


## “Histidine Adem Kit For Recombinant Protein Purification”



- **Recombinant 6xHis-tagged protein purification**

- **High Binding capacity:**

4 mg of recombinant 6xHis-tagged protein/ml of magnetic particles

- **High Yield:**

-From bacterial lysates  
-Reduce non specific binding

- **Convenient:**

-Magnetic technology  
-Ready to use purification buffers

- **Flexible for optimal performance:**

adjust the amount of magnetic particles to the amount of 6xHis-tagged recombinant protein

### SDS Page analysis showing purification of 6xHis-EGFP from bacterial lysate

Lane 1: Molecular weight markers  
Lane 2: Crude *E. Coli* lysate before capture  
Lane 3: Pure 6xHis-EGFP\*  
Lane 4: Elution

The arrow indicates the location of the 6xHis-EGFP\*

\*Enhanced Green Fluorescent Protein



### Ordering information


| Product            | Size                                 | Code number |
|--------------------|--------------------------------------|-------------|
| Histidine Adem kit | 65 Recombinant protein purifications | 04500       |

# “Bio-Adembeads Antibodies : Very easy selection”



## Indirect IgG/IgM Selection : High flexibility

- **Rapid procedure:** capture of any type of cells from different species in less than 5 min
- **High specificity:**
  - all our coated anti-species antibodies are Fc specific
  - no cross reaction
- **Simple:**
  - all frequently used antibodies can be recognized by Bio-Adembeads Antibodies
  - Compatible with flow cytometry and microscopic analysis



## Ordering Information : Bio-Adembeads Antibodies

| Product              | Size       | Code Number          |
|----------------------|------------|----------------------|
| Goat Anti-Mouse IgG  | 2-10-50 ml | G05030-G05031-G05032 |
| Goat Anti-Rat IgG    | 2-10-50 ml | G05130-G05131-G05132 |
| Goat Anti-Rabbit IgG | 2-10-50 ml | G05230-G05231-G05232 |
| Goat Anti-Human IgG  | 2-10-50 ml | G05330-G05331-G05332 |
| Mouse Anti-biotin    | 2-10-50 ml | G05430-G05431-G05432 |
| Goat Anti-Mouse IgM  | 2-10 ml    | M05030-M05031        |
| Goat Anti-Rat IgM    | 2-10 ml    | M05130-M05131        |
| Goat Anti-Rabbit IgM | 2-10 ml    | M05230-M05231        |
| Goat Anti-Human IgM  | 2-10 ml    | M05330-M05331        |

*Diameter 300 nm*

ademtech



ISO 9001 Certified

### ADEMTECH SA

Parc Scientifique Unitec 1

4 Allée du Doyen Georges Brus - 33600 Pessac - France

Tel: +33 (0)5 57 02 02 01 - Fax: +33 (0)5 57 02 02 06

E-mail: ademtech@ademtech.com - www.ademtech.com