

Ademtech Magnetic Devices

Adem-Mag 96
Catalog # 20106



Product Description

Adem-Mag 96
Catalog # 20106

The Adem-Mag 96 magnet offers a high energy neodymium-iron-boron [Nd-Fe-B] configuration.

The magnet is designed to work with a 96-well plate. A 0.62cm/0.25 inch spacer is available to optimize use of the magnetic device. The spacer allows to pellet the nanoparticles at the side of the well for easy aspiration.

Product should be stored in a cool dry environment and should be cleaned with a damp cloth and mild detergent when exposed to harsh solvents. Do not autoclave.

Warning and Handling

This product generates a magnetic field:

- keep away from other magnetic devices.
- people with pacemakers should avoid direct contact.
- credit cards, tapes and disks can be erased: keep away from this product.

This product is assembled with an injection model plastic: do not try to disassemble it or it may not work properly.

Correct handling is required to avoid any injury.

ADEMTECH SA
Parc Scientifique Unitec 1 ;
4 Allée du Doyen Georges Brus ;
33600 Pessac
France

Tel : +33 (0)5 57 02 02 01
Fax : +33 (0)5 57 02 02 06
E-mail : ademtech@ademtech.com