


Antibody Tools


Immunoprecipitation




Western Blot



Sandwich ELISA




Immunohistochemistry




Immunofluorescence



Immunoprecipitation



Flow Cytometry



Antibodies for Immunoprecipitation	1
Secondary Antibodies	36
Antibody Pairs for IP-WB	38
IP&ChIP Kits	51
Immunoprecipitation and Purification Systems	52

Abnova Corporation

www.abnova.com

Email: sales@abnova.com

Address: 9F, No. 108, Jhouzih St., Neihu, Taipei 114, Taiwan

Tel: + 886 2 8751 1888 Fax: + 886 2 6602 1218

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB5411	A1/A2 monoclonal antibody, clone Z2A	PAB5070	Adipoq polyclonal antibody
MAB0794	A1BG monoclonal antibody, clone 54B12	PAB5124	Adipoq polyclonal antibody
MAB0759	A2M monoclonal antibody, clone 3D1	PAB9125	ADIPOQ polyclonal antibody
MAB0758	A2M monoclonal antibody, clone 9A3	PAB9127	Adipoq polyclonal antibody
H00051166-D01	AADAT MaxPab rabbit polyclonal antibody (D01)	PAB9128	Adipoq polyclonal antibody
H00000016-D01	AARS MaxPab rabbit polyclonal antibody (D01)	PAB9129	Adipoq polyclonal antibody
MAB8772	ABCA1 monoclonal antibody, clone AB.H10	PAB9130	Adipoq polyclonal antibody
MAB8291	ABCA1 monoclonal antibody, clone AB1.G6	PAB13780	Adipor2 polyclonal antibody
MAB3345	ABCB1 monoclonal antibody, clone MRK16	PAB4972	ADM polyclonal antibody
H00089845-D01	ABCC10 MaxPab rabbit polyclonal antibody (D01)	PAB5087	Adm polyclonal antibody
H00085320-D01	ABCC11 MaxPab rabbit polyclonal antibody (D01)	H00000135-M17	ADORA2A monoclonal antibody (M17), clone 4E4
PAB12561	ABCE1 polyclonal antibody	H00000154-D01	ADRB2 MaxPab rabbit polyclonal antibody (D01)
H00010061-D01	ABCF2 MaxPab rabbit polyclonal antibody (D01)	H00000166-D01	AES MaxPab rabbit polyclonal antibody (D01)
H00000034-D01	ACADM MaxPab rabbit polyclonal antibody (D01)	MAB3912	AF monoclonal antibody, clone B1822
PAB12498	ACAT1 polyclonal antibody	H00000174-D01	AFP MaxPab rabbit polyclonal antibody (D01)
H00000039-D01	ACAT2 MaxPab rabbit polyclonal antibody (D01)	H00000174-M01	AFP monoclonal antibody (M01), clone 1G7
H00084680-D01	ACCS MaxPab rabbit polyclonal antibody (D01)	MAB2685	AFP monoclonal antibody, clone 3B1
H00065057-D01	ACD MaxPab rabbit polyclonal antibody (D01)	MAB2686	AFP monoclonal antibody, clone 8H2
H00065057-M01	ACD monoclonal antibody (M01), clone 1C11-1A7	MAB0906	AFP monoclonal antibody, clone AFP-01
H00065057-M02	ACD monoclonal antibody (M02), clone 1D8-1B6	MAB6776	AFP monoclonal antibody, clone AFP-Y1
PAB10348	Acetylated lysine polyclonal antibody	H00060509-D01	AGBL5 MaxPab rabbit polyclonal antibody (D01)
PAB16623	Acetylated Lysine polyclonal antibody	H00000177-M02	AGER monoclonal antibody (M02), clone 1C1
PAB5854	Acetylated Lysine polyclonal antibody	H00010554-D01	AGPAT1 MaxPab rabbit polyclonal antibody (D01)
MAB0661	Acetylated-TP53 Lys120 monoclonal antibody, clone 10E5	PAB4973	AGT polyclonal antibody
PAB11344	ACHE polyclonal antibody	PAB5052	AGT polyclonal antibody
MAB9261	Acinar monoclonal antibody, clone 16x2g	PAB5053	AGT polyclonal antibody
H00000048-M01	ACO1 monoclonal antibody (M01), clone 2C1	PAB5054	AGT polyclonal antibody
H00000051-D01	ACOX1 MaxPab rabbit polyclonal antibody (D01)	H00057085-D01	AGTRAP MaxPab rabbit polyclonal antibody (D01)
H00008309-D01	ACOX2 MaxPab rabbit polyclonal antibody (D01)	H00054806-D01	AH1 MaxPab rabbit polyclonal antibody (D01)
H00051703-M01	ACSL5 monoclonal antibody (M01), clone 5H8	MAB8962	Ahi1 monoclonal antibody, clone 645s3
PAB12649	ACTA2 polyclonal antibody	H00079026-M01	AHNAK monoclonal antibody (M01), clone 3G7
PAB12496	ACTB polyclonal antibody	MAB3607	AHNAK monoclonal antibody, clone EM-09
PAB12510	ACTG1 polyclonal antibody	H00000196-M02	AHR monoclonal antibody (M02), clone 3B12
MAB5271	Activated NK Cells monoclonal antibody, clone ANK44	PAB12490	AHR polyclonal antibody
H00000086-D01	ACTL6A MaxPab rabbit polyclonal antibody (D01)	H00010598-D01	AHSA1 MaxPab rabbit polyclonal antibody (D01)
H00010880-D01	ACTL7B MaxPab rabbit polyclonal antibody (D01)	H00010598-M02	AHSA1 monoclonal antibody (M02), clone 4F3-B6
PAB9840	ACVR1 polyclonal antibody	MAB6669	Ahsa1 monoclonal antibody, clone 25F2.D10
H00000093-D01	ACVR2B MaxPab rabbit polyclonal antibody (D01)	MAB6668	Ahsa1 monoclonal antibody, clone 25F2.D9
H00000095-M01	ACY1 monoclonal antibody (M01), clone 4F1-B7	PAB15426	Ahsa1/Ahsa2 polyclonal antibody
H00000095-M02	ACY1 monoclonal antibody (M02), clone 1A11-F6	H00000197-D01	AHSG MaxPab rabbit polyclonal antibody (D01)
H00000100-D01	ADA MaxPab rabbit polyclonal antibody (D01)	H00000197-M02	AHSG monoclonal antibody (M02), clone 1D8
PAB11345	ADA polyclonal antibody	H00064853-D01	AIDA MaxPab rabbit polyclonal antibody (D01)
PAB13735	ADAM17 polyclonal antibody	PAB0324	AIFM1 polyclonal antibody
H00009507-D01	ADAMTS4 MaxPab rabbit polyclonal antibody (D01)	H00084883-D01	AIFM2 MaxPab rabbit polyclonal antibody (D01)
H00009507-M01	ADAMTS4 monoclonal antibody (M01), clone 2A6	H00009447-D01	AIM2 MaxPab rabbit polyclonal antibody (D01)
PAB18057	ADD1 polyclonal antibody	MAB2382	AIP monoclonal antibody, clone 35-2
MAB4645	Adeno-associated virus VP1/VP2/VP3 monoclonal antibody, clone B1	H00023746-M23	AIPL1 monoclonal antibody (M23), clone 3A3
MAB5396	Adenovirus C11 monoclonal antibody, clone 1-N-14	H00000203-M06	AK1 monoclonal antibody (M06), clone 3A6-1F5
MAB5395	Adenovirus E11 monoclonal antibody, clone 1-M-13	H00050808-D01	AK3 MaxPab rabbit polyclonal antibody (D01)
MAB1874	Adenovirus E1A monoclonal antibody, clone M73	H00010270-D01	AKAP8 MaxPab rabbit polyclonal antibody (D01)
MAB3883	Adenovirus hexon monoclonal antibody, clone B111M	H00001645-D01	AKR1C1 MaxPab rabbit polyclonal antibody (D01)
MAB3884	Adenovirus hexon monoclonal antibody, clone B112M	H00001646-D01	AKR1C2 MaxPab rabbit polyclonal antibody (D01)
MAB3885	Adenovirus hexon monoclonal antibody, clone BIOD604	H00001109-D01	AKR1C4 MaxPab rabbit polyclonal antibody (D01)
MAB8757	Adenovirus Hexon monoclonal antibody, clone G097	H00006718-M01	AKR1D1 monoclonal antibody (M01), clone 1A6
MAB3882	Adenovirus type 6 monoclonal antibody, clone 143	H00008574-D01	AKR7A2 MaxPab rabbit polyclonal antibody (D01)
H00000126-D01	ADH1C MaxPab rabbit polyclonal antibody (D01)	H00008574-M15	AKR7A2 monoclonal antibody (M15), clone 2H3
H00000127-D01	ADH4 MaxPab rabbit polyclonal antibody (D01)	MAB0065	AKT1 (phospho S473) monoclonal antibody, clone 104A282
H00000131-D01	ADH7 MaxPab rabbit polyclonal antibody (D01)	PAB12600	AKT1 (phospho S473) polyclonal antibody
MAB3886	ADIPOQ monoclonal antibody, clone Ne.Na	PAB17888	Akt1 (phospho S473) polyclonal antibody
PAB0005	ADIPOQ polyclonal antibody	PAB12601	AKT1 (phospho T308) polyclonal antibody
PAB0006	Adipoq polyclonal antibody	H00000207-M01	AKT1 monoclonal antibody (M01), clone 4C3
PAB13710	Adipoq polyclonal antibody	H00000207-M03	AKT1 monoclonal antibody (M03), clone 2E11
PAB5030	ADIPOQ polyclonal antibody	MAB7913	AKT1 monoclonal antibody, clone 1D1
PAB5031	ADIPOQ polyclonal antibody	PAB12650	AKT1 polyclonal antibody
PAB5069	Adipoq polyclonal antibody	PAB8534	AKT1/AKT2/AKT3 polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB12651	AKT2 polyclonal antibody	MAB9440	Annexins I/II monoclonal antibody, clone 75gd6
PAB12652	AKT3 polyclonal antibody	PAB14756	ANO1 polyclonal antibody
H00000213-M01	ALB monoclonal antibody (M01), clone 1G12-1B3	MAB5147	ANP32A monoclonal antibody, clone RJ1
MAB0742	ALB monoclonal antibody, clone 1G2	MAB9487	ANPEP monoclonal antibody, clone 066j4
PAB15362	ALB polyclonal antibody	MAB6038	ANPEP monoclonal antibody, clone 22A5 (FITC)
PAB5073	Alb polyclonal antibody	MAB6039	ANPEP monoclonal antibody, clone 22A5 (PE)
PAB5079	Alb polyclonal antibody	MAB3827	ANPEP monoclonal antibody, clone WM15
PAB5089	ALB polyclonal antibody	MAB4354	ANPEP monoclonal antibody, clone WM15 (APC)
PAB5090	ALB polyclonal antibody	MAB4355	ANPEP monoclonal antibody, clone WM15 (FITC)
MAB1205	Albumin monoclonal antibody	MAB4356	ANPEP monoclonal antibody, clone WM15 (PE)
PAB11386	<i>Alcaligenes</i> Choline Oxidase polyclonal antibody	MAB4357	ANPEP monoclonal antibody, clone WM15 (PE-Dyomics 590)
MAB6890	ALCAM monoclonal antibody, clone 3A6	H00000301-D01	ANXA1 MaxPab rabbit polyclonal antibody (D01)
PAB10177	Alcohol Dehydrogenase polyclonal antibody	H00000306-M07	ANXA3 monoclonal antibody (M07), clone 1E2
PAB10178	Aldehyde Dehydrogenase polyclonal antibody	H00000307-D01	ANXA4 MaxPab rabbit polyclonal antibody (D01)
H00000216-D01	ALDH1A1 MaxPab rabbit polyclonal antibody (D01)	H00000308-D01	ANXA5 MaxPab rabbit polyclonal antibody (D01)
H00000216-M05	ALDH1A1 monoclonal antibody (M05), clone 1G6	H00000308-M01	ANXA5 monoclonal antibody (M01), clone 1F4-1A5
H00000218-D01	ALDH3A1 MaxPab rabbit polyclonal antibody (D01)	PAB8537	ANXA5 polyclonal antibody
H00000221-M02	ALDH3B1 monoclonal antibody (M02), clone 1F2	H00023028-D01	AOF2 MaxPab rabbit polyclonal antibody (D01)
H00008659-D01	ALDH4A1 MaxPab rabbit polyclonal antibody (D01)	PAB14107	AOF2 polyclonal antibody
H00007915-D01	ALDH5A1 MaxPab rabbit polyclonal antibody (D01)	PAB14134	AOF2 polyclonal antibody
PAB11351	ALDOA polyclonal antibody	PAB8049	Ap1b1 polyclonal antibody
MAB1213	Aldosterone monoclonal antibody	PAB10221	AP1G1 polyclonal antibody
PAB5141	Aldosterone polyclonal antibody	PAB9935	AP2A1 polyclonal antibody
PAB5142	Aldosterone polyclonal antibody	MAB5594	AP-2-alpha monoclonal antibody, clone A6/2/2
PAB5143	Aldosterone polyclonal antibody	PAB9965	AP3D1 polyclonal antibody
H00011036-M01	ALF monoclonal antibody (M01), clone 4F5	MAB0818	APAF1 monoclonal antibody, clone 14E3
H00011036-M02	ALF monoclonal antibody (M02), clone 2E3	MAB3406	APAF1 monoclonal antibody, clone 2E12
H00011036-M03	ALF monoclonal antibody (M03), clone 6F5	PAB0325	APAF1 polyclonal antibody
H00221120-D01	ALKB3 MaxPab rabbit polyclonal antibody (D01)	PAB12116	Apbb1 polyclonal antibody
PAB5173	Allatostatin 1 polyclonal antibody	MAB0778	APCS monoclonal antibody, clone 6E6
PAB5174	Allatostatin 2 polyclonal antibody	H00000328-D02	APEX1 MaxPab rabbit polyclonal antibody (D02)
PAB5175	Allatostatin 3 polyclonal antibody	H00000328-M05	APEX1 monoclonal antibody (M05), clone 3E12
PAB5176	Allatostatin 4 polyclonal antibody	MAB3422	APEX1 monoclonal antibody, clone 13B 8E5C2
H00055821-M07	ALLC monoclonal antibody (M07), clone 3D3	MAB2372	APEX1 monoclonal antibody, clone 13B8E5C2
PAB5144	Allopregnanolone polyclonal antibody	PAB12487	APEX1 polyclonal antibody
PAB13608	ALOX15B polyclonal antibody	H00000335-M01	APOA1 monoclonal antibody (M01), clone 3A11-1A9
PAB13607	Alox8 polyclonal antibody	PAB10089	Apoa1 polyclonal antibody
PAB5137	Alpha-1 acid glycoprotein polyclonal antibody	PAB4999	APOA1 polyclonal antibody
PAB9908	Alpha-Amylase polyclonal antibody	PAB9966	APOA1 polyclonal antibody
MAB9425	alpha-Cat monoclonal antibody, clone 656ds	PAB9967	APOA2 polyclonal antibody
H00000249-D01	ALPL MaxPab rabbit polyclonal antibody (D01)	PAB14895	APOB polyclonal antibody
H00000250-D01	ALPP MaxPab rabbit polyclonal antibody (D01)	PAB5049	APOB polyclonal antibody
H00000251-D01	ALPPL2 MaxPab rabbit polyclonal antibody (D01)	PAB8549	APOB polyclonal antibody
H00000251-M07	ALPPL2 monoclonal antibody (M07), clone 2B3	PAB9968	APOB polyclonal antibody
MAB3788	AMBP/Ambp monoclonal antibody, clone BN11.10	H00027350-D01	APOBEC3C MaxPab rabbit polyclonal antibody (D01)
H00000265-D01	AMELX MaxPab rabbit polyclonal antibody (D01)	H00027350-M01	APOBEC3C monoclonal antibody (M01), clone 3E6
H00084883-M09	AMID monoclonal antibody (M09), clone 3A11	PAB9969	APOC1 polyclonal antibody
H00084883-M13	AMID monoclonal antibody (M13), clone 2C6	PAB9970	APOC2 polyclonal antibody
H00000271-D01	AMPD2 MaxPab rabbit polyclonal antibody (D01)	PAB5048	APOC3 polyclonal antibody
MAB0708	AMPH monoclonal antibody, clone C14-23	PAB9971	APOC3 polyclonal antibody
MAB0699	AMPH monoclonal antibody, clone N8-2	H00000346-M01	APOC4 monoclonal antibody (M01), clone 3D10
H00000275-D01	AMT MaxPab rabbit polyclonal antibody (D01)	MAB4653	APOE monoclonal antibody, clone E6D7
H00000279-D01	AMY2A MaxPab rabbit polyclonal antibody (D01)	PAB9972	APOE polyclonal antibody
MAB5669	Amyloid-beta monoclonal antibody, clone McSA1	H00008542-D01	APOL1 MaxPab rabbit polyclonal antibody (D01)
PAB10017	ANAPC1 (phospho S355) polyclonal antibody	H00023780-D01	APOL2 MaxPab rabbit polyclonal antibody (D01)
PAB10018	ANAPC1 polyclonal antibody	PAB12573	APP polyclonal antibody
H00010393-D01	ANAPC10 MaxPab rabbit polyclonal antibody (D01)	H00000353-D01	APRT MaxPab rabbit polyclonal antibody (D01)
PAB10019	ANAPC10 polyclonal antibody	MAB5324	Aqp1 monoclonal antibody, clone 1/A5F6
PAB10092	Angptl1 polyclonal antibody	PAB12654	AR (phospho S210) polyclonal antibody
PAB10091	Angptl2 polyclonal antibody	PAB12653	AR polyclonal antibody
H00027329-M10	ANGPTL3 monoclonal antibody (M10), clone 3F11	PAB12655	AR polyclonal antibody
H00051129-M02	ANGPTL4 monoclonal antibody (M02), clone 1A12	MAB3599	ARAP1 monoclonal antibody, clone ARAP1-2
H00253935-M01	ANGPTL5 monoclonal antibody (M01), clone 1B2	H00000383-D01	ARG1 MaxPab rabbit polyclonal antibody (D01)
H00000286-M01	ANK1 monoclonal antibody (M01), clone 3C3	H00000384-D01	ARG2 MaxPab rabbit polyclonal antibody (D01)
H00057763-D01	ANKRA2 MaxPab rabbit polyclonal antibody (D01)	H00000384-M05	ARG2 monoclonal antibody (M05), clone 3D10

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00000396-D01	ARHGDI A MaxPab rabbit polyclonal antibody (D01)	H00009529-D01	BAG5 MaxPab rabbit polyclonal antibody (D01)
H00000396-D03	ARHGDI A MaxPab rabbit polyclonal antibody (D03)	PAB0332	BAG5 polyclonal antibody
MAB0059	ARHGDI B monoclonal antibody, clone 97A1015	H00010458-M01	BAIAP2 monoclonal antibody (M01), clone 3F3-2D3
H00007984-M02	ARHGEF5 monoclonal antibody (M02), clone 4D10-D11	H00000578-M01A	BAK1 monoclonal antibody (M01A), clone 1E4
H00115761-D01	ARL11 MaxPab rabbit polyclonal antibody (D01)	PAB0333	BAK1 polyclonal antibody
H00000379-D01	ARL4D MaxPab rabbit polyclonal antibody (D01)	PAB9805	BAPTA polyclonal antibody
PAB12484	ARNTL polyclonal antibody	H00010538-M01	BATF monoclonal antibody (M01), clone 8A12
H00000410-D01	ARSA MaxPab rabbit polyclonal antibody (D01)	MAB3284	Bax monoclonal antibody, clone 5B7
H00000419-D01	ART3 MaxPab rabbit polyclonal antibody (D01)	MAB1809	BAX monoclonal antibody, clone 6A7
H00000419-M09	ART3 monoclonal antibody (M09), clone 1D2	PAB12656	BAX polyclonal antibody
H00000420-D01	ART4 MaxPab rabbit polyclonal antibody (D01)	PAB14668	Bax polyclonal antibody
H00140462-M01	ASB9 monoclonal antibody (M01), clone 1D8	PAB8564	BAX polyclonal antibody
PAB9918	Ascorbate Oxidase polyclonal antibody	PAB9790	BAX polyclonal antibody
PAB8554	ASH2L polyclonal antibody	H00000585-D01	BBS4 MaxPab rabbit polyclonal antibody (D01)
H00008623-D01	ASMTL MaxPab rabbit polyclonal antibody (D01)	H00055973-D01	BCAP29 MaxPab rabbit polyclonal antibody (D01)
H00000440-D01	ASNS MaxPab rabbit polyclonal antibody (D01)	MAB5149	BCAP31 monoclonal antibody, clone A1/182
H00000440-M02	ASNS monoclonal antibody (M02), clone 2B3	PAB7934	Bcar1 (phospho Y762) polyclonal antibody
H00000443-D01	ASPA MaxPab rabbit polyclonal antibody (D01)	MAB1386	Bcar1 monoclonal antibody, clone M144
H00000444-D03	ASPH MaxPab rabbit polyclonal antibody (D03)	H00000587-D01	BCAT2 MaxPab rabbit polyclonal antibody (D01)
H00080150-D01	ASRGL1 MaxPab rabbit polyclonal antibody (D01)	H00053335-D01	BCL11A MaxPab rabbit polyclonal antibody (D01)
H00080150-M03	ASRGL1 monoclonal antibody (M03), clone 5C8	PAB12609	BCL2 (phospho S70) polyclonal antibody
PAB7517	ASS1 polyclonal antibody	PAB12608	BCL2 (phospho T129) polyclonal antibody
H00055252-D01	ASXL2 MaxPab rabbit polyclonal antibody (D01)	PAB12610	BCL2 (phospho T56) polyclonal antibody
H00055210-D01	ATAD3A MaxPab rabbit polyclonal antibody (D01)	MAB0806	BCL2 monoclonal antibody, clone 11C5
PAB17895	ATF2 (phospho T71) polyclonal antibody	MAB0805	BCL2 monoclonal antibody, clone 2B11
H00022926-D01	ATF6 MaxPab rabbit polyclonal antibody (D01)	MAB1745	BCL2 monoclonal antibody, clone Bcl-2/100
MAB6762	ATF6 monoclonal antibody, clone 1-7	PAB0236	BCL2 polyclonal antibody
H00055054-D01	ATG16L1 MaxPab rabbit polyclonal antibody (D01)	PAB12662	BCL2 polyclonal antibody
H00000471-D01	ATIC MaxPab rabbit polyclonal antibody (D01)	PAB2440	BCL2 polyclonal antibody
MAB6132	ATM monoclonal antibody	PAB8567	BCL2 polyclonal antibody
MAB0829	ATM monoclonal antibody, clone 1B10	PAB9708	BCL2 polyclonal antibody
MAB7981	ATM monoclonal antibody, clone 2C1	PAB8528	BCL2A1 polyclonal antibody
MAB7982	ATM monoclonal antibody, clone 3E8	MAB3283	BCL2L1 monoclonal antibody, clone 7B2.5
MAB0998	ATM monoclonal antibody, clone 9A2	MAB5860	BCL2L1 monoclonal antibody, clone 7B2.5 (Biotin)
H00000483-D01	ATP1B3 MaxPab rabbit polyclonal antibody (D01)	MAB5859	BCL2L1 monoclonal antibody, clone 7B2.5 (FITC)
MAB1723	ATP1B3 monoclonal antibody, clone 5G10	MAB5861	BCL2L1 monoclonal antibody, clone 7B2.5 (PE)
H00027032-M01	ATP2C1 monoclonal antibody (M01), clone 2G1	PAB0237	BCL2L1 polyclonal antibody
H00245973-D01	ATP6V1C2 MaxPab rabbit polyclonal antibody (D01)	PAB8573	BCL2L1 polyclonal antibody
H00090423-D01	ATP6V1E2 MaxPab rabbit polyclonal antibody (D01)	MAB3704	Bcl2l11 monoclonal antibody, clone 14A8
H00051606-D01	ATP6V1H MaxPab rabbit polyclonal antibody (D01)	MAB3705	Bcl2l11 monoclonal antibody, clone 5E5
MAB7983	ATR monoclonal antibody, clone 2B5	MAB7553	Bcl2l11 monoclonal antibody, clone Ham151-149
PAB5177	Atriopeptin 2 polyclonal antibody	MAB3449	BCL2L13 monoclonal antibody, clone Rocky-1
PAB5178	Atriopeptin 3 polyclonal antibody	MAB3450	BCL2L13 monoclonal antibody, clone Rocky-2
PAB12476	ATRIP polyclonal antibody	H00079370-M02	BCL2L14 monoclonal antibody (M02), clone 1D11
MAB0766	ATRNL1 monoclonal antibody, clone 9H8	MAB3487	BCL2L2 monoclonal antibody, clone 13F9
PAB9919	Aub polyclonal antibody	H00000602-D01	BCL3 MaxPab rabbit polyclonal antibody (D01)
PAB5179	Auriculin A polyclonal antibody	PAB8568	BCL3 polyclonal antibody
PAB5180	Auriculin B polyclonal antibody	PAB14151	BCL6 polyclonal antibody
MAB9253	<i>Avian myoblastosis virus</i> monoclonal antibody, clone 6w46	PAB8569	BCL6 polyclonal antibody
H00000558-M01	AXL monoclonal antibody (M01), clone 6C8	PAB8600	BCLAF1 polyclonal antibody
MAB8056	B2M monoclonal antibody, clone 6D585	H00000617-M01	BCS1L monoclonal antibody (M01), clone 5F3
MAB0899	B2M monoclonal antibody, clone B2M-01	H00000622-D01	BDH1 MaxPab rabbit polyclonal antibody (D01)
MAB7270	B2M monoclonal antibody, clone BBM.1	MAB9435	beta-tubulin monoclonal antibody, clone 9w431
H00008703-D01	B4GALT3 MaxPab rabbit polyclonal antibody (D01)	PAB0235	BFAR polyclonal antibody
H00008702-D01	B4GALT4 MaxPab rabbit polyclonal antibody (D01)	PAB18014	BFP polyclonal antibody
PAB11390	<i>Bacillus</i> Uricase polyclonal antibody	H00000633-M01	BGN monoclonal antibody (M01), clone 4E1-1G7
MAB7331	Bacteria Cdc13 monoclonal antibody, clone 6F11/2	PAB12512	BHLHE40 polyclonal antibody
PAB17889	Bad (phospho S112) polyclonal antibody	H00023743-D01	BHMT2 MaxPab rabbit polyclonal antibody (D01)
PAB8560	Bad polyclonal antibody	MAB0809	BID monoclonal antibody, clone 8D2
PAB0328	Bag1 polyclonal antibody	PAB0238	Bid polyclonal antibody
H00009532-D01	BAG2 MaxPab rabbit polyclonal antibody (D01)	PAB9131	BID polyclonal antibody
PAB0329	BAG2 polyclonal antibody	PAB9132	Bid polyclonal antibody
PAB0330	BAG3 polyclonal antibody	PAB0239	BIK polyclonal antibody
H00009530-D01	BAG4 MaxPab rabbit polyclonal antibody (D01)	PAB11385	Bilirubin Oxidase polyclonal antibody
PAB0331	BAG4 polyclonal antibody	PAB0252	BIRC2 polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00000330-D01	BIRC3 MaxPab rabbit polyclonal antibody (D01)	H00118924-D01	C10orf4 MaxPab rabbit polyclonal antibody (D01)
PAB0253	BIRC3 polyclonal antibody	H00143384-D01	C10orf46 MaxPab rabbit polyclonal antibody (D01)
PAB0282	BIRC4 polyclonal antibody	H00056672-D01	C11orf17 MaxPab rabbit polyclonal antibody (D01)
H00000332-M01	BIRC5 monoclonal antibody (M01), clone 5B10	H00029902-M01	C12orf24 monoclonal antibody (M01), clone 2G4-1A1
PAB0272	BIRC5 polyclonal antibody	H00080183-D01	C13orf18 MaxPab rabbit polyclonal antibody (D01)
PAB14670	BIRC7 polyclonal antibody	H00084299-D01	C17orf37 MaxPab rabbit polyclonal antibody (D01)
MAB7893	BKbeta3a monoclonal antibody, clone S40B-18	H00000712-D01	C1QA MaxPab rabbit polyclonal antibody (D01)
MAB7313	Blood Group A Antigen monoclonal antibody, clone 9A	H00000714-D01	C1QC MaxPab rabbit polyclonal antibody (D01)
H00000644-M03	BLVRA monoclonal antibody (M03), clone S2	PAB13786	C1qtnf7 polyclonal antibody
MAB3703	Bmf monoclonal antibody, clone 12E10	H00000716-D01	C1S MaxPab rabbit polyclonal antibody (D01)
H00000655-D01	BMP7 MaxPab rabbit polyclonal antibody (D01)	H00128434-M02	C20orf102 monoclonal antibody (M02), clone 1A8
H00000658-D01	BMPR1B MaxPab rabbit polyclonal antibody (D01)	H00055741-M01	C20orf31 monoclonal antibody (M01), clone 2E4
PAB8489	BNAS2 polyclonal antibody	H00079019-M01	C22orf18 monoclonal antibody (M01), clone 4C12-2C8
PAB0264	BNIP3 polyclonal antibody	PAB15380	C3 polyclonal antibody
H00000665-M02	BNIP3L monoclonal antibody (M02), clone 2E11	PAB5002	C3 polyclonal antibody
H00066037-D01	BOLL MaxPab rabbit polyclonal antibody (D01)	MAB1031	C4 monoclonal antibody, clone 49H9
H00066037-M39	BOLL monoclonal antibody (M39), clone 2G4	MAB1032	C4-beta monoclonal antibody, clone 52H10
MAB7863	Bovine papillomavirus type 1 E2 monoclonal antibody, clone 1E2	H00000725-D01	C4BPB MaxPab rabbit polyclonal antibody (D01)
MAB7862	Bovine papillomavirus type 1 E2 monoclonal antibody, clone 3F12	H00000725-M01	C4BPB monoclonal antibody (M01), clone 1G9
MAB7861	Bovine papillomavirus type 1 E2 monoclonal antibody, clone 5E11	H00055319-D01	C4orf43 MaxPab rabbit polyclonal antibody (D01)
MAB2170	BRCA1 monoclonal antibody	MAB0824	C6 monoclonal antibody, clone 13H5
MAB6119	BRCA1 monoclonal antibody	MAB0790	C8B monoclonal antibody, clone 14A10
MAB2363	BRCA1 monoclonal antibody, clone KEN	H00054981-D01	C9orf95 MaxPab rabbit polyclonal antibody (D01)
MAB2364	BRCA1 monoclonal antibody, clone MU	H00054981-M01	C9orf95 monoclonal antibody (M01), clone 2A10
MAB2362	BRCA1 monoclonal antibody, clone RAY	H00000759-M05	CA1 monoclonal antibody (M05), clone M1
PAB8595	BRCA1 polyclonal antibody	H00000771-D01	CA12 MaxPab rabbit polyclonal antibody (D01)
H00006046-D01	BRD2 MaxPab rabbit polyclonal antibody (D01)	MAB6647	CACH2C monoclonal antibody, clone S57-46
PAB8597	BRD8 polyclonal antibody	MAB6649	Cacna1d monoclonal antibody, clone S38-8
PAB13898	BrdU polyclonal antibody	MAB6648	Cacna1d monoclonal antibody, clone S48A-9
PAB9791	BrdU polyclonal antibody	MAB6650	CACNA1H monoclonal antibody, clone S55-10
PAB12468	BRIP1 polyclonal antibody	H00000782-D01	CACNB1 MaxPab rabbit polyclonal antibody (D01)
H00025855-M01	BRMS1 monoclonal antibody (M01), clone 2D4-2G11	MAB6664	Cacnb1 monoclonal antibody, clone S7-18
H00056853-D01	BRUNOL4 MaxPab rabbit polyclonal antibody (D01)	MAB7895	Cacnb2 monoclonal antibody, clone S8b-1
H00000682-D01	BSG MaxPab rabbit polyclonal antibody (D01)	H00000784-D01	CACNB3 MaxPab rabbit polyclonal antibody (D01)
MAB1722	BSG monoclonal antibody, clone 2G2	PAB9803	Cacnb3 polyclonal antibody
MAB3830	BSG monoclonal antibody, clone MEM-M6/1	MAB6665	Cacnb4 monoclonal antibody, clone S10-7
MAB6893	BSG monoclonal antibody, clone UM-8D6	PAB14881	Cacnb4 polyclonal antibody
MAB7040	Bsn monoclonal antibody, clone SAP7F407	H00059285-D01	CACNG6 MaxPab rabbit polyclonal antibody (D01)
H00007809-D01	BSND MaxPab rabbit polyclonal antibody (D01)	H00059284-D01	CACNG7 MaxPab rabbit polyclonal antibody (D01)
H00000683-M08	BST1 monoclonal antibody (M08), clone 4C2	H00027101-D01	CACYBP MaxPab rabbit polyclonal antibody (D01)
H00000684-D01	BST2 MaxPab rabbit polyclonal antibody (D01)	PAB8603	Cadherin polyclonal antibody
H00000684-D03	BST2 MaxPab rabbit polyclonal antibody (D03)	PAB5084	Calca polyclonal antibody
H00000684-M15	BST2 monoclonal antibody (M15), clone 2E6	PAB4979	CALCA/CALCB polyclonal antibody
H00114781-D01	BTBD9 MaxPab rabbit polyclonal antibody (D01)	H00010241-D01	CALCOCO2 MaxPab rabbit polyclonal antibody (D01)
H00000685-D01	BTC MaxPab rabbit polyclonal antibody (D01)	H00000810-D01	CALML3 MaxPab rabbit polyclonal antibody (D01)
H00000686-M01	BTD monoclonal antibody (M01), clone 3B10-2B3	PAB12147	Calr polyclonal antibody
MAB3279	Bu-1 monoclonal antibody, clone 5-11G2	PAB15418	CALR polyclonal antibody
MAB6785	Bu-1 monoclonal antibody, clone AV20	H00125972-D01	CALR3 MaxPab rabbit polyclonal antibody (D01)
MAB5854	Bu-1 monoclonal antibody, clone AV20 (Biotin)	H00057118-M01A	CAMK1D monoclonal antibody (M01A), clone 4G12-1G4
MAB5853	Bu-1 monoclonal antibody, clone AV20 (FITC)	MAB6626	Camk2a monoclonal antibody, clone 6G9
MAB5855	Bu-1 monoclonal antibody, clone AV20 (PE)	H00000816-D01	CAMK2B MaxPab rabbit polyclonal antibody (D01)
MAB3275	Bu-1a monoclonal antibody, clone 21-1A4	H00000818-D01	CAMK2G MaxPab rabbit polyclonal antibody (D01)
MAB5843	Bu-1a monoclonal antibody, clone 21-1A4 (Biotin)	MAB6627	CaMKII monoclonal antibody, clone 22B1
MAB5842	Bu-1a monoclonal antibody, clone 21-1A4 (FITC)	H00010645-D01	CAMKK2 MaxPab rabbit polyclonal antibody (D01)
MAB5844	Bu-1a monoclonal antibody, clone 21-1A4 (PE)	H00010645-M01	CAMKK2 monoclonal antibody (M01), clone 1A11
MAB3276	Bu-1b monoclonal antibody, clone 5-11G2	PAB13902	<i>Candida albicans</i> polyclonal antibody
MAB5845	Bu-1b monoclonal antibody, clone 5-11G2 (Biotin)	H00124583-D01	CANT1 MaxPab rabbit polyclonal antibody (D01)
MAB5846	Bu-1b monoclonal antibody, clone 5-11G2 (PE)	H00124583-M02	CANT1 monoclonal antibody (M02), clone 1A1
H00000699-M02	BUB1 monoclonal antibody (M02), clone 2F9	PAB15419	CANX polyclonal antibody
MAB8963	BUB1 monoclonal antibody, clone P8.14H5.1H1	PAB15425	CANX polyclonal antibody
MAB1929	BUB1B monoclonal antibody, clone 8G1	H00010487-D01	CAP1 MaxPab rabbit polyclonal antibody (D01)
MAB8964	BUB1B monoclonal antibody, clone P10.8G1.G11.H12	MAB1837	CAPN1 monoclonal antibody, clone 15C10
PAB9782	Bub1b polyclonal antibody	H00000828-D01	CAPS MaxPab rabbit polyclonal antibody (D01)
PAB9783	Bub1b polyclonal antibody	H00093661-D01	CAPZA3 MaxPab rabbit polyclonal antibody (D01)
PAB9784	Bub1b polyclonal antibody	MAB2680	Carcinoembryonic antigen monoclonal antibody, clone 36H1

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB2681	Carcinoembryonic antigen monoclonal antibody, clone 37B6	PAB9944	CCND1 polyclonal antibody
PAB0240	CARD6 polyclonal antibody	H00000894-D01	CCND2 MaxPab rabbit polyclonal antibody (D01)
PAB0241	CARD6 polyclonal antibody	MAB1852	CCND2 monoclonal antibody, clone DCS3.1
H00022900-D01	CARD8 MaxPab rabbit polyclonal antibody (D01)	MAB1853	CCND3 monoclonal antibody, clone DCS22
PAB0281	CARD8 polyclonal antibody	MAB1856	CCNE1 monoclonal antibody, clone 168
H00055602-M01	CARF monoclonal antibody (M01), clone 4D11	PAB10197	Ccne2 polyclonal antibody
MAB8582	Cask monoclonal antibody, clone S56A-50	PAB9945	CCNL1 polyclonal antibody
PAB5091	Casomorphin-beta polyclonal antibody	PAB9999	CCNT1 polyclonal antibody
PAB8614	CASP1 polyclonal antibody	PAB11292	CCNT2 polyclonal antibody
PAB0249	CASP10 polyclonal antibody	PAB8939	CCR6 polyclonal antibody
PAB0250	CASP14 polyclonal antibody	H00023552-D01	CCRK MaxPab rabbit polyclonal antibody (D01)
MAB3485	CASP2 monoclonal antibody, clone 11B4	H00009973-D01	CCS MaxPab rabbit polyclonal antibody (D01)
PAB1202	CASP2 polyclonal antibody	MAB0679	CCS monoclonal antibody, clone 2A1
H00000836-D01	CASP3 MaxPab rabbit polyclonal antibody (D01)	MAB0695	CCS monoclonal antibody, clone 3A1
PAB0242	CASP3 polyclonal antibody	PAB1197	CCS polyclonal antibody
PAB0243	CASP3 polyclonal antibody	H00010693-D01	CCT6B MaxPab rabbit polyclonal antibody (D01)
PAB0244	CASP6 polyclonal antibody	H00010574-D01	CCT7 MaxPab rabbit polyclonal antibody (D01)
H00000840-D01	CASP7 MaxPab rabbit polyclonal antibody (D01)	H00150160-D01	CCT8L2 MaxPab rabbit polyclonal antibody (D01)
MAB8485	CASP7 monoclonal antibody, clone B4-G2	MAB6216	CD1 monoclonal antibody, clone 76-7-4
PAB0245	CASP7 polyclonal antibody	MAB4054	CD11b monoclonal antibody, clone 198
PAB8625	CASP7 polyclonal antibody	H00000929-D01	CD14 MaxPab rabbit polyclonal antibody (D01)
MAB8501	CASP8 monoclonal antibody, clone C502S	MAB5621	CD14 monoclonal antibody, clone 61D3
PAB0246	CASP8 polyclonal antibody	MAB6043	CD14 monoclonal antibody, clone 61D3 (APC)
PAB0247	CASP8 polyclonal antibody	MAB6047	CD14 monoclonal antibody, clone 61D3 (APC/Cy7)
PAB0248	CASP9 polyclonal antibody	MAB6041	CD14 monoclonal antibody, clone 61D3 (Biotin)
PAB14669	CASP9 polyclonal antibody	MAB6040	CD14 monoclonal antibody, clone 61D3 (FITC)
PAB9489	CASP9 polyclonal antibody	MAB6042	CD14 monoclonal antibody, clone 61D3 (PE)
H00000847-D01	CAT MaxPab rabbit polyclonal antibody (D01)	MAB6045	CD14 monoclonal antibody, clone 61D3 (PE/Cy5.5)
MAB0671	CAT monoclonal antibody, clone 1A1	MAB6046	CD14 monoclonal antibody, clone 61D3 (PE/Cy7)
MAB0672	CAT monoclonal antibody, clone 2A1	MAB6044	CD14 monoclonal antibody, clone 61D3 (SPRD)
PAB1211	CAT polyclonal antibody	MAB5656	CD14 monoclonal antibody, clone biG 2 (RoMo-1)
MAB2408	CAV1 monoclonal antibody, clone 7C8	MAB3831	CD14 monoclonal antibody, clone MEM-15
PAB12113	Cav1 polyclonal antibody	MAB2531	CD14 monoclonal antibody, clone MEM-18
MAB1674	Cbfa2t2 monoclonal antibody, clone KT42	PAB14381	CD14 polyclonal antibody
PAB9985	CBLC polyclonal antibody	PAB14382	Cd14 polyclonal antibody
H00000873-D01	CBR1 MaxPab rabbit polyclonal antibody (D01)	H00011126-D01	CD160 MaxPab rabbit polyclonal antibody (D01)
H00000874-D01	CBR3 MaxPab rabbit polyclonal antibody (D01)	MAB6978	CD17 monoclonal antibody, clone Huly-m13/H018.3G-6.F5
H00000875-M01	CBS monoclonal antibody (M01), clone 3E1	H00057126-M01	CD177 monoclonal antibody (M01), clone 4C4
H00000875-M02	CBS monoclonal antibody (M02), clone 3D10	MAB0926	CD177 monoclonal antibody, clone MEM-166
PAB15581	CBX1 polyclonal antibody	MAB4365	CD177 monoclonal antibody, clone MEM-166 (FITC)
PAB15582	CBX3 polyclonal antibody	MAB0986	CD18 monoclonal antibody, clone IVA35
PAB10228	CBX4 polyclonal antibody	MAB5256	CD180 monoclonal antibody, clone MHR73
H00054862-D01	CC2D1A MaxPab rabbit polyclonal antibody (D01)	H00000930-M01	CD19 monoclonal antibody (M01), clone 1G3
MAB9069	C-cadherin monoclonal antibody, clone 565s	MAB6499	Cd19 monoclonal antibody, clone 1D3
H00147372-D01	CCBE1 MaxPab rabbit polyclonal antibody (D01)	MAB6500	Cd19 monoclonal antibody, clone 1D3 (PE)
H00006359-D01	CCL15 MaxPab rabbit polyclonal antibody (D01)	MAB3226	Cd19 monoclonal antibody, clone 6D5
H00006362-D01	CCL18 MaxPab rabbit polyclonal antibody (D01)	MAB5690	Cd19 monoclonal antibody, clone 6D5 (APC)
H00006362-M03	CCL18 monoclonal antibody (M03), clone 2C6	MAB5688	Cd19 monoclonal antibody, clone 6D5 (Biotin)
H00006366-D01	CCL21 MaxPab rabbit polyclonal antibody (D01)	MAB5692	Cd19 monoclonal antibody, clone 6D5 (Cy5)
H00006368-D01	CCL23 MaxPab rabbit polyclonal antibody (D01)	MAB5687	Cd19 monoclonal antibody, clone 6D5 (FITC)
PAB8837	Ccl3 polyclonal antibody	MAB5689	Cd19 monoclonal antibody, clone 6D5 (PE)
PAB8838	Ccl4 polyclonal antibody	MAB5693	Cd19 monoclonal antibody, clone 6D5 (PE/Cy5.5)
H00000890-D01	CCNA2 MaxPab rabbit polyclonal antibody (D01)	MAB5694	Cd19 monoclonal antibody, clone 6D5 (PE/Cy7)
MAB7282	CCNA2 monoclonal antibody, clone AT10.2	MAB5691	Cd19 monoclonal antibody, clone 6D5 (SPRD)
MAB7216	CCNA2 monoclonal antibody, clone E67.1	MAB0921	CD19 monoclonal antibody, clone LT19
PAB9996	CCNA2 polyclonal antibody	MAB4369	CD19 monoclonal antibody, clone LT19 (APC)
PAB10042	CCNB1 (phospho S126) polyclonal antibody	MAB4370	CD19 monoclonal antibody, clone LT19 (Biotin)
MAB6997	CCNB1 monoclonal antibody, clone V92.1	MAB4371	CD19 monoclonal antibody, clone LT19 (FITC)
PAB9943	CCNB1 polyclonal antibody	MAB4372	CD19 monoclonal antibody, clone LT19 (PE)
MAB6992	Ccnb2 monoclonal antibody, clone X121.10	MAB4373	CD19 monoclonal antibody, clone LT19 (PE-Dyomics 590)
H00000595-D01	CCND1 MaxPab rabbit polyclonal antibody (D01)	MAB4374	CD19 monoclonal antibody, clone LT19 (PE-Dyomics 647)
MAB7314	CCND1 monoclonal antibody, clone 94B	MAB4375	CD19 monoclonal antibody, clone LT19 (PerCP)
MAB2721	CCND1 monoclonal antibody, clone AF2G5	MAB3227	Cd19 monoclonal antibody, clone MB19-1
MAB6501	CCND1 monoclonal antibody, clone CD1.1	MAB5698	Cd19 monoclonal antibody, clone MB19-1 (APC)
MAB1851	CCND1 monoclonal antibody, clone DCS6	MAB5701	Cd19 monoclonal antibody, clone MB19-1 (APC/Cy5.5)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB5702	Cd19 monoclonal antibody, clone MB19-1 (APC/CY7)	MAB5911	Cd3e monoclonal antibody, clone 145-2C11 (SPRD)
MAB5696	Cd19 monoclonal antibody, clone MB19-1 (Biotin)	MAB3221	Cd3e monoclonal antibody, clone 500-A2
MAB5695	Cd19 monoclonal antibody, clone MB19-1 (FITC)	MAB5915	Cd3e monoclonal antibody, clone 500-A2 (APC)
MAB5697	Cd19 monoclonal antibody, clone MB19-1 (PE)	MAB5913	Cd3e monoclonal antibody, clone 500-A2 (Biotin)
MAB5700	Cd19 monoclonal antibody, clone MB19-1 (PE/Cy5.5)	MAB5917	Cd3e monoclonal antibody, clone 500-A2 (Cy5)
MAB5699	Cd19 monoclonal antibody, clone MB19-1 (SPRD)	MAB5912	Cd3e monoclonal antibody, clone 500-A2 (FITC)
MAB5981	CD19 monoclonal antibody, clone SJ25-C1 (APC)	MAB5914	Cd3e monoclonal antibody, clone 500-A2 (PE)
MAB5986	CD19 monoclonal antibody, clone SJ25-C1 (APC/Cy5.5)	MAB5916	Cd3e monoclonal antibody, clone 500-A2 (SPRD)
MAB5987	CD19 monoclonal antibody, clone SJ25-C1 (APC/CY7)	MAB5604	Cd3e monoclonal antibody, clone C363.29B
MAB5978	CD19 monoclonal antibody, clone SJ25-C1 (Biotin)	MAB6103	Cd3e monoclonal antibody, clone C363.29B (APC)
MAB5983	CD19 monoclonal antibody, clone SJ25-C1 (Cy5)	MAB6101	Cd3e monoclonal antibody, clone C363.29B (Biotin)
MAB5977	CD19 monoclonal antibody, clone SJ25-C1 (FITC)	MAB6100	Cd3e monoclonal antibody, clone C363.29B (FITC)
MAB5979	CD19 monoclonal antibody, clone SJ25-C1 (PE)	MAB6102	Cd3e monoclonal antibody, clone C363.29B (PE)
MAB5984	CD19 monoclonal antibody, clone SJ25-C1 (PE/Cy5.5)	MAB6104	Cd3e monoclonal antibody, clone C363.29B (SPRD)
MAB5985	CD19 monoclonal antibody, clone SJ25-C1 (PE/Cy7)	MAB0870	CD3E monoclonal antibody, clone MEM-57
MAB5980	CD19 monoclonal antibody, clone SJ25-C1 (PE/TXRD)	MAB0952	CD3E monoclonal antibody, clone MEM-92
MAB5982	CD19 monoclonal antibody, clone SJ25-C1 (SPRD)	H00000919-M01	CD3Z monoclonal antibody (M01), clone 4A12-F6
MAB0665	CD1A monoclonal antibody, clone CBT6	MAB4425	CD4 monoclonal antibody, clone EDU-2 (FITC)
H00000910-D01	CD1B MaxPab rabbit polyclonal antibody (D01)	MAB3222	Cd4 monoclonal antibody, clone GK1.5
MAB6894	CD1C monoclonal antibody, clone M241	MAB5925	Cd4 monoclonal antibody, clone GK1.5 (APC/Cy5.5)
MAB4379	CD2 monoclonal antibody, clone LT2 (APC)	MAB5919	Cd4 monoclonal antibody, clone GK1.5 (Biotin)
MAB4381	CD2 monoclonal antibody, clone LT2 (Biotin)	MAB5922	Cd4 monoclonal antibody, clone GK1.5 (Cy5)
MAB4382	CD2 monoclonal antibody, clone LT2 (FITC)	MAB5918	Cd4 monoclonal antibody, clone GK1.5 (FITC)
MAB4380	CD2 monoclonal antibody, clone LT2 (PE)	MAB5923	Cd4 monoclonal antibody, clone GK1.5 (PE/Cy5.5)
MAB0869	CD2 monoclonal antibody, clone MEM-65	MAB5924	Cd4 monoclonal antibody, clone GK1.5 (PE/Cy7)
H00131450-D01	CD200R1 MaxPab rabbit polyclonal antibody (D01)	MAB5920	Cd4 monoclonal antibody, clone GK1.5 (PE/TXRD)
H00050489-D01	CD207 MaxPab rabbit polyclonal antibody (D01)	MAB5921	Cd4 monoclonal antibody, clone GK1.5 (SPRD)
H00030835-D01	CD209 MaxPab rabbit polyclonal antibody (D01)	MAB1319	CD4 monoclonal antibody, clone HCD4-A6
H00030835-M01	CD209 monoclonal antibody (M01), clone 3B10	MAB0871	CD4 monoclonal antibody, clone MEM-115
MAB0908	CD22 monoclonal antibody, clone MEM-01	MAB3627	CD4 monoclonal antibody, clone MEM-16
PAB8584	CD22 polyclonal antibody	MAB6896	CD4 monoclonal antibody, clone QS4120
MAB5270	CD226 monoclonal antibody, clone DX11	MAB5617	CD4 monoclonal antibody, clone RFT-4g
MAB6977	CD24 monoclonal antibody, clone BA-1	MAB6015	CD4 monoclonal antibody, clone RFT-4g (APC)
PAB17803	CD244 (phospho Y266) polyclonal antibody	MAB6020	CD4 monoclonal antibody, clone RFT-4g (APC/CY7)
MAB5125	Cd247 monoclonal antibody, clone H146-968	MAB6012	CD4 monoclonal antibody, clone RFT-4g (Biotin)
PAB12667	CD247 polyclonal antibody	MAB6017	CD4 monoclonal antibody, clone RFT-4g (Cy5)
MAB6214	Cd24a monoclonal antibody, clone 30-F1	MAB6011	CD4 monoclonal antibody, clone RFT-4g (FITC)
MAB6206	Cd28 monoclonal antibody, clone 37.51	MAB6013	CD4 monoclonal antibody, clone RFT-4g (PE)
MAB6904	CD28 monoclonal antibody, clone ANC28.1/5D10	MAB6018	CD4 monoclonal antibody, clone RFT-4g (PE/Cy5.5)
MAB2547	Cd28 monoclonal antibody, clone CD37.51.1	MAB6019	CD4 monoclonal antibody, clone RFT-4g (PE/Cy7)
MAB6944	CD29 monoclonal antibody, clone ANC28.1/5D10	MAB6014	CD4 monoclonal antibody, clone RFT-4g (PE/TXRD)
MAB6010	CD3 monoclonal antibody, clone UCHT1 (APC/CY7)	MAB6016	CD4 monoclonal antibody, clone RFT-4g (SPRD)
MAB6004	CD3 monoclonal antibody, clone UCHT1 (Biotin)	MAB5607	Cd40 monoclonal antibody, clone 1C10
MAB6007	CD3 monoclonal antibody, clone UCHT1 (Cy5)	MAB5934	Cd40 monoclonal antibody, clone 1C10 (APC)
MAB4402	CD3 monoclonal antibody, clone UCHT1 (PE)	MAB5932	Cd40 monoclonal antibody, clone 1C10 (Biotin)
MAB6008	CD3 monoclonal antibody, clone UCHT1 (PE/Cy5.5)	MAB5931	Cd40 monoclonal antibody, clone 1C10 (FITC)
MAB6009	CD3 monoclonal antibody, clone UCHT1 (PE/Cy7)	MAB5933	Cd40 monoclonal antibody, clone 1C10 (PE)
MAB4403	CD3 monoclonal antibody, clone UCHT1 (PE-Dyomics 647)	MAB5935	Cd40 monoclonal antibody, clone 1C10 (SPRD)
MAB6006	CD3 monoclonal antibody, clone UCHT1 (SPRD)	MAB6911	CD40 monoclonal antibody, clone BE-1
MAB8055	Cd34 monoclonal antibody, clone 3H1240	MAB8075	CD40 monoclonal antibody, clone LOB7/6
MAB4038	CD34 monoclonal antibody, clone Immu-133	PAB0251	CD40 polyclonal antibody
MAB7100	Cd34 monoclonal antibody, clone MEC 14.7	MAB3231	Cd40lg monoclonal antibody, clone MR1
MAB2325	Cd36 monoclonal antibody, clone D-2712	MAB5718	Cd40lg monoclonal antibody, clone MR1 (Biotin)
MAB3111	CD37 monoclonal antibody, clone IPO-24	MAB5719	Cd40lg monoclonal antibody, clone MR1 (PE)
MAB5606	Cd38 monoclonal antibody, clone 90	MAB0987	CD41/CD61 monoclonal antibody, clone IVA30
MAB5929	Cd38 monoclonal antibody, clone 90 (APC)	MAB4192	CD44 monoclonal antibody, clone 289-16211
MAB5927	Cd38 monoclonal antibody, clone 90 (Biotin)	MAB6914	CD44 monoclonal antibody, clone BU52
MAB5926	Cd38 monoclonal antibody, clone 90 (FITC)	MAB3936	CD44 monoclonal antibody, clone F10-44-2
MAB5928	Cd38 monoclonal antibody, clone 90 (PE)	MAB4105	Cd44 monoclonal antibody, clone IM7
MAB5930	Cd38 monoclonal antibody, clone 90 (SPRD)	MAB4193	CD44 monoclonal antibody, clone J.173
MAB6909	CD38 monoclonal antibody, clone AT1	MAB6204	Cd44 monoclonal antibody, clone KM201
MAB3220	Cd3e monoclonal antibody, clone 145-2C11	MAB3838	CD44 monoclonal antibody, clone MEM-263
MAB5908	Cd3e monoclonal antibody, clone 145-2C11 (Biotin)	MAB0885	CD44 monoclonal antibody, clone MEM-85
MAB5907	Cd3e monoclonal antibody, clone 145-2C11 (FITC)	MAB1354	Cd44 monoclonal antibody, clone T2-F4
MAB5909	Cd3e monoclonal antibody, clone 145-2C11 (PE)	PAB0912	CD44 polyclonal antibody

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB8746	CD44 polyclonal antibody	MAB4465	CD74 monoclonal antibody, clone LN2 (PE)
MAB6954	CD45R monoclonal antibody, clone 351C5	MAB6922	CD74 monoclonal antibody, clone M-B741
MAB3619	CD45R monoclonal antibody, clone IVA103	MAB6623	CD74 monoclonal antibody, clone PIN.1
MAB6953	CD45RA monoclonal antibody, clone 158.4D3	MAB6923	CD79A monoclonal antibody, clone ZL7-4
MAB5433	CD45RB monoclonal antibody, clone 16A (PE)	MAB6687	Cd79a/Cd79b monoclonal antibody, clone HM79
MAB6948	CD46 monoclonal antibody, clone 169-1-E4.3	H00000974-D01	CD79B MaxPab rabbit polyclonal antibody (D01)
MAB3937	CD46 monoclonal antibody, clone 236-11842 (MEM-258)	H00000974-M01	CD79B monoclonal antibody (M01), clone 4E10-2A10
MAB0930	CD46 monoclonal antibody, clone MEM-258	MAB5628	CD79B monoclonal antibody, clone CB3-1
MAB4439	CD46 monoclonal antibody, clone MEM-258 (APC)	MAB6215	Cd79b monoclonal antibody, clone HM79-12
MAB4440	CD46 monoclonal antibody, clone MEM-258 (Biotin)	MAB3844	CD8 monoclonal antibody, clone MEM-31
MAB4441	CD46 monoclonal antibody, clone MEM-258 (FITC)	MAB0936	CD8 monoclonal antibody, clone MEM-87
MAB4442	CD46 monoclonal antibody, clone MEM-258 (PE)	MAB4466	CD8 monoclonal antibody, clone MEM-87 (Biotin)
MAB5279	CD47 monoclonal antibody, clone BRIC126	MAB1360	CD8.1 monoclonal antibody, clone CD8-E1
MAB0889	CD48 monoclonal antibody, clone MEM-102	MAB1361	CD8.2 monoclonal antibody, clone CD8-F4
MAB8358	CD49e monoclonal antibody, clone SAM-1	H00000941-D01	CD80 MaxPab rabbit polyclonal antibody (D01)
MAB6217	CD4A monoclonal antibody, clone 74-12-4	H00000941-M02A	CD80 monoclonal antibody (M02A), clone 3E10
MAB3840	CD5 monoclonal antibody, clone CRIS1	H00000941-M03A	CD80 monoclonal antibody (M03A), clone 4A4
MAB4443	CD5 monoclonal antibody, clone CRIS1 (APC)	MAB5610	Cd80 monoclonal antibody, clone 1G10
MAB4444	CD5 monoclonal antibody, clone CRIS1 (Biotin)	MAB5951	Cd80 monoclonal antibody, clone 1G10 (APC)
MAB4445	CD5 monoclonal antibody, clone CRIS1 (FITC)	MAB5949	Cd80 monoclonal antibody, clone 1G10 (Biotin)
MAB4446	CD5 monoclonal antibody, clone CRIS1 (PE)	MAB5948	Cd80 monoclonal antibody, clone 1G10 (FITC)
MAB7188	CD5 monoclonal antibody, clone H65	MAB5950	Cd80 monoclonal antibody, clone 1G10 (PE)
MAB0872	CD5 monoclonal antibody, clone MEM-32	MAB5952	Cd80 monoclonal antibody, clone 1G10 (SPRD)
MAB6916	CD53 monoclonal antibody, clone 63.5A3	MAB3845	CD80 monoclonal antibody, clone MEM-233
MAB0890	CD53 monoclonal antibody, clone MEM-53	MAB5306	CD81 monoclonal antibody, clone 1.3.3.22
H00001604-D01	CD55 MaxPab rabbit polyclonal antibody (D01)	MAB3248	Cd81 monoclonal antibody, clone 2F7
MAB0893	CD55 monoclonal antibody, clone MEM-118	MAB5792	Cd81 monoclonal antibody, clone 2F7 (Biotin)
MAB0895	CD58 monoclonal antibody, clone MEM-63	MAB5791	Cd81 monoclonal antibody, clone 2F7 (FITC)
MAB6918	CD59 monoclonal antibody, clone BRA-10G	MAB5793	Cd81 monoclonal antibody, clone 2F7 (PE)
MAB0896	CD59 monoclonal antibody, clone MEM-43	MAB6435	CD81 monoclonal antibody, clone M38
MAB3628	CD59 monoclonal antibody, clone MEM-43/5	H00009308-M01	CD83 monoclonal antibody (M01), clone 3G10-1F4
MAB1719	CD59 monoclonal antibody, clone YUK1	MAB6970	CD84 monoclonal antibody, clone 152.1D5
H00000922-M01	CD5L monoclonal antibody (M01), clone 1C8	MAB5203	CD84 monoclonal antibody, clone 152-1D5
MAB3841	CD6 monoclonal antibody, clone MEM-98	MAB5307	CD84 monoclonal antibody, clone 153-4D9
MAB5280	CD60 monoclonal antibody, clone M-T6004	H00000942-M01	CD86 monoclonal antibody (M01), clone 3D7
MAB6982	CD60b monoclonal antibody, clone UM4D4	MAB3846	CD86 monoclonal antibody, clone BU63
MAB0988	CD62L monoclonal antibody, clone IVA94	MAB4467	CD86 monoclonal antibody, clone BU63 (APC)
MAB6919	CD63 monoclonal antibody, clone AHN16.1/46-4-5	MAB4468	CD86 monoclonal antibody, clone BU63 (FITC)
MAB0931	CD63 monoclonal antibody, clone MEM-259	MAB4469	CD86 monoclonal antibody, clone BU63 (PE)
MAB4456	CD63 monoclonal antibody, clone MEM-259 (APC)	MAB3238	Cd86 monoclonal antibody, clone GL1
MAB4457	CD63 monoclonal antibody, clone MEM-259 (FITC)	MAB5754	Cd86 monoclonal antibody, clone GL1 (APC)
MAB4458	CD63 monoclonal antibody, clone MEM-259 (PE)	MAB5752	Cd86 monoclonal antibody, clone GL1 (Biotin)
MAB4459	CD63 monoclonal antibody, clone MEM-259 (PE-Dyomics 647)	MAB5753	Cd86 monoclonal antibody, clone GL1 (PE)
H00000969-M07	CD69 monoclonal antibody (M07), clone 4H3	MAB5756	Cd86 monoclonal antibody, clone GL1 (PE/Cy7)
MAB3236	Cd69 monoclonal antibody, clone H1.2F3	MAB5755	Cd86 monoclonal antibody, clone GL1 (SPRD)
MAB5745	Cd69 monoclonal antibody, clone H1.2F3 (APC)	H00000925-M10	CD8A monoclonal antibody (M10), clone 4B9
MAB5743	Cd69 monoclonal antibody, clone H1.2F3 (Biotin)	MAB6899	CD8A monoclonal antibody, clone 14
MAB5742	Cd69 monoclonal antibody, clone H1.2F3 (FITC)	MAB5605	Cd8a monoclonal antibody, clone 53-6.7
MAB5744	Cd69 monoclonal antibody, clone H1.2F3 (PE)	MAB4109	Cd8a monoclonal antibody, clone 53-6.72
MAB5747	Cd69 monoclonal antibody, clone H1.2F3 (PE/Cy7)	MAB6218	CD8A monoclonal antibody, clone 76-2-11
MAB5746	Cd69 monoclonal antibody, clone H1.2F3 (SPRD)	MAB7147	Cd8a monoclonal antibody, clone IBL-3/25
MAB6920	CD69 monoclonal antibody, clone HP-4B3	MAB7110	CD8A monoclonal antibody, clone M211
MAB4667	CD69 monoclonal antibody, clone TP1.55.3	MAB7109	CD8A monoclonal antibody, clone M236
H00000924-M04	CD7 monoclonal antibody (M04), clone 1B8	MAB6898	CD8A monoclonal antibody, clone UCHT4
MAB6897	CD7 monoclonal antibody, clone 3A1e	H00000926-D01	CD8B MaxPab rabbit polyclonal antibody (D01)
MAB0873	CD7 monoclonal antibody, clone MEM-186	MAB4669	CD8B monoclonal antibody, clone 2S78.5H7
H00000970-D01	CD70 MaxPab rabbit polyclonal antibody (D01)	MAB6462	Cd9 monoclonal antibody, clone EM-04
MAB6921	CD70 monoclonal antibody, clone BU69	MAB0985	CD9 monoclonal antibody, clone IVA50
MAB0925	CD72 monoclonal antibody, clone 3F3	MAB6900	CD9 monoclonal antibody, clone SN4/C3-3A2
MAB4462	CD72 monoclonal antibody, clone 3F3 (Biotin)	H00010225-D01	CD96 MaxPab rabbit polyclonal antibody (D01)
MAB4463	CD72 monoclonal antibody, clone 3F3 (PE)	MAB0910	CD97 monoclonal antibody, clone MEM-180
H00000972-M01	CD74 monoclonal antibody (M01), clone 1D1	H00000978-D01	CDA MaxPab rabbit polyclonal antibody (D01)
H00000972-M01A	CD74 monoclonal antibody (M01A), clone 1D1	H00081602-D01	CDADC1 MaxPab rabbit polyclonal antibody (D01)
MAB0983	CD74 monoclonal antibody, clone LN2	H00008872-D01	CDC123 MaxPab rabbit polyclonal antibody (D01)
MAB4464	CD74 monoclonal antibody, clone LN2 (APC)	MAB7234	CDC14A monoclonal antibody, clone DCS-291

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00008555-D01	CDC14B MaxPab rabbit polyclonal antibody (D01)	PAB12373	CENPF polyclonal antibody
PAB12614	CDC2 (phospho Y15) polyclonal antibody	H00064946-D01	CENPH MaxPab rabbit polyclonal antibody (D01)
MAB1001	CDC2 monoclonal antibody, clone 26E11	H00055835-D01	CENPJ MaxPab rabbit polyclonal antibody (D01)
MAB1000	CDC2 monoclonal antibody, clone 6C9	H00055835-M02	CENPJ monoclonal antibody (M02), clone 1A5
MAB1921	CDC2 monoclonal antibody, clone POH-1	H00079019-D01	CENPM MaxPab rabbit polyclonal antibody (D01)
PAB12668	CDC2 polyclonal antibody	H00055839-D01	CENPN MaxPab rabbit polyclonal antibody (D01)
PAB12757	CDC2 polyclonal antibody	H00055165-D01	CEP55 MaxPab rabbit polyclonal antibody (D01)
PAB1300	CDC2 polyclonal antibody	H00009702-D01	CEP57 MaxPab rabbit polyclonal antibody (D01)
PAB9939	CDC2 polyclonal antibody	H00009350-D01	CER1 MaxPab rabbit polyclonal antibody (D01)
H00000993-D01	CDC25A MaxPab rabbit polyclonal antibody (D01)	H00009350-M04	CER1 monoclonal antibody (M04), clone 3E11
PAB14480	CDC25A polyclonal antibody	PAB9956	CES1 polyclonal antibody
H00000995-D01	CDC25C MaxPab rabbit polyclonal antibody (D01)	H00008824-D01	CES2 MaxPab rabbit polyclonal antibody (D01)
H00000995-M01	CDC25C monoclonal antibody (M01), clone 3B11	H00008824-M02	CES2 monoclonal antibody (M02), clone 4F12
PAB9987	CDC27 (phospho S427) polyclonal antibody	H00001068-D01	CETN1 MaxPab rabbit polyclonal antibody (D01)
PAB9988	CDC27 (phospho T244) polyclonal antibody	H00055997-D01	CFC1 MaxPab rabbit polyclonal antibody (D01)
MAB7304	CDC27 monoclonal antibody, clone AF3.1	H00055997-M09	CFC1 monoclonal antibody (M09), clone 2G4
H00023097-D01	CDC2L6 MaxPab rabbit polyclonal antibody (D01)	MAB0768	CFH monoclonal antibody, clone 63G5
MAB8265	CDC34 monoclonal antibody	H00003078-D01	CFHR1 MaxPab rabbit polyclonal antibody (D01)
PAB15584	CDC37 polyclonal antibody	H00003080-D01	CFHR2 MaxPab rabbit polyclonal antibody (D01)
H00008318-D01	CDC45L MaxPab rabbit polyclonal antibody (D01)	H00081494-D01	CFHR5 MaxPab rabbit polyclonal antibody (D01)
MAB7254	cdc6 monoclonal antibody, clone Cdc6 9H8/5	MAB7963	CFI monoclonal antibody, clone OX-21
MAB7233	CDC7 monoclonal antibody, clone DCS-341	H00008837-D01	CFLAR MaxPab rabbit polyclonal antibody (D01)
MAB7030	CDC7 monoclonal antibody, clone SPM171	H00008837-M03	CFLAR monoclonal antibody (M03), clone 1G1
PAB14104	CDC73 polyclonal antibody	MAB3465	CFLAR monoclonal antibody, clone Dave-2
H00055143-D01	CDCA8 MaxPab rabbit polyclonal antibody (D01)	PAB0257	CFLAR polyclonal antibody
MAB9084	CDH1 monoclonal antibody, clone 61a3	MAB8760	CFP monoclonal antibody, clone 98s55
MAB1388	Cdh1 monoclonal antibody, clone M168	MAB6384	CGA monoclonal antibody, clone CH5
MAB1381	CDH2 monoclonal antibody, clone M170	MAB6385	CGA monoclonal antibody, clone CH7
PAB7891	CDH2 polyclonal antibody	MAB5316	CGA monoclonal antibody, clone INN-hFSH-132
MAB4657	CDH5 monoclonal antibody, clone TEA1/31	H00001082-D01	CGB MaxPab rabbit polyclonal antibody (D01)
MAB7197	CDK16 monoclonal antibody, clone G6.1	MAB9411	Cha monoclonal antibody, clone 54d9w
H00001017-M01A	CDK2 monoclonal antibody (M01A), clone 3A2-3G6	MAB2356	CHAF1A monoclonal antibody, clone SS 1 1-13
MAB8066	CDK2 monoclonal antibody, clone 4i33	MAB2354	CHAF1B monoclonal antibody, clone SS 24 1-68
PAB9940	CDK2 polyclonal antibody	MAB2355	CHAF1B monoclonal antibody, clone SS 53
H00001020-D01	CDK5 MaxPab rabbit polyclonal antibody (D01)	MAB2720	CHEK1 monoclonal antibody, clone AF7H4
MAB0748	CDK5 monoclonal antibody, clone 2G2	PAB9841	CHEK1 polyclonal antibody
MAB0747	CDK5 monoclonal antibody, clone 8A1	PAB10034	CHEK2 (phospho T68) polyclonal antibody
PAB9990	CDK5 polyclonal antibody	H00011200-D01	CHEK2 MaxPab rabbit polyclonal antibody (D01)
H00001021-D03	CDK6 MaxPab rabbit polyclonal antibody (D03)	PAB11360	CHK1 polyclonal antibody
MAB0830	CDK7 monoclonal antibody, clone 3E3	MAB7112	<i>Chlamydia trachomatis</i> HSPE1 monoclonal antibody, clone M1.2
MAB0831	CDK7 monoclonal antibody, clone 45H1	MAB7111	<i>Chlamydia trachomatis</i> HSPE1 monoclonal antibody, clone M1.4
MAB1855	CDK7 monoclonal antibody, clone MO-1	PAB13910	Cholera Toxin beta subunit polyclonal antibody
PAB10000	CDK9 (phospho T29) polyclonal antibody	PAB11403	Cholesterol Oxidase polyclonal antibody
PAB9941	CDK9 polyclonal antibody	MAB9062	Chromosome 12 monoclonal antibody, clone 653t6
MAB2028	CDKN1A monoclonal antibody, clone DCS-60	PAB17111	ch-TOG polyclonal antibody
MAB2029	Cdkn1b monoclonal antibody, clone DCS-72	MAB0052	CHUK monoclonal antibody, clone 14A231
PAB10300	CDKN1B polyclonal antibody	MAB2676	CHUK monoclonal antibody, clone 3A11
MAB2416	Cdkn2a monoclonal antibody, clone 12-A1-1	PAB8780	CHUK polyclonal antibody
MAB2417	Cdkn2a monoclonal antibody, clone 5-C3-1	PAB5138	Chylomicrons polyclonal antibody
MAB2031	CDKN2A monoclonal antibody, clone DCS50	PAB12457	CIC polyclonal antibody
MAB2030	CDKN2B monoclonal antibody, clone DCS114	PAB12458	CIC polyclonal antibody
H00001032-M08	CDKN2D monoclonal antibody (M08), clone 2E10	H00001149-D01	CIDEA MaxPab rabbit polyclonal antibody (D01)
MAB2027	CDKN2D monoclonal antibody, clone DCS-100	H00001149-M01	CIDEA monoclonal antibody (M01), clone 4B9
H00090273-D01	CEACAM21 MaxPab rabbit polyclonal antibody (D01)	PAB0658	CIITA polyclonal antibody
MAB3848	CEACAM5 monoclonal antibody, clone CB30	PAB12758	CISH polyclonal antibody
MAB3456	CEACAM5 monoclonal antibody, clone COL-1	H00004435-D01	CITED1 MaxPab rabbit polyclonal antibody (D01)
MAB7986	CEBPB monoclonal antibody, clone 21A1	MAB2352	CITED4 monoclonal antibody, clone HT13-2D6.3
PAB10182	Cebpd polyclonal antibody	PAB12456	CKAP5 polyclonal antibody
H00001054-D01	CEBPG MaxPab rabbit polyclonal antibody (D01)	H00001158-M01	CKM monoclonal antibody (M01), clone 3E1-F3
H00001054-M03	CEBPG monoclonal antibody (M03), clone S2	H00001158-M06	CKM monoclonal antibody (M06), clone 2F3-B11
H00051816-M01	CECR1 monoclonal antibody (M01), clone 3D11	MAB9316	CKM monoclonal antibody, clone 79w8
H00001059-D01	CENPB MaxPab rabbit polyclonal antibody (D01)	PAB5007	CKM polyclonal antibody
MAB1924	CENPE monoclonal antibody, clone 1H12	H00001159-D01	CKMT1B MaxPab rabbit polyclonal antibody (D01)
MAB8966	CENPE monoclonal antibody, clone CENP-E.1H12	H00009076-D01	CLDN1 MaxPab rabbit polyclonal antibody (D01)
MAB6121	CENPF monoclonal antibody	H00009076-M01	CLDN1 monoclonal antibody (M01), clone 1C5-D9

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00009076-M02	CLDN1 monoclonal antibody (M02), clone 2E2-H5	MAB6926	CR2 monoclonal antibody, clone BU33
H00023562-D01	CLDN14 MaxPab rabbit polyclonal antibody (D01)	PAB13793	CRADD polyclonal antibody
H00023562-M01	CLDN14 monoclonal antibody (M01), clone 3D11	MAB9317	Creatine kinase monoclonal antibody, clone 4s8c9
H00001364-D01	CLDN4 MaxPab rabbit polyclonal antibody (D01)	MAB9315	Creatine kinase monoclonal antibody, clone 8d79
H00051266-M02	CLEC1B monoclonal antibody (M02), clone 1C10	PAB12522	CREB1 (phospho S121) polyclonal antibody
H00009022-D01	CLIC3 MaxPab rabbit polyclonal antibody (D01)	PAB12617	CREB1 (phospho S133) polyclonal antibody
H00009022-M02	CLIC3 monoclonal antibody (M02), clone 3F8	PAB12673	CREB1 (phospho S133) polyclonal antibody
PAB12145	Clock polyclonal antibody	H00009586-D01	CREB5 MaxPab rabbit polyclonal antibody (D01)
H00001208-D01	CLPS MaxPab rabbit polyclonal antibody (D01)	H00009586-M01	CREB5 monoclonal antibody (M01), clone 1E2
H00001208-M01	CLPS monoclonal antibody (M01), clone 8A8	H00001393-D01	CRHBP MaxPab rabbit polyclonal antibody (D01)
PAB8655	CLSPN polyclonal antibody	H00001393-M01	CRHBP monoclonal antibody (M01), clone 3D9
H00001212-D01	CLTB MaxPab rabbit polyclonal antibody (D01)	H00009244-M01	CRLF1 monoclonal antibody (M01), clone 4F4
MAB3482	CLU monoclonal antibody, clone CLI-9	H00051379-D01	CRLF3 MaxPab rabbit polyclonal antibody (D01)
PAB9973	CLU polyclonal antibody	H00001421-D01	CRYGD MaxPab rabbit polyclonal antibody (D01)
PAB13920	c-Myc tag polyclonal antibody	H00001428-D01	CRYM MaxPab rabbit polyclonal antibody (D01)
H00084735-D01	CNDP1 MaxPab rabbit polyclonal antibody (D01)	H00001428-M03	CRYM monoclonal antibody (M03), clone 6B3
H00154043-D01	CNKSR3 MaxPab rabbit polyclonal antibody (D01)	H00001429-D01	CRYZ MaxPab rabbit polyclonal antibody (D01)
H00154043-M01	CNKSR3 monoclonal antibody (M01), clone 4A11	PAB5024	CS-100 polyclonal antibody
H00001267-D01	CNP MaxPab rabbit polyclonal antibody (D01)	H00008531-D01	CSDA MaxPab rabbit polyclonal antibody (D01)
H00001270-D01	CNTF MaxPab rabbit polyclonal antibody (D01)	H00030818-M01	CSEN monoclonal antibody (M01), clone 2B3-A1
H00080347-D01	COASY MaxPab rabbit polyclonal antibody (D01)	PAB12618	CSF1R (phospho Y723) polyclonal antibody
H00091522-D01	COL23A1 MaxPab rabbit polyclonal antibody (D01)	PAB12674	CSF1R polyclonal antibody
PAB10190	Collagen Type I polyclonal antibody	PAB12675	CSF1R polyclonal antibody
PAB10191	Collagen Type II polyclonal antibody	H00001440-D01	CSF3 MaxPab rabbit polyclonal antibody (D01)
PAB10192	Collagen Type III polyclonal antibody	H00001440-M03	CSF3 monoclonal antibody (M03), clone 2C5
PAB10193	Collagen Type IV polyclonal antibody	MAB4099	CSF3R monoclonal antibody, clone BDI542
PAB10194	Collagen Type V polyclonal antibody	MAB6261	CSF3R monoclonal antibody, clone S-1252
PAB10195	Collagen Type VI polyclonal antibody	MAB6262	CSF3R monoclonal antibody, clone S-1268
MAB9478	Collagenase monoclonal antibody, clone 4h6y	MAB6263	CSF3R monoclonal antibody, clone S-1284
H00150684-M01	COMMD1 monoclonal antibody (M01), clone 2A12	MAB6264	CSF3R monoclonal antibody, clone S-1374
H00150684-M03	COMMD1 monoclonal antibody (M03), clone 2B12	MAB6266	CSF3R monoclonal antibody, clone S-1390
MAB3751	COMP monoclonal antibody, clone 6G5	MAB6265	CSF3R monoclonal antibody, clone S-425
MAB3750	Comp monoclonal antibody, clone 8H7	H00055790-D01	CSGALNACT1 MaxPab rabbit polyclonal antibody (D01)
PAB13707	COMP polyclonal antibody	PAB5015	CSH1 polyclonal antibody
H00001312-D01	COMT MaxPab rabbit polyclonal antibody (D01)	H00001445-M01	CSK monoclonal antibody (M01), clone 3A3
MAB7107	COPB2 monoclonal antibody, clone M3A5	MAB0950	CSK monoclonal antibody, clone CSK-04
MAB8267	COPS5 monoclonal antibody	MAB3766	CSK monoclonal antibody, clone PAG-C1
MAB7987	Cops5 monoclonal antibody, clone 2A10.8	H00001446-D01	CSN1S1 MaxPab rabbit polyclonal antibody (D01)
MAB7988	Cops5 monoclonal antibody, clone 6C3.38	MAB2733	CSNK1D monoclonal antibody, clone AF12G4
MAB9239	cora monoclonal antibody, clone 5df4	MAB2714	CSNK1E monoclonal antibody, clone AF26H1
PAB12760	Corin polyclonal antibody	MAB2713	CSNK1E monoclonal antibody, clone AF6C1
PAB12761	Corin polyclonal antibody	MAB9427	Csp monoclonal antibody, clone 778mk
PAB12762	Corin polyclonal antibody	MAB9426	Csp monoclonal antibody, clone 868f4
PAB4980	CORT polyclonal antibody	MAB7129	CSPG4 monoclonal antibody, clone LHM 2
PAB5132	CORT polyclonal antibody	H00001478-D01	CSTF2 MaxPab rabbit polyclonal antibody (D01)
PAB5147	Cortexolone polyclonal antibody	H00001486-D01	CTBS MaxPab rabbit polyclonal antibody (D01)
PAB5148	Cortexolone polyclonal antibody	H00001491-D01	CTH MaxPab rabbit polyclonal antibody (D01)
PAB5149	Corticosterone 21-acetate polyclonal antibody	MAB3241	Ctla4 monoclonal antibody, clone 9H10
PAB5151	Cortisol polyclonal antibody	MAB6927	CTLA4 monoclonal antibody, clone ANC152.2/8H5
PAB5152	Cortisone polyclonal antibody	MAB3923	CTNNA1 monoclonal antibody, clone BDI950
MAB9413	cos monoclonal antibody, clone 763gh	H00029119-D01	CTNNA3 MaxPab rabbit polyclonal antibody (D01)
PAB8657	COX1 polyclonal antibody	PAB12603	CTNNB1 (phospho S33) polyclonal antibody
H00084701-M01	COX4I2 monoclonal antibody (M01), clone 1F2	PAB12604	CTNNB1 (phospho S33/37) polyclonal antibody
MAB0777	CP monoclonal antibody, clone 3B11	PAB12661	CTNNB1 (phospho S37) polyclonal antibody
PAB5005	CP polyclonal antibody	PAB12605	CTNNB1 (phospho Y228) polyclonal antibody
H00001357-D01	CPA1 MaxPab rabbit polyclonal antibody (D01)	PAB12606	CTNNB1 (phospho Y280) polyclonal antibody
H00001357-M01	CPA1 monoclonal antibody (M01), clone 3F11	PAB7925	CTNNB1 (phospho Y86) polyclonal antibody
H00001358-D01	CPA2 MaxPab rabbit polyclonal antibody (D01)	PAB12607	CTNNB1 (phospho Y96) polyclonal antibody
H00001358-M02	CPA2 monoclonal antibody (M02), clone 2E11	MAB0739	CTNNB1 monoclonal antibody, clone 1C10
H00001360-D01	CPB1 MaxPab rabbit polyclonal antibody (D01)	MAB3927	CTNNB1 monoclonal antibody, clone BDI080
H00001360-M01	CPB1 monoclonal antibody (M01), clone 4D5	MAB3926	CTNNB1 monoclonal antibody, clone BDI109
H00001361-D01	CPB2 MaxPab rabbit polyclonal antibody (D01)	MAB3930	CTNNB1 monoclonal antibody, clone BDI710
H00008904-D01	CPNE1 MaxPab rabbit polyclonal antibody (D01)	MAB3929	CTNNB1 monoclonal antibody, clone BDI770
MAB6925	CR1 monoclonal antibody, clone E11	MAB3928	CTNNB1 monoclonal antibody, clone BDI870
MAB9482	CR2 monoclonal antibody, clone 0y65	MAB3924	CTNNB1 monoclonal antibody, clone BDI920

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB6470	CTNNB1 monoclonal antibody, clone EM-22	H00115265-D01	DDIT4L MaxPab rabbit polyclonal antibody (D01)
MAB1369	CTNNB1 monoclonal antibody, clone M118	H00001652-M04	DDT monoclonal antibody (M04), clone 1D5
PAB1228	CTNNB1 polyclonal antibody	H00168400-D01	DDX53 MaxPab rabbit polyclonal antibody (D01)
PAB1233	CTNNB1 polyclonal antibody	H00079039-D01	DDX54 MaxPab rabbit polyclonal antibody (D01)
PAB12658	CTNNB1 polyclonal antibody	H00010522-M01	DEAF1 monoclonal antibody (M01), clone 6F10
PAB12659	CTNNB1 polyclonal antibody	H00026063-D01	DECR2 MaxPab rabbit polyclonal antibody (D01)
PAB12660	CTNNB1 polyclonal antibody	PAB0255	DEDD2 polyclonal antibody
PAB7883	CTNNB1 polyclonal antibody	H00091614-D01	DEPDC7 MaxPab rabbit polyclonal antibody (D01)
MAB7329	Ctnnd1 monoclonal antibody, clone 6H11	PAB12619	DES (phospho T76/T77) polyclonal antibody
H00001509-D01	CTSD MaxPab rabbit polyclonal antibody (D01)	PAB12677	DES polyclonal antibody
H00001509-M01	CTSD monoclonal antibody (M01), clone 3F12-1B9	PAB12678	DES polyclonal antibody
H00001511-D01	CTSG MaxPab rabbit polyclonal antibody (D01)	PAB5150	Desoxycorticosterone polyclonal antibody
PAB8634	CTSH polyclonal antibody	PAB11387	Dextranase polyclonal antibody
PAB8635	CTSK polyclonal antibody	H00001676-D01	DFFA MaxPab rabbit polyclonal antibody (D01)
PAB8636	Ctsk polyclonal antibody	MAB2702	DFFA monoclonal antibody, clone AF24G5
H00001520-D01	CTSS MaxPab rabbit polyclonal antibody (D01)	MAB2703	DFFA monoclonal antibody, clone AF41E3
H00001523-M01	CUTL1 monoclonal antibody (M01), clone 2A10	H00001718-D01	DHCR24 MaxPab rabbit polyclonal antibody (D01)
H00001523-D01	CUX1 MaxPab rabbit polyclonal antibody (D01)	H00079947-D01	DHDDS MaxPab rabbit polyclonal antibody (D01)
PAB8729	CX3CL1 polyclonal antibody	H00027294-D01	DHDH MaxPab rabbit polyclonal antibody (D01)
PAB8798	Cxcl1 polyclonal antibody	H00010170-D01	DHRS9 MaxPab rabbit polyclonal antibody (D01)
H00006373-M02A	CXCL11 monoclonal antibody (M02A), clone 3C1	H00010170-M06	DHRS9 monoclonal antibody (M06), clone 3E3
PAB9504	Cxcl12 polyclonal antibody	H00079132-D01	DHX58 MaxPab rabbit polyclonal antibody (D01)
H00010563-D01	CXCL13 MaxPab rabbit polyclonal antibody (D01)	MAB3402	Diablo monoclonal antibody, clone 9H10
PAB9505	Cxcl2 polyclonal antibody	PAB11261	DIABLO polyclonal antibody
PAB8943	Cxcr3 polyclonal antibody	PAB14673	DIABLO polyclonal antibody
PAB8944	CXCR3 polyclonal antibody	PAB10201	DIAPH2 polyclonal antibody
PAB9506	CXCR4 polyclonal antibody	PAB5165	Dihydrotestosterone polyclonal antibody
PAB9850	CXCR4 polyclonal antibody	PAB5166	Dihydrotestosterone polyclonal antibody
PAB8945	CXCR6 polyclonal antibody	H00022943-M08	DKK1 monoclonal antibody (M08), clone 2B12
H00001528-M06	CYB5A monoclonal antibody (M06), clone 1A8	H00022943-M11	DKK1 monoclonal antibody (M11), clone 2A5
MAB9115	CycA monoclonal antibody, clone 1a35	H00010395-D01	DLC1 MaxPab rabbit polyclonal antibody (D01)
MAB6240	Cyca monoclonal antibody, clone 6H2.B4	PAB10259	DLD polyclonal antibody
PAB8666	Cyca polyclonal antibody	MAB8585	Dlg1 monoclonal antibody, clone S64-15
H00001543-D01	CYP1A1 MaxPab rabbit polyclonal antibody (D01)	MAB7888	Dlg2 monoclonal antibody, clone S18-30
PAB13660	CYP1A1/CYP1A2 polyclonal antibody	MAB6630	Dlg3 monoclonal antibody, clone S19-2
H00001593-D01	CYP27A1 MaxPab rabbit polyclonal antibody (D01)	MAB8575	DLG4 monoclonal antibody, clone S28-43
H00001571-D01	CYP2E1 MaxPab rabbit polyclonal antibody (D01)	MAB7900	Dlgap1 monoclonal antibody, clone S127-31
H00001576-D01	CYP3A4 MaxPab rabbit polyclonal antibody (D01)	MAB8576	DLGAP1 monoclonal antibody, clone S28-86
H00064816-D01	CYP3A43 MaxPab rabbit polyclonal antibody (D01)	H00001743-D01	DLST MaxPab rabbit polyclonal antibody (D01)
H00010858-M02	CYP46A1 monoclonal antibody (M02), clone 2B5	H00001745-M01	DLX1 monoclonal antibody (M01), clone 2H3
H00001580-D01	CYP4B1 MaxPab rabbit polyclonal antibody (D01)	H00001745-M03	DLX1 monoclonal antibody (M03), clone 3F9
H00009266-D01	CYTH2 MaxPab rabbit polyclonal antibody (D01)	H00001745-M07	DLX1 monoclonal antibody (M07), clone 1D11
PAB9852	Cytochrome c polyclonal antibody	H00001745-M14	DLX1 monoclonal antibody (M14), clone 4H7
PAB15500	Cytosine (5-Methyl) polyclonal antibody	H00001748-D01	DLX4 MaxPab rabbit polyclonal antibody (D01)
PAB14892	Dab1 polyclonal antibody	H00001748-M01	DLX4 monoclonal antibody (M01), clone 1F11
PAB9856	Dab1 polyclonal antibody	H00001756-M07	DMD monoclonal antibody (M07), clone 3A11
PAB9857	Dab1 polyclonal antibody	H00001760-D01	DMPK MaxPab rabbit polyclonal antibody (D01)
PAB12557	Daf-3 polyclonal antibody	H00064446-D01	DNAI2 MaxPab rabbit polyclonal antibody (D01)
PAB10199	DAO1 polyclonal antibody	MAB3738	DNAJA1 monoclonal antibody, clone KA2A5.6
PAB9853	DAP12 polyclonal antibody	H00003337-D01	DNAJB1 MaxPab rabbit polyclonal antibody (D01)
MAB7223	DAPK1 (phospho S308) monoclonal antibody, clone DKPS308	MAB2185	DNAJB1 monoclonal antibody, clone 3B9.E6
H00002532-D01	DARC MaxPab rabbit polyclonal antibody (D01)	PAB15411	DNAJB1 polyclonal antibody
PAB0254	DAXX polyclonal antibody	H00003300-D01	DNAJB2 MaxPab rabbit polyclonal antibody (D01)
H00001620-D01	DBC1 MaxPab rabbit polyclonal antibody (D01)	H00011080-M02	DNAJB4 monoclonal antibody (M02), clone S13
MAB0720	DBI monoclonal antibody, clone 20B10	H00010049-D01	DNAJB6 MaxPab rabbit polyclonal antibody (D01)
MAB0725	DBI monoclonal antibody, clone 27C9	H00010049-M02	DNAJB6 monoclonal antibody (M02), clone 2D12-B9
PAB1215	DBI polyclonal antibody	H00010049-M03	DNAJB6 monoclonal antibody (M03), clone 1A5
H00079877-D01	DCAKD MaxPab rabbit polyclonal antibody (D01)	H00007266-D01	DNAJC7 MaxPab rabbit polyclonal antibody (D01)
H00001633-D03	DCK MaxPab rabbit polyclonal antibody (D03)	H00007802-D01	DNALI1 MaxPab rabbit polyclonal antibody (D01)
PAB2020	DCK polyclonal antibody	H00007802-M04A	DNALI1 monoclonal antibody (M04A), clone 2H3
MAB9245	DCN monoclonal antibody, clone 75h6	PAB9913	DNASE1 polyclonal antibody
H00055802-D01	DCP1A MaxPab rabbit polyclonal antibody (D01)	H00001777-D01	DNASE2 MaxPab rabbit polyclonal antibody (D01)
H00054165-D01	DCUN1D1 MaxPab rabbit polyclonal antibody (D01)	MAB0079	DNMT1 monoclonal antibody, clone 60B1220.1
H00051181-D01	DCXR MaxPab rabbit polyclonal antibody (D01)	PAB12571	DNMT1 polyclonal antibody
PAB8734	DDIT3 polyclonal antibody	PAB15590	Dnmt1 polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB15591	Dnmt1 polyclonal antibody	MAB1914	EGFR monoclonal antibody, clone H11
H00001788-D01	DNMT3A MaxPab rabbit polyclonal antibody (D01)	PAB10223	EGFR polyclonal antibody
MAB0080	Dnmt3a monoclonal antibody, clone 64B1446	PAB12122	Egln1 polyclonal antibody
PAB14129	DNMT3A polyclonal antibody	PAB12440	EGLN1 polyclonal antibody
PAB14130	DNMT3A polyclonal antibody	PAB9915	ego-1 polyclonal antibody
PAB14131	DNMT3A polyclonal antibody	PAB9863	EIA polyclonal antibody
MAB0078	Dnmt3b monoclonal antibody, clone 52A1018	PAB10175	Eif2ak3 polyclonal antibody
PAB12444	DNMT3B polyclonal antibody	H00008892-D01	EIF2B2 MaxPab rabbit polyclonal antibody (D01)
PAB14125	Dnmt3b polyclonal antibody	H00008672-D01	EIF4G3 MaxPab rabbit polyclonal antibody (D01)
PAB15592	Dnmt3b polyclonal antibody	H00056648-D01	EIF5A2 MaxPab rabbit polyclonal antibody (D01)
H00029947-D01	DNMT3L MaxPab rabbit polyclonal antibody (D01)	H00010136-D01	ELA3A MaxPab rabbit polyclonal antibody (D01)
MAB7902	DNMT3L monoclonal antibody, clone S117-9	H00055520-D01	ELAC1 MaxPab rabbit polyclonal antibody (D01)
PAB10043	DNMT3L polyclonal antibody	H00001994-D01	ELAVL1 MaxPab rabbit polyclonal antibody (D01)
PAB8995	DNTT polyclonal antibody	H00001997-D01	ELF1 MaxPab rabbit polyclonal antibody (D01)
H00001796-D01	DOK1 MaxPab rabbit polyclonal antibody (D01)	PAB18470	ELK1 polyclonal antibody
MAB6231	DPP4 monoclonal antibody, clone 134-2C2	PAB12345	EMD polyclonal antibody
MAB5299	DPP4 monoclonal antibody, clone 202.36	H00011117-D01	EMILIN1 MaxPab rabbit polyclonal antibody (D01)
MAB5649	Dpp4 monoclonal antibody, clone 5E8	MAB9305	ena monoclonal antibody, clone 65k5
MAB9489	DPP4 monoclonal antibody, clone 6j995	PAB5181	Endomorphin-1 polyclonal antibody
H00056896-D01	DPYSL5 MaxPab rabbit polyclonal antibody (D01)	PAB5182	Endomorphin-2 polyclonal antibody
PAB8946	Drd1a polyclonal antibody	MAB9476	Endoplasmic reticulum monoclonal antibody, clone 8h46
H00055686-M01	DSU monoclonal antibody (M01), clone 8F9-1B2	PAB5183	Endorphin-alpha polyclonal antibody
H00011266-D01	DUSP12 MaxPab rabbit polyclonal antibody (D01)	PAB5184	Endorphin-gamma polyclonal antibody
H00056940-D01	DUSP22 MaxPab rabbit polyclonal antibody (D01)	H00002022-D01	ENG MaxPab rabbit polyclonal antibody (D01)
MAB2738	DUSP28 monoclonal antibody, clone AF16B7	H00002022-M01	ENG monoclonal antibody (M01), clone 4C11
MAB0835	DUSP3 monoclonal antibody, clone ARM-1	MAB3850	ENG monoclonal antibody, clone MEM-226
H00001854-D01	DUT MaxPab rabbit polyclonal antibody (D01)	MAB4501	ENG monoclonal antibody, clone MEM-226 (FITC)
H00001856-D01	DVL2 MaxPab rabbit polyclonal antibody (D01)	MAB3905	ENG monoclonal antibody, clone SN6
H00001857-D01	DVL3 MaxPab rabbit polyclonal antibody (D01)	MAB6928	ENG monoclonal antibody, clone SN6/N1-3A1
MAB6514	DYKDDDDK tag monoclonal antibody, clone B486M	H00002023-M01	ENO1 monoclonal antibody (M01), clone 8G8
MAB8183	DYKDDDDK tag monoclonal antibody, clone G191	MAB0736	ENO2 monoclonal antibody, clone 18H2
MAB7373	DYKDDDDK tag monoclonal antibody, clone L5	MAB5786	ENO2 monoclonal antibody, clone 2F1 (Biotin)
PAB12767	DYKDDDDK tag polyclonal antibody	MAB5785	ENO2 monoclonal antibody, clone 2F1 (FITC)
H00001783-D01	DYNC1L12 MaxPab rabbit polyclonal antibody (D01)	MAB5787	ENO2 monoclonal antibody, clone 2F1 (PE)
MAB3404	Dynll1/Dynll2 monoclonal antibody, clone 11F7	PAB13697	ep-1/ep-2 polyclonal antibody
PAB9862	Dyrk1a polyclonal antibody	MAB2385	EP300 monoclonal antibody, clone RW105
H00161582-D01	DYX1C1 MaxPab rabbit polyclonal antibody (D01)	MAB2386	EP300 monoclonal antibody, clone RW109
H00161582-M09	DYX1C1 monoclonal antibody (M09), clone 6G1	MAB2384	EP300 monoclonal antibody, clone RW128
PAB12765	E tag polyclonal antibody	PAB15716	Epb4.111 polyclonal antibody
MAB6761	E2F1 (phospho S364) monoclonal antibody, clone #2	H00023136-D01	EPB41L3 MaxPab rabbit polyclonal antibody (D01)
H00001876-D01	E2F6 MaxPab rabbit polyclonal antibody (D01)	H00004072-D01	EPCAM MaxPab rabbit polyclonal antibody (D01)
PAB12726	EBAG9 polyclonal antibody	MAB9471	Epcam monoclonal antibody, clone 4w944
MAB9055	E-cadherin monoclonal antibody, clone 57j5	H00054749-M01	EPDR1 monoclonal antibody (M01), clone 1C1
MAB9095	E-cadherin monoclonal antibody, clone 5c462	H00002045-D01	EPHA7 MaxPab rabbit polyclonal antibody (D01)
H00001891-D01	ECH1 MaxPab rabbit polyclonal antibody (D01)	PAB9865	EPHB1 polyclonal antibody
H00001892-M01A	ECBS1 monoclonal antibody (M01A), clone 1G9	PAB9866	EPHB1 polyclonal antibody
H00001906-M01	EDN1 monoclonal antibody (M01), clone 3D6	H00002053-D01	EPHX2 MaxPab rabbit polyclonal antibody (D01)
PAB4982	EDN1 polyclonal antibody	H00007957-D01	EPM2A MaxPab rabbit polyclonal antibody (D01)
PAB4983	EDN2 polyclonal antibody	H00002057-M01	EPOR monoclonal antibody (M01), clone 3D10
H00001908-M01	EDN3 monoclonal antibody (M01), clone 2A6-2A4	PAB12436	EPS15 polyclonal antibody
PAB4984	EDN3 polyclonal antibody	PAB12437	EPS15 polyclonal antibody
PAB15499	Ednra polyclonal antibody	MAB4694	Epstein-Barr virus early antigen (BAM HI-Z) zeta protein monoclonal antibody, clone BDI506
H00001937-M03	EEF1G monoclonal antibody (M03), clone S1	MAB4690	Epstein-Barr virus early antigen-R Homologue Bcl-2 monoclonal antibody, clone BDI502
H00030008-D01	EFEMP2 MaxPab rabbit polyclonal antibody (D01)	MAB4692	Epstein-Barr virus nuclear antigen monoclonal antibody, clone 0211
H00030008-M01	EFEMP2 monoclonal antibody (M01), clone 2C8	MAB1913	ERBB2 monoclonal antibody, clone N24
H00080303-D01	EFHD1 MaxPab rabbit polyclonal antibody (D01)	MAB7035	ERBB3 monoclonal antibody, clone SGP1
H00001944-M02	EFNA3 monoclonal antibody (M02), clone 1C12	PAB9951	ERBB3 polyclonal antibody
H00001944-M10	EFNA3 monoclonal antibody (M10), clone 2H3	MAB7170	ERBB4 monoclonal antibody, clone HFR1
H00051162-D01	EGFL7 MaxPab rabbit polyclonal antibody (D01)	PAB14786	ERBB4 polyclonal antibody
MAB1765	EGFP monoclonal antibody, clone BV-F4	MAB1877	ERCC1 monoclonal antibody, clone 8F1
PAB8931	EGFP polyclonal antibody	H00054821-D01	ERCC6L MaxPab rabbit polyclonal antibody (D01)
MAB6475	EGFR (phospho Y1173) monoclonal antibody, clone EM-13	H00001161-D01	ERCC8 MaxPab rabbit polyclonal antibody (D01)
MAB6476	EGFR (phospho Y992) monoclonal antibody, clone EM-12	H00114625-D01	ERMAP MaxPab rabbit polyclonal antibody (D01)
MAB3535	EGFR monoclonal antibody, clone 111.6	H00090952-D01	ESAM MaxPab rabbit polyclonal antibody (D01)
MAB7369	EGFR monoclonal antibody, clone 29.1	PAB9922	<i>Escherichia coli</i> Fructosyl-Amino Acid Oxidase polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB5126	<i>Escherichia coli</i> GroEL polyclonal antibody	H00002191-M01	FAP monoclonal antibody (M01), clone 1E5
PAB15567	<i>Escherichia coli</i> LexA polyclonal antibody	PAB8519	Farp2 polyclonal antibody
MAB4258	<i>Escherichia coli</i> Maltose Binding Protein monoclonal antibody, clone 289-10829	H00010667-D01	FARS2 MaxPab rabbit polyclonal antibody (D01)
PAB11397	<i>Escherichia coli</i> Maltose Phosphorylase polyclonal antibody	H00000355-D01	FAS MaxPab rabbit polyclonal antibody (D01)
PAB10310	<i>Escherichia coli</i> Penicillinase polyclonal antibody	H00000355-M01	FAS monoclonal antibody (M01), clone 5B4
PAB11398	<i>Escherichia coli</i> Pyranose Oxidase polyclonal antibody	H00000355-M11	FAS monoclonal antibody (M11), clone 1B6
PAB15568	<i>Escherichia coli</i> RecA polyclonal antibody	MAB6891	FAS monoclonal antibody, clone ANC95.1/5E2
PAB11399	<i>Escherichia coli</i> Sucrose Phosphorylase polyclonal antibody	MAB7220	FAS monoclonal antibody, clone DX2
H00002098-D01	ESD MaxPab rabbit polyclonal antibody (D01)	MAB7219	FAS monoclonal antibody, clone DX3
H00011082-D01	ESM1 MaxPab rabbit polyclonal antibody (D01)	H00000356-D01	FASLG MaxPab rabbit polyclonal antibody (D01)
PAB12621	ESR1 (phospho S167) polyclonal antibody	MAB6892	FASLG monoclonal antibody, clone Alf-2.1a
PAB12622	ESR1 (phospho S305) polyclonal antibody	MAB1756	FASLG monoclonal antibody, clone I-6060
MAB1879	ESR1 monoclonal antibody, clone AER311	MAB8220	FASLG monoclonal antibody, clone SB93a1
MAB1878	ESR1 monoclonal antibody, clone AER-314	H00002194-D01	FASN MaxPab rabbit polyclonal antibody (D01)
MAB0118	ESR1 monoclonal antibody, clone mAN1	H00002194-M01	FASN monoclonal antibody (M01), clone 3F2-1F3
MAB1880	ESR1 monoclonal antibody, clone TE111.5D11	PAB12142	Fasn polyclonal antibody
PAB12620	ESR1 polyclonal antibody	H00002192-D01	FBLN1 MaxPab rabbit polyclonal antibody (D01)
PAB12681	ESR1 polyclonal antibody	H00002192-M01	FBLN1 monoclonal antibody (M01), clone 4C9
PAB12682	ESR2 polyclonal antibody	PAB11346	FBP1 polyclonal antibody
H00002101-D01	ESRRA MaxPab rabbit polyclonal antibody (D01)	PAB5012	FBP1 polyclonal antibody
H00002104-D01	ESRRG MaxPab rabbit polyclonal antibody (D01)	H00055336-D01	FBXL8 MaxPab rabbit polyclonal antibody (D01)
MAB2409	Ethnoadenosine monoclonal antibody, clone 1G4	H00023219-D01	FBXO28 MaxPab rabbit polyclonal antibody (D01)
PAB5157	Ethynylestradiol polyclonal antibody	H00026272-D01	FBXO4 MaxPab rabbit polyclonal antibody (D01)
H00051513-D01	ETV7 MaxPab rabbit polyclonal antibody (D01)	H00026272-M02	FBXO4 monoclonal antibody (M02), clone 2F2
PAB5185	Exendin-3 polyclonal antibody	PAB10206	FBXO43 polyclonal antibody
PAB5186	Exendin-4 polyclonal antibody	H00026270-M02A	FBXO6 monoclonal antibody (M02A), clone 7D1
H00009156-D01	EXO1 MaxPab rabbit polyclonal antibody (D01)	MAB6584	Fcamr monoclonal antibody, clone TX57
H00051010-D01	EXOSC3 MaxPab rabbit polyclonal antibody (D01)	MAB6585	Fcamr monoclonal antibody, clone TX61
H00051010-M11	EXOSC3 monoclonal antibody (M11), clone 3E5	H00002208-M03	FCER2 monoclonal antibody (M03), clone S52
H00056915-D01	EXOSC5 MaxPab rabbit polyclonal antibody (D01)	H00002208-M04A	FCER2 monoclonal antibody (M04A), clone S51
PAB0648	Ezh2 polyclonal antibody	MAB7753	Fcer2a monoclonal antibody, clone 2G8
H00007430-D01	EZR MaxPab rabbit polyclonal antibody (D01)	H00002213-D01	FCGR2B MaxPab rabbit polyclonal antibody (D01)
MAB8812	EZR monoclonal antibody, clone 45j8	MAB6930	FCGR2B monoclonal antibody, clone 7.3
PAB5040	F10 polyclonal antibody	MAB0959	FCGR3A monoclonal antibody, clone LNK16
H00002160-D01	F11 MaxPab rabbit polyclonal antibody (D01)	MAB4358	FCGR3A monoclonal antibody, clone LNK16 (APC)
H00002160-M01	F11 monoclonal antibody (M01), clone 2H8	MAB4359	FCGR3A monoclonal antibody, clone LNK16 (Biotin)
H00050848-D01	F11R MaxPab rabbit polyclonal antibody (D01)	MAB4360	FCGR3A monoclonal antibody, clone LNK16 (FITC)
H00002147-D01	F2 MaxPab rabbit polyclonal antibody (D01)	MAB4361	FCGR3A monoclonal antibody, clone LNK16 (PE)
H00002147-M01	F2 monoclonal antibody (M01), clone 1B7	MAB4362	FCGR3A monoclonal antibody, clone LNK16 (PE-Dyomics 647)
PAB5051	F2 polyclonal antibody	MAB4363	FCGR3A monoclonal antibody, clone LNK16 (PerCP)
PAB5076	F2 polyclonal antibody	MAB5623	FCGR3A/FCGR3B monoclonal antibody, clone GRM1
H00002149-D01	F2R MaxPab rabbit polyclonal antibody (D01)	MAB0879	FCGR3B monoclonal antibody, clone MEM-154
H00002152-D01	F3 MaxPab rabbit polyclonal antibody (D01)	H00002219-D01	FCN1 MaxPab rabbit polyclonal antibody (D01)
H00002152-M01	F3 monoclonal antibody (M01), clone 4G4	H00002220-D01	FCN2 MaxPab rabbit polyclonal antibody (D01)
PAB5008	F5 polyclonal antibody	H00008547-M01A	FCN3 monoclonal antibody (M01A), clone 4B4
H00002155-D01	F7 MaxPab rabbit polyclonal antibody (D01)	H00002224-D01	FDPS MaxPab rabbit polyclonal antibody (D01)
H00002157-D01	F8 MaxPab rabbit polyclonal antibody (D01)	PAB15743	Fdps polyclonal antibody
H00002158-M01	F9 monoclonal antibody (M01), clone 2C9	PAB1186	FDXR polyclonal antibody
PAB5041	F9 polyclonal antibody	H00055612-D01	FERMT1 MaxPab rabbit polyclonal antibody (D01)
PAB5063	FABP polyclonal antibody	H00002243-M01	FGA monoclonal antibody (M01), clone 1D12
H00002168-M01	FABP1 monoclonal antibody (M01), clone 4G5	PAB4985	FGA polyclonal antibody
H00002167-M01	FABP4 monoclonal antibody (M01), clone 2H3-1G10	H00002246-D01	FGF1 MaxPab rabbit polyclonal antibody (D01)
H00002173-D02	FABP7 MaxPab rabbit polyclonal antibody (D02)	H00002246-M02	FGF1 monoclonal antibody (M02), clone 2E12
PAB8706	FADD polyclonal antibody	H00009965-D01	FGF19 MaxPab rabbit polyclonal antibody (D01)
PAB8707	FAF1 polyclonal antibody	H00026291-M01	FGF21 monoclonal antibody (M01), clone 2F11
PAB12766	Fah polyclonal antibody	H00026291-M02	FGF21 monoclonal antibody (M02), clone 1A8
H00054097-M07	FAM3B monoclonal antibody (M07), clone 1E7	H00026291-M03	FGF21 monoclonal antibody (M03), clone 3G10
MAB7338	FANCA monoclonal antibody, clone 5G9	H00002250-M01	FGF5 monoclonal antibody (M01), clone 1B4
H00002176-D01	FANCC MaxPab rabbit polyclonal antibody (D01)	MAB4707	FGFR1 monoclonal antibody, clone M5G10
MAB7315	FANCC monoclonal antibody, clone 8F3	H00002267-M01	FGL1 monoclonal antibody (M01), clone 2A4
MAB2398	FANCD2 monoclonal antibody, clone FI-17	H00010875-M01	FGL2 monoclonal antibody (M01), clone 6D9
PAB12428	FANCD2 polyclonal antibody	H00002268-D01	FGR MaxPab rabbit polyclonal antibody (D01)
H00002189-D01	FANCG MaxPab rabbit polyclonal antibody (D01)	H00009457-D01	FHL5 MaxPab rabbit polyclonal antibody (D01)
H00002189-M01	FANCG monoclonal antibody (M01), clone 2C8	PAB14913	Fibrinogen fragment D polyclonal antibody
H00055120-D01	FANCL MaxPab rabbit polyclonal antibody (D01)	PAB14914	Fibrinogen fragment E polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB1212	Fibrinogen monoclonal antibody	H00007855-M01	FZD5 monoclonal antibody (M01), clone 6A3
PAB5134	Fibrinogen polyclonal antibody	H00050486-D01	G0S2 MaxPab rabbit polyclonal antibody (D01)
PAB5139	Fibrinogen polyclonal antibody	H00010146-D01	G3BP1 MaxPab rabbit polyclonal antibody (D01)
PAB5140	Fibrinogen polyclonal antibody	MAB1139	GAB1 monoclonal antibody, clone 1AT1979
H00002277-D01	FIGF MaxPab rabbit polyclonal antibody (D01)	H00002553-D01	GABPB1 MaxPab rabbit polyclonal antibody (D01)
PAB8718	FIS1 polyclonal antibody	PAB9679	Gabrd polyclonal antibody
H00051303-D01	FKBP11 MaxPab rabbit polyclonal antibody (D01)	MAB9472	GAD monoclonal antibody, clone 3s65d
PAB8720	Fkbp1a polyclonal antibody	H00002572-D01	GAD2 MaxPab rabbit polyclonal antibody (D01)
H00002288-D01	FKBP4 MaxPab rabbit polyclonal antibody (D01)	PAB10212	Gadd45g polyclonal antibody
H00002288-M01	FKBP4 monoclonal antibody (M01), clone 5C11	H00002580-D01	GAK MaxPab rabbit polyclonal antibody (D01)
MAB9464	fl(2)d monoclonal antibody, clone 64w95	PAB5055	GAL polyclonal antibody
H00055011-M01	FLJ20643 monoclonal antibody (M01), clone 4F5-1E10	PAB5088	GAL polyclonal antibody
PAB13531	FLT1 polyclonal antibody	MAB1883	GAL4 monoclonal antibody, clone 5C8
MAB8212	FLT3 monoclonal antibody, clone BV10A4	MAB1884	GAL4 monoclonal antibody, clone 8C-1
H00002328-D01	FMO3 MaxPab rabbit polyclonal antibody (D01)	PAB5187	Galantide polyclonal antibody
MAB9086	Fmr1 monoclonal antibody, clone 46a6	H00002582-D01	GALE MaxPab rabbit polyclonal antibody (D01)
MAB2440	Fmr1 monoclonal antibody, clone 6A15	H00002588-D01	GALNS MaxPab rabbit polyclonal antibody (D01)
MAB0770	FN1 monoclonal antibody, clone 16E5	H00011226-D01	GALNT6 MaxPab rabbit polyclonal antibody (D01)
MAB0762	FN1 monoclonal antibody, clone 23C9	H00050614-D01	GALNT9 MaxPab rabbit polyclonal antibody (D01)
MAB9134	FN1 monoclonal antibody, clone 976h3	H00002592-D01	GALT MaxPab rabbit polyclonal antibody (D01)
MAB7201	FN1 monoclonal antibody, clone Fn-3	H00023193-D01	GANAB MaxPab rabbit polyclonal antibody (D01)
MAB7200	FN1 monoclonal antibody, clone Fn-4	H00002597-M01	GAPDH monoclonal antibody (M01), clone 3C2
PAB10210	FN1 polyclonal antibody	MAB2549	GAPDH monoclonal antibody, clone 6C5
PAB5071	Fn1 polyclonal antibody	MAB0694	GAPDH monoclonal antibody, clone 7B
PAB5085	Fn1 polyclonal antibody	MAB0687	GAPDH monoclonal antibody, clone 8C2
H00064122-D01	FN3K MaxPab rabbit polyclonal antibody (D01)	PAB1195	GAPDH polyclonal antibody
H00064838-D01	FNDC4 MaxPab rabbit polyclonal antibody (D01)	H00026330-D01	GAPDHS MaxPab rabbit polyclonal antibody (D01)
PAB12367	FNIP1 polyclonal antibody	H00002617-D01	GARS MaxPab rabbit polyclonal antibody (D01)
MAB3611	FOLH1 monoclonal antibody, clone GCP-05	H00002618-D01	GART MaxPab rabbit polyclonal antibody (D01)
H00002354-D01	FOSB MaxPab rabbit polyclonal antibody (D01)	H00002618-M01	GART monoclonal antibody (M01), clone 4D6-1D5
H00002355-D01	FOSL2 MaxPab rabbit polyclonal antibody (D01)	H00283431-D01	GAS2L3 MaxPab rabbit polyclonal antibody (D01)
H00003170-M01	FOXA2 monoclonal antibody (M01), clone 7E6	H00002621-D01	GAS6 MaxPab rabbit polyclonal antibody (D01)
H00003170-M10	FOXA2 monoclonal antibody (M10), clone 1C7	H00002623-D01	GATA1 MaxPab rabbit polyclonal antibody (D01)
H00003171-D01	FOXA3 MaxPab rabbit polyclonal antibody (D01)	H00002629-D01	GBA MaxPab rabbit polyclonal antibody (D01)
PAB12683	FOXP1 polyclonal antibody	H00002632-D01	GBE1 MaxPab rabbit polyclonal antibody (D01)
PAB12684	FOXP2 polyclonal antibody	H00064223-D01	GBL MaxPab rabbit polyclonal antibody (D01)
H00050943-D01	FOXP3 MaxPab rabbit polyclonal antibody (D01)	H00002633-D01	GBP1 MaxPab rabbit polyclonal antibody (D01)
PAB12141	Foxp3 polyclonal antibody	H00002634-D01	GBP2 MaxPab rabbit polyclonal antibody (D01)
PAB12685	FOXP3 polyclonal antibody	H00002634-M03	GBP2 monoclonal antibody (M03), clone 2A10
PAB8952	Foxp3 polyclonal antibody	PAB5066	GC polyclonal antibody
PAB12686	FOXP4 polyclonal antibody	PAB15688	Gcc2 polyclonal antibody
H00002358-D01	FPR2 MaxPab rabbit polyclonal antibody (D01)	PAB4991	GCG polyclonal antibody
H00002483-M01	FRG1 monoclonal antibody (M01), clone 4A5	PAB5056	GCG polyclonal antibody
H00002483-M02	FRG1 monoclonal antibody (M02), clone 2B9	PAB5057	GCG polyclonal antibody
H00010818-D01	FRS2 MaxPab rabbit polyclonal antibody (D01)	PAB5077	Gcg polyclonal antibody
MAB7754	FSBA monoclonal antibody, clone 2F5/5	H00002643-M01	GCH1 monoclonal antibody (M01), clone 4A12
H00079187-D01	FSD1 MaxPab rabbit polyclonal antibody (D01)	H00002730-D01	GCLM MaxPab rabbit polyclonal antibody (D01)
H00079187-M01	FSD1 monoclonal antibody (M01), clone 1G3	PAB8740	Gdnf polyclonal antibody
H00002512-M16	FTL monoclonal antibody (M16), clone X1	H00079760-M09	GEMIN7 monoclonal antibody (M09), clone 4H6
H00002512-M11	FTL monoclonal antibody, clone 4C10	H00002670-D01	GFAP MaxPab rabbit polyclonal antibody (D01)
H00002512-M08	FTL monoclonal antibody, clone 4G10	MAB8070	GFAP monoclonal antibody, clone 0.N.308
PAB8732	FUS polyclonal antibody	MAB5132	GFAP monoclonal antibody, clone GA-5
H00002523-D01	FUT1 MaxPab rabbit polyclonal antibody (D01)	MAB3612	GFAP monoclonal antibody, clone GF-01
H00002524-D01	FUT2 MaxPab rabbit polyclonal antibody (D01)	MAB8184	GFP monoclonal antibody
MAB5622	FUT4 monoclonal antibody, clone 28	MAB6746	GFP monoclonal antibody, clone 1A5
MAB6048	FUT4 monoclonal antibody, clone 28 (FITC)	PAB1264	GFP polyclonal antibody
MAB6049	FUT4 monoclonal antibody, clone 28 (PE)	PAB12768	GFP polyclonal antibody
MAB6931	FUT4 monoclonal antibody, clone AHN1.1	PAB13744	GFP polyclonal antibody
MAB0878	FUT4 monoclonal antibody, clone MEM-158	MAB0849	GFP tag monoclonal antibody, clone 18A11
H00002395-D01	FXN MaxPab rabbit polyclonal antibody (D01)	MAB8758	GFP tag monoclonal antibody, clone 946as
H00002395-M01	FXN monoclonal antibody (M01), clone 1D9	PAB1253	GFP tag polyclonal antibody
MAB2493	FXR2 monoclonal antibody, clone A42	PAB10213	GGA1 polyclonal antibody
MAB1882	Fyn monoclonal antibody, clone 1S	PAB10050	GGA3 polyclonal antibody
MAB0947	FYN monoclonal antibody, clone FYN-01	H00009453-D01	GGPS1 MaxPab rabbit polyclonal antibody (D01)
MAB7969	FYN monoclonal antibody, clone FYN-59	H00002678-M01	GGT1 monoclonal antibody (M01), clone 1F9

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00002688-D01	GH1 MaxPab rabbit polyclonal antibody (D01)	MAB6237	Grin1 monoclonal antibody, clone R1JHL
PAB5099	GH1 polyclonal antibody	PAB9699	Grin2a polyclonal antibody
PAB4992	GHRH polyclonal antibody	MAB8566	Grin2b monoclonal antibody, clone S59-20
PAB5068	Ghrh polyclonal antibody	MAB7899	Grin2b monoclonal antibody, clone S59-36
PAB5078	Ghrh polyclonal antibody	PAB9700	Grin2b polyclonal antibody
PAB4987	GHRL polyclonal antibody	H00002905-D01	GRIN2C MaxPab rabbit polyclonal antibody (D01)
PAB4988	GHRL polyclonal antibody	PAB9705	Grin2c polyclonal antibody
PAB5059	GHRL polyclonal antibody	MAB8285	GRIN3A/GRIN3B monoclonal antibody
H00002694-M03	GIF monoclonal antibody (M03), clone 1D9	H00002870-D01	GRK6 MaxPab rabbit polyclonal antibody (D01)
H00055340-D01	GIMAP5 MaxPab rabbit polyclonal antibody (D01)	PAB4986	GRP polyclonal antibody
H00055340-M10	GIMAP5 monoclonal antibody (M10), clone 1E10	PAB4996	GRP polyclonal antibody
H00055340-M14A	GIMAP5 monoclonal antibody (M14A), clone 3F11	MAB6629	Grp75 monoclonal antibody, clone S52A-42
H00009837-D01	GINS1 MaxPab rabbit polyclonal antibody (D01)	H00145258-D01	GSC MaxPab rabbit polyclonal antibody (D01)
PAB4990	GIP polyclonal antibody	H00145258-M01	GSC monoclonal antibody (M01), clone 4H7
H00009815-M03	GIT2 monoclonal antibody (M03), clone M1	H00145258-M18	GSC monoclonal antibody (M18), clone 4D6
PAB12759	GJA1 polyclonal antibody	PAB12623	GSK3A (phospho S21) polyclonal antibody
PAB9916	GK polyclonal antibody	MAB5634	GSK3A monoclonal antibody, clone SBGSK3a4
H00002717-D01	GLA MaxPab rabbit polyclonal antibody (D01)	MAB6088	GSK3A monoclonal antibody, clone SBGSK3a4 (SEPH)
MAB9107	Glass Drosophila protein monoclonal antibody, clone 4162h	PAB12687	GSK3A polyclonal antibody
H00002720-D01	GLB1 MaxPab rabbit polyclonal antibody (D01)	PAB17915	GSK3A polyclonal antibody
H00002720-M01	GLB1 monoclonal antibody (M01), clone 6E7	PAB12625	GSK3B (phospho S9) polyclonal antibody
H00026035-D01	GLCE MaxPab rabbit polyclonal antibody (D01)	PAB12624	GSK3B (phospho Y216) polyclonal antibody
H00011146-D01	GLMN MaxPab rabbit polyclonal antibody (D01)	PAB12688	GSK3B polyclonal antibody
H00011146-M01	GLMN monoclonal antibody (M01), clone 1C12	PAB12689	GSK3B polyclonal antibody
H00002739-D01	GLO1 MaxPab rabbit polyclonal antibody (D01)	PAB12690	GSK3B polyclonal antibody
MAB0717	GLRX monoclonal antibody, clone 28C3	MAB4712	GSK3-beta monoclonal antibody, clone 21A
MAB0726	GLRX monoclonal antibody, clone 30A1	H00002934-M01	GSN monoclonal antibody (M01), clone 3G5
H00002746-M01	GLUD1 monoclonal antibody (M01), clone 3C2	MAB0786	GSN monoclonal antibody, clone 35B2
PAB9914	GLUD1 polyclonal antibody	H00002936-D01	GSR MaxPab rabbit polyclonal antibody (D01)
MAB0735	GLUL monoclonal antibody, clone 27D10	MAB0693	GSR monoclonal antibody, clone 1A7
MAB0734	GLUL monoclonal antibody, clone 8G9	MAB0685	GSR monoclonal antibody, clone 2B3
PAB11391	Glutamine Synthetase polyclonal antibody	MAB8236	GST monoclonal antibody, clone 1B3
MAB6676	Glutathione monoclonal antibody, clone L4H	MAB3650	GST monoclonal antibody, clone S-tag-02
MAB2343	GMNN monoclonal antibody, clone 1B9	MAB3651	GST monoclonal antibody, clone S-tag-05
H00027232-D01	GNMT MaxPab rabbit polyclonal antibody (D01)	PAB10346	GST polyclonal antibody
H00008443-D01	GNPAT MaxPab rabbit polyclonal antibody (D01)	PAB15371	GST polyclonal antibody
H00051280-M04	GOLPH2 monoclonal antibody (M04), clone 3B10	PAB15372	GST polyclonal antibody
H00009527-D01	GOSR1 MaxPab rabbit polyclonal antibody (D01)	PAB17015	GST polyclonal antibody
MAB7171	Gosr1 monoclonal antibody, clone HFD9-5	MAB0851	GST tag monoclonal antibody, clone LF4G2
H00002805-D01	GOT1 MaxPab rabbit polyclonal antibody (D01)	H00002938-D01	GSTA1 MaxPab rabbit polyclonal antibody (D01)
H00002811-M02	GP1BA monoclonal antibody (M02), clone 1C6	H00373156-D01	GSTK1 MaxPab rabbit polyclonal antibody (D01)
MAB4144	GP1BA monoclonal antibody, clone SZ2	H00002948-D01	GSTM4 MaxPab rabbit polyclonal antibody (D01)
MAB5454	GP9 monoclonal antibody, clone SZ1 (FITC)	H00009446-M01	GSTO1 monoclonal antibody (M01), clone 3D9
H00002262-D01	GPC5 MaxPab rabbit polyclonal antibody (D01)	H00119391-D01	GSTO2 MaxPab rabbit polyclonal antibody (D01)
H00002262-M01	GPC5 monoclonal antibody (M01), clone 1C9	H00002950-M01	GSTP1 monoclonal antibody (M01), clone 2G6-F6
H00002262-M02	GPC5 monoclonal antibody (M02), clone 1B8	H00002953-M01	GSTT2 monoclonal antibody (M01), clone 1C12
H00023171-D01	GPD1L MaxPab rabbit polyclonal antibody (D01)	H00002954-D01	GSTZ1 MaxPab rabbit polyclonal antibody (D01)
PAB15686	Gpd1l polyclonal antibody	H00002960-D01	GTF2E1 MaxPab rabbit polyclonal antibody (D01)
H00002821-D01	GPI MaxPab rabbit polyclonal antibody (D01)	H00002965-D01	GTF2H1 MaxPab rabbit polyclonal antibody (D01)
H00011321-D01	GPN1 MaxPab rabbit polyclonal antibody (D01)	H00002966-D01	GTF2H2 MaxPab rabbit polyclonal antibody (D01)
H00131601-D01	GPR175 MaxPab rabbit polyclonal antibody (D01)	H00002967-D01	GTF2H3 MaxPab rabbit polyclonal antibody (D01)
H00084706-D01	GPT2 MaxPab rabbit polyclonal antibody (D01)	H00002968-D01	GTF2H4 MaxPab rabbit polyclonal antibody (D01)
MAB0845	Gpx1 monoclonal antibody, clone 13B2AF	H00002969-D01	GTF2I MaxPab rabbit polyclonal antibody (D01)
MAB0730	GPX1 monoclonal antibody, clone 2A10	H00009569-D01	GTF2IRD1 MaxPab rabbit polyclonal antibody (D01)
MAB0731	GPX1 monoclonal antibody, clone 42C9	H00051512-D01	GTSE1 MaxPab rabbit polyclonal antibody (D01)
PAB8035	GPX1 polyclonal antibody	H00009626-D01	GUCA1C MaxPab rabbit polyclonal antibody (D01)
MAB0765	GPX3 monoclonal antibody, clone 55A	H00009626-M03	GUCA1C monoclonal antibody (M03), clone 3F10
MAB0729	GPX4 monoclonal antibody, clone 1H11	PAB4993	GUCA2A polyclonal antibody
MAB0707	GPX4 monoclonal antibody, clone 7A4	H00002982-D01	GUCY1A3 MaxPab rabbit polyclonal antibody (D01)
H00002880-M02	GPX5 monoclonal antibody (M02), clone 3B9	H00002987-D01	GUK1 MaxPab rabbit polyclonal antibody (D01)
H00009402-D01	GRAP2 MaxPab rabbit polyclonal antibody (D01)	H00002987-M01	GUK1 monoclonal antibody (M01), clone 4C3-1A7
PAB11330	Grasp polyclonal antibody	H00002987-M02A	GUK1 monoclonal antibody (M02A), clone 1H3
H00002886-D01	GRB7 MaxPab rabbit polyclonal antibody (D01)	H00003002-D01	GZMB MaxPab rabbit polyclonal antibody (D01)
H00002886-M03	GRB7 monoclonal antibody (M03), clone 3C12	H00003002-M01	GZMB monoclonal antibody (M01), clone 4F5
H00064388-D01	GREM2 MaxPab rabbit polyclonal antibody (D01)	MAB1329	H2-Ab1 monoclonal antibody, clone IA-B2

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB1330	H2-Ab1 monoclonal antibody, clone IA-C6	H00055869-D01	HDAC8 MaxPab rabbit polyclonal antibody (D01)
MAB1332	H2-Ab1 monoclonal antibody, clone IA-E6	H00055869-M07	HDAC8 monoclonal antibody (M07), clone 2F4
MAB1331	H2-Ab1 monoclonal antibody, clone NR-D11	PAB8754	Hdac8 polyclonal antibody
PAB8764	H2AFX polyclonal antibody	PAB2338	HDAC9 polyclonal antibody
PAB15497	H2B polyclonal antibody	PAB2339	HDAC9 polyclonal antibody
MAB1359	H2-D1 monoclonal antibody, clone KD-B4	PAB8755	HDAC9 polyclonal antibody
MAB5821	H2-D4 monoclonal antibody, clone 34-5-8S (Biotin)	PAB5014	HDL polyclonal antibody
MAB5820	H2-D4 monoclonal antibody, clone 34-5-8S (FITC)	MAB1014	Heat Shock Protein 70 monoclonal antibody, clone 16A12
MAB5822	H2-D4 monoclonal antibody, clone 34-5-8S (PE)	MAB5458	Hepatitis B virus core antigen (ayw) monoclonal antibody, clone 13A9 (Biotin)
PAB14096	H3 polyclonal antibody	MAB5403	Hepatitis B virus core antigen monoclonal antibody, clone 10C6
PAB14115	H3K27me1 polyclonal antibody	MAB5402	Hepatitis B virus core antigen monoclonal antibody, clone 10-E-11
PAB14116	H3K27me2 polyclonal antibody	MAB4774	Hepatitis B virus core antigen monoclonal antibody, clone 13A9
MAB4947	H3K27me2/3 monoclonal antibody	MAB5404	Hepatitis B virus core antigen monoclonal antibody, clone 14-E-11
PAB14051	H3K27me3 polyclonal antibody	MAB4775	Hepatitis B virus core Antigen monoclonal antibody, clone C1
PAB14058	H3K27me3S28p polyclonal antibody	MAB5401	Hepatitis B virus core antigen monoclonal antibody, clone I-5
PAB14132	H3K36me1 polyclonal antibody	MAB0463	Hepatitis B virus monoclonal antibody, clone C 1
PAB14090	H3K36me2 polyclonal antibody	MAB5399	Hepatitis B virus pre-surface antigen monoclonal antibody, clone S 26
PAB14112	H3K4me1 polyclonal antibody	MAB4776	Hepatitis B virus pre-surface antigen monoclonal antibody, clone S26
PAB14111	H3K4me2 polyclonal antibody	PAB13967	Hepatitis B virus surface antigen (ad/ay) polyclonal antibody
PAB14056	H3K79me1 polyclonal antibody	MAB5400	Hepatitis B virus surface antigen surface monoclonal antibody, clone S 38
PAB14047	H3K79me2 polyclonal antibody	MAB4787	Hepatitis C virus core antigen monoclonal antibody, clone 6A1
PAB14122	H3K79me3 polyclonal antibody	MAB5512	Hepatitis C virus core antigen monoclonal antibody, clone 6A1 (FITC)
PAB14069	H3K9acS10p polyclonal antibody	MAB4791	Hepatitis C virus core antigen monoclonal antibody, clone 7A1
PAB14049	H3K9me1 polyclonal antibody	MAB5513	Hepatitis C virus core antigen monoclonal antibody, clone 7A1 (Biotin)
PAB14099	H3K9me3 polyclonal antibody	MAB5514	Hepatitis C virus core antigen monoclonal antibody, clone 7A1 (FITC)
PAB14091	H3K9me3S10p polyclonal antibody	MAB5393	Hepatitis C virus core monoclonal antibody, clone 1F6
PAB14059	H3R17me2 polyclonal antibody	MAB0487	Hepatitis C virus monoclonal antibody, clone 1-E-5
PAB14082	H3S10p polyclonal antibody	MAB0490	Hepatitis C virus monoclonal antibody, clone 2-H1
PAB14070	H4K8ac polyclonal antibody	MAB0489	Hepatitis C virus monoclonal antibody, clone 4-F5
MAB8504	HA tag monoclonal antibody	MAB0491	Hepatitis C virus monoclonal antibody, clone 7-D4
MAB8505	HA tag monoclonal antibody, clone 2C16 (Biotin)	MAB0488	Hepatitis C virus monoclonal antibody, clone 9-G2
MAB7911	HA tag monoclonal antibody, clone LF7H5	MAB4797	Hepatitis C virus NS-4 antigen monoclonal antibody, clone 8A1
H00003026-M01	HABP2 monoclonal antibody (M01), clone 1H4	MAB5394	Hepatitis C virus NS4 monoclonal antibody, clone 1G7
H00003033-D01	HADH MaxPab rabbit polyclonal antibody (D01)	MAB5515	Hepatitis C virus NS4 monoclonal antibody, clone 8A1 (Biotin)
H00003033-M01	HADHSC monoclonal antibody (M01), clone 4B5	MAB5516	Hepatitis C virus NS4 monoclonal antibody, clone 8A1 (FITC)
MAB9494	Half pint protein monoclonal antibody, clone 4h95j	PAB14171	Hepatitis C virus NS5b polyclonal antibody
H00051179-D01	HAO2 MaxPab rabbit polyclonal antibody (D01)	MAB7798	Herpes simplex virus type 1 ICP8 Major DNA binding protein monoclonal antibody, clone 11E2
MAB2312	Hap1 monoclonal antibody, clone 1B6	MAB7749	Herpes simplex virus type 1 UL42 monoclonal antibody, clone 2H4
H00003035-D01	HARS MaxPab rabbit polyclonal antibody (D01)	MAB4801	Herpes simplex virus type 2 Glycoprotein B monoclonal antibody, clone 1-G-10
MAB1912	HA-tag monoclonal antibody, clone 12CA5 (Agarose)	H00009709-D01	HERPUD1 MaxPab rabbit polyclonal antibody (D01)
H00026762-M03	HAVCR1 monoclonal antibody (M03), clone 2G11	H00009709-M04	HERPUD1 monoclonal antibody (M04), clone 2G7
MAB6756	HBEGF monoclonal antibody, clone 4G10	MAB0793	HES1 monoclonal antibody, clone 7H11
H00010767-D01	HBS1L MaxPab rabbit polyclonal antibody (D01)	H00003073-D01	HEXA MaxPab rabbit polyclonal antibody (D01)
H00003052-M01	HCCS monoclonal antibody (M01), clone 3C7	H00003081-M09	HGD monoclonal antibody (M09), clone 2C10
H00003055-D01	HCK MaxPab rabbit polyclonal antibody (D01)	H00003081-M11	HGD monoclonal antibody (M11), clone 1F1
H00003055-M01	HCK monoclonal antibody (M01), clone 1D9-1A9	H00064399-M01	HHIP monoclonal antibody (M01), clone 5D11
H00003059-D01	HCLS1 MaxPab rabbit polyclonal antibody (D01)	H00003091-D01	HIF1A MaxPab rabbit polyclonal antibody (D01)
MAB6651	Hcn1 monoclonal antibody, clone S70-28	MAB2342	HIF1A monoclonal antibody, clone ESEE122
MAB6652	Hcn2 monoclonal antibody, clone S71-37	MAB1892	HIF1A monoclonal antibody, clone H1alpha67
H00003062-D01	HCRTR2 MaxPab rabbit polyclonal antibody (D01)	MAB1893	HIF1A monoclonal antibody, clone OZ12
H00003065-D01	HDAC1 MaxPab rabbit polyclonal antibody (D01)	PAB12418	HIF1A polyclonal antibody
PAB0647	HDAC1 polyclonal antibody	PAB12417	HIF1AN polyclonal antibody
PAB8747	HDAC1 polyclonal antibody	PAB14562	HIF3A polyclonal antibody
H00079885-D01	HDAC11 MaxPab rabbit polyclonal antibody (D01)	H00003092-M01	HIP1 monoclonal antibody (M01), clone 1F12
PAB0001	HDAC1 polyclonal antibody	H00147746-D01	HIPK4 MaxPab rabbit polyclonal antibody (D01)
PAB0002	HDAC2 polyclonal antibody	MAB0852	His tag monoclonal antibody, clone 6A7
H00008841-D01	HDAC3 MaxPab rabbit polyclonal antibody (D01)	MAB8503	His tag monoclonal antibody, clone BV-G020
H00008841-M01	HDAC3 monoclonal antibody (M01), clone 2F9-4B7	MAB3417	His tag monoclonal antibody, clone HIS.H8
PAB0171	HDAC3 polyclonal antibody	MAB2785	His tag monoclonal antibody, clone His17
PAB0181	HDAC4 polyclonal antibody	MAB1274	His tag monoclonal antibody, clone HT84
PAB0182	HDAC5 polyclonal antibody	PAB17049	His tag polyclonal antibody
PAB8752	HDAC5 polyclonal antibody	PAB18029	His tag polyclonal antibody
PAB0183	HDAC6 polyclonal antibody	PAB8929	His tag polyclonal antibody
H00051564-D01	HDAC7 MaxPab rabbit polyclonal antibody (D01)	PAB15498	HIST1 polyclonal antibody
PAB8928	Hdac7 polyclonal antibody	PAB12679	HIST2H3C polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB12691	HIST2H3C polyclonal antibody	MAB3855	HLA-G monoclonal antibody, clone MEM-G/9
PAB8766	HIST3H2BB polyclonal antibody	H00003141-D01	HLCS MaxPab rabbit polyclonal antibody (D01)
PAB15451	Histon H3 polyclonal antibody	PAB12415	HLTF polyclonal antibody
MAB5685	Histone H2B (ubiquitinated) monoclonal antibody, clone NRO3	H00003142-D01	HLX MaxPab rabbit polyclonal antibody (D01)
PAB0656	Histone H3 (K36me3) polyclonal antibody	H00003142-M05	HLX1 monoclonal antibody (M05), clone 1B9
PAB15453	Histone H3 (K4me2) polyclonal antibody	H00010362-D01	HMG20B MaxPab rabbit polyclonal antibody (D01)
PAB0650	Histone H3 (K4me3) polyclonal antibody	H00003146-D01	HMG20B MaxPab rabbit polyclonal antibody (D01)
PAB15454	Histone H3 (K4me3) polyclonal antibody	H00003155-D01	HMGCL MaxPab rabbit polyclonal antibody (D01)
PAB0652	Histone H3 (K9/14ac) polyclonal antibody	H00003156-D01	HMGCR MaxPab rabbit polyclonal antibody (D01)
PAB0651	Histone H3 (K9ac) polyclonal antibody	MAB6628	HMOX1 monoclonal antibody, clone 1F12-A6
PAB0653	Histone H3 (K9me2) polyclonal antibody	PAB8760	HMOX1 polyclonal antibody
PAB0654	Histone H3 (K9me3) polyclonal antibody	PAB8761	HMOX2 polyclonal antibody
PAB0655	Histone H3 (K9me3) polyclonal antibody	H00003172-D01	HNF4A MaxPab rabbit polyclonal antibody (D01)
PAB0657	Histone H4 (K20me3) polyclonal antibody	MAB2483	HNRNPA1 monoclonal antibody, clone 4B10
PAB8768	Histone H4 polyclonal antibody	MAB2492	HNRNPA1 monoclonal antibody, clone 9H10
PAB15496	Histone Nucleosome Core polyclonal antibody	MAB2497	HNRNPA2B1 monoclonal antibody, clone DP3B3
H00003098-M01	HK1 monoclonal antibody (M01), clone 5G9	MAB2486	HNRNPC monoclonal antibody, clone 4F4
MAB3854	HLA Class I monoclonal antibody, clone MEM-123	MAB2478	HNRNPK monoclonal antibody, clone 3C2
MAB4206	HLA Class II (DQ), Monomorphic monoclonal antibody, clone SPV-L3	MAB2485	HNRNPL monoclonal antibody, clone 4D11
H00003105-D01	HLA-A MaxPab rabbit polyclonal antibody (D01)	MAB2339	HNRNPM monoclonal antibody, clone 2A6
MAB6449	HLA-A monoclonal antibody, clone W6/32	MAB1758	Holocytochrome c monoclonal antibody, clone 2CYZ-702
MAB6434	HLA-A2 monoclonal antibody, clone BB7.2	H00009456-D01	HOMER1 MaxPab rabbit polyclonal antibody (D01)
H00003107-D01	HLA-C MaxPab rabbit polyclonal antibody (D01)	H00009455-D01	HOMER2 MaxPab rabbit polyclonal antibody (D01)
MAB0903	HLA-Class I monoclonal antibody, clone MEM-147	H00003198-D01	HOXA1 MaxPab rabbit polyclonal antibody (D01)
H00003108-D01	HLA-DMA MaxPab rabbit polyclonal antibody (D01)	H00003211-D01	HOXB1 MaxPab rabbit polyclonal antibody (D01)
H00003109-M06	HLA-DMB monoclonal antibody (M06), clone 5G11	H00003211-M06	HOXB1 monoclonal antibody (M06), clone 2F12
H00003109-M07	HLA-DMB monoclonal antibody (M07), clone 7D11	H00003211-M13	HOXB1 monoclonal antibody (M13), clone 2E5
H00003109-M08	HLA-DMB monoclonal antibody (M08), clone 4D12	H00003219-D03	HOXB9 MaxPab rabbit polyclonal antibody (D03)
H00003109-M09	HLA-DMB monoclonal antibody (M09), clone 6E10	H00003219-M03	HOXB9 monoclonal antibody (M03), clone 4C3
H00003111-D01	HLA-DOA MaxPab rabbit polyclonal antibody (D01)	H00003219-M05	HOXB9 monoclonal antibody (M05), clone 4C11
H00003112-D01	HLA-DOB MaxPab rabbit polyclonal antibody (D01)	H00003219-M29	HOXB9 monoclonal antibody (M29), clone 2E8
MAB7221	HLA-DPA1 monoclonal antibody, clone DP 11.1	H00003224-M02	HOXC8 monoclonal antibody (M02), clone 1H2
MAB3615	HLA-DQ1/DQ3 monoclonal antibody, clone HL-37	H00003231-D01	HOXD1 MaxPab rabbit polyclonal antibody (D01)
MAB7748	HLA-DQ3 monoclonal antibody, clone 2HB6	H00003231-M01	HOXD1 monoclonal antibody (M01), clone 4F4
H00003119-D01	HLA-DQB1 MaxPab rabbit polyclonal antibody (D01)	H00003234-D01	HOXD8 MaxPab rabbit polyclonal antibody (D01)
MAB5237	HLA-DR monoclonal antibody, clone HL39	H00003248-M01	HPGD monoclonal antibody (M01), clone 1D8
MAB0902	HLA-DR monoclonal antibody, clone HL-39	H00003263-M01	HPX monoclonal antibody (M01), clone 3A9-1A9
MAB0978	HLA-DR monoclonal antibody, clone MEM-12	H00003273-D01	HRG MaxPab rabbit polyclonal antibody (D01)
MAB4477	HLA-DR monoclonal antibody, clone MEM-12 (APC)	H00009653-D01	HS2ST1 MaxPab rabbit polyclonal antibody (D01)
MAB4478	HLA-DR monoclonal antibody, clone MEM-12 (FITC)	MAB2191	Hsc70 monoclonal antibody, clone 1B5
MAB4479	HLA-DR monoclonal antibody, clone MEM-12 (PE)	H00003291-D01	HSD11B2 MaxPab rabbit polyclonal antibody (D01)
MAB4481	HLA-DR monoclonal antibody, clone MEM-12 (PE-Dyomics 647)	H00003294-D01	HSD17B2 MaxPab rabbit polyclonal antibody (D01)
MAB4482	HLA-DR monoclonal antibody, clone MEM-12 (PerCP)	H00051478-D01	HSD17B7 MaxPab rabbit polyclonal antibody (D01)
MAB0900	HLA-DR/HLA-DP monoclonal antibody, clone HL-38	H00007923-D01	HSD17B8 MaxPab rabbit polyclonal antibody (D01)
MAB0901	HLA-DR/HLA-DP monoclonal antibody, clone HL-40	H00003283-D01	HSD3B1 MaxPab rabbit polyclonal antibody (D01)
MAB0917	HLA-DR/HLA-DP monoclonal antibody, clone MEM-136	MAB2285	Hsf1 monoclonal antibody, clone 10H8
MAB4483	HLA-DR/HLA-DP monoclonal antibody, clone MEM-136 (APC)	MAB3736	Hsf1 monoclonal antibody, clone Cocktail
MAB4484	HLA-DR/HLA-DP monoclonal antibody, clone MEM-136 (FITC)	MAB2199	Hsf2 monoclonal antibody, clone 3E2
MAB4485	HLA-DR/HLA-DP monoclonal antibody, clone MEM-136 (PE)	PAB5125	HSP 104 polyclonal antibody
H00003122-D01	HLA-DRA MaxPab rabbit polyclonal antibody (D01)	MAB2968	HSP25 monoclonal antibody, clone HSP25-19
MAB7134	HLA-DRA monoclonal antibody, clone L243	MAB2969	HSP25 monoclonal antibody, clone HSP25-26
MAB7128	HLA-DRA monoclonal antibody, clone LHM 4	MAB2970	HSP25 monoclonal antibody, clone HSP25-31
MAB7125	HLA-DRA monoclonal antibody, clone LN3	PAB13575	HSP70 polyclonal antibody
H00003127-D01	HLA-DRB5 MaxPab rabbit polyclonal antibody (D01)	PAB13580	HSP70 polyclonal antibody
H00003133-D01	HLA-E MaxPab rabbit polyclonal antibody (D01)	MAB6620	HSP90 monoclonal antibody, clone 8D3
MAB3630	HLA-E monoclonal antibody, clone MEM-E/06	MAB8494	Hsp90 monoclonal antibody, clone S88
MAB3631	HLA-E monoclonal antibody, clone MEM-E/07	MAB2186	HSP90AA1 monoclonal antibody, clone 2G5.G3
MAB3632	HLA-E monoclonal antibody, clone MEM-E/08	MAB2187	HSP90AA1 monoclonal antibody, clone 4F3.E8
MAB6463	HLA-G monoclonal antibody, clone 01G	MAB2197	HSP90AA1 monoclonal antibody, clone D7alpha
MAB3858	HLA-G monoclonal antibody, clone 4H84	PAB13569	Hsp90aa1 polyclonal antibody
MAB4487	HLA-G monoclonal antibody, clone 4H84 (Biotin)	PAB8776	HSP90AA1 polyclonal antibody
MAB6225	HLA-G monoclonal antibody, clone F21-21	MAB6457	HSP90AA1/HSP90AB1 monoclonal antibody, clone MBH90AB
MAB3609	HLA-G monoclonal antibody, clone G233	MAB6458	HSP90AB1 monoclonal antibody, clone MBH90B
MAB3856	HLA-G monoclonal antibody, clone MEM-G/11	PAB18096	HSP90AB1 polyclonal antibody
MAB4520	HLA-G monoclonal antibody, clone MEM-G/11 (FITC)	MAB2195	HSP90B1 monoclonal antibody, clone 9G10

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name
MAB6618	HSP90B1 monoclonal antibody, clone H9010
PAB13570	Hsp90b1 polyclonal antibody
H00003305-D01	HSPA1L MaxPab rabbit polyclonal antibody (D01)
H00003306-M04	HSPA2 monoclonal antibody (M04), clone 3A11-1A2
MAB3737	HSPA4 monoclonal antibody, clone W27
PAB15700	Hspa4 polyclonal antibody
MAB6130	HSPA4/HSPA8 monoclonal antibody, clone N27F3-4
PAB15446	HSPA7 polyclonal antibody
MAB6636	HSPA8 monoclonal antibody, clone 1F2-H5
PAB15412	HSPA8 polyclonal antibody
PAB13598	HSPA9 polyclonal antibody
H00003315-M04	HSPB1 monoclonal antibody (M04), clone 3G3
MAB6642	HSPB1 monoclonal antibody, clone 5D12-A12
MAB1889	HSPB1 monoclonal antibody, clone G3.1
MAB5557	HSPB2 (phospho S28) monoclonal antibody, clone 41A8
PAB15424	HSPB8 polyclonal antibody
H00023640-D01	HSPBP1 MaxPab rabbit polyclonal antibody (D01)
H00051504-D01	HSPC152 MaxPab rabbit polyclonal antibody (D01)
H00003329-D01	HSPD1 MaxPab rabbit polyclonal antibody (D01)
MAB6621	HSPD1 monoclonal antibody, clone LK1
PAB15431	HSPE1 polyclonal antibody
PAB8773	Hsph1 polyclonal antibody
PAB16672	htb1 polyclonal antibody
MAB4819	HTT monoclonal antibody, clone HDB4E10
MAB4805	Human herpes virus type 6 gp60/110 monoclonal antibody, clone 200
MAB4806	Human herpes virus type 6 gp90 monoclonal antibody, clone 300
MAB4807	Human herpes virus type 6 p150 monoclonal antibody, clone 20
MAB1894	Human immunodeficiency virus reverse transcriptase monoclonal antibody, clone AS34
MAB5392	Human immunodeficiency virus type 1 gp41 monoclonal antibody, clone 1D4
MAB2453	Human immunodeficiency virus type 1 integrase monoclonal antibody, clone IN-2
MAB7835	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 01-002 (3D12)
MAB7836	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 2F2 (01-003)
MAB7839	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 2H12 (01-007)
MAB7838	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 3A2 (01-006)
MAB7837	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 3E6 (01-004)
MAB7834	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone 3F2 (01-001)
MAB7840	Human immunodeficiency virus type 1 Nef monoclonal antibody, clone JR6 (01-009)
MAB2434	Human immunodeficiency virus type 1 P17 Gag monoclonal antibody, clone 17-1
MAB5390	Human immunodeficiency virus type 1 p17 monoclonal antibody, clone 1C3
MAB4811	Human immunodeficiency virus type 1 p17 monoclonal antibody, clone 4A1
MAB5389	Human immunodeficiency virus type 1 p24 monoclonal antibody, clone 1A1
MAB7849	Human immunodeficiency virus type 1 p24 monoclonal antibody, clone 5 (05-001)
MAB7848	Human immunodeficiency virus type 1 Rev monoclonal antibody, clone 10 (03-002)
MAB7842	Human immunodeficiency virus type 1 Tat monoclonal antibody, clone 1 (02-001)
MAB5391	Human immunodeficiency virus type 2 gp36 monoclonal antibody, clone 3B10/G2
MAB4817	Human immunodeficiency virus type 2 gp36 monoclonal antibody, clone BD1411
MAB2438	Human immunodeficiency virus type 2 Nef monoclonal antibody, clone 4B5
MAB4227	Human papilloma virus type 16 E7 monoclonal antibody, clone 289-17013
MAB4225	Human papilloma virus type 16 L1 monoclonal antibody, clone 289-16981
MAB4229	Human papilloma virus type 18 E6 monoclonal antibody, clone 289-13965
H00003364-D01	HUS1 MaxPab rabbit polyclonal antibody (D01)
H00003364-M01	HUS1 monoclonal antibody (M01), clone 2B12
PAB10059	HUS1B polyclonal antibody
H00003373-M01	HYAL1 monoclonal antibody (M01), clone 2H7
H00008372-D01	HYAL3 MaxPab rabbit polyclonal antibody (D01)
PAB5153	Hydrocortisone polyclonal antibody
MAB1333	I-Ak monoclonal antibody, clone IA-G3
PAB5033	IAPP polyclonal antibody
PAB5080	Iapp polyclonal antibody
MAB1334	I-As monoclonal antibody, clone IA-G1
H00003383-D01	ICAM1 MaxPab rabbit polyclonal antibody (D01)
H00003383-M01	ICAM1 monoclonal antibody (M01), clone 3H8-2G6
MAB4195	ICAM1 monoclonal antibody, clone 289-10152
MAB8078	ICAM1 monoclonal antibody, clone 3H1547
MAB6935	ICAM1 monoclonal antibody, clone 8.4A6

Cat. Num.	Product Name
MAB9131	ICAM1 monoclonal antibody, clone 978w4
MAB4107	Icam1 monoclonal antibody, clone BE29G1
MAB5292	ICAM1 monoclonal antibody, clone MEM 112
MAB0892	ICAM1 monoclonal antibody, clone MEM-112
MAB5609	Icam1 monoclonal antibody, clone YN1/1.7.4
MAB5946	Icam1 monoclonal antibody, clone YN1/1.7.4 (Biotin)
MAB5945	Icam1 monoclonal antibody, clone YN1/1.7.4 (FITC)
MAB5947	Icam1 monoclonal antibody, clone YN1/1.7.4 (PE)
MAB1718	ICAM1 monoclonal antibody, clone YUK11
MAB3267	Icam2 monoclonal antibody, clone 3C4
MAB5840	Icam2 monoclonal antibody, clone 3C4 (Biotin)
MAB5839	Icam2 monoclonal antibody, clone 3C4 (FITC)
MAB5841	Icam2 monoclonal antibody, clone 3C4 (PE)
MAB0909	ICAM3 monoclonal antibody, clone MEM-171
H00003386-D01	ICAM4 MaxPab rabbit polyclonal antibody (D01)
MAB8028	ICOSLG monoclonal antibody, clone 2D3
MAB8266	ID3 monoclonal antibody
H00003418-D01	IDH2 MaxPab rabbit polyclonal antibody (D01)
H00003423-D01	IDS MaxPab rabbit polyclonal antibody (D01)
MAB3254	I-Ek monoclonal antibody, clone 14-4-4S
MAB5802	I-Ek monoclonal antibody, clone 14-4-4S (Biotin)
MAB5801	I-Ek monoclonal antibody, clone 14-4-4S (FITC)
MAB5803	I-Ek monoclonal antibody, clone 14-4-4S (PE)
MAB1335	I-Ek monoclonal antibody, clone IE-D6
MAB9307	if monoclonal antibody, clone 46d5
H00010437-D01	IFI30 MaxPab rabbit polyclonal antibody (D01)
H00064135-D01	IFIH1 MaxPab rabbit polyclonal antibody (D01)
H00003434-D01	IFIT1 MaxPab rabbit polyclonal antibody (D01)
H00003433-D01	IFIT2 MaxPab rabbit polyclonal antibody (D01)
H00003437-D01	IFIT3 MaxPab rabbit polyclonal antibody (D01)
H00010410-D01	IFITM3 MaxPab rabbit polyclonal antibody (D01)
H00010410-D03	IFITM3 MaxPab rabbit polyclonal antibody (D03)
H00010410-M02	IFITM3 monoclonal antibody (M02), clone 2H4-1D5
H00003439-D01	IFNA1 MaxPab rabbit polyclonal antibody (D01)
H00003446-D01	IFNA10 MaxPab rabbit polyclonal antibody (D01)
H00003451-D01	IFNA17 MaxPab rabbit polyclonal antibody (D01)
H00003452-D01	IFNA21 MaxPab rabbit polyclonal antibody (D01)
H00003441-D01	IFNA4 MaxPab rabbit polyclonal antibody (D01)
H00003442-D01	IFNA5 MaxPab rabbit polyclonal antibody (D01)
H00003443-M01	IFNA6 monoclonal antibody (M01), clone 3C9
H00003445-D01	IFNA8 MaxPab rabbit polyclonal antibody (D01)
PAB9607	IFNAR1 (phospho S535/S539) polyclonal antibody
H00003456-D01	IFNB1 MaxPab rabbit polyclonal antibody (D01)
H00338376-D01	IFNE MaxPab rabbit polyclonal antibody (D01)
H00338376-M01	IFNE1 monoclonal antibody (M01), clone 3B8
MAB1628	IFNG monoclonal antibody, clone K14-B
PAB5110	IFNG polyclonal antibody
H00003459-D01	IFNGR1 MaxPab rabbit polyclonal antibody (D01)
H00003460-D01	IFNGR2 MaxPab rabbit polyclonal antibody (D01)
H00007866-D01	IFRD2 MaxPab rabbit polyclonal antibody (D01)
MAB5556	IgA monoclonal antibody, clone 47C12
MAB8071	IgA monoclonal antibody, clone Mc24-2E11
MAB4834	IGF1R monoclonal antibody, clone 3B7
H00010644-D01	IGF2BP2 MaxPab rabbit polyclonal antibody (D01)
H00003489-M01	IGFBP6 monoclonal antibody (M01), clone 1H2-2D10
MAB5554	IgG monoclonal antibody, clone 3E8
MAB8514	IgG1 monoclonal antibody, clone MOPC-21
H00003500-D01	IGHG1 MaxPab rabbit polyclonal antibody (D01)
H00003514-D01	IGKC MaxPab rabbit polyclonal antibody (D01)
H00003514-M05	IGKC monoclonal antibody (M05), clone 2E7
H00003543-D01	IGLL1 MaxPab rabbit polyclonal antibody (D01)
MAB5555	IgM monoclonal antibody, clone 13A11
MAB6221	IgM monoclonal antibody, clone M-1
MAB5217	IGSF2 monoclonal antibody, clone BB27
H00010261-D01	IGSF6 MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB0075	IKBKB monoclonal antibody, clone 10AG2	H00003628-D01	INPP1 MaxPab rabbit polyclonal antibody (D01)
MAB0799	IKBKB monoclonal antibody, clone 42D1	H00003628-M14	INPP1 monoclonal antibody (M14), clone 1A6
H00064375-D01	IKZF4 MaxPab rabbit polyclonal antibody (D01)	H00003632-M05	INPP5A monoclonal antibody (M05), clone 3D8
MAB2539	IL10 monoclonal antibody, clone SPM464	MAB9475	INS monoclonal antibody, clone 86s5
H00003588-D01	IL10RB MaxPab rabbit polyclonal antibody (D01)	MAB9474	INS monoclonal antibody, clone 96s79
MAB7259	Il12 monoclonal antibody, clone C17.8	PAB5029	INS polyclonal antibody
H00003594-D01	IL12RB1 MaxPab rabbit polyclonal antibody (D01)	PAB5058	INS polyclonal antibody
H00003598-D01	IL13RA2 MaxPab rabbit polyclonal antibody (D01)	PAB5118	INS polyclonal antibody
H00003600-D01	IL15 MaxPab rabbit polyclonal antibody (D01)	MAB4835	INSR monoclonal antibody, clone 18-44
H00027189-D01	IL17C MaxPab rabbit polyclonal antibody (D01)	MAB4836	INSR monoclonal antibody, clone 83-14
H00053342-D01	IL17D MaxPab rabbit polyclonal antibody (D01)	MAB4837	INSR monoclonal antibody, clone 83-7
H00112744-D01	IL17F MaxPab rabbit polyclonal antibody (D01)	MAB2389	INSR monoclonal antibody, clone MA-20
PAB9533	Il18 polyclonal antibody	MAB8375	Insulin receptor monoclonal antibody, clone 47-9
H00008809-D01	IL18R1 MaxPab rabbit polyclonal antibody (D01)	MAB8374	Insulin receptor monoclonal antibody, clone 83-7
H00003553-D01	IL1B MaxPab rabbit polyclonal antibody (D01)	MAB9085	Integrin beta-1 monoclonal antibody, clone 67a9
H00003553-M01	IL1B monoclonal antibody (M01), clone 2A8	H00128239-M01	IQGAP3 monoclonal antibody (M01), clone 6E3
PAB10321	Il1b polyclonal antibody	H00003654-D01	IRAK1 MaxPab rabbit polyclonal antibody (D01)
PAB11294	IL1B polyclonal antibody	H00003654-M01	IRAK1 monoclonal antibody (M01), clone 3F7
PAB16928	Il1b polyclonal antibody	H00003654-M02	IRAK1 monoclonal antibody (M02), clone 3A9
PAB10060	IL1R1 polyclonal antibody	PAB0158	IRAK1 polyclonal antibody
PAB9952	IL1R2 polyclonal antibody	PAB10249	IRAK1 polyclonal antibody
H00003556-M06	IL1RAP monoclonal antibody (M06), clone 2A3	MAB2538	IRAK4 (phospho T345) monoclonal antibody, clone A8A8
H00003557-M03	IL1RN monoclonal antibody (M03), clone 1H5	H00003659-D01	IRF1 MaxPab rabbit polyclonal antibody (D01)
MAB1902	Il2 monoclonal antibody, clone 7H11AB7	H00003661-D01	IRF3 MaxPab rabbit polyclonal antibody (D01)
MAB1899	Il2 monoclonal antibody, clone JES6-5H4	H00003661-M10	IRF3 monoclonal antibody (M10), clone 3C8
H00050604-M01	IL20 monoclonal antibody (M01), clone 2H8	H00003663-M01	IRF5 monoclonal antibody (M01), clone 2E3-1A11
H00051561-D01	IL23A MaxPab rabbit polyclonal antibody (D01)	MAB1068	IRF5 monoclonal antibody, clone 10T1
H00011009-M01	IL24 monoclonal antibody (M01), clone 4B6	H00010379-D01	IRF9 MaxPab rabbit polyclonal antibody (D01)
H00246778-D01	IL27 MaxPab rabbit polyclonal antibody (D01)	PAB12627	IRS1 (phospho S312) polyclonal antibody
H00246778-M01	IL27 monoclonal antibody (M01), clone 3F12	PAB12628	IRS1 (phospho Y612) polyclonal antibody
MAB6205	Il2ra monoclonal antibody, clone 3C7	PAB12694	IRS1 polyclonal antibody
MAB3228	Il2ra monoclonal antibody, clone 7D4	PAB12695	IRS1 polyclonal antibody
MAB5706	Il2ra monoclonal antibody, clone 7D4 (APC)	PAB12629	IRS2 (phospho S731) polyclonal antibody
MAB5704	Il2ra monoclonal antibody, clone 7D4 (Biotin)	PAB12696	IRS2 polyclonal antibody
MAB5703	Il2ra monoclonal antibody, clone 7D4 (FITC)	PAB12697	IRS2 polyclonal antibody
MAB5705	Il2ra monoclonal antibody, clone 7D4 (PE)	H00009636-D01	ISG15 MaxPab rabbit polyclonal antibody (D01)
MAB5708	Il2ra monoclonal antibody, clone 7D4 (PE/Cy7)	H00003669-D01	ISG20 MaxPab rabbit polyclonal antibody (D01)
MAB5707	Il2ra monoclonal antibody, clone 7D4 (SPRD)	H00003669-M05	ISG20 monoclonal antibody (M05), clone S1
MAB0918	IL2RA monoclonal antibody, clone MEM-140	PAB15698	Islr2 polyclonal antibody
MAB4396	IL2RA monoclonal antibody, clone MEM-140 (Biotin)	MAB5224	ITGA1 monoclonal antibody, clone TS2/7
H00003560-D01	IL2RB MaxPab rabbit polyclonal antibody (D01)	MAB7187	ITGA2 monoclonal antibody, clone HAS-3
MAB3256	Il2rb monoclonal antibody, clone 5H4	MAB6939	ITGA2 monoclonal antibody, clone HAS6
MAB5810	Il2rb monoclonal antibody, clone 5H4 (APC)	MAB7098	Itga3 monoclonal antibody, clone MHF4
MAB5808	Il2rb monoclonal antibody, clone 5H4 (Biotin)	MAB5290	ITGA4 monoclonal antibody, clone 44H6
MAB5807	Il2rb monoclonal antibody, clone 5H4 (FITC)	MAB6941	ITGA4 monoclonal antibody, clone BU49
MAB5809	Il2rb monoclonal antibody, clone 5H4 (PE)	MAB5603	Itga4 monoclonal antibody, clone PS/2
MAB6937	IL2RB monoclonal antibody, clone 9A2 (DU-2)	MAB9128	ITGA5 monoclonal antibody, clone 598h5
H00009235-D01	IL32 MaxPab rabbit polyclonal antibody (D01)	MAB9130	ITGA5 monoclonal antibody, clone 978a
MAB7328	IL3RA monoclonal antibody, clone 6H6	MAB7087	ITGA5 monoclonal antibody, clone NKI-SAM-1
PAB9953	Il6 polyclonal antibody	MAB9129	ITGA5/CD49e monoclonal antibody, clone 59f6
H00003576-M12	IL8 monoclonal antibody (M12), clone 4G11	MAB9302	ITGA6 monoclonal antibody, clone 6g56
H00003576-M16	IL8 monoclonal antibody (M16), clone 1B4	MAB7096	ITGA6 monoclonal antibody, clone MP 4F10
H00003576-M17	IL8 monoclonal antibody (M17), clone 1A3	MAB6213	Itgae monoclonal antibody, clone 2E7
H00003576-M18	IL8 monoclonal antibody (M18), clone 4E3	MAB4124	ITGAE monoclonal antibody, clone 2G5.1
H00003576-M21	IL8 monoclonal antibody (M21), clone 4F1	MAB5465	ITGAE monoclonal antibody, clone 2G5.1 (Biotin)
PAB10317	IL8 polyclonal antibody	MAB5260	ITGAE monoclonal antibody, clone LF61
PAB8941	IL8RA polyclonal antibody	MAB7916	Itgal monoclonal antibody, clone M17/4
PAB8942	Il8rb polyclonal antibody	MAB0876	ITGAL monoclonal antibody, clone MEM-25
H00003608-D01	ILF2 MaxPab rabbit polyclonal antibody (D01)	MAB0951	ITGAL monoclonal antibody, clone MEM-83
H00003609-D01	ILF3 MaxPab rabbit polyclonal antibody (D01)	MAB7371	ITGAM monoclonal antibody, clone 24
H00003614-D01	IMPDH1 MaxPab rabbit polyclonal antibody (D01)	MAB3223	Itgam monoclonal antibody, clone 3A33
H00003615-M01	IMPDH2 monoclonal antibody (M01), clone 1E12-B6	MAB7255	ITGAM monoclonal antibody, clone CBRM1/5
H00003621-D01	ING1 MaxPab rabbit polyclonal antibody (D01)	MAB3224	Itgam monoclonal antibody, clone M1/70
H00051147-D01	ING4 MaxPab rabbit polyclonal antibody (D01)	MAB0877	ITGAM monoclonal antibody, clone MEM-174
H00083729-D01	INHBE MaxPab rabbit polyclonal antibody (D01)	MAB7080	Itgam monoclonal antibody, clone OX 42

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB6943	ITGAV monoclonal antibody, clone P2W7	MAB6655	KCNQ2 monoclonal antibody, clone S26A-23
MAB3864	ITGAX monoclonal antibody, clone BU15	MAB6656	KCNQ4 monoclonal antibody, clone S43-6
MAB4351	ITGAX monoclonal antibody, clone BU15 (APC)	H00386618-D01	KCTD4 MaxPab rabbit polyclonal antibody (D01)
MAB4350	ITGAX monoclonal antibody, clone BU15 (FITC)	MAB2194	KDEL1 monoclonal antibody, clone KR-10
MAB4352	ITGAX monoclonal antibody, clone BU15 (PE)	PAB12648	KDR (phospho Y951) polyclonal antibody
MAB4353	ITGAX monoclonal antibody, clone BU15 (PerCP)	PAB12154	Kdr polyclonal antibody
MAB9250	ITGB1 monoclonal antibody, clone 646s	PAB12647	KDR polyclonal antibody
MAB0883	ITGB1 monoclonal antibody, clone MEM-101A	PAB12746	KDR polyclonal antibody
MAB7010	ITGB1 monoclonal antibody, clone TS2/16	MAB5250	KEL monoclonal antibody, clone BRIC203
MAB9098	itgb1-a/itgb1-b monoclonal antibody, clone 589c	MAB1991	Keratin monoclonal antibody, clone C11
MAB9122	ITGB2 monoclonal antibody, clone 33d5	H00010657-D01	KHDRBS1 MaxPab rabbit polyclonal antibody (D01)
MAB9123	Itgb2 monoclonal antibody, clone 56f62	H00003795-D01	KHK MaxPab rabbit polyclonal antibody (D01)
MAB3225	Itgb2 monoclonal antibody, clone C71/16	H00009933-D01	KIAA0020 MaxPab rabbit polyclonal antibody (D01)
MAB7917	Itgb2 monoclonal antibody, clone M18/2	H00023392-D01	KIAA0368 MaxPab rabbit polyclonal antibody (D01)
MAB0929	ITGB2 monoclonal antibody, clone MEM-148	H00057214-D01	KIAA1199 MaxPab rabbit polyclonal antibody (D01)
MAB0880	ITGB2 monoclonal antibody, clone MEM-48	H00057214-M01	KIAA1199 monoclonal antibody (M01), clone 3C12
MAB6226	ITGB3 monoclonal antibody, clone PM6/13	MAB7078	Kidins220 monoclonal antibody, clone P220 1F8/3
MAB3249	Itgb4 monoclonal antibody, clone 346-11A	PAB15692	Kidins220 polyclonal antibody
MAB5239	ITGB4 monoclonal antibody, clone 450-9D	PAB12408	KIF11 polyclonal antibody
MAB7722	ITGB4 monoclonal antibody, clone 58XB4	PAB15697	Kif17 polyclonal antibody
MAB6946	ITGB4 monoclonal antibody, clone UMA9	PAB15699	Kif21a polyclonal antibody
MAB7204	ITGB7 monoclonal antibody, clone FIB27	H00026153-D01	KIF26A MaxPab rabbit polyclonal antibody (D01)
H00003696-M01	ITGB8 monoclonal antibody (M01), clone 2B4	PAB12407	KIF2A polyclonal antibody
H00003697-D01	ITIH1 MaxPab rabbit polyclonal antibody (D01)	H00011004-M01	KIF2C monoclonal antibody (M01), clone 1G2
H00003702-D01	ITK MaxPab rabbit polyclonal antibody (D01)	PAB15721	Kif5c polyclonal antibody
MAB1905	Itk monoclonal antibody, clone 2F12	MAB1078	KIR2DL1 monoclonal antibody, clone 2F9
H00003705-D01	ITPK1 MaxPab rabbit polyclonal antibody (D01)	MAB4187	KIR2DL2 monoclonal antibody, clone GL183
MAB7026	IVL monoclonal antibody, clone SY5	H00084063-M01	KIRREL2 monoclonal antibody (M01), clone 2B9-1D3
MAB7025	IVL monoclonal antibody, clone SY8	MAB3412	KIT monoclonal antibody, clone 57A5
PAB12749	JAK2 (phospho Y1007/Y1008) polyclonal antibody	PAB6148	KLF1 polyclonal antibody
PAB12769	JAK2 polyclonal antibody	H00028999-M02	KLF15 monoclonal antibody (M02), clone 1F3
PAB12770	JAK2 polyclonal antibody	H00011279-D01	KLF8 MaxPab rabbit polyclonal antibody (D01)
H00008681-D01	JMJD7-PLA2G4B MaxPab rabbit polyclonal antibody (D01)	PAB15377	KLH polyclonal antibody
PAB12568	JNK1 polyclonal antibody	H00054758-D01	KLHDC4 MaxPab rabbit polyclonal antibody (D01)
PAB8477	JUB polyclonal antibody	PAB15695	Klhl9 polyclonal antibody
PAB9930	JUB polyclonal antibody	H00005655-D01	KLK10 MaxPab rabbit polyclonal antibody (D01)
PAB12616	JUN (phospho S63) polyclonal antibody	H00003817-D01	KLK2 MaxPab rabbit polyclonal antibody (D01)
PAB12615	JUN (phospho T239) polyclonal antibody	H00003817-M03	KLK2 monoclonal antibody (M03), clone 3C5
H00003725-D01	JUN MaxPab rabbit polyclonal antibody (D01)	H00005650-D01	KLK7 MaxPab rabbit polyclonal antibody (D01)
PAB12669	JUN polyclonal antibody	MAB3244	Klr1c monoclonal antibody, clone PK136
H00003726-D01	JUNB MaxPab rabbit polyclonal antibody (D01)	MAB5778	Klr1c monoclonal antibody, clone PK136 (APC)
MAB1383	Jup monoclonal antibody, clone M111	MAB5776	Klr1c monoclonal antibody, clone PK136 (Biotin)
MAB2328	Kaposi's sarcoma-associated herpesvirus LANA2 monoclonal antibody, clone CM-A807	MAB5775	Klr1c monoclonal antibody, clone PK136 (FITC)
MAB6445	Kappa light chains monoclonal antibody, clone MEM-09	MAB5777	Klr1c monoclonal antibody, clone PK136 (PE)
H00010524-D01	KAT5 MaxPab rabbit polyclonal antibody (D01)	MAB5779	Klr1c monoclonal antibody, clone PK136 (PE/Cy7)
PAB8948	kay polyclonal antibody	H00003821-M01	KLRC1 monoclonal antibody (M01), clone 2C3
MAB8584	Kchip2 monoclonal antibody, clone S60-73	H00003821-M01A	KLRC1 monoclonal antibody (M01A), clone 2C3
MAB8574	Kcna1 monoclonal antibody, clone S20-78	H00010219-D01	KLRG1 MaxPab rabbit polyclonal antibody (D01)
MAB8572	Kcna2 monoclonal antibody, clone S14-16	MAB3246	Klrg1 monoclonal antibody, clone 2F1
MAB8571	Kcna4 monoclonal antibody, clone S13-31	H00003827-D01	KNG1 MaxPab rabbit polyclonal antibody (D01)
MAB8569	Kcna5 monoclonal antibody, clone S7-45	H00003827-M02	KNG1 monoclonal antibody (M02), clone 4A1
MAB8570	Kcnab1 monoclonal antibody, clone S9-40	PAB4976	KNG1 polyclonal antibody
MAB8573	Kcnab2 monoclonal antibody, clone S17-70	MAB3663	KRT10 monoclonal antibody, clone VIK-10
MAB7897	Kcnc4 monoclonal antibody, clone S72-16	H00003866-D01	KRT15 MaxPab rabbit polyclonal antibody (D01)
MAB8583	Kcnd2 monoclonal antibody, clone S57-1	H00003872-M01	KRT17 monoclonal antibody (M01), clone 2D4-1G9
MAB8588	Kcnd3 monoclonal antibody, clone S75-41	H00003875-M01	KRT18 monoclonal antibody (M01), clone 2F8
H00003753-D01	KCNE1 MaxPab rabbit polyclonal antibody (D01)	MAB0857	KRT18 monoclonal antibody, clone C-04
PAB15524	KCNH2 polyclonal antibody	MAB6504	KRT18 monoclonal antibody, clone C-04 (Biotin)
MAB8580	Kcnp1 monoclonal antibody, clone S55-7	MAB0860	KRT18 monoclonal antibody, clone DA-7
MAB8581	Kcnp1 monoclonal antibody, clone S55-82	MAB6506	KRT18 monoclonal antibody, clone DA-7 (Biotin)
H00030818-D01	KCNIP3 MaxPab rabbit polyclonal antibody (D01)	MAB6505	KRT18 monoclonal antibody, clone DC-10 (Biotin)
MAB8586	Kcnp3 monoclonal antibody, clone S66-38	PAB18098	KRT18 polyclonal antibody
MAB8589	Kcnma1 monoclonal antibody, clone S6-60	MAB3966	KRT19 monoclonal antibody, clone 236-11221
MAB7894	Kcnmb2 monoclonal antibody, clone S53-32	MAB3967	KRT19 monoclonal antibody, clone 236-11791
MAB6654	KCNQ1 monoclonal antibody, clone S37A-10	MAB0862	KRT19 monoclonal antibody, clone A53-B/A2

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB6502	KRT19 monoclonal antibody, clone A53-B/A2 (Biotin)	MAB9483	LCT monoclonal antibody, clone 9h058
MAB0861	KRT19 monoclonal antibody, clone BA-17	MAB6677	L-Cysteine monoclonal antibody, clone F2D
PAB12676	KRT19 polyclonal antibody	H00003939-D01	LDHA MaxPab rabbit polyclonal antibody (D01)
H00054474-D01	KRT20 MaxPab rabbit polyclonal antibody (D01)	MAB2737	LDHA monoclonal antibody, clone AF14A11
H00003885-D01	KRT34 MaxPab rabbit polyclonal antibody (D01)	MAB2732	LDHB monoclonal antibody, clone 60H11
H00003886-D01	KRT35 MaxPab rabbit polyclonal antibody (D01)	H00003948-D01	LDHC MaxPab rabbit polyclonal antibody (D01)
H00003852-M08	KRT5 monoclonal antibody (M08), clone 1A12	PAB11350	LDH-M polyclonal antibody
MAB0853	KRT5/KRT18 monoclonal antibody, clone C-50	PAB5135	LDL polyclonal antibody
MAB0859	KRT7/KRT17 monoclonal antibody, clone C-46	H00011061-D01	LECT1 MaxPab rabbit polyclonal antibody (D01)
H00003856-D01	KRT8 MaxPab rabbit polyclonal antibody (D01)	H00010637-M01	LEFTY1 monoclonal antibody (M01), clone 4D7
MAB0856	KRT8 monoclonal antibody, clone C-43	H00007044-D01	LEFTY2 MaxPab rabbit polyclonal antibody (D01)
H00003889-D01	KRT83 MaxPab rabbit polyclonal antibody (D01)	MAB5675	<i>Leishmania</i> A2 monoclonal antibody, clone C9
H00008942-D01	KYNU MaxPab rabbit polyclonal antibody (D01)	PAB14102	LEO1 polyclonal antibody
PAB14105	L3MBTL polyclonal antibody	PAB5093	LEP polyclonal antibody
H00091133-D01	L3MBTL4 MaxPab rabbit polyclonal antibody (D01)	PAB5111	Lep polyclonal antibody
MAB8235	La Crosse Virus monoclonal antibody, clone 10G5.4	PAB5113	Lep polyclonal antibody
MAB8233	La Crosse Virus monoclonal antibody, clone 4A5.5	MAB9304	LEP-100 monoclonal antibody, clone 7y95a
MAB8342	La Crosse Virus monoclonal antibody, clone 8C2.2	H00003954-M03	LETM1 monoclonal antibody (M03), clone 6F7
MAB8234	La Crosse Virus monoclonal antibody, clone 9E7.2	H00003955-M08	LFNG monoclonal antibody (M08), clone 4F8
H00003903-M09	LAIR1 monoclonal antibody (M09), clone 3B2	H00003956-M01	LGALS1 monoclonal antibody (M01), clone 1E8-1B2
MAB7869	LAMA4 monoclonal antibody, clone 3H2	H00003958-D01	LGALS3 MaxPab rabbit polyclonal antibody (D01)
MAB7870	LAMA4 monoclonal antibody, clone 6C3	H00003958-M01	LGALS3 monoclonal antibody (M01), clone 3C2-2A3
MAB7871	LAMA4 monoclonal antibody, clone 9B2	H00003959-D01	LGALS3BP MaxPab rabbit polyclonal antibody (D01)
MAB7872	LAMA5 monoclonal antibody, clone 4B12	H00005641-D01	LGMN MaxPab rabbit polyclonal antibody (D01)
MAB9420	LAMB2 monoclonal antibody, clone 135ta	H00005641-M03	LGMN monoclonal antibody (M03), clone M2
MAB9417	LAMB2 monoclonal antibody, clone 946s5	MAB7994	LIG1 monoclonal antibody, clone 10H5
MAB6979	Lambda chain monoclonal antibody, clone ICO106	H00010859-D01	LILRB1 MaxPab rabbit polyclonal antibody (D01)
MAB2519	LAMC1 monoclonal antibody, clone BC17	MAB0960	LIME1 monoclonal antibody, clone LIME-06
MAB9232	Lamin (Dm0) monoclonal antibody, clone 4n95	MAB8509	Lime1 monoclonal antibody, clone mLIME-05
MAB9233	Lamin (Dm0) monoclonal antibody, clone 5n98	PAB13550	Lime1 polyclonal antibody
PAB10312	Laminin polyclonal antibody	H00079727-D01	LIN28 MaxPab rabbit polyclonal antibody (D01)
MAB9468	LAMP1 monoclonal antibody, clone 116vd	H00286826-D01	LIN9 MaxPab rabbit polyclonal antibody (D01)
MAB3265	Lamp1 monoclonal antibody, clone 1D4B	MAB9105	Link protein monoclonal antibody, clone 5m795
MAB5836	Lamp1 monoclonal antibody, clone 1D4B (APC)	H00003988-M01	LIPA monoclonal antibody (M01), clone 1F9
MAB5835	Lamp1 monoclonal antibody, clone 1D4B (FITC)	H00009388-M01	LIPG monoclonal antibody (M01), clone 4A9
MAB6634	Lamp1 monoclonal antibody, clone Ly1C6	H00009516-D01	LITAF MaxPab rabbit polyclonal antibody (D01)
MAB9226	Lamp2 monoclonal antibody, clone 1y653	H00004000-D01	LMNA MaxPab rabbit polyclonal antibody (D01)
MAB9473	Lamp2 monoclonal antibody, clone 34d9	MAB2499	LMNA monoclonal antibody, clone JOL2
MAB3266	Lamp2 monoclonal antibody, clone ABL-93	MAB2501	LMNA monoclonal antibody, clone JOL4
MAB5838	Lamp2 monoclonal antibody, clone ABL-93 (APC)	H00004001-M01	LMNB1 monoclonal antibody (M01), clone 4B10
MAB5837	Lamp2 monoclonal antibody, clone ABL-93 (FITC)	MAB1761	LMNB1 monoclonal antibody, clone ZL-5
MAB6635	Lamp2 monoclonal antibody, clone GL2A7	MAB2336	LMO2 monoclonal antibody, clone 1A9-3B11
H00051056-D01	LAP3 MaxPab rabbit polyclonal antibody (D01)	H00493829-M04	LOC493829 monoclonal antibody (M04), clone 2B8
H00055323-D01	LARP6 MaxPab rabbit polyclonal antibody (D01)	H00493829-M05	LOC493829 monoclonal antibody (M05), clone 3B5
MAB3865	LAT monoclonal antibody, clone LAT-01	H00091614-M01	LOC91614 monoclonal antibody (M01), clone 4D9
MAB8029	LAT monoclonal antibody, clone LAT1111	H00004016-D01	LOXL1 MaxPab rabbit polyclonal antibody (D01)
MAB3623	Lat monoclonal antibody, clone LAT-1D1	H00004023-D01	LPL MaxPab rabbit polyclonal antibody (D01)
MAB0945	LAT2 monoclonal antibody, clone NAP-03	PAB15495	LPO polyclonal antibody
MAB0946	LAT2 monoclonal antibody, clone NAP-07	PAB13792	LRDD polyclonal antibody
MAB4522	LAT2 monoclonal antibody, clone NAP-07 (PE)	MAB5197	LRP1 monoclonal antibody, clone A2Mr alpha-2
MAB0019	LATS2 monoclonal antibody, clone ST-3D10	PAB12631	LRP5 (phospho T1492) polyclonal antibody
MAB5666	Lbp monoclonal antibody, clone biG 311	PAB12699	LRP5 polyclonal antibody
PAB15583	Lbr polyclonal antibody	PAB12632	LRP6 (phospho T1479) polyclonal antibody
PAB12750	LCK (phospho Y505) polyclonal antibody	PAB12700	LRP6 polyclonal antibody
PAB8909	LCK (phospho Y505) polyclonal antibody	H00079705-D01	LRRK1 MaxPab rabbit polyclonal antibody (D01)
H00003932-M01	LCK monoclonal antibody (M01), clone 3F7-F5	H00079705-M03	LRRK1 monoclonal antibody (M03), clone 3G8
MAB1939	Lck monoclonal antibody, clone 3A5	PAB12701	LRRK2 polyclonal antibody
MAB2724	LCK monoclonal antibody, clone AF4D5	H00090678-D01	LRSAM1 MaxPab rabbit polyclonal antibody (D01)
MAB0912	LCK monoclonal antibody, clone LCK-01	PAB15450	LSD1 polyclonal antibody
PAB12774	LCK polyclonal antibody	PAB12633	LSP1 (phospho S252) polyclonal antibody
PAB12775	LCK polyclonal antibody	H00004046-D01	LSP1 MaxPab rabbit polyclonal antibody (D01)
PAB8803	LCK polyclonal antibody	PAB12702	LSP1 polyclonal antibody
H00051451-D01	LCMT1 MaxPab rabbit polyclonal antibody (D01)	MAB6480	LST1 monoclonal antibody, clone LST1/02
H00003934-D01	LCN2 MaxPab rabbit polyclonal antibody (D01)	MAB7746	LTBR monoclonal antibody, clone 31G4D8
MAB8067	LCP1 monoclonal antibody, clone SPM414	MAB3721	LTBR monoclonal antibody, clone Carlos-1

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB0745	LTF monoclonal antibody, clone 1B8	PAB0327	MAP3K5 polyclonal antibody
PAB5050	LTF polyclonal antibody	H00006885-D01	MAP3K7 MaxPab rabbit polyclonal antibody (D01)
MAB3239	Ly6a monoclonal antibody, clone D7	H00010454-D01	MAP3K7IP1 MaxPab rabbit polyclonal antibody (D01)
MAB5760	Ly6a monoclonal antibody, clone D7 (APC)	H00023118-D01	MAP3K7IP2 MaxPab rabbit polyclonal antibody (D01)
MAB5758	Ly6a monoclonal antibody, clone D7 (Biotin)	MAB0813	MAP4K1 monoclonal antibody, clone 2A1
MAB5757	Ly6a monoclonal antibody, clone D7 (FITC)	MAB1636	MAPK1 monoclonal antibody, clone B20-U
MAB5759	Ly6a monoclonal antibody, clone D7 (PE)	PAB14286	MAPK1 polyclonal antibody
MAB5762	Ly6a monoclonal antibody, clone D7 (PE/Cy5.5)	PAB17866	Mapk1 polyclonal antibody
MAB5761	Ly6a monoclonal antibody, clone D7 (SPRD)	PAB9958	MAPK1 polyclonal antibody
MAB5613	Ly6g monoclonal antibody, clone RB6-8C5	MAB3337	MAPK1/3 monoclonal antibody, clone SB46b
MAB5971	Ly6g monoclonal antibody, clone RB6-8C5 (APC)	MAB5901	MAPK1/3 monoclonal antibody, clone SB46b (Biotin)
MAB5976	Ly6g monoclonal antibody, clone RB6-8C5 (APC/CY7)	MAB5900	MAPK1/3 monoclonal antibody, clone SB46b (HRP)
MAB5969	Ly6g monoclonal antibody, clone RB6-8C5 (Biotin)	MAB1637	MAPK1/MAPK3 (phospho T202/204) monoclonal antibody, clone G15-B
MAB5973	Ly6g monoclonal antibody, clone RB6-8C5 (Cy5)	PAB12748	MAPK1/MAPK3 (phospho T202/Y204) polyclonal antibody
MAB5968	Ly6g monoclonal antibody, clone RB6-8C5 (FITC)	PAB12763	MAPK1/MAPK3 polyclonal antibody
MAB5970	Ly6g monoclonal antibody, clone RB6-8C5 (PE)	PAB12764	MAPK1/MAPK3 polyclonal antibody
MAB5974	Ly6g monoclonal antibody, clone RB6-8C5 (PE/Cy5.5)	H00005600-D01	MAPK11 MaxPab rabbit polyclonal antibody (D01)
MAB5975	Ly6g monoclonal antibody, clone RB6-8C5 (PE/Cy7)	H00006300-D01	MAPK12 MaxPab rabbit polyclonal antibody (D01)
MAB5972	Ly6g monoclonal antibody, clone RB6-8C5 (SPRD)	H00005603-M06A	MAPK13 monoclonal antibody (M06A), clone 3D8
PAB0184	LY86 polyclonal antibody	H00005603-M09	MAPK13 monoclonal antibody (M09), clone 2B6
PAB10260	LYG2 polyclonal antibody	H00005603-M13	MAPK13 monoclonal antibody (M13), clone 3H7
MAB1675	Lyl1 monoclonal antibody, clone KT43	PAB8911	MAPK14 (phospho T180/Y182) polyclonal antibody
MAB1003	LYN monoclonal antibody, clone 11A7	MAB6637	MAPK14 monoclonal antibody, clone 9F12
MAB0949	LYN monoclonal antibody, clone LYN-01	PAB8882	MAPK14 polyclonal antibody
PAB5044	LYZ polyclonal antibody	PAB8906	MAPK3 (phospho T202/Y204) polyclonal antibody
H00054585-D01	LZTFL1 MaxPab rabbit polyclonal antibody (D01)	PAB10301	MAPK3 polyclonal antibody
H00054585-M09	LZTFL1 monoclonal antibody (M09), clone 1B1	H00005597-D01	MAPK6 MaxPab rabbit polyclonal antibody (D01)
MAB1948	M13 Bacteriophage (gp8) monoclonal antibody, clone RL-ph2	PAB12630	MAPK8 (phospho T183/Y185) polyclonal antibody
MAB1930	MAD1L1 monoclonal antibody, clone 9B10	H00005599-M01	MAPK8 monoclonal antibody (M01), clone 3H2
MAB8826	MAD1L1 monoclonal antibody, clone P22.9B10.C12	H00005599-M03	MAPK8 monoclonal antibody (M03), clone 4F10
PAB8815	MAD1L1 polyclonal antibody	H00005599-M04	MAPK8 monoclonal antibody (M04), clone 2F3
MAB1926	MAD2L1 monoclonal antibody, clone 17D10	H00005599-M05	MAPK8 monoclonal antibody (M05), clone 4H5
MAB8965	MAD2L1 monoclonal antibody, clone KD6.17D10.2G9	H00005599-M06	MAPK8 monoclonal antibody (M06), clone 4A11
H00009587-D01	MAD2L1BP MaxPab rabbit polyclonal antibody (D01)	H00005599-M07	MAPK8 monoclonal antibody (M07), clone 4H6
H00009587-M13	MAD2L1BP monoclonal antibody (M13), clone 4B9	PAB8795	MAPK8 polyclonal antibody
H00010459-D01	MAD2L2 MaxPab rabbit polyclonal antibody (D01)	PAB17934	MAPK9 (phospho T183/Y185) polyclonal antibody
MAB3251	Madcam1 monoclonal antibody, clone MECA-367	H00009261-D01	MAPKAPK2 MaxPab rabbit polyclonal antibody (D01)
MAB5795	Madcam1 monoclonal antibody, clone MECA-367 (Biotin)	H00009261-M01	MAPKAPK2 monoclonal antibody (M01), clone 2B3
MAB5794	Madcam1 monoclonal antibody, clone MECA-367 (FITC)	H00009261-M06	MAPKAPK2 monoclonal antibody (M06), clone 1F9
H00004110-D01	MAGEA11 MaxPab rabbit polyclonal antibody (D01)	PAB8818	Mapkapk2 polyclonal antibody
H00004102-D01	MAGEA3 MaxPab rabbit polyclonal antibody (D01)	H00008550-D01	MAPKAPK5 MaxPab rabbit polyclonal antibody (D01)
H00004103-D01	MAGEA4 MaxPab rabbit polyclonal antibody (D01)	MAB1678	Mapre1 monoclonal antibody, clone KT51
H00004107-D01	MAGEA8 MaxPab rabbit polyclonal antibody (D01)	H00010982-D01	MAPRE2 MaxPab rabbit polyclonal antibody (D01)
H00028986-D01	MAGEH1 MaxPab rabbit polyclonal antibody (D01)	MAB1681	Mapre2 monoclonal antibody, clone KT52
H00028986-M08	MAGEH1 monoclonal antibody (M08), clone 2B3	MAB9139	MAPT monoclonal antibody, clone 59m6
PAB0270	MALT1 polyclonal antibody	MAB5429	MAPT monoclonal antibody, clone 7B10.F2.E5
H00004128-D01	MAOA MaxPab rabbit polyclonal antibody (D01)	MAB1773	MAPT monoclonal antibody, clone B11E8
H00004129-D01	MAOB MaxPab rabbit polyclonal antibody (D01)	MAB1832	MAPT monoclonal antibody, clone TAU-5
MAB3866	MAP2 monoclonal antibody, clone MT-01	PAB14198	MAPT polyclonal antibody
MAB3636	MAP2 monoclonal antibody, clone MT-07	PAB14199	MAPT polyclonal antibody
MAB3637	MAP2 monoclonal antibody, clone MT-08	PAB9049	MAPT polyclonal antibody
PAB8829	Map2k1 polyclonal antibody	PAB9846	MARK3 polyclonal antibody
PAB12635	MAP2K1/MAP2K2 (phospho S218/S222) polyclonal antibody	H00004141-D01	MARS MaxPab rabbit polyclonal antibody (D01)
PAB12706	MAP2K1/MAP2K2 polyclonal antibody	PAB14126	MBD1 polyclonal antibody
PAB12707	MAP2K1/MAP2K2 polyclonal antibody	H00053615-D01	MBD3 MaxPab rabbit polyclonal antibody (D01)
H00005605-D01	MAP2K2 MaxPab rabbit polyclonal antibody (D01)	MAB3278	MCAM monoclonal antibody, clone C264
MAB0800	MAP2K2 monoclonal antibody, clone 1A1	MAB5851	MCAM monoclonal antibody, clone C264 (Biotin)
H00006416-D01	MAP2K4 MaxPab rabbit polyclonal antibody (D01)	MAB5850	MCAM monoclonal antibody, clone C264 (FITC)
MAB2716	MAP2K4 monoclonal antibody, clone AF23A1	MAB5852	MCAM monoclonal antibody, clone C264 (PE)
PAB12634	MAP2K5 (phospho S311/T315) polyclonal antibody	H00027349-D01	MCAT MaxPab rabbit polyclonal antibody (D01)
H00005607-D01	MAP2K5 MaxPab rabbit polyclonal antibody (D01)	PAB0260	MCL1 polyclonal antibody
PAB12704	MAP2K5 polyclonal antibody	PAB8824	MCL1 polyclonal antibody
PAB12705	MAP2K5 polyclonal antibody	H00004171-D01	MCM2 MaxPab rabbit polyclonal antibody (D01)
MAB2695	MAP2K6 monoclonal antibody, clone AF4G6	H00004171-M01	MCM2 monoclonal antibody (M01), clone 6A8
MAB2696	MAP2K6 monoclonal antibody, clone AF8B6	H00010445-D01	MCRS1 MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB10262	MDH1/MDH2 polyclonal antibody	H00004314-D01	MMP3 MaxPab rabbit polyclonal antibody (D01)
H00004192-D01	MDK MaxPab rabbit polyclonal antibody (D01)	H00004314-M02	MMP3 monoclonal antibody (M02), clone 4C11
PAB10334	Mdm2 (phospho S185) polyclonal antibody	PAB9548	MMP3 polyclonal antibody
PAB10274	Mdm2 polyclonal antibody	H00004316-D01	MMP7 MaxPab rabbit polyclonal antibody (D01)
MAB6760	MDM4 (phospho S367) monoclonal antibody, clone #15	PAB12712	MMP7 polyclonal antibody
H00004194-D01	MDM4 MaxPab rabbit polyclonal antibody (D01)	H00004317-D01	MMP8 MaxPab rabbit polyclonal antibody (D01)
PAB8758	MDM4 polyclonal antibody	PAB12713	MMP8 polyclonal antibody
PAB0646	Mecp2 polyclonal antibody	PAB9550	MMP8 polyclonal antibody
H00051102-D01	MECR MaxPab rabbit polyclonal antibody (D01)	PAB12714	MMP9 polyclonal antibody
PAB8661	MED1 polyclonal antibody	PAB13701	MMP9 polyclonal antibody
H00054797-D01	MED18 MaxPab rabbit polyclonal antibody (D01)	PAB8850	Mmp9 polyclonal antibody
H00004221-D01	MEN1 MaxPab rabbit polyclonal antibody (D01)	H00027304-D01	MOCS3 MaxPab rabbit polyclonal antibody (D01)
H00004223-D01	MEOX2 MaxPab rabbit polyclonal antibody (D01)	H00004340-D01	MOG MaxPab rabbit polyclonal antibody (D01)
MAB9488	MEP1A monoclonal antibody, clone 88154	MAB6219	Monocyte/Granulocyte monoclonal antibody, clone 74-22-15
PAB11308	MER2 (phospho S30) polyclonal antibody	H00004342-D01	MOS MaxPab rabbit polyclonal antibody (D01)
PAB11356	MER2 polyclonal antibody	MAB7045	Mos monoclonal antibody, clone S3.1
PAB10350	Methylated Lysine polyclonal antibody	H00004352-D01	MPL MaxPab rabbit polyclonal antibody (D01)
PAB16626	Methylated Lysine polyclonal antibody	MAB0772	MPO monoclonal antibody, clone 1A1
H00004240-D01	MFGE8 MaxPab rabbit polyclonal antibody (D01)	H00004355-D01	MPP2 MaxPab rabbit polyclonal antibody (D01)
H00004240-M03	MFGE8 monoclonal antibody (M03), clone 1E7	H00004355-M07	MPP2 monoclonal antibody (M07), clone 2B1
MAB3420	Mfge8 monoclonal antibody, clone 2422	H00004356-D01	MPP3 MaxPab rabbit polyclonal antibody (D01)
MAB9484	MGAM monoclonal antibody, clone 0h9j4	H00004357-M01	MPST monoclonal antibody (M01), clone 1B5
H00004248-M01	MGAT3 monoclonal antibody (M01), clone 2G4	H00009019-D01	MPZL1 MaxPab rabbit polyclonal antibody (D01)
H00011282-D01	MGAT4B MaxPab rabbit polyclonal antibody (D01)	H00010205-D01	MPZL2 MaxPab rabbit polyclonal antibody (D01)
H00010724-D01	MGEA5 MaxPab rabbit polyclonal antibody (D01)	PAB10277	MRE11 polyclonal antibody
H00004255-D01	MGMT MaxPab rabbit polyclonal antibody (D01)	MAB7996	MRE11A monoclonal antibody, clone 12D7
H00004255-M01	MGMT monoclonal antibody (M01), clone 2G4	PAB12381	MRE11A polyclonal antibody
H00004258-M01A	MGST2 monoclonal antibody (M01A), clone 1D11	H00055686-D01	MREG MaxPab rabbit polyclonal antibody (D01)
H00004258-M04	MGST2 monoclonal antibody (M04), clone 4E1	H00065993-D01	MRPS34 MaxPab rabbit polyclonal antibody (D01)
MAB9241	MHC class I monoclonal antibody, clone 54n5	H00051154-D01	MRTO4 MaxPab rabbit polyclonal antibody (D01)
MAB1363	MHC Class I monoclonal antibody, clone KD-H3	H00000931-D01	MS4A1 MaxPab rabbit polyclonal antibody (D01)
MAB6934	MHC Class II monoclonal antibody, clone DF1588	MAB6902	MS4A1 monoclonal antibody, clone 2H7
MAB6981	MHC Class II monoclonal antibody, clone TDR31.1	MAB5614	MS4A1 monoclonal antibody, clone B-Ly1
MAB7335	MICA monoclonal antibody, clone 6D4	MAB5993	MS4A1 monoclonal antibody, clone B-Ly1 (APC)
H00004277-D01	MICB MaxPab rabbit polyclonal antibody (D01)	MAB5996	MS4A1 monoclonal antibody, clone B-Ly1 (APC/CY7)
H00004282-M01	MIF monoclonal antibody (M01), clone 2A10-4D3	MAB5991	MS4A1 monoclonal antibody, clone B-Ly1 (Biotin)
PAB8836	MIF polyclonal antibody	MAB5990	MS4A1 monoclonal antibody, clone B-Ly1 (FITC)
PAB9001	mis11 polyclonal antibody	MAB5992	MS4A1 monoclonal antibody, clone B-Ly1 (PE)
H00004286-D01	MITF MaxPab rabbit polyclonal antibody (D01)	MAB5995	MS4A1 monoclonal antibody, clone B-Ly1 (PE/Cy7)
MAB1955	MITF monoclonal antibody, clone C5	MAB5994	MS4A1 monoclonal antibody, clone B-Ly1 (SPRD)
PAB12698	MKI67 polyclonal antibody	MAB0881	MS4A1 monoclonal antibody, clone MEM-97
H00008195-D01	MKKS MaxPab rabbit polyclonal antibody (D01)	PAB12666	MS4A1 polyclonal antibody
H00023609-D01	MKRN2 MaxPab rabbit polyclonal antibody (D01)	H00004440-D01	MSI1 MaxPab rabbit polyclonal antibody (D01)
H00004291-D01	MLF1 MaxPab rabbit polyclonal antibody (D01)	H00124540-D01	MSI2 MaxPab rabbit polyclonal antibody (D01)
H00055904-D01	MLL5 MaxPab rabbit polyclonal antibody (D01)	H00010232-D01	MSLN MaxPab rabbit polyclonal antibody (D01)
PAB12520	MLXIPL polyclonal antibody	H00010232-M02	MSLN monoclonal antibody (M02), clone 1A10
MAB8549	MME monoclonal antibody, clone CB-Calla	MAB0727	MSRA monoclonal antibody, clone 1C8
MAB0875	MME monoclonal antibody, clone MEM-78	MAB0728	MSRA monoclonal antibody, clone 5B5
MAB6949	MME monoclonal antibody, clone SN5c/L4-1A1	PAB12611	MST1R (phospho Y1238) polyclonal antibody
H00004312-D01	MMP1 MaxPab rabbit polyclonal antibody (D01)	PAB12612	MST1R (phospho Y1238/S1239) polyclonal antibody
MAB9060	Mmp1 monoclonal antibody, clone 57o5h	PAB12613	MST1R (phospho Y1239) polyclonal antibody
PAB12708	MMP1 polyclonal antibody	PAB12664	MST1R polyclonal antibody
H00004319-D01	MMP10 MaxPab rabbit polyclonal antibody (D01)	PAB12665	MST1R polyclonal antibody
H00004321-D01	MMP12 MaxPab rabbit polyclonal antibody (D01)	H00004487-D01	MSX1 MaxPab rabbit polyclonal antibody (D01)
H00004322-D01	MMP13 MaxPab rabbit polyclonal antibody (D01)	H00004487-M03	MSX1 monoclonal antibody (M03), clone 2B12
H00004322-M07	MMP13 monoclonal antibody (M07), clone 3B11	H00004487-M04	MSX1 monoclonal antibody (M04), clone 3A8
PAB12710	MMP13 polyclonal antibody	H00004487-M05	MSX1 monoclonal antibody (M05), clone 1D2
PAB8846	Mmp13 polyclonal antibody	H00004487-M11	MSX1 monoclonal antibody (M11), clone 1E2
PAB17925	MMP17 polyclonal antibody	H00004488-M10	MSX2 monoclonal antibody (M10), clone 1E8
MAB2334	MMP2 monoclonal antibody, clone MMP2/2C1	H00092140-D01	MTDH MaxPab rabbit polyclonal antibody (D01)
MAB2335	MMP2 monoclonal antibody, clone MMP2/8B4	H00004524-D01	MTHFR MaxPab rabbit polyclonal antibody (D01)
MAB2541	MMP2 monoclonal antibody, clone SPM346	PAB12715	MTOR polyclonal antibody
PAB12711	MMP2 polyclonal antibody	H00051537-D01	MTP18 MaxPab rabbit polyclonal antibody (D01)
PAB8848	Mmp2 polyclonal antibody	MAB2474	MUS81 monoclonal antibody, clone 2G10/3
PAB17924	MMP24 polyclonal antibody	PAB10279	MUT polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00004595-M01	MUTYH monoclonal antibody (M01), clone 4D10	H00004692-M06	NDN monoclonal antibody (M06), clone 3C1
H00004597-D01	MVD MaxPab rabbit polyclonal antibody (D01)	H00004692-M15	NDN monoclonal antibody (M15), clone 2G8
H00004598-D01	MVK MaxPab rabbit polyclonal antibody (D01)	H00051079-D01	NDUFA13 MaxPab rabbit polyclonal antibody (D01)
MAB1945	MVP monoclonal antibody, clone 1032	H00004716-D01	NDUFB10 MaxPab rabbit polyclonal antibody (D01)
H00004602-D01	MYB MaxPab rabbit polyclonal antibody (D01)	H00004712-D01	NDUFB6 MaxPab rabbit polyclonal antibody (D01)
H00010514-D01	MYBBP1A MaxPab rabbit polyclonal antibody (D01)	H00004720-D01	NDUFS2 MaxPab rabbit polyclonal antibody (D01)
PAB8903	MYC (phospho T58/S62) polyclonal antibody	H00004722-D01	NDUFS3 MaxPab rabbit polyclonal antibody (D01)
H00004609-D01	MYC MaxPab rabbit polyclonal antibody (D01)	H00004722-M02	NDUFS3 monoclonal antibody (M02), clone 1D6
MAB0967	MYC monoclonal antibody, clone 9E10	MAB2091	NEDD9 monoclonal antibody, clone 14A11
MAB4518	Myc monoclonal antibody, clone 9E10 (FITC)	MAB2476	NEDD9 monoclonal antibody, clone 2G9
MAB4519	Myc monoclonal antibody, clone 9E10 (HRP)	MAB8969	NEDD9 monoclonal antibody, clone C7.16D5.2G9.3C6
MAB7910	MYC monoclonal antibody, clone LF10E3	H00004747-D01	NEFL MaxPab rabbit polyclonal antibody (D01)
MAB5543	MYC monoclonal antibody, clone Myc.A7	MAB7757	NEFL monoclonal antibody, clone 2F11
PAB8856	MYC polyclonal antibody	H00252969-D01	NEIL2 MaxPab rabbit polyclonal antibody (D01)
PAB12716	MYC tag polyclonal antibody	H00004751-D01	NEK2 MaxPab rabbit polyclonal antibody (D01)
H00004615-D01	MYD88 MaxPab rabbit polyclonal antibody (D01)	H00004751-M01	NEK2 monoclonal antibody (M01), clone 2F6
MAB7264	Myeloid Marker monoclonal antibody, clone BM-1	H00004752-D01	NEK3 MaxPab rabbit polyclonal antibody (D01)
PAB15685	Myh9 polyclonal antibody	H00091754-D01	NEK9 MaxPab rabbit polyclonal antibody (D01)
H00004633-M01	MYL2 monoclonal antibody (M01), clone 3B9-B4	H00004745-M01	NELL1 monoclonal antibody (M01), clone 6A8
PAB9961	MYL6 (phospho S19/20) polyclonal antibody	PAB10408	Nephronophthisis (pic) polyclonal antibody
PAB11255	MYL9 polyclonal antibody	MAB2473	NES monoclonal antibody, clone 2C1.3A11
PAB10281	MYO1G polyclonal antibody	MAB9117	NES monoclonal antibody, clone 499s
H00004653-M01	MYOC monoclonal antibody (M01), clone 4F8	MAB9436	NES monoclonal antibody, clone 5h456
PAB8971	MYOC polyclonal antibody	MAB9463	Neural precursor cells monoclonal antibody, clone 31a4w
MAB1962	Myod1 monoclonal antibody, clone 5.2F	MAB9081	Neurocan receptor monoclonal antibody, clone 9f75
MAB0077	Myod1 monoclonal antibody, clone 5.8A	H00004760-D01	NEUROD1 MaxPab rabbit polyclonal antibody (D01)
MAB9458	Myog monoclonal antibody, clone 2t35e	H00004761-D01	NEUROD2 MaxPab rabbit polyclonal antibody (D01)
MAB1961	Myog monoclonal antibody, clone F12B	H00050674-M01	NEUROG3 monoclonal antibody (M01), clone 4B7
MAB3557	Myog monoclonal antibody, clone F5D	H00050674-M02	NEUROG3 monoclonal antibody (M02), clone 4E9
MAB9308	mys monoclonal antibody, clone 97t9	PAB12637	NF1 (phospho S2741) polyclonal antibody
PAB15516	MYST1 polyclonal antibody	MAB2306	NF1 monoclonal antibody, clone McFn27
H00011143-D01	MYST2 MaxPab rabbit polyclonal antibody (D01)	MAB2305	NF1 monoclonal antibody, clone McFn27a
H00029104-D01	N6AMT1 MaxPab rabbit polyclonal antibody (D01)	PAB12636	NF1 polyclonal antibody
H00027163-M02	NAAA monoclonal antibody (M02), clone 3F4	PAB12717	NF1 polyclonal antibody
H00004665-D01	NAB2 MaxPab rabbit polyclonal antibody (D01)	H00004771-D01	NF2 MaxPab rabbit polyclonal antibody (D01)
MAB2333	NACC1 monoclonal antibody, clone 3	H00084901-D01	NFATC2IP MaxPab rabbit polyclonal antibody (D01)
H00008883-D01	NAE1 MaxPab rabbit polyclonal antibody (D01)	H00084901-M02	NFATC2IP monoclonal antibody (M02), clone 3E9-B7
MAB8231	NAMPT monoclonal antibody, clone AF-1E12	MAB0056	NFKBIA (phospho S32/S36) monoclonal antibody, clone 39A1413
PAB13693	Nanog polyclonal antibody	PAB12626	NFKBIA (phospho S32/S36) polyclonal antibody
H00004674-D01	NAP1L2 MaxPab rabbit polyclonal antibody (D01)	PAB7922	NFKBIA (phospho Y42) polyclonal antibody
H00004675-D01	NAP1L3 MaxPab rabbit polyclonal antibody (D01)	MAB1760	NFKBIA monoclonal antibody, clone 6A055
MAB1795	NAPA monoclonal antibody, clone 15D4	MAB0057	NFKBIA monoclonal antibody, clone 6A920
MAB1807	NAPA monoclonal antibody, clone 4E4	PAB12692	NFKBIA polyclonal antibody
MAB9433	NAPA-73 monoclonal antibody, clone 34g55	PAB12693	NFKBIA polyclonal antibody
H00093100-D01	NAPRT1 MaxPab rabbit polyclonal antibody (D01)	PAB0659	NFYB polyclonal antibody
H00004677-D01	NARS MaxPab rabbit polyclonal antibody (D01)	MAB3446	NGFR monoclonal antibody, clone Nagy-1
PAB15597	NAS6 polyclonal antibody	MAB2324	NHE-3 monoclonal antibody, clone 19F5
H00000009-D01	NAT1 MaxPab rabbit polyclonal antibody (D01)	MAB1679	Nipbl monoclonal antibody, clone KT54
H00000010-D01	NAT2 MaxPab rabbit polyclonal antibody (D01)	MAB1026	Nitrotyrosine monoclonal antibody, clone 1F1
MAB5198	Natural Killer Cell Activation Structures monoclonal antibody, clone ANK61	H00085409-D01	NKD2 MaxPab rabbit polyclonal antibody (D01)
MAB7997	NBN monoclonal antibody, clone 1D7	MAB1964	NKR-P1 monoclonal antibody, clone 1C10
MAB9068	N-cadherin monoclonal antibody, clone 5j46h	H00001482-M04	NKX2-5 monoclonal antibody (M04), clone 3C1
H00004684-D01	NCAM1 MaxPab rabbit polyclonal antibody (D01)	H00001482-M05	NKX2-5 monoclonal antibody (M05), clone 4B11
MAB9132	Ncam1 monoclonal antibody, clone 465c4	H00022829-M02	NLGN4Y monoclonal antibody (M02), clone 2F7
MAB9061	NCAM1 monoclonal antibody, clone 85j46	PAB0258	NLRC4 polyclonal antibody
MAB0894	NCAM1 monoclonal antibody, clone MEM-188	MAB3372	NLRP1 monoclonal antibody, clone Nalpy1-4
MAB2727	NCK2 monoclonal antibody, clone AF3F11	PAB0261	NLRP1 polyclonal antibody
MAB2726	NCK2 monoclonal antibody, clone AF6D7	PAB0262	NLRP1 polyclonal antibody
PAB12541	NCL polyclonal antibody	PAB0269	NLRP12 polyclonal antibody
PAB8877	NCL polyclonal antibody	PAB0267	NLRP2 polyclonal antibody
MAB7998	NCOA1 monoclonal antibody, clone 1135	MAB3724	NLRP3 monoclonal antibody, clone Nalpy3-b
PAB10282	NCOA3 polyclonal antibody	PAB0268	NLRP4 polyclonal antibody
MAB8270	NCOR2 monoclonal antibody	H00004830-M02	NME1 monoclonal antibody (M02), clone 1D7
MAB8000	NDC80 monoclonal antibody, clone 9G3.23	H00004831-D01	NME2 MaxPab rabbit polyclonal antibody (D01)
H00004692-M05	NDN monoclonal antibody (M05), clone 3B9	H00008382-D01	NME5 MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00008382-M04	NME5 monoclonal antibody (M04), clone 4F5	H00011162-D01	NUDT6 MaxPab rabbit polyclonal antibody (D01)
H00009111-D01	NMI MaxPab rabbit polyclonal antibody (D01)	H00004926-M01	NUMA1 monoclonal antibody (M01), clone 1C5
H00004837-M03	NNMT monoclonal antibody (M03), clone 2F2	PAB12361	NUP153 polyclonal antibody
PAB15596	NOB1 polyclonal antibody	H00004928-D01	NUP98 MaxPab rabbit polyclonal antibody (D01)
PAB0266	NOD1 polyclonal antibody	H00115861-D01	NXNL1 MaxPab rabbit polyclonal antibody (D01)
H00064127-D01	NOD2 MaxPab rabbit polyclonal antibody (D01)	H00004938-M01	OAS1 monoclonal antibody (M01), clone 2D4
MAB2303	NOD2 monoclonal antibody, clone 2D9	H00004950-D01	OCLN MaxPab rabbit polyclonal antibody (D01)
PAB0265	NOD2 polyclonal antibody	H00004950-M01	OCLN monoclonal antibody (M01), clone 1G7
H00009241-D01	NOG MaxPab rabbit polyclonal antibody (D01)	H00004953-M01	ODC1 monoclonal antibody (M01), clone 2G5
PAB8871	NONO polyclonal antibody	H00004968-D01	OGG1 MaxPab rabbit polyclonal antibody (D01)
MAB3423	Nonspecific Cytotoxic Cell monoclonal antibody, clone 5C.6	H00004969-D01	OGN MaxPab rabbit polyclonal antibody (D01)
MAB3338	NOP1 monoclonal antibody, clone 28F2	H00011339-D01	OIP5 MaxPab rabbit polyclonal antibody (D01)
H00051491-D01	NOP16 MaxPab rabbit polyclonal antibody (D01)	H00011339-M01	OIP5 monoclonal antibody (M01), clone 3G10-1D5
MAB3347	NOP58 monoclonal antibody, clone 34B12	H00116448-M04	OLIG1 monoclonal antibody (M04), clone 2B11
MAB3346	NOP58/SIK1 monoclonal antibody, clone 37C12	H00116448-M05	OLIG1 monoclonal antibody (M05), clone 3B3
MAB1627	NOS2 monoclonal antibody, clone K13-A	H00116448-M06	OLIG1 monoclonal antibody (M06), clone 3E3
PAB12680	NOS3 polyclonal antibody	H00116448-M08	OLIG1 monoclonal antibody (M08), clone 2A4
MAB9240	NOTCH2 monoclonal antibody, clone 55h6	PAB16939	Olig1 polyclonal antibody
H00004856-D01	NOV MaxPab rabbit polyclonal antibody (D01)	H00010215-D01	OLIG2 MaxPab rabbit polyclonal antibody (D01)
H00004856-M02	NOV monoclonal antibody (M02), clone 2G8	H00010215-M02	OLIG2 monoclonal antibody (M02), clone 3D7
H00010811-D01	NOXA1 MaxPab rabbit polyclonal antibody (D01)	H00010215-M05	OLIG2 monoclonal antibody (M05), clone 2A8
H00004862-D01	NPAS2 MaxPab rabbit polyclonal antibody (D01)	H00010215-M06	OLIG2 monoclonal antibody (M06), clone 2B11
PAB12093	NPC1 polyclonal antibody	PAB16940	Olig2 polyclonal antibody
H00010577-M11	NPC2 monoclonal antibody (M11), clone 1F4	MAB7372	OLLAS tag monoclonal antibody, clone L2
H00004867-D01	NPHP1 MaxPab rabbit polyclonal antibody (D01)	H00004973-D01	OLR1 MaxPab rabbit polyclonal antibody (D01)
MAB1638	NPHS1 monoclonal antibody, clone Y17-R	H00004958-D01	OMD MaxPab rabbit polyclonal antibody (D01)
H00007827-D01	NPHS2 MaxPab rabbit polyclonal antibody (D01)	MAB1742	OPRD1 monoclonal antibody, clone 1E7
H00055666-D01	NPLOC4 MaxPab rabbit polyclonal antibody (D01)	MAB1740	OPRK1 monoclonal antibody, clone 5E3
PAB4974	NPPA polyclonal antibody	H00010133-D01	OPTN MaxPab rabbit polyclonal antibody (D01)
PAB5081	Nppa polyclonal antibody	MAB2489	ORC1L monoclonal antibody, clone 7A7
PAB5131	NPPA polyclonal antibody	MAB2471	ORC3L monoclonal antibody, clone 1D6
PAB8975	NPPA polyclonal antibody	H00023594-D01	ORC6L MaxPab rabbit polyclonal antibody (D01)
H00004879-M01	NPPB monoclonal antibody (M01), clone 2D11	H00023594-D02	ORC6L MaxPab rabbit polyclonal antibody (D02)
PAB4975	NPPB polyclonal antibody	H00023594-M04	ORC6L monoclonal antibody (M04), clone 3C4
PAB5082	Nppb polyclonal antibody	H00005004-M01	ORM1 monoclonal antibody (M01), clone 2F9-1F10
PAB5083	Nppb polyclonal antibody	PAB12082	OS9 polyclonal antibody
PAB8976	NPPB polyclonal antibody	H00005008-D01	OSM MaxPab rabbit polyclonal antibody (D01)
PAB4978	NPPC polyclonal antibody	H00009180-D01	OSMR MaxPab rabbit polyclonal antibody (D01)
H00004885-D01	NPTX2 MaxPab rabbit polyclonal antibody (D01)	H00078990-M14	OTUB2 monoclonal antibody (M14), clone 2E7
H00001728-M01	NQO1 monoclonal antibody (M01), clone 1E3-A6	H00005013-D01	OTX1 MaxPab rabbit polyclonal antibody (D01)
MAB2302	NQO1 monoclonal antibody, clone A180	PAB10298	Ovalbumin polyclonal antibody
H00008856-D01	NR112 MaxPab rabbit polyclonal antibody (D01)	H00005018-D01	OXA1L MaxPab rabbit polyclonal antibody (D01)
H00009970-D01	NR113 MaxPab rabbit polyclonal antibody (D01)	H00009943-M05	OXSR1 monoclonal antibody (M05), clone 5E1
H00007026-D01	NR2F2 MaxPab rabbit polyclonal antibody (D01)	H00009943-M07	OXSR1 monoclonal antibody (M07), clone 2E9
MAB0119	NR3C1 monoclonal antibody, clone m2F8	H00009943-M08	OXSR1 monoclonal antibody (M08), clone 1F6
H00002649-D01	NR6A1 MaxPab rabbit polyclonal antibody (D01)	H00009943-M11	OXSR1 monoclonal antibody (M11), clone 4D12
H00003084-D01	NRG1 MaxPab rabbit polyclonal antibody (D01)	H00009943-M14	OXSR1 monoclonal antibody (M14), clone 1B10
MAB1890	Nrg1 monoclonal antibody, clone 7D5	MAB1997	p13 Myeloid Marker monoclonal antibody, clone BM3
PAB12639	NRP1 (phospho T916) polyclonal antibody	MAB9479	P46 cell surface protein monoclonal antibody, clone 4756h
H00008829-M05	NRP1 monoclonal antibody (M05), clone 1B3	H00005034-D01	P4HB MaxPab rabbit polyclonal antibody (D01)
PAB12638	NRP1 polyclonal antibody	PAB15416	P4hb polyclonal antibody
PAB12718	NRP1 polyclonal antibody	MAB1982	p73 alpha/beta monoclonal antibody, clone ER-15
H00050814-D01	NSDHL MaxPab rabbit polyclonal antibody (D01)	MAB9477	p80 monoclonal antibody, clone 4t56
MAB1966	NSF monoclonal antibody, clone NSF-1	H00005036-D01	PA2G4 MaxPab rabbit polyclonal antibody (D01)
PAB8875	Nsf polyclonal antibody	MAB2465	PABPC1 monoclonal antibody, clone 10E10
H00004909-D01	NTF4 MaxPab rabbit polyclonal antibody (D01)	H00005042-D01	PABPC3 MaxPab rabbit polyclonal antibody (D01)
PAB17918	NTRK1 (phospho Y490) polyclonal antibody	PAB14100	PAF1 polyclonal antibody
PAB9064	NTRK1 polyclonal antibody	PAB14101	Paf1 polyclonal antibody
PAB11334	Ntrk3 polyclonal antibody	MAB0934	PAG1 monoclonal antibody, clone PAG-C1
MAB8178	Nuclei monoclonal antibody, clone 235-1	PAB1257	PAG1 polyclonal antibody
MAB9429	Nucleolar protein monoclonal antibody, clone 090j6	PAB8889	PAK1 polyclonal antibody
MAB9249	Nucleolar protein monoclonal antibody, clone 6a84	PAB12642	PAK4 (phospho S474) polyclonal antibody
H00004521-D01	NUDT1 MaxPab rabbit polyclonal antibody (D01)	PAB12721	PAK4 polyclonal antibody
H00000318-D01	NUDT2 MaxPab rabbit polyclonal antibody (D01)	PAB2303	PAK6 polyclonal antibody
H00011163-D01	NUDT4 MaxPab rabbit polyclonal antibody (D01)	PAB2302	PAK7 polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB7377	PALLD monoclonal antibody, clone 1E6	H00064236-D01	PDLIM2 MaxPab rabbit polyclonal antibody (D01)
PAB10305	Pancreatic islets polyclonal antibody	H00009260-D01	PDLIM7 MaxPab rabbit polyclonal antibody (D01)
PAB15540	Pancreatic Trypsinogen polyclonal antibody	H00009260-M01	PDLIM7 monoclonal antibody (M01), clone 2D6
H00053354-D01	PANK1 MaxPab rabbit polyclonal antibody (D01)	H00005170-D01	PDPK1 MaxPab rabbit polyclonal antibody (D01)
MAB0583	PAPPA monoclonal antibody, clone 10A5	PAB8898	PDPK1 polyclonal antibody
MAB0579	PAPPA monoclonal antibody, clone 10E1	H00005173-M01	PDYN monoclonal antibody (M01), clone 2E12
MAB0580	PAPPA monoclonal antibody, clone 10E2	PAB4981	PDYN polyclonal antibody
MAB0585	PAPPA monoclonal antibody, clone 10H9	PAB5133	PDYN polyclonal antibody
MAB0587	PAPPA monoclonal antibody, clone 11E4	H00079849-D01	PDZD3 MaxPab rabbit polyclonal antibody (D01)
MAB0586	PAPPA monoclonal antibody, clone 18A10	H00005175-M01	PECAM1 monoclonal antibody (M01), clone 1D2-1A5
MAB0584	PAPPA monoclonal antibody, clone 3C8	MAB3230	Pecam1 monoclonal antibody, clone 390
MAB0582	PAPPA monoclonal antibody, clone 4G11	MAB5715	Pecam1 monoclonal antibody, clone 390 (APC)
MAB0581	PAPPA monoclonal antibody, clone 5H9	MAB5713	Pecam1 monoclonal antibody, clone 390 (Biotin)
MAB0588	PAPPA monoclonal antibody, clone 7A6	MAB5712	Pecam1 monoclonal antibody, clone 390 (FITC)
H00009060-D01	PAPSS2 MaxPab rabbit polyclonal antibody (D01)	MAB5714	Pecam1 monoclonal antibody, clone 390 (PE)
H00050855-D01	PARD6A MaxPab rabbit polyclonal antibody (D01)	MAB5716	Pecam1 monoclonal antibody, clone 390 (PE/Cy7)
H00005071-D01	PARK2 MaxPab rabbit polyclonal antibody (D01)	MAB5210	PECAM1 monoclonal antibody, clone LCI-4
PAB1105	PARK2 polyclonal antibody	MAB7099	Pecam1 monoclonal antibody, clone MEC13.3
PAB12080	PARK7 polyclonal antibody	H00010455-D01	PECI MaxPab rabbit polyclonal antibody (D01)
MAB5193	PARP1 monoclonal antibody, clone A6.4.12	PAB5130	PENK polyclonal antibody
PAB13692	Parp2 polyclonal antibody	H00005184-D01	PEPD MaxPab rabbit polyclonal antibody (D01)
PAB8893	PAWR polyclonal antibody	H00005184-M01	PEPD monoclonal antibody (M01), clone 1D5-H3
H00005080-D01	PAX6 MaxPab rabbit polyclonal antibody (D01)	H00008863-D01	PER3 MaxPab rabbit polyclonal antibody (D01)
H00005089-D01	PBX2 MaxPab rabbit polyclonal antibody (D01)	H00008863-M01	PER3 monoclonal antibody (M01), clone 3A7
H00005094-D01	PCBP2 MaxPab rabbit polyclonal antibody (D01)	H00005192-D01	PEX10 MaxPab rabbit polyclonal antibody (D01)
H00005096-D01	PCCB MaxPab rabbit polyclonal antibody (D01)	H00005196-M01	PF4 monoclonal antibody (M01), clone 3F6
H00005097-M03	PCDH1 monoclonal antibody (M03), clone 2E8	H00005208-D01	PFKFB2 MaxPab rabbit polyclonal antibody (D01)
H00005105-D01	PCK1 MaxPab rabbit polyclonal antibody (D01)	H00005209-D01	PFKFB3 MaxPab rabbit polyclonal antibody (D01)
H00005105-M01	PCK1 monoclonal antibody (M01), clone 3E4	H00005209-M08	PFKFB3 monoclonal antibody (M08), clone 3F3
H00005106-D01	PCK2 MaxPab rabbit polyclonal antibody (D01)	H00005223-M01	PGAM1 monoclonal antibody (M01), clone 2G1-A6
H00005111-M02	PCNA monoclonal antibody (M02), clone 1G7	H00005225-D01	PGC MaxPab rabbit polyclonal antibody (D01)
PAB10179	PCNA polyclonal antibody	H00025796-M01	PGLS monoclonal antibody (M01), clone 2D9
PAB8896	PCNA polyclonal antibody	H00005236-M01	PGM1 monoclonal antibody (M01), clone 3B8-H4
MAB8228	PCSK9 (prodomain) monoclonal antibody, clone KS-2C8	H00005238-M01	PGM3 monoclonal antibody (M01), clone 1E2-1B12
MAB8229	PCSK9 monoclonal antibody, clone KS-4H12	MAB2467	PGR (phospho S162) monoclonal antibody, clone 1064/E2
PAB17045	PCSK9 polyclonal antibody	MAB1975	PGR monoclonal antibody, clone hPRa2+hPRa3
H00005127-D01	PCTK1 MaxPab rabbit polyclonal antibody (D01)	PAB12723	PGR polyclonal antibody
H00005127-M02	PCTK1 monoclonal antibody (M02), clone 2E5	PAB9007	PHB polyclonal antibody
H00005134-D01	PDCD2 MaxPab rabbit polyclonal antibody (D01)	H00001911-D01	PHC1 MaxPab rabbit polyclonal antibody (D01)
MAB3890	PDCD6IP monoclonal antibody, clone B234M	H00055274-D01	PHF10 MaxPab rabbit polyclonal antibody (D01)
PAB12362	PDE10A polyclonal antibody	H00026227-M01	PHGDH monoclonal antibody (M01), clone 4A3-1D6
PAB12366	PDE11A polyclonal antibody	H00005261-D01	PHKG2 MaxPab rabbit polyclonal antibody (D01)
PAB12350	PDE3A polyclonal antibody	MAB2545	Phosphoserine monoclonal antibody, clone A8G9
PAB8960	PDE3A polyclonal antibody	PAB10351	Phosphoserine polyclonal antibody
PAB12351	PDE3B polyclonal antibody	PAB13796	Phosphoserine polyclonal antibody
PAB12353	PDE4B polyclonal antibody	PAB16619	Phosphoserine polyclonal antibody
PAB12354	PDE4D polyclonal antibody	MAB6640	Phosphothreonine monoclonal antibody, clone 18F9
H00009659-D01	PDE4DIP MaxPab rabbit polyclonal antibody (D01)	PAB10352	Phosphothreonine polyclonal antibody
H00009659-M02	PDE4DIP monoclonal antibody (M02), clone 1C4	PAB16622	Phosphothreonine polyclonal antibody
PAB12359	PDE5A polyclonal antibody	MAB6670	Phosphotyrosine monoclonal antibody, clone G104
PAB8978	Pde7b polyclonal antibody	MAB1971	Phosphotyrosine monoclonal antibody, clone PY20
H00005152-D01	PDE9A MaxPab rabbit polyclonal antibody (D01)	PAB10353	Phosphotyrosine polyclonal antibody
PAB12356	PDE9A polyclonal antibody	PAB13650	Phosphotyrosine polyclonal antibody
H00005154-D01	PDGFA MaxPab rabbit polyclonal antibody (D01)	PAB13767	Phosphotyrosine polyclonal antibody
H00080310-D01	PDGFD MaxPab rabbit polyclonal antibody (D01)	PAB16629	Phosphotyrosine polyclonal antibody
H00080310-M09	PDGFD monoclonal antibody (M09), clone 4H2	H00000401-D01	PHOX2A MaxPab rabbit polyclonal antibody (D01)
H00005159-M08	PDGFRB monoclonal antibody (M08), clone 4C12	H00055361-D01	PI4K2A MaxPab rabbit polyclonal antibody (D01)
MAB7066	PDGFRB monoclonal antibody, clone PDGFR-B2	H00008554-D01	PIAS1 MaxPab rabbit polyclonal antibody (D01)
H00008050-D01	PDHX MaxPab rabbit polyclonal antibody (D01)	H00009063-D01	PIAS2 MaxPab rabbit polyclonal antibody (D01)
PAB11963	Pdia2 polyclonal antibody	MAB2257	PIAS4 polyclonal antibody
PAB12587	PDIA2 polyclonal antibody	H00009463-D01	PICK1 MaxPab rabbit polyclonal antibody (D01)
H00002923-D01	PDIA3 MaxPab rabbit polyclonal antibody (D01)	H00005283-D01	PIGH MaxPab rabbit polyclonal antibody (D01)
MAB6645	PDIA3 monoclonal antibody, clone MaP.ERP57	H00055011-D01	PIH1D1 MaxPab rabbit polyclonal antibody (D01)
H00005165-D01	PDK3 MaxPab rabbit polyclonal antibody (D01)	PAB17948	Pim1 polyclonal antibody
H00009124-D01	PDLIM1 MaxPab rabbit polyclonal antibody (D01)	H00011040-D01	PIM2 MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB2305	PIM2 polyclonal antibody	H00051082-D01	POLR1D MaxPab rabbit polyclonal antibody (D01)
PAB11244	PIN1 polyclonal antibody	MAB0120	POLR2A monoclonal antibody
MAB7762	PIP2 monoclonal antibody, clone 2C11	MAB2439	POLR2A monoclonal antibody, clone 4H8
H00008396-D01	PIP4K2B MaxPab rabbit polyclonal antibody (D01)	H00005434-D01	POLR2E MaxPab rabbit polyclonal antibody (D01)
H00079837-D01	PIP4K2C MaxPab rabbit polyclonal antibody (D01)	H00005437-D01	POLR2H MaxPab rabbit polyclonal antibody (D01)
H00023396-M02	PIP5K1C monoclonal antibody (M02), clone 2E9	MAB7056	Polyoma virus early antigen monoclonal antibody, clone PyC
H00051268-M10	PIPOX monoclonal antibody (M10), clone 3D1	MAB7375	Polyoma virus, Large T antigen monoclonal antibody, clone PyLT
H00005307-M01	PITX1 monoclonal antibody (M01), clone 5G4	MAB7374	Polyoma virus, Medium T antigen monoclonal antibody, clone PyMT
MAB0843	PIWIL2 monoclonal antibody, clone 11C11	H00005443-D01	POMC MaxPab rabbit polyclonal antibody (D01)
MAB5593	PKC monoclonal antibody, clone MC5	MAB8230	PON1 monoclonal antibody, clone KS-5F5
H00005313-D01	PKLR MaxPab rabbit polyclonal antibody (D01)	MAB2743	PON2 monoclonal antibody, clone AF3E6
PAB0271	PKN1 polyclonal antibody	PAB11972	Por polyclonal antibody
H00063876-D01	PKNOX2 MaxPab rabbit polyclonal antibody (D01)	PAB12771	POU2F1 polyclonal antibody
H00051365-D01	PLA1A MaxPab rabbit polyclonal antibody (D01)	PAB8878	POU2F1 polyclonal antibody
H00023659-D01	PLA2G15 MaxPab rabbit polyclonal antibody (D01)	H00005452-D01	POU2F2 MaxPab rabbit polyclonal antibody (D01)
H00011145-D01	PLA2G16 MaxPab rabbit polyclonal antibody (D01)	PAB12772	POU2F2 polyclonal antibody
H00005327-D01	PLAT MaxPab rabbit polyclonal antibody (D01)	H00005459-D01	POU4F3 MaxPab rabbit polyclonal antibody (D01)
PAB5112	Plat polyclonal antibody	PAB10296	POU5F1 polyclonal antibody
PAB14920	PLAU polyclonal antibody	PAB12773	POU5F1 polyclonal antibody
PAB14928	PLAU polyclonal antibody	PAB8880	Pou5f1 polyclonal antibody
PAB14929	PLAU polyclonal antibody	H00005463-D01	POU6F1 MaxPab rabbit polyclonal antibody (D01)
PAB5098	PLAU polyclonal antibody	PAB10354	PP2C polyclonal antibody
PAB5114	Plau polyclonal antibody	PAB10355	PP2C polyclonal antibody
PAB5117	PLAU polyclonal antibody	PAB5095	PPARA polyclonal antibody
PAB5123	Plau polyclonal antibody	H00005471-D01	PPAT MaxPab rabbit polyclonal antibody (D01)
PAB7929	PLCG1 (phospho Y775) polyclonal antibody	H00005475-D01	PPEF1 MaxPab rabbit polyclonal antibody (D01)
MAB0700	PLCG1 monoclonal antibody, clone 1F1	PAB15689	Ppfia3 polyclonal antibody
MAB0701	PLCG1 monoclonal antibody, clone 2B1	H00005481-D01	PPID MaxPab rabbit polyclonal antibody (D01)
MAB1377	PLCG1 monoclonal antibody, clone M156	H00009360-D01	PPIG MaxPab rabbit polyclonal antibody (D01)
H00059338-D01	PLEKHA1 MaxPab rabbit polyclonal antibody (D01)	MAB1050	PPM1A monoclonal antibody, clone p6c7
MAB0797	PLG monoclonal antibody, clone 44A1	PAB13655	PPM1A/PPM1B polyclonal antibody
MAB0796	PLG monoclonal antibody, clone 9A1	H00009647-D01	PPM1F MaxPab rabbit polyclonal antibody (D01)
PAB5072	Plg polyclonal antibody	H00009647-M02	PPM1F monoclonal antibody (M02), clone 1G10
PAB5086	Plg polyclonal antibody	H00333926-D01	PPM1J MaxPab rabbit polyclonal antibody (D01)
H00005347-D01	PLK1 MaxPab rabbit polyclonal antibody (D01)	PAB10410	PPP1CA polyclonal antibody
H00005347-M01	PLK1 monoclonal antibody (M01), clone 2G12	PAB10412	PPP1CB polyclonal antibody
MAB0801	PLK1 monoclonal antibody, clone 18D4	PAB10414	PPP1CC polyclonal antibody
H00010733-D01	PLK4 MaxPab rabbit polyclonal antibody (D01)	PAB13653	PPP1CC polyclonal antibody
H00005352-D01	PLOD2 MaxPab rabbit polyclonal antibody (D01)	H00153743-D01	PPP1R2P3 MaxPab rabbit polyclonal antibody (D01)
H00005354-D01	PLP1 MaxPab rabbit polyclonal antibody (D01)	H00005507-D03	PPP1R3C MaxPab rabbit polyclonal antibody (D03)
H00005359-D01	PLSCR1 MaxPab rabbit polyclonal antibody (D01)	PAB10421	PPP2CA polyclonal antibody
H00005359-M02	PLSCR1 monoclonal antibody (M02), clone 3B4	MAB2677	PPP2CA/PPP2CB monoclonal antibody, clone 13D4
MAB6482	Plscr1 monoclonal antibody, clone 13A6 [TEC-23]	MAB2678	PPP2CA/PPP2CB monoclonal antibody, clone 20G6
H00057048-D01	PLSCR3 MaxPab rabbit polyclonal antibody (D01)	PAB10416	PPP2R1A polyclonal antibody
PAB7879	PLXNA1 polyclonal antibody	PAB10417	PPP2R1A polyclonal antibody
MAB3726	Pml monoclonal antibody, clone 36-1-104	PAB13656	PPP2R1A polyclonal antibody
H00005373-D01	PMM2 MaxPab rabbit polyclonal antibody (D01)	H00005519-D01	PPP2R1B MaxPab rabbit polyclonal antibody (D01)
H00010654-D01	PMVK MaxPab rabbit polyclonal antibody (D01)	H00005521-M02	PPP2R2B monoclonal antibody (M02), clone 1F3
H00010654-M07	PMVK monoclonal antibody (M07), clone 2B8	H00005521-M23	PPP2R2B monoclonal antibody (M23), clone 1B10
H00010654-M07A	PMVK monoclonal antibody (M07A), clone 2B8	H00028227-M08	PPP2R3B monoclonal antibody (M08), clone 2E12
H00009240-D01	PNMA1 MaxPab rabbit polyclonal antibody (D01)	H00028227-M14	PPP2R3B monoclonal antibody (M14), clone 1F4
H00005409-D01	PNMT MaxPab rabbit polyclonal antibody (D01)	PAB8471	Ppp2r3c polyclonal antibody
H00005409-M06	PNMT monoclonal antibody (M06), clone 3G4	MAB2390	PPP2R5C monoclonal antibody, clone TQ11-1G6
PAB8557	PNPLA2 polyclonal antibody	PAB8605	PPP3CA polyclonal antibody
H00087178-D01	PNPT1 MaxPab rabbit polyclonal antibody (D01)	H00005535-D01	PPP3R2 MaxPab rabbit polyclonal antibody (D01)
H00057645-D01	POGK MaxPab rabbit polyclonal antibody (D01)	H00023532-D01	PRAME MaxPab rabbit polyclonal antibody (D01)
H00023649-D01	POLA2 MaxPab rabbit polyclonal antibody (D01)	MAB9438	prc monoclonal antibody, clone 1655y
H00023649-M02	POLA2 monoclonal antibody (M02), clone 2A8	H00000639-D01	PRDM1 MaxPab rabbit polyclonal antibody (D01)
H00005423-D01	POLB MaxPab rabbit polyclonal antibody (D01)	MAB1151	PRDM2 monoclonal antibody, clone 33AT1045
MAB3437	Polb monoclonal antibody, clone 18S	PAB2383	Prdm2 polyclonal antibody
PAB15576	Polb polyclonal antibody	MAB0714	PRDX1 monoclonal antibody, clone 13E7
MAB6748	POLD3 monoclonal antibody, clone 2A1C11	MAB0713	PRDX1 monoclonal antibody, clone 2A4
MAB3438	POLE monoclonal antibody, clone 93H3A	MAB0846	PRDX1 monoclonal antibody, clone 3G5
H00011232-D01	POLG2 MaxPab rabbit polyclonal antibody (D01)	MAB0690	PRDX1 monoclonal antibody, clone 9D2
H00051426-M02	POLK monoclonal antibody (M02), clone 4F12	MAB0716	PRDX2 monoclonal antibody, clone 12B1

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB0715	PRDX2 monoclonal antibody, clone 9A1	H00143471-D01	PSMA8 MaxPab rabbit polyclonal antibody (D01)
H00010935-M01	PRDX3 monoclonal antibody (M01), clone 1E3-B2	H00005696-D01	PSMB8 MaxPab rabbit polyclonal antibody (D01)
MAB0696	PRDX3 monoclonal antibody, clone 2B11	H00005696-M01	PSMB8 monoclonal antibody (M01), clone 1B3
MAB0698	PRDX3 monoclonal antibody, clone 4G10	PAB9089	PSMC5 polyclonal antibody
MAB2725	PRDX3 monoclonal antibody, clone AF17D2FF	MAB7301	Psmc4 monoclonal antibody, clone AH1.1
PAB1200	PRDX3 polyclonal antibody	H00005711-D01	PSMD5 MaxPab rabbit polyclonal antibody (D01)
H00010549-D01	PRDX4 MaxPab rabbit polyclonal antibody (D01)	H00005713-D01	PSMD7 MaxPab rabbit polyclonal antibody (D01)
MAB0673	PRDX4 monoclonal antibody, clone 1A1	H00005721-M02	PSME2 monoclonal antibody (M02), clone 1G4
MAB0674	PRDX4 monoclonal antibody, clone 3A1	H00005725-M01	PTBP1 monoclonal antibody (M01), clone 3H8
MAB0677	PRDX4 monoclonal antibody, clone 7A1	PAB12751	PTEN (phospho S385) polyclonal antibody
PAB1189	PRDX4 polyclonal antibody	H00005728-M13	PTEN monoclonal antibody (M13), clone 1C3
H00025824-D01	PRDX5 MaxPab rabbit polyclonal antibody (D01)	PAB11075	PTEN polyclonal antibody
MAB0670	PRDX5 monoclonal antibody, clone 3F11	PAB12724	PTEN polyclonal antibody
MAB0669	PRDX5 monoclonal antibody, clone 4C3	PAB12725	PTEN polyclonal antibody
MAB0676	PRDX6 monoclonal antibody, clone 1A11	PAB13619	Ptgs2 polyclonal antibody
MAB0712	PRDX6 monoclonal antibody, clone 5E1	H00009536-D01	PTGES MaxPab rabbit polyclonal antibody (D01)
MAB0741	PRDX6 monoclonal antibody, clone 6H5	H00009536-M04	PTGES monoclonal antibody (M04), clone 2B9
PAB1191	PRDX6 polyclonal antibody	H00080142-D01	PTGES2 MaxPab rabbit polyclonal antibody (D01)
H00005551-D01	PRF1 MaxPab rabbit polyclonal antibody (D01)	MAB6639	PTGES3 monoclonal antibody, clone JJ6
H00005551-M04	PRF1 monoclonal antibody (M04), clone 3B4	H00005740-D01	PTGIS MaxPab rabbit polyclonal antibody (D01)
MAB3710	Prf1 monoclonal antibody, clone CB5.4	MAB3798	PTGIS monoclonal antibody
MAB3711	Prf1 monoclonal antibody, clone CE2.10	H00005742-D01	PTGS1 MaxPab rabbit polyclonal antibody (D01)
MAB1764	PRF1 monoclonal antibody, clone deltaG9	H00005743-D01	PTGS2 MaxPab rabbit polyclonal antibody (D01)
PAB17882	Prf1 polyclonal antibody	H00005741-M02	PTH monoclonal antibody (M02), clone 2E4
PAB12747	PRKAA1 (phospho S79) polyclonal antibody	H00005741-M04	PTH monoclonal antibody (M04), clone 4C12
PAB12602	PRKAA1 (phospho T172) polyclonal antibody	H00005741-M05	PTH monoclonal antibody (M05), clone 4A2
PAB14475	PRKAA1 polyclonal antibody	H00005741-M06	PTH monoclonal antibody (M06), clone 1E1
PAB8919	PRKACA polyclonal antibody	H00005741-M07	PTH monoclonal antibody (M07), clone 4D7
H00051422-D01	PRKAG2 MaxPab rabbit polyclonal antibody (D01)	H00005741-M08	PTH monoclonal antibody (M08), clone 4E10
PAB7910	PRKCA (phospho S657/Y658) polyclonal antibody	PAB5103	PTH polyclonal antibody
MAB1629	PRKCA monoclonal antibody, clone G13-R	H00005747-D01	PTK2 MaxPab rabbit polyclonal antibody (D01)
MAB1630	PRKCB monoclonal antibody, clone A10-F	H00005747-M01	PTK2 monoclonal antibody (M01), clone 2C3
MAB1631	PRKCB monoclonal antibody, clone N11-S	MAB0839	PTK2 monoclonal antibody, clone 63D5
MAB1633	PRKCD monoclonal antibody, clone K15-K	PAB1237	PTK2 polyclonal antibody
MAB1634	PRKCE monoclonal antibody, clone G13-R	PAB9880	Ptk2 polyclonal antibody
MAB1632	PRKCG monoclonal antibody, clone A12-H	H00005756-M02	PTK9 monoclonal antibody (M02), clone 1E2
MAB1367	PRKCG monoclonal antibody, clone M110	MAB8069	PTMA monoclonal antibody, clone 6D709
H00005589-D01	PRKCSH MaxPab rabbit polyclonal antibody (D01)	H00005770-D01	PTPN1 MaxPab rabbit polyclonal antibody (D01)
MAB3430	PRKDC monoclonal antibody, clone 18-2	MAB0795	PTPN11 monoclonal antibody, clone 5F2
MAB3433	PRKDC monoclonal antibody, clone 18-2, 25-4, 42-psc	MAB5240	PTPN11 monoclonal antibody, clone SIT-01
MAB3431	PRKDC monoclonal antibody, clone 25-4	PAB9023	PTPN11 polyclonal antibody
H00008575-D01	PRKRA MaxPab rabbit polyclonal antibody (D01)	PAB0256	PTPN13 polyclonal antibody
H00008575-M01	PRKRA monoclonal antibody (M01), clone 1B9-1A7	H00005771-D01	PTPN2 MaxPab rabbit polyclonal antibody (D01)
H00005613-D01	PRKX MaxPab rabbit polyclonal antibody (D01)	MAB5684	PTPN2 monoclonal antibody, clone 3E2
H00005613-M01	PRKX monoclonal antibody (M01), clone 1H7	MAB5683	PTPN2 monoclonal antibody, clone 6F3
MAB7060	PRMT1 monoclonal antibody, clone PRMT1-171	H00084867-M03	PTPN5 monoclonal antibody (M03), clone 2H5
H00055170-D01	PRMT6 MaxPab rabbit polyclonal antibody (D01)	MAB2017	PTPN6 monoclonal antibody, clone 1SH01
H00054496-D01	PRMT7 MaxPab rabbit polyclonal antibody (D01)	MAB1379	PTPN6 monoclonal antibody, clone M160
H00010544-D01	PROCR MaxPab rabbit polyclonal antibody (D01)	H00005778-D01	PTPN7 MaxPab rabbit polyclonal antibody (D01)
PAB5159	Progesterone polyclonal antibody	MAB2429	Ptprb monoclonal antibody, clone 122.2
PAB12663	PROM1 polyclonal antibody	MAB3233	Ptprc exon B-dependent form monoclonal antibody, clone C363.16A
H00005626-D01	PROP1 MaxPab rabbit polyclonal antibody (D01)	MAB5732	Ptprc exon B-dependent form monoclonal antibody, clone C363.16A (Biotin)
H00005627-D01	PROS1 MaxPab rabbit polyclonal antibody (D01)	MAB5731	Ptprc exon B-dependent form monoclonal antibody, clone C363.16A (FITC)
MAB9093	Protein tyrosine phosphatase monoclonal antibody, clone 8qw7	MAB5733	Ptprc exon B-dependent form monoclonal antibody, clone C363.16A (PE)
H00005639-D01	PRRG2 MaxPab rabbit polyclonal antibody (D01)	MAB5734	Ptprc exon B-dependent form monoclonal antibody, clone C363.16A (SPRD)
H00051450-D01	PRRX2 MaxPab rabbit polyclonal antibody (D01)	MAB6210	Ptprc monoclonal antibody, clone 104
H00005644-D01	PRSS1 MaxPab rabbit polyclonal antibody (D01)	MAB4106	Ptprc monoclonal antibody, clone 13/2.3
MAB1642	PRSS1 monoclonal antibody, clone G15-B	MAB8062	PTPRC monoclonal antibody, clone 2Q1375
H00010942-M01	PRSS2 monoclonal antibody (M01), clone 2E10	MAB6209	Ptprc monoclonal antibody, clone A20
H00005652-M11A	PRSS8 monoclonal antibody (M11A), clone 3C4	MAB6955	PTPRC monoclonal antibody, clone C11
H00005660-M01	PSAP monoclonal antibody (M01), clone 1D1-C12	MAB6467	Ptprc monoclonal antibody, clone EM-05
PAB12733	PSAP polyclonal antibody	MAB5626	PTPRC monoclonal antibody, clone F10-89-4
PAB11400	<i>Pseudomonas</i> NAM-DH polyclonal antibody	MAB6067	PTPRC monoclonal antibody, clone F10-89-4 (APC)
H00005688-D01	PSMA7 MaxPab rabbit polyclonal antibody (D01)	MAB6065	PTPRC monoclonal antibody, clone F10-89-4 (Biotin)
H00005688-M01	PSMA7 monoclonal antibody (M01), clone 1A10-3G12	MAB6069	PTPRC monoclonal antibody, clone F10-89-4 (Cy5)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB6064	PTPRC monoclonal antibody, clone F10-89-4 (FITC)	PAB9011	RAB5A polyclonal antibody
MAB6066	PTPRC monoclonal antibody, clone F10-89-4 (PE)	H00338382-M01	RAB7B monoclonal antibody (M01), clone 3B3
MAB6070	PTPRC monoclonal antibody, clone F10-89-4 (PE/Cy5.5)	H00023637-D01	RABGAP1 MaxPab rabbit polyclonal antibody (D01)
MAB6068	PTPRC monoclonal antibody, clone F10-89-4 (SPRD)	H00011158-D01	RABL2B MaxPab rabbit polyclonal antibody (D01)
MAB5608	Ptprc monoclonal antibody, clone I3/2.3	PAB15575	RAD18 polyclonal antibody
MAB5940	Ptprc monoclonal antibody, clone I3/2.3 (APC)	PAB1173	rad22 polyclonal antibody
MAB5944	Ptprc monoclonal antibody, clone I3/2.3 (APC/CY7)	MAB8002	RAD50 monoclonal antibody, clone 13B3
MAB5937	Ptprc monoclonal antibody, clone I3/2.3 (Biotin)	PAB12052	RAD50 polyclonal antibody
MAB5936	Ptprc monoclonal antibody, clone I3/2.3 (FITC)	PAB12153	Rad50 polyclonal antibody
MAB5938	Ptprc monoclonal antibody, clone I3/2.3 (PE)	H00005888-M02	RAD51 monoclonal antibody (M02), clone 4C11-B4
MAB5942	Ptprc monoclonal antibody, clone I3/2.3 (PE/Cy5.5)	PAB1172	rad51 polyclonal antibody
MAB5943	Ptprc monoclonal antibody, clone I3/2.3 (PE/Cy7)	PAB1184	RAD51 polyclonal antibody
MAB5939	Ptprc monoclonal antibody, clone I3/2.3 (PE/TXRD)	PAB15573	RAD51 polyclonal antibody
MAB5941	Ptprc monoclonal antibody, clone I3/2.3 (SPRD)	PAB15574	RAD52 polyclonal antibody
MAB0886	PTPRC monoclonal antibody, clone MEM-28	H00025788-D01	RAD54B MaxPab rabbit polyclonal antibody (D01)
MAB0888	PTPRC monoclonal antibody, clone MEM-55	MAB0063	RAD9A monoclonal antibody, clone 93A535
MAB0887	PTPRC monoclonal antibody, clone MEM-56	H00005897-D01	RAG2 MaxPab rabbit polyclonal antibody (D01)
MAB3232	Ptprc monoclonal antibody, clone RA3-6B2	H00005899-M04	RALB monoclonal antibody (M04), clone 4D1
MAB5724	Ptprc monoclonal antibody, clone RA3-6B2 (APC)	H00010928-M02	RALBP1 monoclonal antibody (M02), clone 2A1
MAB5729	Ptprc monoclonal antibody, clone RA3-6B2 (APC/Cy5.5)	H00055103-D01	RALGPS2 MaxPab rabbit polyclonal antibody (D01)
MAB5730	Ptprc monoclonal antibody, clone RA3-6B2 (APC/CY7)	MAB7286	RAN monoclonal antibody, clone ARAN1
MAB5721	Ptprc monoclonal antibody, clone RA3-6B2 (Biotin)	PAB12357	RAN polyclonal antibody
MAB5726	Ptprc monoclonal antibody, clone RA3-6B2 (Cy5)	PAB15577	RANBP2 polyclonal antibody
MAB5720	Ptprc monoclonal antibody, clone RA3-6B2 (FITC)	PAB16671	RANBP9 polyclonal antibody
MAB5722	Ptprc monoclonal antibody, clone RA3-6B2 (PE)	H00010411-D01	RAPGEF3 MaxPab rabbit polyclonal antibody (D01)
MAB5727	Ptprc monoclonal antibody, clone RA3-6B2 (PE/Cy5.5)	H00005919-D01	RARRES2 MaxPab rabbit polyclonal antibody (D01)
MAB5728	Ptprc monoclonal antibody, clone RA3-6B2 (PE/Cy7)	H00005917-D01	RARS MaxPab rabbit polyclonal antibody (D01)
MAB5723	Ptprc monoclonal antibody, clone RA3-6B2 (PE/TXRD)	H00005921-D01	RASA1 MaxPab rabbit polyclonal antibody (D01)
MAB5725	Ptprc monoclonal antibody, clone RA3-6B2 (SPRD)	PAB9014	RASA1 polyclonal antibody
MAB3873	PTPRC monoclonal antibody, clone UCHL1	H00023551-M02	RASD2 monoclonal antibody (M02), clone 1C7
MAB4434	PTPRC monoclonal antibody, clone UCHL1 (APC)	H00011228-D01	RASSF8 MaxPab rabbit polyclonal antibody (D01)
MAB4435	PTPRC monoclonal antibody, clone UCHL1 (FITC)	H00011228-M01	RASSF8 monoclonal antibody (M01), clone 2G1
MAB4436	PTPRC monoclonal antibody, clone UCHL1 (PE)	MAB0003	RB1 (phospho S612) monoclonal antibody, clone 3C11
MAB6540	PTPRM monoclonal antibody, clone T10-AF1A8	MAB0002	RB1 (phospho S612) monoclonal antibody, clone 4E4
MAB5645	PTPRR monoclonal antibody, clone 6A6	MAB0005	RB1 (phospho S807) monoclonal antibody, clone 5H12
MAB6538	PTPRR monoclonal antibody, clone T10-AF4H11	MAB0004	RB1 (phospho T356) monoclonal antibody, clone 4E3
MAB5681	Ptprs monoclonal antibody, clone 17G7.2	MAB0295	RB1 monoclonal antibody
MAB2740	PTPRT monoclonal antibody, clone T20-AF3C7	MAB2168	RB1 monoclonal antibody, clone 1D1
H00284119-D01	PTRF MaxPab rabbit polyclonal antibody (D01)	MAB2167	RB1 monoclonal antibody, clone 1E5
H00005805-D01	PTS MaxPab rabbit polyclonal antibody (D01)	MAB2166	RB1 monoclonal antibody, clone 1G5
H00080324-D01	PUS1 MaxPab rabbit polyclonal antibody (D01)	MAB1029	RB1 monoclonal antibody, clone 32C8
H00080324-M01	PUS1 monoclonal antibody (M01), clone 3A3	MAB2169	RB1 monoclonal antibody, clone 3C8
H00005818-D01	PVRL1 MaxPab rabbit polyclonal antibody (D01)	MAB2163	RB1 monoclonal antibody, clone 4G2
H00005819-M01	PVRL2 monoclonal antibody (M01), clone 2A6-2C1	MAB2164	RB1 monoclonal antibody, clone 4H6
MAB3387	Pvrl2 monoclonal antibody, clone 502-57	MAB0783	RB1 monoclonal antibody, clone 7E6
H00025945-D01	PVRL3 MaxPab rabbit polyclonal antibody (D01)	MAB2165	RB1 monoclonal antibody, clone 7F12
MAB5184	Pvrl3 monoclonal antibody, clone 103-A1	MAB2008	RB1 monoclonal antibody, clone C36
H00005822-D01	PWP2 MaxPab rabbit polyclonal antibody (D01)	MAB2006	RB1 monoclonal antibody, clone XZ104
H00005829-D01	PXN MaxPab rabbit polyclonal antibody (D01)	MAB2007	RB1 monoclonal antibody, clone XZ133
MAB1366	PXN monoclonal antibody, clone M107	MAB2010	RB1 monoclonal antibody, clone XZ55
PAB13764	PYCARD polyclonal antibody	MAB8004	RBBP4 monoclonal antibody, clone 11G10
H00005831-D01	PYCR1 MaxPab rabbit polyclonal antibody (D01)	H00010616-D01	RBCK1 MaxPab rabbit polyclonal antibody (D01)
H00065263-D01	PYCR1 MaxPab rabbit polyclonal antibody (D01)	PAB11266	RBL2 polyclonal antibody
PAB0326	PYDC1 polyclonal antibody	PAB8611	RBM39 polyclonal antibody
H00005836-D01	PYGL MaxPab rabbit polyclonal antibody (D01)	H00010181-D01	RBM5 MaxPab rabbit polyclonal antibody (D01)
PAB13592	PYGM polyclonal antibody	MAB2484	RBM8A monoclonal antibody, clone 4C4
H00005859-M01	QARS monoclonal antibody (M01), clone 5F5	H00023543-D01	RBM9 MaxPab rabbit polyclonal antibody (D01)
H00025797-D01	QPCT MaxPab rabbit polyclonal antibody (D01)	H00005937-D01	RBMS1 MaxPab rabbit polyclonal antibody (D01)
H00057403-D01	RAB22A MaxPab rabbit polyclonal antibody (D01)	H00001104-D01	RCC1 MaxPab rabbit polyclonal antibody (D01)
H00051715-D01	RAB23 MaxPab rabbit polyclonal antibody (D01)	H00001104-M02	RCC1 monoclonal antibody (M02), clone 1C1
H00010981-D01	RAB32 MaxPab rabbit polyclonal antibody (D01)	PAB12339	RCC1 polyclonal antibody
H00010981-M01	RAB32 monoclonal antibody (M01), clone 1C7	H00005957-M02A	RCV1 monoclonal antibody (M02A), clone 3G10
H00023682-D01	RAB38 MaxPab rabbit polyclonal antibody (D01)	H00005957-M03A	RCV1 monoclonal antibody (M03A), clone 1D3
H00116442-D01	RAB39B MaxPab rabbit polyclonal antibody (D01)	H00005957-M42A	RCV1 monoclonal antibody (M42A), clone 3H1
H00005865-M02	RAB3B monoclonal antibody (M02), clone 1A7	H00005957-D01	RCVRN MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB5104	RCVRN polyclonal antibody	H00009045-D01	RPL14 MaxPab rabbit polyclonal antibody (D01)
H00005962-D01	RDX MaxPab rabbit polyclonal antibody (D01)	PAB15593	RPN3 polyclonal antibody
H00005971-D01	RELB MaxPab rabbit polyclonal antibody (D01)	PAB15594	RPN5 polyclonal antibody
MAB9470	Reovirus T3D monoclonal antibody, clone 3356s	PAB15595	RPN7 polyclonal antibody
MAB4329	Respiratory syncytial virus monoclonal antibody, clone RSV3216(B016)	H00055197-D01	RPRD1A MaxPab rabbit polyclonal antibody (D01)
MAB4330	Respiratory syncytial virus nucleoprotein monoclonal antibody, clone B023	MAB0792	RPS6KA1 monoclonal antibody, clone 8C12
PAB13584	RFA1/RFA2/RFA3 polyclonal antibody	PAB12645	RPS6KB1 (phospho S434) polyclonal antibody
PAB11311	RFA2 (phospho S122) polyclonal antibody	H00006198-D01	RPS6KB1 MaxPab rabbit polyclonal antibody (D01)
PAB11364	RFA2 polyclonal antibody	PAB10302	RPS6KB1 polyclonal antibody
H00005984-D01	RFC4 MaxPab rabbit polyclonal antibody (D01)	PAB12720	RPS6KB1 polyclonal antibody
PAB9104	RFP polyclonal antibody	PAB14283	RPS6KB1/RPS6KB2 polyclonal antibody
MAB8763	RFP tag monoclonal antibody, clone 465d3	H00006199-D02	RPS6KB2 MaxPab rabbit polyclonal antibody (D01)
H00005993-D01	RFX5 MaxPab rabbit polyclonal antibody (D01)	H00057521-D01	RPTOR MaxPab rabbit polyclonal antibody (D01)
PAB0661	RFXAP polyclonal antibody	H00010325-D01	RRAGB MaxPab rabbit polyclonal antibody (D01)
H00008601-D01	RGS20 MaxPab rabbit polyclonal antibody (D01)	H00006237-M01	RRAS monoclonal antibody (M01), clone 2E12
PAB9098	RGS3 polyclonal antibody	H00006237-M02	RRAS monoclonal antibody (M02), clone 2A6
H00008490-D01	RGS5 MaxPab rabbit polyclonal antibody (D01)	PAB15696	Rrbp1 polyclonal antibody
H00006009-D01	RHEB MaxPab rabbit polyclonal antibody (D01)	H00006241-D01	RRM2 MaxPab rabbit polyclonal antibody (D01)
PAB12644	RHO (phospho S188) polyclonal antibody	H00006241-M01	RRM2 monoclonal antibody (M01), clone 1E1
MAB8559	RHO monoclonal antibody, clone 1D4	MAB2412	rRNA monoclonal antibody
MAB8558	RHO monoclonal antibody, clone 4D2	H00023212-D01	RRS1 MaxPab rabbit polyclonal antibody (D01)
PAB12727	RHO polyclonal antibody	MAB4949	RSF1 monoclonal antibody
PAB8982	RIC8A polyclonal antibody	PAB10306	RSL1D1 polyclonal antibody
H00253260-M01	RICTOR monoclonal antibody (M01), clone 1F3	MAB2786	RT1-B monoclonal antibody, clone 3D6
PAB14880	Rim1 polyclonal antibody	H00065078-D01	RTN4R MaxPab rabbit polyclonal antibody (D01)
H00083732-D01	RIOK1 MaxPab rabbit polyclonal antibody (D01)	H00000864-D01	RUNX3 MaxPab rabbit polyclonal antibody (D01)
H00055781-D01	RIOK2 MaxPab rabbit polyclonal antibody (D01)	H00010856-D01	RUVBL2 MaxPab rabbit polyclonal antibody (D01)
H00008780-D01	RIOK3 MaxPab rabbit polyclonal antibody (D01)	H00006256-M03	RXRA monoclonal antibody (M03), clone 4B5
H00008780-M01	RIOK3 monoclonal antibody (M01), clone 4A12-1D2	H00006256-M17	RXRA monoclonal antibody (M17), clone 1G1
MAB0836	RIPK1 monoclonal antibody, clone 7H10	H00023429-D01	RYBP MaxPab rabbit polyclonal antibody (D01)
MAB2692	RIPK2 monoclonal antibody, clone AF28D3	H00023429-D03	RYBP MaxPab rabbit polyclonal antibody (D03)
MAB3270	RLA-DR monoclonal antibody, clone 2C4	H00023429-M08	RYBP monoclonal antibody (M08), clone 1F4
MAB3269	RLA-DQ monoclonal antibody, clone RDR34	MAB2546	S100A1 monoclonal antibody, clone SPM354
H00006017-D01	RLBP1 MaxPab rabbit polyclonal antibody (D01)	PAB12732	S100A1 polyclonal antibody
H00006013-M01	RLN1 monoclonal antibody (M01), clone 1H6	H00006278-D03	S100A7 MaxPab rabbit polyclonal antibody (D03)
PAB11340	RNASE1 polyclonal antibody	H00006278-M07A	S100A7 monoclonal antibody (M07A), clone 3E11
H00140432-D01	RNF113B MaxPab rabbit polyclonal antibody (D01)	PAB8984	S100B polyclonal antibody
H00140432-M02	RNF113B monoclonal antibody (M02), clone 2B11	H00006291-D01	SAA4 MaxPab rabbit polyclonal antibody (D01)
H00051132-D01	RNF12 MaxPab rabbit polyclonal antibody (D01)	H00006291-M08	SAA4 monoclonal antibody (M08), clone 3C11
H00051132-M01	RNF12 monoclonal antibody (M01), clone 1G10	PAB8989	sad polyclonal antibody
H00168433-D01	RNF133 MaxPab rabbit polyclonal antibody (D01)	PAB8990	sad polyclonal antibody
H00255488-D01	RNF144B MaxPab rabbit polyclonal antibody (D01)	MAB8006	SAFB monoclonal antibody, clone 6F7
H00056254-D01	RNF20 MaxPab rabbit polyclonal antibody (D01)	H00025939-D01	SAMHD1 MaxPab rabbit polyclonal antibody (D01)
H00064320-D01	RNF25 MaxPab rabbit polyclonal antibody (D01)	H00010284-D01	SAP18 MaxPab rabbit polyclonal antibody (D01)
PAB11256	RNF25 polyclonal antibody	PAB11392	Sarcosine Oxidase polyclonal antibody
H00140545-M01	RNF32 monoclonal antibody (M01), clone 5C1	MAB3435	SARM1 monoclonal antibody, clone Sarmy-1
H00006049-D01	RNF6 MaxPab rabbit polyclonal antibody (D01)	PAB11369	SARS orf1ab polyclonal antibody
H00006050-D01	RNH1 MaxPab rabbit polyclonal antibody (D01)	PAB11370	SARS orf1ab polyclonal antibody
H00006051-M02	RNPEP monoclonal antibody (M02), clone 4E1	PAB11371	SARS orf1ab polyclonal antibody
MAB9162	robo monoclonal antibody, clone 7a746	PAB11372	SARS orf1ab polyclonal antibody
PAB12731	ROCK2 (phospho T249) polyclonal antibody	H00054938-D01	SARS2 MaxPab rabbit polyclonal antibody (D01)
PAB12643	ROCK2 (phospho T396) polyclonal antibody	H00060485-D01	SAV1 MaxPab rabbit polyclonal antibody (D01)
PAB12728	ROCK2 polyclonal antibody	H00060485-M02	SAV1 monoclonal antibody (M02), clone 3B2
PAB12729	ROCK2 polyclonal antibody	PAB8985	sax polyclonal antibody
PAB12730	ROCK2 polyclonal antibody	PAB8986	sax polyclonal antibody
PAB9075	ROCK2 polyclonal antibody	H00010609-D01	SC65 MaxPab rabbit polyclonal antibody (D01)
MAB5397	Rotavirus p42 monoclonal antibody, clone 1-O-15	PAB12042	Scarb1 polyclonal antibody
H00051765-D01	RP6-213H19.1 MaxPab rabbit polyclonal antibody (D01)	PAB12597	Scarb1 polyclonal antibody
H00006117-D01	RPA1 MaxPab rabbit polyclonal antibody (D01)	PAB12157	Scarb2 polyclonal antibody
MAB3410	RPA1 monoclonal antibody, clone 2H10, 70C or alphaSSB70C	MAB5546	SCD monoclonal antibody, clone CD.E10
MAB8269	RPA2 monoclonal antibody	PAB17107	SCG2/Scg2 polyclonal antibody
MAB8005	RPA3 monoclonal antibody, clone 11.1	PAB12703	SCGB2A2 polyclonal antibody
H00029935-D01	RPA4 MaxPab rabbit polyclonal antibody (D01)	MAB7310	<i>Schizosaccharomyces pombe</i> CDK monoclonal antibody, clone AA10.2
MAB3413	RPA4 monoclonal antibody, clone alphaSSB34A	MAB7252	<i>Schizosaccharomyces pombe</i> Cig2 monoclonal antibody, clone CIG 3A11/5
H00079657-D01	RPAP3 MaxPab rabbit polyclonal antibody (D01)	H00051540-D01	SCLY MaxPab rabbit polyclonal antibody (D01)

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00006322-D01	SCML1 MaxPab rabbit polyclonal antibody (D01)	H00006317-M01	SERPINB3 monoclonal antibody (M01), clone 2F5
MAB8587	Scn2a1 monoclonal antibody, clone S69-3	MAB1972	SERPINB6 monoclonal antibody, clone 3A
MAB6661	SCN9A monoclonal antibody, clone S68-6	MAB1774	SERPINB6 monoclonal antibody, clone Pl6-3A
H00006342-M01	SCP2 monoclonal antibody (M01), clone 2E9-1B3	H00008710-D01	SERPINB7 MaxPab rabbit polyclonal antibody (D01)
H00057147-D01	SCYL3 MaxPab rabbit polyclonal antibody (D01)	MAB1204	Serpinc1 monoclonal antibody
MAB3455	SDC1 monoclonal antibody, clone MI15	PAB5074	Serpinc1 polyclonal antibody
H00006385-D01	SDC4 MaxPab rabbit polyclonal antibody (D01)	PAB5075	Serpinc1 polyclonal antibody
PAB9045	SDC4 polyclonal antibody	H00003053-D01	SERPIND1 MaxPab rabbit polyclonal antibody (D01)
H00006386-D01	SDCBP MaxPab rabbit polyclonal antibody (D01)	H00005054-D01	SERPINE1 MaxPab rabbit polyclonal antibody (D01)
H00006386-M01	SDCBP monoclonal antibody (M01), clone 2C12	H00005054-M01	SERPINE1 monoclonal antibody (M01), clone 3F2
H00195814-D01	SDR16C5 MaxPab rabbit polyclonal antibody (D01)	MAB2661	SERPINE1 monoclonal antibody, clone I-2
H00113675-D01	SDSL MaxPab rabbit polyclonal antibody (D01)	PAB14897	SERPINE1 polyclonal antibody
H00113675-M01A	SDSL monoclonal antibody (M01A), clone 1D6-6B8	PAB5094	SERPINE1 polyclonal antibody
H00113675-M03	SDSL monoclonal antibody (M03), clone 2F8-6C9	H00005345-D01	SERPINF2 MaxPab rabbit polyclonal antibody (D01)
H00006396-M02	SEC13 monoclonal antibody (M02), clone 1G7	PAB5042	SERPINF2 polyclonal antibody
PAB7505	SECISBP2 polyclonal antibody	H00000710-M01	SERPING1 monoclonal antibody (M01), clone 3F4-1D9
MAB3644	Secretory component monoclonal antibody, clone SC-05	MAB0756	SERPING1 monoclonal antibody, clone 4G12
MAB6227	SELE monoclonal antibody, clone 1.2B6	PAB8774	SERPINH1 polyclonal antibody
MAB9303	SELE monoclonal antibody, clone 56c6	H00027244-D01	SESN1 MaxPab rabbit polyclonal antibody (D01)
MAB6956	SELE monoclonal antibody, clone HAE-1f	H00083667-M03	SESN2 monoclonal antibody (M03), clone 3B8
MAB1915	Sele monoclonal antibody, clone UZ5	PAB13763	SET polyclonal antibody
PAB8702	SELE polyclonal antibody	H00002810-M01	SFN monoclonal antibody (M01), clone 3C3
MAB8213	SELL monoclonal antibody, clone DREG56	MAB7036	SFRS2 monoclonal antibody, clone SC-35
MAB6957	SELL monoclonal antibody, clone LAM1-116	H00006439-D01	SFTPB MaxPab rabbit polyclonal antibody (D01)
MAB3235	Sell monoclonal antibody, clone MEL-14	H00006444-D01	SGCD MaxPab rabbit polyclonal antibody (D01)
MAB5738	Sell monoclonal antibody, clone MEL-14 (APC)	H00006445-D01	SGCG MaxPab rabbit polyclonal antibody (D01)
MAB5741	Sell monoclonal antibody, clone MEL-14 (APC/CY7)	H00006445-M02	SGCG monoclonal antibody (M02), clone 3C5
MAB5736	Sell monoclonal antibody, clone MEL-14 (Biotin)	PAB9020	SGK1 polyclonal antibody
MAB5735	Sell monoclonal antibody, clone MEL-14 (FITC)	H00084197-D01	SGK196 MaxPab rabbit polyclonal antibody (D01)
MAB5737	Sell monoclonal antibody, clone MEL-14 (PE)	H00010110-D01	SGK2 MaxPab rabbit polyclonal antibody (D01)
MAB5740	Sell monoclonal antibody, clone MEL-14 (PE/Cy7)	H00010110-M17	SGK2 monoclonal antibody (M17), clone 2F6
MAB5739	Sell monoclonal antibody, clone MEL-14 (SPRD)	H00151648-D01	SGOL1 MaxPab rabbit polyclonal antibody (D01)
PAB8813	SELL polyclonal antibody	H00151648-M01	SGOL1 monoclonal antibody (M01), clone 3C11
MAB0706	SELM monoclonal antibody, clone 10C1	H00151648-M02	SGOL1 monoclonal antibody (M02), clone 4G6
MAB0705	SELM monoclonal antibody, clone 4C2	H00006448-D01	SGSH MaxPab rabbit polyclonal antibody (D01)
MAB6958	SELP monoclonal antibody, clone G1/G1-4	H00006449-D01	SGTA MaxPab rabbit polyclonal antibody (D01)
MAB7346	SELPLG monoclonal antibody, clone KPL-1	H00006449-M06	SGTA monoclonal antibody (M06), clone 2E11
H00055829-D01	SELS MaxPab rabbit polyclonal antibody (D01)	PAB10367	SH2D1A polyclonal antibody
PAB8494	Sema3b polyclonal antibody	H00117157-M04	SH2D1B monoclonal antibody (M04), clone 2B4
MAB6581	Sema4a monoclonal antibody, clone HIAT-2	H00009047-D01	SH2D2A MaxPab rabbit polyclonal antibody (D01)
MAB5207	SEMA4D monoclonal antibody, clone 133.1G6	H00063898-D01	SH2D4A MaxPab rabbit polyclonal antibody (D01)
MAB6973	SEMA4D monoclonal antibody, clone 133-1C6	PAB10066	SH3BP2 (phospho S427) polyclonal antibody
MAB6583	Sema4d monoclonal antibody, clone SK-3	PAB12776	Shaker polyclonal antibody
MAB0914	SEMA7A monoclonal antibody, clone MEM-150	MAB7890	Shank2 monoclonal antibody, clone S23B-49
MAB6582	Sema7a monoclonal antibody, clone SKK-7	MAB7891	Shank2 monoclonal antibody, clone S23b-6
PAB14672	Semg1 polyclonal antibody	MAB7898	Shank3 monoclonal antibody, clone S69-46
PAB14671	Semg2 polyclonal antibody	MAB0807	SHC1 monoclonal antibody, clone 47F4
H00123228-D01	SENP8 MaxPab rabbit polyclonal antibody (D01)	H00053358-D01	SHC3 MaxPab rabbit polyclonal antibody (D01)
PAB11323	Sephs2 polyclonal antibody	MAB9424	shg monoclonal antibody, clone 018nf
H00001731-D01	SEPT1 MaxPab rabbit polyclonal antibody (D01)	H00006474-D01	SHOX2 MaxPab rabbit polyclonal antibody (D01)
H00001731-M03	SEPT1 monoclonal antibody (M03), clone 1F12	MAB9485	SI monoclonal antibody, clone 10n65
H00151011-D01	SEPT10 MaxPab rabbit polyclonal antibody (D01)	MAB2299	SIAH1/SIAH2 monoclonal antibody, clone 8G7H12
H00005413-D01	SEPT5 MaxPab rabbit polyclonal antibody (D01)	PAB12777	Sialyltransferase polyclonal antibody
PAB15710	Sept8 polyclonal antibody	H00059307-D01	SIGIRR MaxPab rabbit polyclonal antibody (D01)
PAB11324	Sepw1 polyclonal antibody	H00000946-D01	SIGLEC6 MaxPab rabbit polyclonal antibody (D01)
H00051156-M02	SERPINA10 monoclonal antibody (M02), clone 1E11	H00000946-M02	SIGLEC6 monoclonal antibody (M02), clone 2G6
H00145264-D01	SERPINA12 MaxPab rabbit polyclonal antibody (D01)	MAB8014	SIK2 monoclonal antibody, clone S15G10
H00145264-M01	SERPINA12 monoclonal antibody (M01), clone 3B3	H00006490-D01	SILV MaxPab rabbit polyclonal antibody (D01)
H00000012-D01	SERPINA3 MaxPab rabbit polyclonal antibody (D01)	MAB7856	Simian immunodeficiency virus type 1 Nef monoclonal antibody, clone 17 (04-001)
H00000012-M01	SERPINA3 monoclonal antibody (M01), clone 1E6	MAB2466	SIP1 monoclonal antibody, clone 2E17
PAB5027	SERPINA7 polyclonal antibody	PAB0004	SIRT1 polyclonal antibody
H00001992-D01	SERPINB1 MaxPab rabbit polyclonal antibody (D01)	PAB0875	SIRT1 polyclonal antibody
H00005055-M08	SERPINB2 monoclonal antibody (M08), clone 3A9	H00051547-D01	SIRT7 MaxPab rabbit polyclonal antibody (D01)
H00006317-D01	SERPINB3 MaxPab rabbit polyclonal antibody (D01)	PAB12032	SIRT7 polyclonal antibody
H00006317-D03	SERPINB3 MaxPab rabbit polyclonal antibody (D03)	PAB12033	SIRT7 polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB12034	SIRT7 polyclonal antibody	MAB0684	SOD1 monoclonal antibody, clone 72B1
MAB0933	SIT1 monoclonal antibody, clone SIT-01	MAB0688	SOD1 monoclonal antibody, clone 8A1
MAB4525	SIT1 monoclonal antibody, clone SIT-01 (PE)	PAB1193	SOD1 polyclonal antibody
MAB4908	Ski gene protein monoclonal antibody, clone B548M	PAB11970	Sod1 polyclonal antibody
H00006502-D01	SKP2 MaxPab rabbit polyclonal antibody (D01)	PAB14492	SOD1 polyclonal antibody
H00084174-D01	SLA2 MaxPab rabbit polyclonal antibody (D01)	PAB9030	SOD1 polyclonal antibody
MAB5264	SLAMF1 monoclonal antibody, clone A12	MAB0692	SOD2 monoclonal antibody, clone 1E8
MAB3908	SLAMF1 monoclonal antibody, clone IPO-3	MAB0711	SOD2 monoclonal antibody, clone 23G5
MAB7082	SLAMF6 monoclonal antibody, clone NT-7	MAB0710	SOD2 monoclonal antibody, clone 4F10
H00006563-D01	SLC14A1 MaxPab rabbit polyclonal antibody (D01)	PAB1196	SOD2 polyclonal antibody
H00000292-D03	SLC25A5 MaxPab rabbit polyclonal antibody (D03)	PAB11971	Sod2 polyclonal antibody
H00376497-D01	SLC27A1 MaxPab rabbit polyclonal antibody (D01)	PAB14493	SOD2 polyclonal antibody
PAB8954	Sic2a10 polyclonal antibody	MAB0749	SOD3 monoclonal antibody, clone 1H12
MAB0911	SLC3A2 monoclonal antibody, clone MEM-108	H00054937-D01	SOHLH2 MaxPab rabbit polyclonal antibody (D01)
MAB6960	SLC3A2 monoclonal antibody, clone UM7F8	H00006652-M03	SORD monoclonal antibody (M03), clone 1C2
MAB1721	SLC3A2 monoclonal antibody, clone WK4	H00009580-D01	SOX13 MaxPab rabbit polyclonal antibody (D01)
H00006523-D01	SLC5A1 MaxPab rabbit polyclonal antibody (D01)	PAB12735	SOX2 polyclonal antibody
PAB12026	Sic7a11 polyclonal antibody	H00011063-D01	SOX30 MaxPab rabbit polyclonal antibody (D01)
MAB2294	Sic9a3 (phospho S552) monoclonal antibody, clone 14D5	PAB12646	SOX9 (phospho S181) polyclonal antibody
MAB2293	Sic9a3 (phospho S605) monoclonal antibody, clone 10A8	H00006662-M02	SOX9 monoclonal antibody (M02), clone 3C10
PAB9079	SLK polyclonal antibody	PAB12736	SOX9 polyclonal antibody
H00057152-M06	SLURP1 monoclonal antibody (M06), clone 4D1	PAB12737	SOX9 polyclonal antibody
MAB1006	SMAD2 monoclonal antibody, clone 7A5	H00006672-D01	SP100 MaxPab rabbit polyclonal antibody (D01)
MAB0996	SMAD4 monoclonal antibody, clone 31G4	H00003431-D01	SP110 MaxPab rabbit polyclonal antibody (D01)
H00004091-D01	SMAD6 MaxPab rabbit polyclonal antibody (D01)	H00003431-M01	SP110 monoclonal antibody (M01), clone 8C8
PAB8486	SMARCA2 polyclonal antibody	PAB8794	SPAG9 polyclonal antibody
PAB8598	SMARCA4 polyclonal antibody	MAB9257	SPARC monoclonal antibody, clone 97a2
MAB4948	SMARCA5 monoclonal antibody	H00008877-D01	SPHK1 MaxPab rabbit polyclonal antibody (D01)
H00056916-D01	SMARCAD1 MaxPab rabbit polyclonal antibody (D01)	H00008877-M02	SPHK1 monoclonal antibody (M02), clone 2A8
PAB15737	Smarcad1 polyclonal antibody	PAB15559	SPHK1 polyclonal antibody
H00006598-D01	SMARCB1 MaxPab rabbit polyclonal antibody (D01)	PAB2801	SPHK1 polyclonal antibody
H00006605-M01	SMARCE1 monoclonal antibody (M01), clone 6G11	H00056848-D01	SPHK2 MaxPab rabbit polyclonal antibody (D01)
H00009126-D01	SMC3 MaxPab rabbit polyclonal antibody (D01)	H00006690-M01	SPINK1 monoclonal antibody (M01), clone 4D4
PAB12734	SMN1 polyclonal antibody	H00006693-D01	SPN MaxPab rabbit polyclonal antibody (D01)
H00006607-D01	SMN2 MaxPab rabbit polyclonal antibody (D01)	MAB0884	SPN monoclonal antibody, clone MEM-59
H00055512-D01	SMPD3 MaxPab rabbit polyclonal antibody (D01)	H00006697-D01	SPR MaxPab rabbit polyclonal antibody (D01)
H00027293-M01	SMPDL3B monoclonal antibody (M01), clone 5E12	H00006697-M01	SPR monoclonal antibody (M01), clone 4F2
H00057154-M01	SMURF1 monoclonal antibody (M01), clone 1D7	H00010252-D01	SPRY1 MaxPab rabbit polyclonal antibody (D01)
H00056950-D01	SMYD2 MaxPab rabbit polyclonal antibody (D01)	MAB2515	SPTA1 monoclonal antibody, clone AF10
H00006615-D01	SNAI1 MaxPab rabbit polyclonal antibody (D01)	MAB2516	SPTB monoclonal antibody, clone DB2
H00006615-M41	SNAI1 monoclonal antibody (M41), clone 1A5	PAB12022	SPTBN2 polyclonal antibody
H00006591-D01	SNAI2 MaxPab rabbit polyclonal antibody (D01)	PAB12023	SPTBN2 polyclonal antibody
H00008773-M01	SNAP23 monoclonal antibody (M01), clone 2F5-3D4	PAB10405	SR-alpha polyclonal antibody
MAB1796	SNAP25 monoclonal antibody, clone 16D1	PAB10406	SR-beta polyclonal antibody
H00006622-M01	SNCA monoclonal antibody (M01), clone 2E4	MAB1155	SRC monoclonal antibody, clone 17AT28
PAB5100	SNCA polyclonal antibody	H00006728-D01	SRP19 MaxPab rabbit polyclonal antibody (D01)
PAB5101	SNCB polyclonal antibody	H00006732-D01	SRPK1 MaxPab rabbit polyclonal antibody (D01)
PAB5102	SNCG polyclonal antibody	H00063826-D01	SRR MaxPab rabbit polyclonal antibody (D01)
H00079753-D01	SNIP1 MaxPab rabbit polyclonal antibody (D01)	MAB8814	SRY monoclonal antibody, clone SRY.G12
H00006626-D01	SNRPA MaxPab rabbit polyclonal antibody (D01)	MAB8488	SSEA-1 monoclonal antibody, clone MC-480
MAB6749	Snrpc monoclonal antibody, clone 4H12	MAB8489	SSEA-3 monoclonal antibody, clone MC-631
H00010073-D01	SNUPN MaxPab rabbit polyclonal antibody (D01)	MAB8490	SSEA-4 monoclonal antibody, clone MC-813-70
H00006642-D01	SNX1 MaxPab rabbit polyclonal antibody (D01)	PAB10358	SSeCKS polyclonal antibody
H00029916-D01	SNX11 MaxPab rabbit polyclonal antibody (D01)	H00085464-D01	SSH2 MaxPab rabbit polyclonal antibody (D01)
H00029916-M05	SNX11 monoclonal antibody (M05), clone 2G1	H00008082-D01	SSPN MaxPab rabbit polyclonal antibody (D01)
H00257364-D01	SNX33 MaxPab rabbit polyclonal antibody (D01)	PAB10372	SSRP1 polyclonal antibody
H00051429-D01	SNX9 MaxPab rabbit polyclonal antibody (D01)	PAB10373	SSRP1 polyclonal antibody
H00051429-M03	SNX9 monoclonal antibody (M03), clone 3C6	PAB12019	Sstr1 polyclonal antibody
H00006646-D01	SOAT1 MaxPab rabbit polyclonal antibody (D01)	H00006757-D01	SSX2 MaxPab rabbit polyclonal antibody (D01)
PAB12778	SOCS1 polyclonal antibody	H00117178-D01	SSX2IP MaxPab rabbit polyclonal antibody (D01)
PAB12779	SOCS2 polyclonal antibody	H00006487-D01	ST3GAL3 MaxPab rabbit polyclonal antibody (D01)
PAB12780	SOCS3 polyclonal antibody	H00006480-M01	ST6GAL1 monoclonal antibody (M01), clone 2E12
PAB12781	SOCS4 polyclonal antibody	MAB6959	ST6GAL1 monoclonal antibody, clone LN1
H00009655-D01	SOCS5 MaxPab rabbit polyclonal antibody (D01)	PAB10067	<i>Staphylococcus aureus</i> Protein A polyclonal antibody (Agarose)
H00009306-D01	SOCS6 MaxPab rabbit polyclonal antibody (D01)	MAB3641	STAT1 (phospho S727) monoclonal antibody, clone PSM1

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB12752	STAT1 (phospho Y701) polyclonal antibody	MAB9143	SUMO2 monoclonal antibody, clone 6j4k9
MAB3648	STAT1 monoclonal antibody, clone SM1	H00006821-D01	SUOX MaxPab rabbit polyclonal antibody (D01)
MAB0962	STAT1 monoclonal antibody, clone SM2	PAB10370	SUPT16H polyclonal antibody
PAB12782	STAT1 polyclonal antibody	MAB7368	SUV39H1 monoclonal antibody, clone 44.1
PAB12783	STAT1 polyclonal antibody	PAB14108	SUZ12 polyclonal antibody
PAB12820	STAT1 polyclonal antibody	MAB8761	SV40 tag monoclonal antibody, clone 8f7s9
PAB13772	STAT2 polyclonal antibody	H00023075-D01	SWAP70 MaxPab rabbit polyclonal antibody (D01)
PAB12753	STAT3 (phospho S727) polyclonal antibody	PAB11378	swi6 polyclonal antibody
PAB17909	Stat3 (phospho S727) polyclonal antibody	PAB15585	SWI6 polyclonal antibody
PAB12754	STAT3 (phospho Y705) polyclonal antibody	MAB0838	SYK monoclonal antibody, clone 12E3
PAB12784	STAT3 polyclonal antibody	MAB1806	SYK monoclonal antibody, clone 4D10.1
PAB13773	STAT3 polyclonal antibody	MAB0837	SYK monoclonal antibody, clone 8H5
PAB12755	STAT4 (phospho Y693) polyclonal antibody	MAB0948	SYK monoclonal antibody, clone SYK-01
H00006775-M01	STAT4 monoclonal antibody (M01), clone 1C2-1C12	PAB9702	SYN1 polyclonal antibody
PAB12756	STAT5A (phospho Y694) polyclonal antibody	H00008224-D01	SYN3 MaxPab rabbit polyclonal antibody (D01)
H00006776-D01	STAT5A MaxPab rabbit polyclonal antibody (D01)	MAB2498	SYNCRIP monoclonal antibody, clone I8E4
PAB12785	STAT5A polyclonal antibody	MAB2315	Synj1 monoclonal antibody, clone AC1
PAB12786	STAT5A polyclonal antibody	PAB4994	TAC1 polyclonal antibody
PAB15867	STAT5A polyclonal antibody	PAB4995	TAC3 polyclonal antibody
PAB9095	STAT5B polyclonal antibody	H00006881-D01	TAF10 MaxPab rabbit polyclonal antibody (D01)
MAB0100	STAT6 (phospho Y641) monoclonal antibody, clone 177C322	H00008148-D01	TAF15 MaxPab rabbit polyclonal antibody (D01)
PAB15557	STAT6 (phospho Y641) polyclonal antibody	H00009014-D01	TAF1B MaxPab rabbit polyclonal antibody (D01)
H00006780-D01	STAU1 MaxPab rabbit polyclonal antibody (D01)	H00010629-D01	TAF6L MaxPab rabbit polyclonal antibody (D01)
H00008987-D01	STBD1 MaxPab rabbit polyclonal antibody (D01)	H00006888-D01	TALDO1 MaxPab rabbit polyclonal antibody (D01)
MAB3602	STIM1 monoclonal antibody, clone CDN3H4	PAB0259	TANK polyclonal antibody
H00010963-M01	STIP1 monoclonal antibody (M01), clone 4B6	PAB9048	TANK polyclonal antibody
H00008576-D01	STK16 MaxPab rabbit polyclonal antibody (D01)	PAB17946	Taok1 polyclonal antibody
H00009263-D01	STK17A MaxPab rabbit polyclonal antibody (D01)	PAB9085	TAOK1 polyclonal antibody
H00008859-D01	STK19 MaxPab rabbit polyclonal antibody (D01)	H00051347-D01	TAOK3 MaxPab rabbit polyclonal antibody (D01)
H00008859-M01	STK19 monoclonal antibody (M01), clone 4E11	PAB8791	TAOK3 polyclonal antibody
H00202374-D01	STK32A MaxPab rabbit polyclonal antibody (D01)	H00006895-D01	TARBP2 MaxPab rabbit polyclonal antibody (D01)
H00065975-D01	STK33 MaxPab rabbit polyclonal antibody (D01)	H00006895-M01	TARBP2 monoclonal antibody (M01), clone 1B1
H00065975-M02	STK33 monoclonal antibody (M02), clone 6G1	H00023435-M01	TARDBP monoclonal antibody (M01), clone 2E2-D3
H00065975-M03	STK33 monoclonal antibody (M03), clone 4A12	H00006905-D01	TBCE MaxPab rabbit polyclonal antibody (D01)
H00065975-M06	STK33 monoclonal antibody (M06), clone 4F7	H00006907-D01	TBL1X MaxPab rabbit polyclonal antibody (D01)
H00065975-M09	STK33 monoclonal antibody (M09), clone 1F10	H00026608-D01	TBL2 MaxPab rabbit polyclonal antibody (D01)
H00011329-M01	STK38 monoclonal antibody (M01), clone 2G8-1F3	MAB0121	TBP monoclonal antibody
PAB3184	STK39 polyclonal antibody	MAB4734	TBP monoclonal antibody, clone 1TBP18
H00083931-D01	STK40 MaxPab rabbit polyclonal antibody (D01)	H00006919-D01	TCEA2 MaxPab rabbit polyclonal antibody (D01)
H00003925-M01A	STMN1 monoclonal antibody (M01A), clone 3A9	H00006920-M06	TCEA3 monoclonal antibody (M06), clone 2C9
H00030968-D01	STOML2 MaxPab rabbit polyclonal antibody (D01)	H00009338-D01	TCEAL1 MaxPab rabbit polyclonal antibody (D01)
H00092335-D01	STRADA MaxPab rabbit polyclonal antibody (D01)	H00009338-M01	TCEAL1 monoclonal antibody (M01), clone 3B9
PAB15556	Streptavidin polyclonal antibody	MAB4922	T-Cell Receptor Pan gamma/delta monoclonal antibody, clone Immu510
MAB4915	<i>Streptococcus</i> Group A (carbohydrate) monoclonal antibody, clone 3401	MAB4923	T-Cell receptor V alpha-24 monoclonal antibody, clone C15
MAB4916	<i>Streptococcus</i> Group B monoclonal antibody, clone 072	MAB1812	TCF4 monoclonal antibody, clone 6H5-3
PAB15599	<i>Streptococcus</i> NADase polyclonal antibody	PAB13547	TCF7L2 polyclonal antibody
PAB11375	<i>Streptococcus</i> Protein G polyclonal antibody	H00008115-M02A	TCL1A monoclonal antibody (M02A), clone 3G5
PAB12016	STRN polyclonal antibody	H00008115-M10A	TCL1A monoclonal antibody (M10A), clone 2D7
MAB2291	STRN3 monoclonal antibody, clone S68	H00006947-D01	TCN1 MaxPab rabbit polyclonal antibody (D01)
MAB9089	Stromal monoclonal antibody, clone 87a95	H00006947-M12	TCN1 monoclonal antibody (M12), clone 3F10
H00000412-D01	STS MaxPab rabbit polyclonal antibody (D01)	H00006949-D01	TCOF1 MaxPab rabbit polyclonal antibody (D01)
H00008676-D01	STX11 MaxPab rabbit polyclonal antibody (D01)	MAB2457	TCR-alpha/beta monoclonal antibody, clone KJ12-98
H00023673-D01	STX12 MaxPab rabbit polyclonal antibody (D01)	MAB3240	Tcrb monoclonal antibody, clone H57-597
H00023673-M01	STX12 monoclonal antibody (M01), clone 3B9	MAB4735	Tcrg monoclonal antibody, clone 536
H00006804-D01	STX1A MaxPab rabbit polyclonal antibody (D01)	MAB7194	Tcrg monoclonal antibody, clone GL3
H00002054-D01	STX2 MaxPab rabbit polyclonal antibody (D01)	H00055775-D01	TDP1 MaxPab rabbit polyclonal antibody (D01)
H00006810-D01	STX4 MaxPab rabbit polyclonal antibody (D01)	H00007004-M01	TEAD4 monoclonal antibody (M01), clone 5H3
H00008417-D01	STX7 MaxPab rabbit polyclonal antibody (D01)	H00054937-M04	TEB1 monoclonal antibody (M04), clone 1E1
H00055359-D01	STYK1 MaxPab rabbit polyclonal antibody (D01)	MAB1831	Tenascin monoclonal antibody, clone T2H5
H00006817-D01	SULT1A1 MaxPab rabbit polyclonal antibody (D01)	MAB0074	TERF2 monoclonal antibody, clone 4A794.15
MAB0722	SULT1A3 monoclonal antibody, clone 3F10	PAB12598	TERF2 polyclonal antibody
H00027284-M04A	SULT1B1 monoclonal antibody (M04A), clone 4C9	PAB15121	TERF2 polyclonal antibody (Biotin)
H00006783-D01	SULT1E1 MaxPab rabbit polyclonal antibody (D01)	PAB15120	TERF2 polyclonal antibody (DyLight 488)
H00006820-D01	SULT2B1 MaxPab rabbit polyclonal antibody (D01)	PAB15119	TERF2 polyclonal antibody (DyLight 549)
MAB9142	SUMO1 monoclonal antibody, clone 85j96	PAB15118	TERF2 polyclonal antibody (DyLight 649)

Antibody Tools for Immunoprecipitation

Antibodies for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB15117	TERF2 polyclonal antibody (HRP)	H00008795-M01	TNFRSF10B monoclonal antibody (M01), clone 2D6
PAB5164	Testosterone polyclonal antibody	MAB3478	TNFRSF10B monoclonal antibody, clone HS201
H00056157-D01	TEX13A MaxPab rabbit polyclonal antibody (D01)	MAB3481	TNFRSF10D monoclonal antibody, clone HS402
H00051368-M01	TEX264 monoclonal antibody (M01), clone 2A3-1A10	MAB0091	TNFRSF11A monoclonal antibody, clone 64C1385.1
PAB14522	TF polyclonal antibody	H00004982-D01	TNFRSF11B MaxPab rabbit polyclonal antibody (D01)
H00007019-D01	TFAM MaxPab rabbit polyclonal antibody (D01)	MAB3414	TNFRSF11B monoclonal antibody, clone Bony-1
MAB9057	TFAP2A monoclonal antibody, clone 84r69	MAB3720	TNFRSF14 monoclonal antibody, clone Laury-1
H00007024-D01	TFCP2 MaxPab rabbit polyclonal antibody (D01)	H0000608-D01	TNFRSF17 MaxPab rabbit polyclonal antibody (D01)
MAB0938	TFRC monoclonal antibody, clone MEM-189	MAB3374	Tnfrsf17 monoclonal antibody, clone Vicky-2
MAB5273	TFRC monoclonal antibody, clone MEM75	H00008784-D01	TNFRSF18 MaxPab rabbit polyclonal antibody (D01)
MAB0897	TFRC monoclonal antibody, clone MEM-75	MAB6244	Tnfrsf18 monoclonal antibody, clone DTA-1
MAB4108	Tfrc monoclonal antibody, clone R17217.1.4	MAB3877	TNFRSF21 monoclonal antibody, clone DR-6-04-EC
MAB1720	TFRC monoclonal antibody, clone YUK9	MAB4474	TNFRSF21 monoclonal antibody, clone DR-6-04-EC (PE)
PAB5026	TG polyclonal antibody	MAB6965	TNFRSF4 monoclonal antibody, clone Ber Act 35
PAB12738	TGFB1 polyclonal antibody	H00000939-M01	TNFRSF7 monoclonal antibody (M01), clone 1E2-A3
PAB12739	TGFB2 polyclonal antibody	MAB6906	TNFRSF8 monoclonal antibody, clone AC10
H00007045-D01	TGFB1 MaxPab rabbit polyclonal antibody (D01)	MAB7144	Tnfsf11 monoclonal antibody, clone IK22/5
H00007052-D01	TGM2 MaxPab rabbit polyclonal antibody (D01)	MAB3366	TNFSF13 monoclonal antibody, clone Aprily-8
MAB0733	TGM2 monoclonal antibody, clone 14G2	H00010673-M04	TNFSF13B monoclonal antibody (M04), clone 2C11
MAB0732	TGM2 monoclonal antibody, clone 2A1	H00010673-M19	TNFSF13B monoclonal antibody (M19), clone 3G6
PAB1225	TGM2 polyclonal antibody	H00008740-M01	TNFSF14 monoclonal antibody (M01), clone 4E3
H00055145-D01	THAP1 MaxPab rabbit polyclonal antibody (D01)	H00008995-D01	TNFSF18 MaxPab rabbit polyclonal antibody (D01)
MAB7009	THBS1 monoclonal antibody, clone TSP-B7	MAB3725	Tnfsf4 monoclonal antibody, clone MPG34
MAB8007	THOC1 monoclonal antibody, clone 5E10	H00010188-D01	TNK2 MaxPab rabbit polyclonal antibody (D01)
MAB2468	THOC4 monoclonal antibody, clone 11G5	MAB3065	TNNI3 monoclonal antibody, clone 16A11
PAB11276	THRA polyclonal antibody	MAB3063	TNNI3 monoclonal antibody, clone 581
PAB11277	THRB polyclonal antibody	PAB12743	TNNI3 polyclonal antibody
PAB5065	Thrombin polyclonal antibody	H00007139-D01	TNNT2 MaxPab rabbit polyclonal antibody (D01)
H00079178-D01	THTPA MaxPab rabbit polyclonal antibody (D01)	PAB12744	TNNT2 polyclonal antibody
MAB4084	Thy1 monoclonal antibody, clone 30-H12	MAB2496	TNPO1 monoclonal antibody, clone D45
MAB5520	Thy1 monoclonal antibody, clone 30-H12 (Biotin)	H00007148-M01	TNXB monoclonal antibody (M01), clone 2D3-1B10
MAB5521	Thy1 monoclonal antibody, clone 30-H12 (FITC)	H00054472-D01	TOLLIP MaxPab rabbit polyclonal antibody (D01)
MAB5522	Thy1 monoclonal antibody, clone 30-H12 (PE)	MAB3709	Tollip monoclonal antibody, clone Kimmy-1
MAB7281	THY1 monoclonal antibody, clone aTHy-1A1	H00010953-D01	TOMM34 MaxPab rabbit polyclonal antibody (D01)
MAB2525	THY1 monoclonal antibody, clone F15-42-1	H00010953-M02	TOMM34 monoclonal antibody (M02), clone 1D2
MAB1356	Thy1 monoclonal antibody, clone T1-A3	PAB13698	TOP2A polyclonal antibody
MAB1357	Thy1 monoclonal antibody, clone T1-H7	PAB12640	TP53 (phospho S15) polyclonal antibody
H00007072-D01	TIA1 MaxPab rabbit polyclonal antibody (D01)	PAB8912	TP53 (phospho S15) polyclonal antibody
H00026230-D01	TIAM2 MaxPab rabbit polyclonal antibody (D01)	MAB0660	TP53 (phospho S315) monoclonal antibody, clone 18
PAB13795	TICAM1 polyclonal antibody	PAB12641	TP53 (phospho S37) polyclonal antibody
H00353376-D01	TICAM2 MaxPab rabbit polyclonal antibody (D01)	H00007157-D01	TP53 MaxPab rabbit polyclonal antibody (D01)
H00007076-M01	TIMP1 monoclonal antibody (M01), clone 4D12	H00007157-M01	TP53 monoclonal antibody (M01), clone 2C3
H00114609-D01	TIRAP MaxPab rabbit polyclonal antibody (D01)	MAB2711	TP53 monoclonal antibody, clone AF2C1
MAB3447	TIRAP monoclonal antibody, clone Sebi-1	MAB1969	TP53 monoclonal antibody, clone HR231
H00007083-M02	TK1 monoclonal antibody (M02), clone F12	MAB7073	TP53 monoclonal antibody, clone PAb 1103
H00007086-D01	TKT MaxPab rabbit polyclonal antibody (D01)	MAB7072	TP53 monoclonal antibody, clone PAb 1104
PAB15690	Tln1 polyclonal antibody	MAB1626	TP53 monoclonal antibody, clone S11-K
MAB0066	TLR2 monoclonal antibody, clone 1030A5.138	PAB12719	TP53 polyclonal antibody
MAB7833	Tlr2 monoclonal antibody, clone T2.5	PAB14128	TP53 polyclonal antibody
MAB0067	TLR2 monoclonal antibody, clone TL2.1	H00009540-D02	TP53I3 MaxPab rabbit polyclonal antibody (D02)
MAB3706	TLR3 monoclonal antibody, clone 40C1285	PAB0660	TP63 polyclonal antibody
MAB0088	TLR3 monoclonal antibody, clone 40C1285.6	MAB8038	TP73 monoclonal antibody, clone 38C674.2
MAB0071	TLR4 monoclonal antibody, clone HTA125	MAB8039	TP73 monoclonal antibody, clone 5B429
MAB0072	TLR4 monoclonal antibody, clone HTA125 (NaN3 Free)	PAB10146	TP73 polyclonal antibody
PAB0374	TLR4 polyclonal antibody	MAB7327	TPO monoclonal antibody, clone 6H7
H00030012-D01	TLX3 MaxPab rabbit polyclonal antibody (D01)	H00001200-D01	TPP1 MaxPab rabbit polyclonal antibody (D01)
H00011070-D01	TMEM115 MaxPab rabbit polyclonal antibody (D01)	PAB11993	TPX2 polyclonal antibody
H00055863-D01	TMEM126B MaxPab rabbit polyclonal antibody (D01)	PAB9088	TPX2 polyclonal antibody
H00029767-D01	TMOD2 MaxPab rabbit polyclonal antibody (D01)	PAB9051	TRA@ polyclonal antibody
MAB7093	Tnc monoclonal antibody, clone MTn-12	MAB2391	TRA-1-60 monoclonal antibody, clone TRA-1-60(R)
MAB7014	Tnf monoclonal antibody, clone TN3-19.12	MAB8491	TRA-1-60 monoclonal antibody, clone TRA-1-60R
PAB15549	TNF polyclonal antibody	MAB2392	TRA-1-81 monoclonal antibody, clone TRA-1-81
PAB5105	TNF polyclonal antibody	MAB8492	TRA-1-81 monoclonal antibody, clone TRA-1-81
MAB3473	TNFRSF10A monoclonal antibody, clone HS101	H00008717-M03	TRADD monoclonal antibody (M03), clone 3G3
H00008795-D01	TNFRSF10B MaxPab rabbit polyclonal antibody (D01)	PAB0273	TRADD polyclonal antibody

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB9061	TRADD polyclonal antibody	H00081629-M03	TSSK3 monoclonal antibody (M03), clone 3F9
H00007185-D01	TRAF1 MaxPab rabbit polyclonal antibody (D01)	H00081629-M06	TSSK3 monoclonal antibody (M06), clone 6H4
H00007185-M03	TRAF1 monoclonal antibody (M03), clone 1A6	H00083983-D01	TSSK6 MaxPab rabbit polyclonal antibody (D01)
PAB0274	TRAF1 polyclonal antibody	H00007265-D01	TTC1 MaxPab rabbit polyclonal antibody (D01)
PAB0275	TRAF2 polyclonal antibody	H00123016-D01	TTC8 MaxPab rabbit polyclonal antibody (D01)
PAB0276	TRAF3 polyclonal antibody	H00007272-D01	TTK MaxPab rabbit polyclonal antibody (D01)
H00010758-D01	TRAF3IP2 MaxPab rabbit polyclonal antibody (D01)	MAB0784	TTR monoclonal antibody, clone 10E1
PAB0277	TRAF4 polyclonal antibody	MAB0822	TTR monoclonal antibody, clone 9G6
H00007188-D01	TRAF5 MaxPab rabbit polyclonal antibody (D01)	PAB1221	TTR polyclonal antibody
MAB0062	TRAF5 monoclonal antibody, clone 55A219	PAB4998	TTR polyclonal antibody
PAB0278	Traf5 polyclonal antibody	PAB12569	Tub-1 polyclonal antibody
H00007189-D01	TRAF6 MaxPab rabbit polyclonal antibody (D01)	H00079861-D01	TUBAL3 MaxPab rabbit polyclonal antibody (D01)
PAB0279	TRAF6 polyclonal antibody	H00081027-D01	TUBB1 MaxPab rabbit polyclonal antibody (D01)
PAB0280	TRAIIP polyclonal antibody	H00056604-D01	TUBB4Q MaxPab rabbit polyclonal antibody (D01)
MAB5398	Transmissible Gastroenteritis virus 1E11 monoclonal antibody, clone 1Q17	H00051174-D01	TUBD1 MaxPab rabbit polyclonal antibody (D01)
H00010131-D01	TRAP1 MaxPab rabbit polyclonal antibody (D01)	H00051174-M05	TUBD1 monoclonal antibody (M05), clone 3C7
MAB0932	TRAT1 monoclonal antibody, clone TRIM-04	H00010844-D01	TUBGCP2 MaxPab rabbit polyclonal antibody (D01)
MAB4526	TRAT1 monoclonal antibody, clone TRIM-04 (PE)	H00027229-D01	TUBGCP4 MaxPab rabbit polyclonal antibody (D01)
PAB9052	TRB@ polyclonal antibody	MAB9224	Tubulin-alpha monoclonal antibody, clone 46q6
MAB3879	TRBV5-3 monoclonal antibody, clone MEM-262	MAB3660	Tubulin-alpha monoclonal antibody, clone TU-16
H00001787-D01	TRDMT1 MaxPab rabbit polyclonal antibody (D01)	MAB9223	Tubulin-beta monoclonal antibody, clone 465s
PAB5170	Trenbolone polyclonal antibody	H00007284-D01	TUFM MaxPab rabbit polyclonal antibody (D01)
H00007200-D01	TRH MaxPab rabbit polyclonal antibody (D01)	H00007284-M03	TUFM monoclonal antibody (M03), clone 1A2
H00028951-M11	TRIB2 monoclonal antibody (M11), clone 1D11	H00007288-D01	TULP2 MaxPab rabbit polyclonal antibody (D01)
H00057761-M07	TRIB3 monoclonal antibody (M07), clone 3D7	H00005756-D01	TWF1 MaxPab rabbit polyclonal antibody (D01)
H00010626-M01	TRIM16 monoclonal antibody (M01), clone 5F4	H00011344-D01	TWF2 MaxPab rabbit polyclonal antibody (D01)
H00010346-D01	TRIM22 MaxPab rabbit polyclonal antibody (D01)	H00007291-M01	TWIST1 monoclonal antibody (M01), clone 3E11
H00008805-M01	TRIM24 monoclonal antibody (M01), clone 2F2	H00007291-M13	TWIST1 monoclonal antibody (M13), clone 2G12
H00007706-D01	TRIM25 MaxPab rabbit polyclonal antibody (D01)	H00007291-M19	TWIST1 monoclonal antibody (M19), clone 1D3
H00022954-D01	TRIM32 MaxPab rabbit polyclonal antibody (D01)	H00007291-M20	TWIST1 monoclonal antibody (M20), clone 3H5
PAB15691	Trim33 polyclonal antibody	H00117581-M06	TWIST2 monoclonal antibody (M06), clone 1E2
H00084675-D01	TRIM55 MaxPab rabbit polyclonal antibody (D01)	H00117581-M12	TWIST2 monoclonal antibody (M12), clone 3C12
H00055223-D01	TRIM62 MaxPab rabbit polyclonal antibody (D01)	H00007295-M01	TXN monoclonal antibody (M01), clone 2A7
H00084676-D01	TRIM63 MaxPab rabbit polyclonal antibody (D01)	MAB0704	TXN monoclonal antibody, clone 3A1
H00493829-D01	TRIM72 MaxPab rabbit polyclonal antibody (D01)	MAB0721	TXN monoclonal antibody, clone 8A1
H00493829-M01	TRIM72 monoclonal antibody (M01), clone 2G1	MAB0723	TXN2 monoclonal antibody, clone 4C5
PAB12740	Trimethyl-Histone H3 polyclonal antibody	MAB0724	TXN2 monoclonal antibody, clone 71G4
PAB12741	Trimethyl-Histone H3 polyclonal antibody	H00023071-M01	TXNDC4 monoclonal antibody (M01), clone 3C7
PAB12742	Trimethyl-Histone H3 polyclonal antibody	H00010190-D01	TXNDC9 MaxPab rabbit polyclonal antibody (D01)
MAB7047	Trinitrophenol monoclonal antibody, clone RTK2758	MAB0678	TXNRD1 monoclonal antibody, clone 19A1
H00009322-M01	TRIP10 monoclonal antibody (M01), clone 1A9	MAB0683	TXNRD1 monoclonal antibody, clone 5A5
MAB2172	TRIP11 monoclonal antibody	MAB0703	TXNRD2 monoclonal antibody, clone 25B3
H00007205-M07	TRIP6 monoclonal antibody (M07), clone 3D12	MAB0686	TXNRD2 monoclonal antibody, clone 7B2
MAB9311	Tropomyosin monoclonal antibody, clone 6m45	H00007297-M01	TYK2 monoclonal antibody (M01), clone 6G12
MAB9310	Tropomyosin monoclonal antibody, clone 87u65	H00001890-D01	TYMP MaxPab rabbit polyclonal antibody (D01)
MAB1027	Troponin I monoclonal antibody, clone 3D10	H00007298-D01	TYMS MaxPab rabbit polyclonal antibody (D01)
MAB3714	Trp53 monoclonal antibody, clone PAb 122	H00007298-M02	TYMS monoclonal antibody (M02), clone 2B2
MAB3712	Trp53 monoclonal antibody, clone PAb 240	MAB3411	TYMS monoclonal antibody, clone TS 106
MAB1970	Trp53 monoclonal antibody, clone Pab122	MAB1833	TYMS monoclonal antibody, clone TS106
MAB6963	Trp53 monoclonal antibody, clone PAb240	H00007299-D01	TYR MaxPab rabbit polyclonal antibody (D01)
PAB8504	Trpa1 polyclonal antibody	H00007306-D01	TYRP1 MaxPab rabbit polyclonal antibody (D01)
MAB6662	Trpc4 monoclonal antibody, clone S77-15	H00053347-M01	UBASH3A monoclonal antibody (M01), clone 4D2
MAB8568	TRPC5 monoclonal antibody, clone S67-15	PAB12745	UBB polyclonal antibody
MAB7904	TRPC7 monoclonal antibody, clone S64A-36	PAB11288	UBE2J1 polyclonal antibody
MAB6663	Trpm7 monoclonal antibody, clone S74-25	H00009246-D01	UBE2L6 MaxPab rabbit polyclonal antibody (D01)
PAB12118	Trpm7 polyclonal antibody	PAB0191	UBE2N polyclonal antibody
MAB7896	Trpv3 monoclonal antibody, clone S15-39	H00027338-D01	UBE2S MaxPab rabbit polyclonal antibody (D01)
MAB6666	Trpv3 monoclonal antibody, clone S15-4	H00029089-M01	UBE2T monoclonal antibody (M01), clone 1E12-4A3
PAB15687	Tsc1 polyclonal antibody	H00029089-M02	UBE2T monoclonal antibody (M02), clone 4G1-4C2
H00008848-D01	TSC22D1 MaxPab rabbit polyclonal antibody (D01)	H00029089-M03	UBE2T monoclonal antibody (M03), clone 2A12-4F11
H00001831-M01	TSC22D3 monoclonal antibody (M01), clone 3A5	PAB8693	UBE3A polyclonal antibody
H00116461-D01	TSEN15 MaxPab rabbit polyclonal antibody (D01)	MAB6459	Ubinuclein 1 monoclonal antibody, clone UBN1-02
MAB8273	TSG1 monoclonal antibody, clone 4A10	MAB1016	Ubiquitin monoclonal antibody, clone 2A1
H00007251-D01	TSG101 MaxPab rabbit polyclonal antibody (D01)	MAB1017	Ubiquitin monoclonal antibody, clone 4A1
H00081629-D01	TSSK3 MaxPab rabbit polyclonal antibody (D01)	MAB1012	Ubiquitin+1 monoclonal antibody, clone 40B3

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
H00023304-D01	UBR2 MaxPab rabbit polyclonal antibody (D01)	MAB5857	VSIG1 monoclonal antibody, clone CT1 (Biotin)
H00007371-M03	UCK2 monoclonal antibody (M03), clone 1G12	MAB5856	VSIG1 monoclonal antibody, clone CT1 (FITC)
H00054963-D01	UCKL1 MaxPab rabbit polyclonal antibody (D01)	MAB5858	VSIG1 monoclonal antibody, clone CT1 (PE)
H00054600-D01	UGT1A9 MaxPab rabbit polyclonal antibody (D01)	MAB8350	v-SKI monoclonal antibody, clone G8
H00054600-M04	UGT1A9 monoclonal antibody (M04), clone 1A2	H00128434-D01	VSTM2L MaxPab rabbit polyclonal antibody (D01)
H00051720-D01	UIMC1 MaxPab rabbit polyclonal antibody (D01)	H00007448-D01	VTN MaxPab rabbit polyclonal antibody (D01)
H00007369-D01	UMOD MaxPab rabbit polyclonal antibody (D01)	MAB8077	VWF monoclonal antibody, clone 2Q2134
H00007369-M03	UMOD monoclonal antibody (M03), clone 3F10	H00023063-D01	WAPAL MaxPab rabbit polyclonal antibody (D01)
H00009094-D01	UNC119 MaxPab rabbit polyclonal antibody (D01)	H00007453-D01	WARS MaxPab rabbit polyclonal antibody (D01)
PAB11382	Urease polyclonal antibody	H00007453-M02	WARS monoclonal antibody (M02), clone 3A12
H00007390-D01	UROS MaxPab rabbit polyclonal antibody (D01)	H00010352-D01	WARS2 MaxPab rabbit polyclonal antibody (D01)
H00011274-D01	USP18 MaxPab rabbit polyclonal antibody (D01)	H00010810-D01	WASF3 MaxPab rabbit polyclonal antibody (D01)
H00009099-D01	USP2 MaxPab rabbit polyclonal antibody (D01)	H00114049-M11	WBSCR22 monoclonal antibody (M11), clone 2E12
PAB3916	USP2 polyclonal antibody	H00057728-D01	WDR19 MaxPab rabbit polyclonal antibody (D01)
H00010868-D01	USP20 MaxPab rabbit polyclonal antibody (D01)	H00091833-D01	WDR20 MaxPab rabbit polyclonal antibody (D01)
H00023032-M01	USP33 monoclonal antibody (M01), clone 5B5	H00010785-D01	WDR4 MaxPab rabbit polyclonal antibody (D01)
PAB9090	USP7 polyclonal antibody	PAB9581	WEE1 polyclonal antibody
PAB4997	UTS2 polyclonal antibody	PAB11985	WFS1 polyclonal antibody
PAB9804	Vamp polyclonal antibody	H00008839-D01	WISP2 MaxPab rabbit polyclonal antibody (D01)
PAB5109	VAMP2 polyclonal antibody	H00008839-M09	WISP2 monoclonal antibody (M09), clone 3D10
H00081839-D01	VANGL1 MaxPab rabbit polyclonal antibody (D01)	H00089780-D01	WNT3A MaxPab rabbit polyclonal antibody (D01)
H00009218-D01	VAPA MaxPab rabbit polyclonal antibody (D01)	H00081029-D01	WNT5B MaxPab rabbit polyclonal antibody (D01)
H00009217-D01	VAPB MaxPab rabbit polyclonal antibody (D01)	H00007478-D01	WNT8A MaxPab rabbit polyclonal antibody (D01)
MAB6698	Varicella zoster virus (gpI/IV) monoclonal antibody, clone SG1	MAB2289	WT1 monoclonal antibody, clone 6F-H2
MAB6699	Varicella zoster virus (gpI/IV) monoclonal antibody, clone SG1-1	H00011059-M01	WWP1 monoclonal antibody (M01), clone 1A7
MAB6700	Varicella zoster virus (gpII) monoclonal antibody, clone 2E6	H00011060-D01	WWP2 MaxPab rabbit polyclonal antibody (D01)
MAB6701	Varicella zoster virus (gpIII) monoclonal antibody, clone SG3	PAB11979	WWTR1 polyclonal antibody
MAB6702	Varicella zoster virus (gpIV) monoclonal antibody, clone SG4	PAB15161	WWTR1 polyclonal antibody (Biotin)
MAB6703	Varicella zoster virus monoclonal antibody, clone NCP-1	PAB15160	WWTR1 polyclonal antibody (DyLight 488)
H00114990-D01	VASN MaxPab rabbit polyclonal antibody (D01)	PAB15159	WWTR1 polyclonal antibody (DyLight 549)
PAB5188	Vasonatin polyclonal antibody	PAB15158	WWTR1 polyclonal antibody (DyLight 649)
PAB5189	Vasopressin polyclonal antibody	PAB15157	WWTR1 polyclonal antibody (HRP)
PAB9710	VAV2 polyclonal antibody	PAB15693	X99384 polyclonal antibody
PAB9711	VAV2 polyclonal antibody	PAB11978	XAGE1 polyclonal antibody
H00007411-D01	VBP1 MaxPab rabbit polyclonal antibody (D01)	H00007494-D01	XBP1 MaxPab rabbit polyclonal antibody (D01)
MAB1797	VCAM1 monoclonal antibody, clone 1G11.B1	PAB9587	Xcl1 polyclonal antibody
MAB9121	VCAM1 monoclonal antibody, clone 49s56	PAB9588	Xcl1 polyclonal antibody
MAB5602	Vcam1 monoclonal antibody, clone M/K-2	PAB11343	XDH polyclonal antibody
MAB8211	VCAM1 monoclonal antibody, clone STA	MAB5405	xgsk-3 monoclonal antibody, clone 1H8
MAB9156	VCAN monoclonal antibody, clone 135s5	MAB3755	xgsk-3 monoclonal antibody, clone G3
MAB8019	VCP monoclonal antibody, clone 4G9	H00000331-D01	XIAP MaxPab rabbit polyclonal antibody (D01)
PAB9091	VCP polyclonal antibody	PAB9589	XIAP polyclonal antibody
H00007421-D01	VDR MaxPab rabbit polyclonal antibody (D01)	H00007507-D01	XPA MaxPab rabbit polyclonal antibody (D01)
MAB2544	VEGFA monoclonal antibody, clone SPM225	MAB3418	XPA monoclonal antibody, clone 12F5
H00007423-D01	VEGFB MaxPab rabbit polyclonal antibody (D01)	H00057510-M01	XPO5 monoclonal antibody (M01), clone 2C5-1B3
H00007423-M03	VEGFB monoclonal antibody (M03), clone 4H5	PAB15694	Xpo5 polyclonal antibody
PAB13530	Vegfc polyclonal antibody	MAB3419	XRCC1 monoclonal antibody, clone 33-2-5
H00389136-D01	VGLL3 MaxPab rabbit polyclonal antibody (D01)	H00007517-D01	XRCC3 MaxPab rabbit polyclonal antibody (D01)
MAB1825	VHL monoclonal antibody, clone Ig33	H00007518-D01	XRCC4 MaxPab rabbit polyclonal antibody (D01)
H00007429-D01	VIL1 MaxPab rabbit polyclonal antibody (D01)	MAB3426	XRCC5 monoclonal antibody, clone 111
H00007431-M01	VIM monoclonal antibody (M01), clone 3E9	MAB6468	XRCC5 monoclonal antibody, clone MEM-54
MAB1008	VIM monoclonal antibody, clone 33C7	H00002547-D01	XRCC6 MaxPab rabbit polyclonal antibody (D01)
MAB9141	VIM monoclonal antibody, clone 679k	H00010413-M01	YAP1 monoclonal antibody (M01), clone 2F12
MAB3880	VIM monoclonal antibody, clone VI-01	PAB12570	YAP1 polyclonal antibody
MAB8013	VIM monoclonal antibody, clone VI-10	PAB15065	YAP1 polyclonal antibody (Biotin)
PAB1302	VIM polyclonal antibody	PAB15064	YAP1 polyclonal antibody (DyLight 488)
PAB9580	VIM polyclonal antibody	PAB15063	YAP1 polyclonal antibody (DyLight 549)
H00007433-D01	VIPR1 MaxPab rabbit polyclonal antibody (D01)	PAB15062	YAP1 polyclonal antibody (DyLight 649)
PAB5136	VLDL polyclonal antibody	PAB15061	YAP1 polyclonal antibody (HRP)
H00008876-D01	VNN1 MaxPab rabbit polyclonal antibody (D01)	H00051087-D01	YBX2 MaxPab rabbit polyclonal antibody (D01)
H00008876-M05	VNN1 monoclonal antibody (M05), clone 2B10	MAB2200	YDJ1 monoclonal antibody, clone 1G10.H8
H00057617-M04	VPS18 monoclonal antibody (M04), clone 4E9	MAB2201	YDJ1 monoclonal antibody, clone 2A7.H6
H00023339-D01	VPS39 MaxPab rabbit polyclonal antibody (D01)	H00007525-D01	YES1 MaxPab rabbit polyclonal antibody (D01)
H00007443-D01	VRK1 MaxPab rabbit polyclonal antibody (D01)	MAB8759	YFP monoclonal antibody, clone 64a59
MAB3280	VSIG1 monoclonal antibody, clone CT1	MAB6224	YFV monoclonal antibody, clone F21-2

Antibody Tools for Immunoprecipitation

Cat. Num.	Product Name	Cat. Num.	Product Name
MAB1804	YWHAQ monoclonal antibody, clone 3B9	H00084307-D01	ZNF397 MaxPab rabbit polyclonal antibody (D01)
PAB14303	YWHAQ polyclonal antibody	H00084838-D01	ZNF496 MaxPab rabbit polyclonal antibody (D01)
PAB18088	YWHAZ polyclonal antibody	H00090441-D01	ZNF622 MaxPab rabbit polyclonal antibody (D01)
MAB0804	ZAP70 monoclonal antibody, clone 49B4	H00007579-D01	ZSCAN20 MaxPab rabbit polyclonal antibody (D01)
PAB9592	ZBTB7A polyclonal antibody	H00007589-D01	ZSCAN21 MaxPab rabbit polyclonal antibody (D01)
PAB0263	ZDHHC23 polyclonal antibody	H00011130-D01	ZWINT MaxPab rabbit polyclonal antibody (D01)
PAB15492	Z-DNA polyclonal antibody	H00007791-M01	ZYX monoclonal antibody (M01), clone 2C10-4A7
H00084936-D01	ZFYVE19 MaxPab rabbit polyclonal antibody (D01)	PAB5167	11a-hydroxytestosterone polyclonal antibody
H00084936-M02	ZFYVE19 monoclonal antibody (M02), clone 3G4-2B11	PAB5168	11-ketotestosterone polyclonal antibody
H00009372-D01	ZFYVE9 MaxPab rabbit polyclonal antibody (D01)	PAB5162	17-hydroxyprogesterone polyclonal antibody
H00064393-D01	ZMAT3 MaxPab rabbit polyclonal antibody (D01)	PAB5169	19-nortestosterone polyclonal antibody
H00079698-D01	ZMAT4 MaxPab rabbit polyclonal antibody (D01)	PAB5163	20a-hydroxyprogesterone polyclonal antibody
H00079698-M02	ZMAT4 monoclonal antibody (M02), clone 2A4	PAB5172	3-keto-5a, 16-androstene polyclonal antibody
H00007690-D01	ZNF131 MaxPab rabbit polyclonal antibody (D01)	PAB5145	4-androsten-17a-OL-3-one polyclonal antibody
H00007566-D01	ZNF18 MaxPab rabbit polyclonal antibody (D01)	PAB5171	4-pregnen-17, 20-diol-3-one polyclonal antibody
H00007569-D01	ZNF182 MaxPab rabbit polyclonal antibody (D01)	PAB5146	5-androsten-3-OL-7, 17-dione polyclonal antibody
H00007762-D01	ZNF215 MaxPab rabbit polyclonal antibody (D01)	MAB4951	5-mC monoclonal antibody, clone 33D3
H00007775-D01	ZNF232 MaxPab rabbit polyclonal antibody (D01)	MAB4952	5-mC monoclonal antibody, clone b
H00010127-D01	ZNF263 MaxPab rabbit polyclonal antibody (D01)	PAB5161	6-hydroxyprogesterone polyclonal antibody
H00010781-M05	ZNF266 monoclonal antibody (M05), clone 1B9	PAB5154	6-keto-17a-estradiol polyclonal antibody
H00055905-M03	ZNF313 monoclonal antibody (M03), clone S2	PAB5155	6-keto-17beta-estradiol polyclonal antibody
PAB8868	ZNF335 polyclonal antibody	PAB5156	6-ketoestriol polyclonal antibody
H00023567-D01	ZNF346 MaxPab rabbit polyclonal antibody (D01)	PAB5158	6-ketoestrone polyclonal antibody
H00059348-M01	ZNF350 monoclonal antibody (M01), clone 1A9	PAB5160	7a-progesterone polyclonal antibody
H00252884-D01	ZNF396 MaxPab rabbit polyclonal antibody (D01)		

Secondary Antibodies

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB9204	Donkey Anti-Goat IgG (H&L) secondary antibody	PAB9193	Goat Anti-Donkey IgG F(c) secondary antibody
PAB9205	Donkey Anti-Goat IgG (H&L) secondary antibody	PAB9200	Goat Anti-Ferret IgA (alpha chain) secondary antibody
PAB9276	Donkey Anti-Human IgG (H&L) secondary antibody	PAB9201	Goat Anti-Ferret IgG (gamma chain) secondary antibody
PAB9277	Donkey Anti-Human IgG (H&L) secondary antibody	PAB9202	Goat Anti-Ferret IgG IgA IgM (H&L) secondary antibody
PAB9328	Donkey Anti-Mouse IgG (H&L) secondary antibody	PAB9203	Goat Anti-Ferret IgM (mu) secondary antibody
PAB9329	Donkey Anti-Mouse IgG (H&L) secondary antibody	PAB9228	Goat Anti-Golden Syrian hamster IgG (H&L) secondary antibody
PAB9381	Donkey Anti-Rabbit IgG (H&L) secondary antibody	PAB9229	Goat Anti-Golden Syrian hamster IgG (H&L) secondary antibody
PAB9382	Donkey Anti-Rabbit IgG (H&L) secondary antibody	PAB9231	Goat Anti-Golden Syrian hamster IgG (H&L) secondary antibody
PAB9383	Donkey Anti-Rabbit IgG (H&L) secondary antibody	PAB9244	Goat Anti-Guinea pig IgG (H&L) secondary antibody
PAB9413	Donkey Anti-Rat IgG (H&L) secondary antibody	PAB9245	Goat Anti-Guinea pig IgG (H&L) secondary antibody
PAB9414	Donkey Anti-Rat IgG (H&L) secondary antibody	PAB9250	Goat Anti-Guinea pig IgG (H&L) secondary antibody
PAB9445	Donkey Anti-Sheep IgG (H&L) secondary antibody	PAB9246	Goat Anti-Guinea pig IgG F(ab')2 secondary antibody
PAB9446	Donkey Anti-Sheep IgG (H&L) secondary antibody	PAB9251	Goat Anti-Guinea pig IgG F(ab')2 secondary antibody
PAB9449	Donkey Anti-Sheep IgG (H&L) secondary antibody	PAB9247	Goat Anti-Guinea pig IgG F(c) secondary antibody
PAB9133	Goat Anti-Armenian hamster IgG (H&L) secondary antibody	PAB9261	Goat Anti-Horse IgG (H&L) secondary antibody
PAB9134	Goat Anti-Armenian hamster IgG (H&L) secondary antibody	PAB9263	Goat Anti-Horse IgG (H&L) secondary antibody
PAB9135	Goat Anti-Armenian hamster IgG (H&L) secondary antibody	PAB9264	Goat Anti-Horse IgG F(ab')2 secondary antibody
PAB9136	Goat Anti-Bovine IgG (H&L) secondary antibody	PAB9262	Goat Anti-Horse IgG F(c) secondary antibody
PAB9138	Goat Anti-Bovine IgG (H&L) secondary antibody	PAB9278	Goat Anti-Human (kappa chain) secondary antibody
PAB9137	Goat Anti-Bovine IgG F(ab')2 secondary antibody	PAB9279	Goat Anti-Human (lambda chain) secondary antibody
PAB9151	Goat Anti-Cat IgG (H&L) secondary antibody	PAB9280	Goat Anti-Human IgA (alpha chain) secondary antibody
PAB9154	Goat Anti-Cat IgG (H&L) secondary antibody	PAB9281	Goat Anti-Human IgG (gamma chain) secondary antibody
PAB9152	Goat Anti-Cat IgG F(ab')2 secondary antibody	PAB9282	Goat Anti-Human IgG (H&L) secondary antibody
PAB9153	Goat Anti-Cat IgG F(c) secondary antibody	PAB9283	Goat Anti-Human IgG (H&L) secondary antibody
PAB9166	Goat Anti-Chicken IgG (H&L) secondary antibody	PAB9284	Goat Anti-Human IgG (H&L) secondary antibody
PAB9167	Goat Anti-Chicken IgG (H&L) secondary antibody	PAB9300	Goat Anti-Human IgG (H&L) secondary antibody
PAB9168	Goat Anti-Chicken IgM (mu) secondary antibody	PAB9285	Goat Anti-Human IgG F(ab')2 secondary antibody
PAB9179	Goat Anti-Dog IgG (H&L) secondary antibody	PAB9301	Goat Anti-Human IgG F(ab')2 secondary antibody
PAB9180	Goat Anti-Dog IgG F(c) secondary antibody	PAB9286	Goat Anti-Human IgG F(c) secondary antibody
PAB9191	Goat Anti-Donkey IgG (H&L) secondary antibody	PAB9287	Goat Anti-Human IgG F(c) secondary antibody
PAB9192	Goat Anti-Donkey IgG (H&L) secondary antibody	PAB9302	Goat Anti-Human IgG F(c) secondary antibody

Secondary Antibodies

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB9288	Goat Anti-Human IgG IgA IgM (H&L) secondary antibody	PAB9183	Rabbit Anti-Dog IgG F(c) secondary antibody
PAB9289	Goat Anti-Human IgM (mu) secondary antibody	PAB9186	Rabbit Anti-Dog IgM (mu) secondary antibody
PAB9324	Goat Anti-Monkey IgA (alpha chain) secondary antibody	PAB9196	Rabbit Anti-Donkey IgG (H&L) secondary antibody
PAB9325	Goat Anti-Monkey IgG (gamma chain) secondary antibody	PAB9197	Rabbit Anti-Donkey IgG (H&L) secondary antibody
PAB9326	Goat Anti-Monkey IgG IgA IgM (H&L) secondary antibody	PAB9211	Rabbit Anti-Goat IgG (H&L) secondary antibody
PAB9327	Goat Anti-Monkey IgM (mu) secondary antibody	PAB9212	Rabbit Anti-Goat IgG (H&L) secondary antibody
PAB9330	Goat Anti-Mouse IgG (H&L) secondary antibody	PAB9220	Rabbit Anti-Goat IgG (H&L) secondary antibody
PAB9331	Goat Anti-Mouse IgG (H&L) secondary antibody	PAB9213	Rabbit Anti-Goat IgG F(ab') ₂ secondary antibody
PAB9332	Goat Anti-Mouse IgG (H&L) secondary antibody	PAB9214	Rabbit Anti-Goat IgG F(ab') ₂ secondary antibody
PAB9347	Goat Anti-Mouse IgG (H&L) secondary antibody	PAB9221	Rabbit Anti-Goat IgG F(ab') ₂ secondary antibody
PAB9334	Goat Anti-Mouse IgG F(ab') ₂ secondary antibody	PAB9215	Rabbit Anti-Goat IgG F(c) secondary antibody
PAB9348	Goat Anti-Mouse IgG F(ab') ₂ secondary antibody	PAB9222	Rabbit Anti-Goat IgG F(c) secondary antibody
PAB9335	Goat Anti-Mouse IgG F(c) secondary antibody	PAB9223	Rabbit Anti-Goat IgM (mu) secondary antibody
PAB9343	Goat Anti-Mouse IgG F(c) secondary antibody	PAB9233	Rabbit Anti-Golden Syrian hamster IgG (H&L) secondary antibody
PAB9349	Goat Anti-Mouse IgG F(c) secondary antibody	PAB9237	Rabbit Anti-Golden Syrian hamster IgG (H&L) secondary antibody
PAB9336	Goat Anti-Mouse IgG IgA IgM (H&L) secondary antibody	PAB9234	Rabbit Anti-Golden Syrian hamster IgG F(ab') ₂ secondary antibody
PAB9333	Goat Anti-Mouse IgG IgM (H&L) secondary antibody	PAB9238	Rabbit Anti-Golden Syrian hamster IgG F(ab') ₂ secondary antibody
PAB9350	Goat Anti-Mouse IgM (H&L) secondary antibody	PAB9239	Rabbit Anti-Golden Syrian hamster IgM (mu) secondary antibody
PAB9338	Goat Anti-Mouse IgM (mu) secondary antibody	PAB9253	Rabbit Anti-Guinea pig IgG (H&L) secondary antibody
PAB9353	Goat Anti-Mouse serum secondary antibody	PAB9257	Rabbit Anti-Guinea pig IgG (H&L) secondary antibody
PAB9465	Goat Anti-Pig IgG (H&L) secondary antibody	PAB9254	Rabbit Anti-Guinea pig IgG F(ab') ₂ secondary antibody
PAB9385	Goat Anti-Rabbit IgG (H&L) secondary antibody	PAB9255	Rabbit Anti-Guinea pig IgG F(c) secondary antibody
PAB9386	Goat Anti-Rabbit IgG (H&L) secondary antibody	PAB9266	Rabbit Anti-Horse IgG (H&L) secondary antibody
PAB9387	Goat Anti-Rabbit IgG (H&L) secondary antibody	PAB9269	Rabbit Anti-Horse IgG (H&L) secondary antibody
PAB9397	Goat Anti-Rabbit IgG (H&L) secondary antibody	PAB9267	Rabbit Anti-Horse IgG F(ab') ₂ secondary antibody
PAB9388	Goat Anti-Rabbit IgG F(ab') ₂ secondary antibody	PAB9270	Rabbit Anti-Horse IgG F(ab') ₂ secondary antibody
PAB9389	Goat Anti-Rabbit IgG F(ab') ₂ secondary antibody	PAB9268	Rabbit Anti-Horse IgG F(c) secondary antibody
PAB9398	Goat Anti-Rabbit IgG F(ab') ₂ secondary antibody	PAB9271	Rabbit Anti-Horse IgM (mu) secondary antibody
PAB9390	Goat Anti-Rabbit IgG F(c) secondary antibody	PAB9308	Rabbit Anti-Human IgA (alpha chain) secondary antibody
PAB9391	Goat Anti-Rabbit IgG F(c) secondary antibody	PAB9309	Rabbit Anti-Human IgG (gamma chain) secondary antibody
PAB9399	Goat Anti-Rabbit IgG F(c) secondary antibody	PAB9310	Rabbit Anti-Human IgG (H&L) secondary antibody
PAB9400	Goat Anti-Rabbit IgM (mu) secondary antibody	PAB9311	Rabbit Anti-Human IgG (H&L) secondary antibody
PAB9402	Goat Anti-Rabbit serum secondary antibody	PAB9316	Rabbit Anti-Human IgG (H&L) secondary antibody
PAB9415	Goat Anti-Rat IgA (alpha chain) secondary antibody	PAB9312	Rabbit Anti-Human IgG F(ab') ₂ secondary antibody
PAB9416	Goat Anti-Rat IgG (H&L) secondary antibody	PAB9313	Rabbit Anti-Human IgG F(ab') ₂ secondary antibody
PAB9417	Goat Anti-Rat IgG (H&L) secondary antibody	PAB9317	Rabbit Anti-Human IgG F(ab') ₂ secondary antibody
PAB9418	Goat Anti-Rat IgG (H&L) secondary antibody	PAB9314	Rabbit Anti-Human IgG F(c) secondary antibody
PAB9428	Goat Anti-Rat IgG (H&L) secondary antibody	PAB9318	Rabbit Anti-Human IgG F(c) secondary antibody
PAB9419	Goat Anti-Rat IgG F(c) secondary antibody	PAB9319	Rabbit Anti-Human IgM (mu) secondary antibody
PAB9420	Goat Anti-Rat IgG IgA IgM (H&L) secondary antibody	PAB9315	Rabbit Anti-Human IgM Fc _{5u} secondary antibody
PAB9421	Goat Anti-Rat IgM (mu) secondary antibody	PAB9358	Rabbit Anti-Mouse (kappa chain) secondary antibody
PAB9407	Guinea pig Anti-Rabbit IgG (H&L) secondary antibody	PAB9359	Rabbit Anti-Mouse (lambda chain) secondary antibody
PAB9209	Mouse Anti-Goat IgG (H&L) secondary antibody	PAB9360	Rabbit Anti-Mouse IgA (alpha chain) secondary antibody
PAB9307	Mouse Anti-Human IgG (H&L) secondary antibody	PAB9361	Rabbit Anti-Mouse IgG (gamma chain) secondary antibody
PAB9408	Mouse Anti-Rabbit IgG (H&L) secondary antibody	PAB9362	Rabbit Anti-Mouse IgG (H&L) secondary antibody
PAB9141	Rabbit Anti-Bovine IgG (H&L) secondary antibody	PAB9363	Rabbit Anti-Mouse IgG F(ab') ₂ secondary antibody
PAB9144	Rabbit Anti-Bovine IgG (H&L) secondary antibody	PAB9364	Rabbit Anti-Mouse IgG F(ab') ₂ secondary antibody
PAB9142	Rabbit Anti-Bovine IgG F(ab') ₂ secondary antibody	PAB9374	Rabbit Anti-Mouse IgG F(ab') ₂ secondary antibody
PAB9145	Rabbit Anti-Bovine IgG F(ab') ₂ secondary antibody	PAB9365	Rabbit Anti-Mouse IgG F(c) secondary antibody
PAB9143	Rabbit Anti-Bovine IgG F(c) secondary antibody	PAB9366	Rabbit Anti-Mouse IgG1 (Gamma 1 chain) secondary antibody
PAB9146	Rabbit Anti-Bovine IgM (mu) secondary antibody	PAB9367	Rabbit Anti-Mouse IgG2a (Gamma 2a chain) secondary antibody
PAB9156	Rabbit Anti-Cat IgG (H&L) secondary antibody	PAB9368	Rabbit Anti-Mouse IgG2b (Gamma 2b chain) secondary antibody
PAB9159	Rabbit Anti-Cat IgG (H&L) secondary antibody	PAB9369	Rabbit Anti-Mouse IgG3 (Gamma 3 chain) secondary antibody
PAB9157	Rabbit Anti-Cat IgG F(ab') ₂ secondary antibody	PAB9375	Rabbit Anti-Mouse IgM (H&L) secondary antibody
PAB9160	Rabbit Anti-Cat IgG F(ab') ₂ secondary antibody	PAB9370	Rabbit Anti-Mouse IgM (mu) secondary antibody
PAB9158	Rabbit Anti-Cat IgG F(c) secondary antibody	PAB9376	Rabbit Anti-Mouse serum secondary antibody
PAB9161	Rabbit Anti-Cat IgM (mu) secondary antibody	PAB9466	Rabbit Anti-Pig IgG (H&L) secondary antibody
PAB9170	Rabbit Anti-Chicken IgG (H&L) secondary antibody	PAB9469	Rabbit Anti-Pig IgG (H&L) secondary antibody
PAB9174	Rabbit Anti-Chicken IgG (H&L) secondary antibody	PAB9467	Rabbit Anti-Pig IgG F(ab') ₂ secondary antibody
PAB9171	Rabbit Anti-Chicken IgG F(ab') ₂ secondary antibody	PAB9470	Rabbit Anti-Pig IgG F(ab') ₂ secondary antibody
PAB9175	Rabbit Anti-Chicken IgG F(ab') ₂ secondary antibody	PAB9468	Rabbit Anti-Pig IgG F(c) secondary antibody
PAB9172	Rabbit Anti-Chicken IgG F(c) secondary antibody	PAB9471	Rabbit Anti-Pig IgM (mu) secondary antibody
PAB9181	Rabbit Anti-Dog IgG (H&L) secondary antibody	PAB9430	Rabbit Anti-Rat IgG (H&L) secondary antibody
PAB9182	Rabbit Anti-Dog IgG F(ab') ₂ secondary antibody	PAB9434	Rabbit Anti-Rat IgG (H&L) secondary antibody
PAB9185	Rabbit Anti-Dog IgG F(ab') ₂ secondary antibody	PAB9431	Rabbit Anti-Rat IgG F(ab') ₂ secondary antibody

Secondary Antibodies

Cat. Num.	Product Name	Cat. Num.	Product Name
PAB9438	Rabbit Anti-Rat IgG F(ab') ₂ secondary antibody	PAB9456	Rabbit Anti-Sheep IgG F(c) secondary antibody
PAB9432	Rabbit Anti-Rat IgG F(c) secondary antibody	PAB9460	Rabbit Anti-Sheep IgM (mu) secondary antibody
PAB9433	Rabbit Anti-Rat IgM (mu) secondary antibody	PAB9476	Rabbit Anti-Turkey IgG (H&L) secondary antibody
PAB9451	Rabbit Anti-Sheep IgG (H&L) secondary antibody	PAB9410	Rat Anti-Rabbit IgG (H&L) secondary antibody
PAB9452	Rabbit Anti-Sheep IgG (H&L) secondary antibody	PAB9323	Sheep Anti-Human IgG (H&L) secondary antibody
PAB9458	Rabbit Anti-Sheep IgG (H&L) secondary antibody	PAB9380	Sheep Anti-Mouse IgG (H&L) secondary antibody
PAB9453	Rabbit Anti-Sheep IgG F(ab') ₂ secondary antibody	PAB9411	Sheep Anti-Rabbit IgG (H&L) secondary antibody
PAB9454	Rabbit Anti-Sheep IgG F(ab') ₂ secondary antibody	PAB9412	Sheep Anti-Rabbit IgG (H&L) secondary antibody
PAB9459	Rabbit Anti-Sheep IgG F(ab') ₂ secondary antibody	PAB9444	Sheep Anti-Rat IgG (H&L) secondary antibody
PAB9455	Rabbit Anti-Sheep IgG F(c) secondary antibody		

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00051166-PW1	AADAT (Human) IP-WB Antibody Pair	H00023746-PW1	AIPL1 (Human) IP-WB Antibody Pair
H00000016-PW1	AARS (Human) IP-WB Antibody Pair	H00000203-PW1	AK1 (Human) IP-WB Antibody Pair
H00089845-PW1	ABCC10 (Human) IP-WB Antibody Pair	H00010270-PW1	AKAP8 (Human) IP-WB Antibody Pair
H00085320-PW1	ABCC11 (Human) IP-WB Antibody Pair	H00001645-PW1	AKR1C1 (Human) IP-WB Antibody Pair
H00010061-PW1	ABCF2 (Human) IP-WB Antibody Pair	H00001646-PW1	AKR1C2 (Human) IP-WB Antibody Pair
H00000034-PW1	ACADM (Human) IP-WB Antibody Pair	H00001109-PW1	AKR1C4 (Human) IP-WB Antibody Pair
H00000039-PW1	ACAT2 (Human) IP-WB Antibody Pair	H00006718-PW1	AKR1D1 (Human) IP-WB Antibody Pair
H00084680-PW1	ACCS (Human) IP-WB Antibody Pair	H00008574-PW1	AKR7A2 (Human) IP-WB Antibody Pair
H00065057-PW1	ACD (Human) IP-WB Antibody Pair	H00008574-PW2	AKR7A2 (Human) IP-WB Antibody Pair
H00065057-PW2	ACD (Human) IP-WB Antibody Pair	H00000207-PW1	AKT1 (Human) IP-WB Antibody Pair
H00065057-PW3	ACD (Human) IP-WB Antibody Pair	H00000207-PW2	AKT1 (Human) IP-WB Antibody Pair
H00000048-PW1	ACO1 (Human) IP-WB Antibody Pair	H00000213-PW1	ALB (Human) IP-WB Antibody Pair
H00000051-PW1	ACOX1 (Human) IP-WB Antibody Pair	H00000216-PW1	ALDH1A1 (Human) IP-WB Antibody Pair
H00008309-PW1	ACOX2 (Human) IP-WB Antibody Pair	H00000216-PW2	ALDH1A1 (Human) IP-WB Antibody Pair
H00051703-PW1	ACSL5 (Human) IP-WB Antibody Pair	H00000218-PW1	ALDH3A1 (Human) IP-WB Antibody Pair
H00000086-PW1	ACTL6A (Human) IP-WB Antibody Pair	H00000221-PW1	ALDH3B1 (Human) IP-WB Antibody Pair
H00010880-PW1	ACTL7B (Human) IP-WB Antibody Pair	H00008659-PW1	ALDH4A1 (Human) IP-WB Antibody Pair
H00000093-PW1	ACVR2B (Human) IP-WB Antibody Pair	H00007915-PW1	ALDH5A1 (Human) IP-WB Antibody Pair
H00000095-PW1	ACY1 (Human) IP-WB Antibody Pair	H00221120-PW1	ALKBH3 (Human) IP-WB Antibody Pair
H00000095-PW2	ACY1 (Human) IP-WB Antibody Pair	H00055821-PW1	ALLC (Human) IP-WB Antibody Pair
H00000100-PW1	ADA (Human) IP-WB Antibody Pair	H00000249-PW1	ALPL (Human) IP-WB Antibody Pair
H00009507-PW1	ADAMT5 (Human) IP-WB Antibody Pair	H00000250-PW1	ALPP (Human) IP-WB Antibody Pair
H00009507-PW2	ADAMT5 (Human) IP-WB Antibody Pair	H00000251-PW1	ALPPL2 (Human) IP-WB Antibody Pair
H00000126-PW1	ADH1C (Human) IP-WB Antibody Pair	H00000251-PW2	ALPPL2 (Human) IP-WB Antibody Pair
H00000127-PW1	ADH4 (Human) IP-WB Antibody Pair	H00000265-PW1	AMELX (Human) IP-WB Antibody Pair
H00000131-PW1	ADH7 (Human) IP-WB Antibody Pair	H00000271-PW1	AMPD2 (Human) IP-WB Antibody Pair
H00000135-PW1	ADORA2A (Human) IP-WB Antibody Pair	H00000275-PW1	AMT (Human) IP-WB Antibody Pair
H00000154-PW1	ADRB2 (Human) IP-WB Antibody Pair	H00000279-PW1	AMY2A (Human) IP-WB Antibody Pair
H00000166-PW1	AES (Human) IP-WB Antibody Pair	H00010393-PW1	ANAPC10 (Human) IP-WB Antibody Pair
H00000174-PW1	AFP (Human) IP-WB Antibody Pair	H00027329-PW1	ANGPTL3 (Human) IP-WB Antibody Pair
H00000174-PW2	AFP (Human) IP-WB Antibody Pair	H00051129-PW1	ANGPTL4 (Human) IP-WB Antibody Pair
H00060509-PW1	AGBL5 (Human) IP-WB Antibody Pair	H00253935-PW1	ANGPTL5 (Human) IP-WB Antibody Pair
H00000177-PW1	AGER (Human) IP-WB Antibody Pair	H00000286-PW1	ANK1 (Human) IP-WB Antibody Pair
H00010554-PW1	AGPAT1 (Human) IP-WB Antibody Pair	H00057763-PW1	ANKRA2 (Human) IP-WB Antibody Pair
H00057085-PW1	AGTRAP (Human) IP-WB Antibody Pair	H00000301-PW1	ANXA1 (Human) IP-WB Antibody Pair
H00054806-PW1	AHI1 (Human) IP-WB Antibody Pair	H00000306-PW1	ANXA3 (Human) IP-WB Antibody Pair
H00079026-PW1	AHNAK (Human) IP-WB Antibody Pair	H00000307-PW1	ANXA4 (Human) IP-WB Antibody Pair
H00000196-PW1	AHR (Human) IP-WB Antibody Pair	H00000308-PW1	ANXA5 (Human) IP-WB Antibody Pair
H00010598-PW1	AHSA1 (Human) IP-WB Antibody Pair	H00000308-PW2	ANXA5 (Human) IP-WB Antibody Pair
H00010598-PW2	AHSA1 (Human) IP-WB Antibody Pair	H00023028-PW1	AOF2 (Human) IP-WB Antibody Pair
H00000197-PW1	AHSG (Human) IP-WB Antibody Pair	H00000328-PW1	APEX1 (Human) IP-WB Antibody Pair
H00000197-PW2	AHSG (Human) IP-WB Antibody Pair	H00000328-PW2	APEX1 (Human) IP-WB Antibody Pair
H00064853-PW1	AIDA (Human) IP-WB Antibody Pair	H00000335-PW1	APOA1 (Human) IP-WB Antibody Pair
H00084883-PW1	AIFM2 (Human) IP-WB Antibody Pair	H00027350-PW1	APOBEC3C (Human) IP-WB Antibody Pair
H00084883-PW2	AIFM2 (Human) IP-WB Antibody Pair	H00027350-PW2	APOBEC3C (Human) IP-WB Antibody Pair
H00084883-PW3	AIFM2 (Human) IP-WB Antibody Pair	H00000346-PW1	APOC4 (Human) IP-WB Antibody Pair
H00009447-PW1	AIM2 (Human) IP-WB Antibody Pair	H00023780-PW1	APOL2 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00000353-PW1	APRT (Human) IP-WB Antibody Pair	H00000684-PW2	BST2 (Human) IP-WB Antibody Pair
H00000383-PW1	ARG1 (Human) IP-WB Antibody Pair	H00000684-PW3	BST2 (Human) IP-WB Antibody Pair
H00000384-PW1	ARG2 (Human) IP-WB Antibody Pair	H00114781-PW1	BTBD9 (Human) IP-WB Antibody Pair
H00000384-PW2	ARG2 (Human) IP-WB Antibody Pair	H00000685-PW1	BTC (Human) IP-WB Antibody Pair
H00000396-PW1	ARHGDI A (Human) IP-WB Antibody Pair	H00000686-PW1	BTD (Human) IP-WB Antibody Pair
H00000396-PW2	ARHGDI A (Human) IP-WB Antibody Pair	H00000699-PW1	BUB1 (Human) IP-WB Antibody Pair
H00007984-PW1	ARHGEF5 (Human) IP-WB Antibody Pair	H00118924-PW1	C10orf4 (Human) IP-WB Antibody Pair
H00115761-PW1	ARL11 (Human) IP-WB Antibody Pair	H00143384-PW1	C10orf46 (Human) IP-WB Antibody Pair
H00000379-PW1	ARL4D (Human) IP-WB Antibody Pair	H00056672-PW1	C11orf17 (Human) IP-WB Antibody Pair
H00000410-PW1	ARSA (Human) IP-WB Antibody Pair	H00029902-PW1	C12orf24 (Human) IP-WB Antibody Pair
H00000419-PW1	ART3 (Human) IP-WB Antibody Pair	H00080183-PW1	C13orf18 (Human) IP-WB Antibody Pair
H00000419-PW2	ART3 (Human) IP-WB Antibody Pair	H00084299-PW1	C17orf37 (Human) IP-WB Antibody Pair
H00000420-PW1	ART4 (Human) IP-WB Antibody Pair	H00000716-PW1	C1S (Human) IP-WB Antibody Pair
H00140462-PW1	ASB9 (Human) IP-WB Antibody Pair	H00000725-PW1	C4BPB (Human) IP-WB Antibody Pair
H00008623-PW1	ASMTL (Human) IP-WB Antibody Pair	H00055319-PW1	C4orf43 (Human) IP-WB Antibody Pair
H00000440-PW1	ASNS (Human) IP-WB Antibody Pair	H00054981-PW1	C9orf95 (Human) IP-WB Antibody Pair
H00000440-PW2	ASNS (Human) IP-WB Antibody Pair	H00054981-PW2	C9orf95 (Human) IP-WB Antibody Pair
H00000443-PW1	ASPA (Human) IP-WB Antibody Pair	H00000759-PW1	CA1 (Human) IP-WB Antibody Pair
H00000444-PW1	ASPH (Human) IP-WB Antibody Pair	H00000771-PW1	CA12 (Human) IP-WB Antibody Pair
H00080150-PW1	ASRGL1 (Human) IP-WB Antibody Pair	H00000782-PW1	CACNB1 (Human) IP-WB Antibody Pair
H00080150-PW2	ASRGL1 (Human) IP-WB Antibody Pair	H00000784-PW1	CACNB3 (Human) IP-WB Antibody Pair
H00055252-PW1	ASXL2 (Human) IP-WB Antibody Pair	H00059285-PW1	CACNG6 (Human) IP-WB Antibody Pair
H00055210-PW1	ATAD3A (Human) IP-WB Antibody Pair	H00059284-PW1	CACNG7 (Human) IP-WB Antibody Pair
H00022926-PW1	ATF6 (Human) IP-WB Antibody Pair	H00010241-PW1	CALCOCO2 (Human) IP-WB Antibody Pair
H00055054-PW1	ATG16L1 (Human) IP-WB Antibody Pair	H00000810-PW1	CALML3 (Human) IP-WB Antibody Pair
H00010533-PW1	ATG7 (Human) IP-WB Antibody Pair	H00125972-PW1	CALR3 (Human) IP-WB Antibody Pair
H00000471-PW1	ATIC (Human) IP-WB Antibody Pair	H00057118-PW1	CAMK1D (Human) IP-WB Antibody Pair
H00000483-PW1	ATP1B3 (Human) IP-WB Antibody Pair	H00000816-PW1	CAMK2B (Human) IP-WB Antibody Pair
H00027032-PW1	ATP2C1 (Human) IP-WB Antibody Pair	H00000818-PW1	CAMK2G (Human) IP-WB Antibody Pair
H00245973-PW1	ATP6V1C2 (Human) IP-WB Antibody Pair	H00010645-PW1	CAMKK2 (Human) IP-WB Antibody Pair
H00090423-PW1	ATP6V1E2 (Human) IP-WB Antibody Pair	H00010645-PW2	CAMKK2 (Human) IP-WB Antibody Pair
H00051606-PW1	ATP6V1H (Human) IP-WB Antibody Pair	H00124583-PW1	CANT1 (Human) IP-WB Antibody Pair
H00000558-PW1	AXL (Human) IP-WB Antibody Pair	H00124583-PW2	CANT1 (Human) IP-WB Antibody Pair
H00008703-PW1	B4GALT3 (Human) IP-WB Antibody Pair	H00010487-PW1	CAP1 (Human) IP-WB Antibody Pair
H00008702-PW1	B4GALT4 (Human) IP-WB Antibody Pair	H00000828-PW1	CAPS (Human) IP-WB Antibody Pair
H00009532-PW1	BAG2 (Human) IP-WB Antibody Pair	H00093661-PW1	CAPZA3 (Human) IP-WB Antibody Pair
H00009530-PW1	BAG4 (Human) IP-WB Antibody Pair	H00022900-PW1	CARD8 (Human) IP-WB Antibody Pair
H00009529-PW1	BAG5 (Human) IP-WB Antibody Pair	H00000836-PW1	CASP3 (Human) IP-WB Antibody Pair
H00010458-PW1	BAIAP2 (Human) IP-WB Antibody Pair	H00000840-PW1	CASP7 (Human) IP-WB Antibody Pair
H00000578-PW1	BAK1 (Human) IP-WB Antibody Pair	H00000847-PW1	CAT (Human) IP-WB Antibody Pair
H00010538-PW1	BATF (Human) IP-WB Antibody Pair	H00000873-PW1	CBR1 (Human) IP-WB Antibody Pair
H00000585-PW1	BBS4 (Human) IP-WB Antibody Pair	H00000874-PW1	CBR3 (Human) IP-WB Antibody Pair
H00000587-PW1	BCAT2 (Human) IP-WB Antibody Pair	H00000875-PW1	CBS (Human) IP-WB Antibody Pair
H00079370-PW1	BCL2L14 (Human) IP-WB Antibody Pair	H00000875-PW2	CBS (Human) IP-WB Antibody Pair
H00000602-PW1	BCL3 (Human) IP-WB Antibody Pair	H00054862-PW1	CC2D1A (Human) IP-WB Antibody Pair
H00000617-PW1	BCS1L (Human) IP-WB Antibody Pair	H00147372-PW1	CCBE1 (Human) IP-WB Antibody Pair
H00000622-PW1	BDH1 (Human) IP-WB Antibody Pair	H00006359-PW1	CCL15 (Human) IP-WB Antibody Pair
H00000633-PW1	BGN (Human) IP-WB Antibody Pair	H00006362-PW1	CCL18 (Human) IP-WB Antibody Pair
H00023743-PW1	BHMT2 (Human) IP-WB Antibody Pair	H00006362-PW2	CCL18 (Human) IP-WB Antibody Pair
H00000330-PW1	BIRC3 (Human) IP-WB Antibody Pair	H00006366-PW1	CCL21 (Human) IP-WB Antibody Pair
H00000332-PW1	BIRC5 (Human) IP-WB Antibody Pair	H00006368-PW1	CCL23 (Human) IP-WB Antibody Pair
H00000332-PW2	BIRC5 (Human) IP-WB Antibody Pair	H00000890-PW1	CCNA2 (Human) IP-WB Antibody Pair
H00000644-PW1	BLVRA (Human) IP-WB Antibody Pair	H00000595-PW1	CCND1 (Human) IP-WB Antibody Pair
H00000655-PW1	BMP7 (Human) IP-WB Antibody Pair	H00000894-PW1	CCND2 (Human) IP-WB Antibody Pair
H00000658-PW1	BMPR1B (Human) IP-WB Antibody Pair	H00023552-PW1	CCRK (Human) IP-WB Antibody Pair
H00000665-PW1	BNIP3L (Human) IP-WB Antibody Pair	H00010693-PW1	CCT6B (Human) IP-WB Antibody Pair
H00066037-PW1	BOLL (Human) IP-WB Antibody Pair	H00010574-PW1	CCT7 (Human) IP-WB Antibody Pair
H00066037-PW2	BOLL (Human) IP-WB Antibody Pair	H00150160-PW1	CCT8L2 (Human) IP-WB Antibody Pair
H00006046-PW1	BRD2 (Human) IP-WB Antibody Pair	H00000929-PW1	CD14 (Human) IP-WB Antibody Pair
H00025855-PW1	BRMS1 (Human) IP-WB Antibody Pair	H00011126-PW1	CD160 (Human) IP-WB Antibody Pair
H00056853-PW1	BRUNOL4 (Human) IP-WB Antibody Pair	H00057126-PW1	CD177 (Human) IP-WB Antibody Pair
H00000682-PW1	BSG (Human) IP-WB Antibody Pair	H00000930-PW1	CD19 (Human) IP-WB Antibody Pair
H00007809-PW1	BSND (Human) IP-WB Antibody Pair	H00000910-PW1	CD1B (Human) IP-WB Antibody Pair
H00000683-PW1	BST1 (Human) IP-WB Antibody Pair	H00131450-PW1	CD200R1 (Human) IP-WB Antibody Pair
H00000684-PW1	BST2 (Human) IP-WB Antibody Pair	H00050489-PW1	CD207 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00030835-PW1	CD209 (Human) IP-WB Antibody Pair	H00001082-PW2	CGB (Human) IP-WB Antibody Pair
H00030835-PW2	CD209 (Human) IP-WB Antibody Pair	H00011200-PW1	CHEK2 (Human) IP-WB Antibody Pair
H00000919-PW1	CD247 (Human) IP-WB Antibody Pair	H00001149-PW1	CIDEA (Human) IP-WB Antibody Pair
H00000939-PW1	CD27 (Human) IP-WB Antibody Pair	H00004435-PW1	CITED1 (Human) IP-WB Antibody Pair
H00001604-PW1	CD55 (Human) IP-WB Antibody Pair	H00001158-PW1	CKM (Human) IP-WB Antibody Pair
H00000922-PW1	CD5L (Human) IP-WB Antibody Pair	H00001158-PW2	CKM (Human) IP-WB Antibody Pair
H00000969-PW1	CD69 (Human) IP-WB Antibody Pair	H00001159-PW1	CKMT1B (Human) IP-WB Antibody Pair
H00000924-PW1	CD7 (Human) IP-WB Antibody Pair	H00009076-PW1	CLDN1 (Human) IP-WB Antibody Pair
H00000970-PW1	CD70 (Human) IP-WB Antibody Pair	H00009076-PW2	CLDN1 (Human) IP-WB Antibody Pair
H00000972-PW1	CD74 (Human) IP-WB Antibody Pair	H00009076-PW3	CLDN1 (Human) IP-WB Antibody Pair
H00000972-PW2	CD74 (Human) IP-WB Antibody Pair	H00023562-PW1	CLDN14 (Human) IP-WB Antibody Pair
H00000974-PW1	CD79B (Human) IP-WB Antibody Pair	H00001364-PW1	CLDN4 (Human) IP-WB Antibody Pair
H00000941-PW1	CD80 (Human) IP-WB Antibody Pair	H00051266-PW1	CLEC1B (Human) IP-WB Antibody Pair
H00000941-PW2	CD80 (Human) IP-WB Antibody Pair	H00009022-PW1	CLIC3 (Human) IP-WB Antibody Pair
H00000941-PW3	CD80 (Human) IP-WB Antibody Pair	H00001208-PW1	CLPS (Human) IP-WB Antibody Pair
H00009308-PW1	CD83 (Human) IP-WB Antibody Pair	H00001208-PW2	CLPS (Human) IP-WB Antibody Pair
H00000942-PW1	CD86 (Human) IP-WB Antibody Pair	H00001212-PW1	CLTB (Human) IP-WB Antibody Pair
H00000925-PW1	CD8A (Human) IP-WB Antibody Pair	H00084735-PW1	CNDP1 (Human) IP-WB Antibody Pair
H00010225-PW1	CD96 (Human) IP-WB Antibody Pair	H00154043-PW1	CNKSR3 (Human) IP-WB Antibody Pair
H00000978-PW1	CDA (Human) IP-WB Antibody Pair	H00154043-PW2	CNKSR3 (Human) IP-WB Antibody Pair
H00081602-PW1	CDADC1 (Human) IP-WB Antibody Pair	H00001267-PW1	CNP (Human) IP-WB Antibody Pair
H00081602-PW2	CDADC1 (Human) IP-WB Antibody Pair	H00001270-PW1	CNTF (Human) IP-WB Antibody Pair
H00008872-PW1	CDC123 (Human) IP-WB Antibody Pair	H00080347-PW1	COASY (Human) IP-WB Antibody Pair
H00008555-PW1	CDC14B (Human) IP-WB Antibody Pair	H00091522-PW1	COL23A1 (Human) IP-WB Antibody Pair
H00000993-PW1	CDC25A (Human) IP-WB Antibody Pair	H00091522-PW2	COL23A1 (Human) IP-WB Antibody Pair
H00000993-PW2	CDC25A (Human) IP-WB Antibody Pair	H00150684-PW1	COMMD1 (Human) IP-WB Antibody Pair
H00000995-PW1	CDC25C (Human) IP-WB Antibody Pair	H00150684-PW2	COMMD1 (Human) IP-WB Antibody Pair
H00000995-PW2	CDC25C (Human) IP-WB Antibody Pair	H00001312-PW1	COMT (Human) IP-WB Antibody Pair
H00023097-PW1	CDC2L6 (Human) IP-WB Antibody Pair	H00084701-PW1	COX4I2 (Human) IP-WB Antibody Pair
H00008318-PW1	CDC45L (Human) IP-WB Antibody Pair	H00001357-PW1	CPA1 (Human) IP-WB Antibody Pair
H00008318-PW2	CDC45L (Human) IP-WB Antibody Pair	H00001357-PW2	CPA1 (Human) IP-WB Antibody Pair
H00055143-PW1	CDCA8 (Human) IP-WB Antibody Pair	H00001358-PW1	CPA2 (Human) IP-WB Antibody Pair
H00001012-PW1	CDH13 (Human) IP-WB Antibody Pair	H00001358-PW2	CPA2 (Human) IP-WB Antibody Pair
H00001017-PW1	CDK2 (Human) IP-WB Antibody Pair	H00001360-PW1	CPB1 (Human) IP-WB Antibody Pair
H00001020-PW1	CDK5 (Human) IP-WB Antibody Pair	H00001360-PW2	CPB1 (Human) IP-WB Antibody Pair
H00001021-PW1	CDK6 (Human) IP-WB Antibody Pair	H00001361-PW1	CPB2 (Human) IP-WB Antibody Pair
H00055602-PW1	CDKN2AIP (Human) IP-WB Antibody Pair	H00008904-PW1	CPNE1 (Human) IP-WB Antibody Pair
H00001032-PW1	CDKN2D (Human) IP-WB Antibody Pair	H00009586-PW1	CREB5 (Human) IP-WB Antibody Pair
H00090273-PW1	CEACAM21 (Human) IP-WB Antibody Pair	H00009586-PW2	CREB5 (Human) IP-WB Antibody Pair
H00001054-PW1	CEBPG (Human) IP-WB Antibody Pair	H00001393-PW1	CRHBP (Human) IP-WB Antibody Pair
H00001054-PW2	CEBPG (Human) IP-WB Antibody Pair	H00001393-PW2	CRHBP (Human) IP-WB Antibody Pair
H00051816-PW1	CECR1 (Human) IP-WB Antibody Pair	H00009244-PW1	CRLF1 (Human) IP-WB Antibody Pair
H00001059-PW1	CENPB (Human) IP-WB Antibody Pair	H00051379-PW1	CRLF3 (Human) IP-WB Antibody Pair
H00055835-PW1	CENPJ (Human) IP-WB Antibody Pair	H00001421-PW1	CRYGD (Human) IP-WB Antibody Pair
H00055835-PW2	CENPJ (Human) IP-WB Antibody Pair	H00001428-PW1	CRYM (Human) IP-WB Antibody Pair
H00079019-PW1	CENPM (Human) IP-WB Antibody Pair	H00001428-PW2	CRYM (Human) IP-WB Antibody Pair
H00079019-PW2	CENPM (Human) IP-WB Antibody Pair	H00001428-PW3	CRYM (Human) IP-WB Antibody Pair
H00079019-PW3	CENPM (Human) IP-WB Antibody Pair	H00001429-PW1	CRYZ (Human) IP-WB Antibody Pair
H00055839-PW1	CENPN (Human) IP-WB Antibody Pair	H00008531-PW1	CSDA (Human) IP-WB Antibody Pair
H00055165-PW1	CEP55 (Human) IP-WB Antibody Pair	H00001440-PW1	CSF3 (Human) IP-WB Antibody Pair
H00009702-PW1	CEP57 (Human) IP-WB Antibody Pair	H00001440-PW2	CSF3 (Human) IP-WB Antibody Pair
H00009350-PW1	CER1 (Human) IP-WB Antibody Pair	H00055790-PW1	CSGALNACT1 (Human) IP-WB Antibody Pair
H00009350-PW2	CER1 (Human) IP-WB Antibody Pair	H00001445-PW1	CSK (Human) IP-WB Antibody Pair
H00008824-PW1	CES2 (Human) IP-WB Antibody Pair	H00001446-PW1	CSN1S1 (Human) IP-WB Antibody Pair
H00008824-PW2	CES2 (Human) IP-WB Antibody Pair	H00001478-PW1	CSTF2 (Human) IP-WB Antibody Pair
H00001068-PW1	CETN1 (Human) IP-WB Antibody Pair	H00001478-PW2	CSTF2 (Human) IP-WB Antibody Pair
H00055997-PW1	CFC1 (Human) IP-WB Antibody Pair	H00001486-PW1	CTBS (Human) IP-WB Antibody Pair
H00055997-PW2	CFC1 (Human) IP-WB Antibody Pair	H00001491-PW1	CTH (Human) IP-WB Antibody Pair
H00003078-PW1	CFHR1 (Human) IP-WB Antibody Pair	H00029119-PW1	CTNNA3 (Human) IP-WB Antibody Pair
H00003080-PW1	CFHR2 (Human) IP-WB Antibody Pair	H00001509-PW1	CTSD (Human) IP-WB Antibody Pair
H00003080-PW2	CFHR2 (Human) IP-WB Antibody Pair	H00001509-PW2	CTSD (Human) IP-WB Antibody Pair
H00081494-PW1	CFHR5 (Human) IP-WB Antibody Pair	H00001520-PW1	CTSS (Human) IP-WB Antibody Pair
H00008837-PW1	CFLAR (Human) IP-WB Antibody Pair	H00001520-PW2	CTSS (Human) IP-WB Antibody Pair
H00008837-PW2	CFLAR (Human) IP-WB Antibody Pair	H00001523-PW1	CUX1 (Human) IP-WB Antibody Pair
H00001082-PW1	CGB (Human) IP-WB Antibody Pair	H00001523-PW2	CUX1 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00006373-PW1	CXCL11 (Human) IP-WB Antibody Pair	H00001892-PW1	ECHS1 (Human) IP-WB Antibody Pair
H00010563-PW1	CXCL13 (Human) IP-WB Antibody Pair	H00055741-PW1	EDEM2 (Human) IP-WB Antibody Pair
H00001528-PW1	CYB5A (Human) IP-WB Antibody Pair	H00001906-PW1	EDN1 (Human) IP-WB Antibody Pair
H00001543-PW1	CYP1A1 (Human) IP-WB Antibody Pair	H00001908-PW1	EDN3 (Human) IP-WB Antibody Pair
H00001593-PW1	CYP27A1 (Human) IP-WB Antibody Pair	H00001937-PW1	EEF1G (Human) IP-WB Antibody Pair
H00001571-PW1	CYP2E1 (Human) IP-WB Antibody Pair	H00030008-PW1	EFEMP2 (Human) IP-WB Antibody Pair
H00001576-PW1	CYP3A4 (Human) IP-WB Antibody Pair	H00030008-PW2	EFEMP2 (Human) IP-WB Antibody Pair
H00064816-PW1	CYP3A43 (Human) IP-WB Antibody Pair	H00001944-PW1	EFNA3 (Human) IP-WB Antibody Pair
H00010858-PW1	CYP46A1 (Human) IP-WB Antibody Pair	H00001944-PW2	EFNA3 (Human) IP-WB Antibody Pair
H00001580-PW1	CYP4B1 (Human) IP-WB Antibody Pair	H00051162-PW1	EGFL7 (Human) IP-WB Antibody Pair
H00009266-PW1	CYTH2 (Human) IP-WB Antibody Pair	H00008892-PW1	EIF2B2 (Human) IP-WB Antibody Pair
H00001620-PW1	DBC1 (Human) IP-WB Antibody Pair	H00008672-PW1	EIF4G3 (Human) IP-WB Antibody Pair
H00001633-PW1	DCK (Human) IP-WB Antibody Pair	H00056648-PW1	EIF5A2 (Human) IP-WB Antibody Pair
H00055802-PW1	DCP1A (Human) IP-WB Antibody Pair	H00010136-PW1	ELA3A (Human) IP-WB Antibody Pair
H00054165-PW1	DCUN1D1 (Human) IP-WB Antibody Pair	H00055520-PW1	ELAC1 (Human) IP-WB Antibody Pair
H00115265-PW1	DDIT4L (Human) IP-WB Antibody Pair	H00001994-PW1	ELAVL1 (Human) IP-WB Antibody Pair
H00001652-PW1	DDT (Human) IP-WB Antibody Pair	H00001994-PW2	ELAVL1 (Human) IP-WB Antibody Pair
H00168400-PW1	DDX53 (Human) IP-WB Antibody Pair	H00001997-PW1	ELF1 (Human) IP-WB Antibody Pair
H00079039-PW1	DDX54 (Human) IP-WB Antibody Pair	H00002006-PW1	ELN (Human) IP-WB Antibody Pair
H00010522-PW1	DEAF1 (Human) IP-WB Antibody Pair	H00011117-PW1	EMILIN1 (Human) IP-WB Antibody Pair
H00026063-PW1	DECR2 (Human) IP-WB Antibody Pair	H00002022-PW1	ENG (Human) IP-WB Antibody Pair
H00091614-PW1	DEPDC7 (Human) IP-WB Antibody Pair	H00002022-PW2	ENG (Human) IP-WB Antibody Pair
H00091614-PW2	DEPDC7 (Human) IP-WB Antibody Pair	H00002023-PW1	ENO1 (Human) IP-WB Antibody Pair
H00001676-PW1	DFFA (Human) IP-WB Antibody Pair	H00023136-PW1	EPB41L3 (Human) IP-WB Antibody Pair
H00001718-PW1	DHCR24 (Human) IP-WB Antibody Pair	H00004072-PW1	EPCAM (Human) IP-WB Antibody Pair
H00079947-PW1	DHDDS (Human) IP-WB Antibody Pair	H00054749-PW1	EPDR1 (Human) IP-WB Antibody Pair
H00027294-PW1	DHDH (Human) IP-WB Antibody Pair	H00002045-PW1	EPHA7 (Human) IP-WB Antibody Pair
H00010170-PW1	DHRS9 (Human) IP-WB Antibody Pair	H00002053-PW1	EPHX2 (Human) IP-WB Antibody Pair
H00010170-PW2	DHRS9 (Human) IP-WB Antibody Pair	H00007957-PW1	EPM2A (Human) IP-WB Antibody Pair
H00079132-PW1	DHX58 (Human) IP-WB Antibody Pair	H00002057-PW1	EPOR (Human) IP-WB Antibody Pair
H00022943-PW1	DKK1 (Human) IP-WB Antibody Pair	H00054821-PW1	ERCC6L (Human) IP-WB Antibody Pair
H00022943-PW2	DKK1 (Human) IP-WB Antibody Pair	H00001161-PW1	ERCC8 (Human) IP-WB Antibody Pair
H00010395-PW1	DLC1 (Human) IP-WB Antibody Pair	H00114625-PW1	ERMAP (Human) IP-WB Antibody Pair
H00001743-PW1	DLST (Human) IP-WB Antibody Pair	H00114625-PW2	ERMAP (Human) IP-WB Antibody Pair
H00001745-PW1	DLX1 (Human) IP-WB Antibody Pair	H00090952-PW1	ESAM (Human) IP-WB Antibody Pair
H00001745-PW2	DLX1 (Human) IP-WB Antibody Pair	H00002098-PW1	ESD (Human) IP-WB Antibody Pair
H00001745-PW3	DLX1 (Human) IP-WB Antibody Pair	H00011082-PW1	ESM1 (Human) IP-WB Antibody Pair
H00001745-PW4	DLX1 (Human) IP-WB Antibody Pair	H00002101-PW1	ESRRA (Human) IP-WB Antibody Pair
H00001748-PW1	DLX4 (Human) IP-WB Antibody Pair	H00002104-PW2	ESRRG (Human) IP-WB Antibody Pair
H00001756-PW1	DMD (Human) IP-WB Antibody Pair	H00051513-PW1	ETV7 (Human) IP-WB Antibody Pair
H00001760-PW1	DMPK (Human) IP-WB Antibody Pair	H00009156-PW1	EXO1 (Human) IP-WB Antibody Pair
H00064446-PW1	DNAI2 (Human) IP-WB Antibody Pair	H00051010-PW1	EXOSC3 (Human) IP-WB Antibody Pair
H00003337-PW1	DNAJB1 (Human) IP-WB Antibody Pair	H00051010-PW2	EXOSC3 (Human) IP-WB Antibody Pair
H00003300-PW1	DNAJB2 (Human) IP-WB Antibody Pair	H00056915-PW1	EXOSC5 (Human) IP-WB Antibody Pair
H00011080-PW1	DNAJB4 (Human) IP-WB Antibody Pair	H00007430-PW1	EZR (Human) IP-WB Antibody Pair
H00010049-PW1	DNAJB6 (Human) IP-WB Antibody Pair	H00002160-PW1	F11 (Human) IP-WB Antibody Pair
H00010049-PW2	DNAJB6 (Human) IP-WB Antibody Pair	H00002160-PW2	F11 (Human) IP-WB Antibody Pair
H00010049-PW3	DNAJB6 (Human) IP-WB Antibody Pair	H00050848-PW1	F11R (Human) IP-WB Antibody Pair
H00007266-PW1	DNAJC7 (Human) IP-WB Antibody Pair	H00002147-PW1	F2 (Human) IP-WB Antibody Pair
H00007802-PW1	DNAL1 (Human) IP-WB Antibody Pair	H00002147-PW2	F2 (Human) IP-WB Antibody Pair
H00007802-PW2	DNAL1 (Human) IP-WB Antibody Pair	H00002149-PW1	F2R (Human) IP-WB Antibody Pair
H00001777-PW1	DNASE2 (Human) IP-WB Antibody Pair	H00002152-PW1	F3 (Human) IP-WB Antibody Pair
H00001788-PW1	DNMT3A (Human) IP-WB Antibody Pair	H00002152-PW2	F3 (Human) IP-WB Antibody Pair
H00029947-PW1	DNMT3L (Human) IP-WB Antibody Pair	H00002155-PW1	F7 (Human) IP-WB Antibody Pair
H00001796-PW1	DOK1 (Human) IP-WB Antibody Pair	H00002157-PW1	F8 (Human) IP-WB Antibody Pair
H00056896-PW1	DPYSL5 (Human) IP-WB Antibody Pair	H00002158-PW1	F9 (Human) IP-WB Antibody Pair
H00011266-PW1	DUSP12 (Human) IP-WB Antibody Pair	H00002168-PW1	FABP1 (Human) IP-WB Antibody Pair
H00056940-PW1	DUSP22 (Human) IP-WB Antibody Pair	H00002167-PW1	FABP4 (Human) IP-WB Antibody Pair
H00001856-PW1	DVL2 (Human) IP-WB Antibody Pair	H00002173-PW1	FABP7 (Human) IP-WB Antibody Pair
H00001857-PW1	DVL3 (Human) IP-WB Antibody Pair	H00054097-PW1	FAM3B (Human) IP-WB Antibody Pair
H00001783-PW1	DYNC1L12 (Human) IP-WB Antibody Pair	H00002176-PW1	FANCC (Human) IP-WB Antibody Pair
H00161582-PW1	DYX1C1 (Human) IP-WB Antibody Pair	H00002189-PW1	FANCG (Human) IP-WB Antibody Pair
H00161582-PW2	DYX1C1 (Human) IP-WB Antibody Pair	H00002189-PW2	FANCG (Human) IP-WB Antibody Pair
H00001876-PW1	E2F6 (Human) IP-WB Antibody Pair	H00055120-PW1	FANCL (Human) IP-WB Antibody Pair
H00001891-PW1	ECH1 (Human) IP-WB Antibody Pair	H00002191-PW1	FAP (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00010667-PW1	FARS2 (Human) IP-WB Antibody Pair	H00002553-PW1	GABPB1 (Human) IP-WB Antibody Pair
H00000355-PW1	FAS (Human) IP-WB Antibody Pair	H00002572-PW1	GAD2 (Human) IP-WB Antibody Pair
H00000355-PW2	FAS (Human) IP-WB Antibody Pair	H00002580-PW1	GAK (Human) IP-WB Antibody Pair
H00000355-PW3	FAS (Human) IP-WB Antibody Pair	H00002582-PW1	GALE (Human) IP-WB Antibody Pair
H00000356-PW1	FASLG (Human) IP-WB Antibody Pair	H00002588-PW1	GALNS (Human) IP-WB Antibody Pair
H00002194-PW1	FASN (Human) IP-WB Antibody Pair	H00011226-PW1	GALNT6 (Human) IP-WB Antibody Pair
H00002194-PW2	FASN (Human) IP-WB Antibody Pair	H00050614-PW1	GALNT9 (Human) IP-WB Antibody Pair
H00002192-PW1	FBLN1 (Human) IP-WB Antibody Pair	H00002592-PW1	GALT (Human) IP-WB Antibody Pair
H00002192-PW2	FBLN1 (Human) IP-WB Antibody Pair	H00023193-PW1	GANAB (Human) IP-WB Antibody Pair
H00055336-PW1	FBXL8 (Human) IP-WB Antibody Pair	H00002597-PW1	GAPDH (Human) IP-WB Antibody Pair
H00023219-PW1	FBXO28 (Human) IP-WB Antibody Pair	H00026330-PW1	GAPDHS (Human) IP-WB Antibody Pair
H00026272-PW1	FBXO4 (Human) IP-WB Antibody Pair	H00002617-PW1	GARS (Human) IP-WB Antibody Pair
H00026272-PW2	FBXO4 (Human) IP-WB Antibody Pair	H00002618-PW1	GART (Human) IP-WB Antibody Pair
H00054455-PW1	FBXO42 (Human) IP-WB Antibody Pair	H00002618-PW2	GART (Human) IP-WB Antibody Pair
H00026270-PW1	FBXO6 (Human) IP-WB Antibody Pair	H00283431-PW1	GAS2L3 (Human) IP-WB Antibody Pair
H00002208-PW1	FCER2 (Human) IP-WB Antibody Pair	H00002621-PW1	GAS6 (Human) IP-WB Antibody Pair
H00002208-PW2	FCER2 (Human) IP-WB Antibody Pair	H00002623-PW1	GATA1 (Human) IP-WB Antibody Pair
H00002213-PW1	FCGR2B (Human) IP-WB Antibody Pair	H00002629-PW1	GBA (Human) IP-WB Antibody Pair
H00002213-PW2	FCGR2B (Human) IP-WB Antibody Pair	H00002632-PW1	GBE1 (Human) IP-WB Antibody Pair
H00002219-PW1	FCN1 (Human) IP-WB Antibody Pair	H00064223-PW1	GBL (Human) IP-WB Antibody Pair
H00002220-PW1	FCN2 (Human) IP-WB Antibody Pair	H00002633-PW1	GBP1 (Human) IP-WB Antibody Pair
H00002220-PW2	FCN2 (Human) IP-WB Antibody Pair	H00002634-PW1	GBP2 (Human) IP-WB Antibody Pair
H00008547-PW1	FCN3 (Human) IP-WB Antibody Pair	H00002634-PW2	GBP2 (Human) IP-WB Antibody Pair
H00002224-PW1	FDPS (Human) IP-WB Antibody Pair	H00002643-PW1	GCH1 (Human) IP-WB Antibody Pair
H00055612-PW1	FERMT1 (Human) IP-WB Antibody Pair	H00002730-PW1	GCLM (Human) IP-WB Antibody Pair
H00002243-PW1	FGA (Human) IP-WB Antibody Pair	H00079760-PW1	GEMIN7 (Human) IP-WB Antibody Pair
H00089846-PW1	FGD3 (Human) IP-WB Antibody Pair	H00002670-PW1	GFAP (Human) IP-WB Antibody Pair
H00002246-PW1	FGF1 (Human) IP-WB Antibody Pair	H00002670-PW2	GFAP (Human) IP-WB Antibody Pair
H00002246-PW2	FGF1 (Human) IP-WB Antibody Pair	H00009453-PW1	GGPS1 (Human) IP-WB Antibody Pair
H00009965-PW1	FGF19 (Human) IP-WB Antibody Pair	H00002678-PW1	GGT1 (Human) IP-WB Antibody Pair
H00026291-PW1	FGF21 (Human) IP-WB Antibody Pair	H00002688-PW1	GH1 (Human) IP-WB Antibody Pair
H00026291-PW2	FGF21 (Human) IP-WB Antibody Pair	H00002694-PW1	GIF (Human) IP-WB Antibody Pair
H00026291-PW3	FGF21 (Human) IP-WB Antibody Pair	H00055340-PW1	GIMAP5 (Human) IP-WB Antibody Pair
H00002250-PW1	FGF5 (Human) IP-WB Antibody Pair	H00055340-PW2	GIMAP5 (Human) IP-WB Antibody Pair
H00002267-PW1	FGL1 (Human) IP-WB Antibody Pair	H00055340-PW3	GIMAP5 (Human) IP-WB Antibody Pair
H00010875-PW1	FGL2 (Human) IP-WB Antibody Pair	H00009837-PW1	GINS1 (Human) IP-WB Antibody Pair
H00002268-PW1	FGR (Human) IP-WB Antibody Pair	H00009815-PW1	GIT2 (Human) IP-WB Antibody Pair
H00009457-PW1	FHL5 (Human) IP-WB Antibody Pair	H00002717-PW1	GLA (Human) IP-WB Antibody Pair
H00002277-PW1	FIGF (Human) IP-WB Antibody Pair	H00002720-PW1	GLB1 (Human) IP-WB Antibody Pair
H00051303-PW1	FKBP11 (Human) IP-WB Antibody Pair	H00002720-PW2	GLB1 (Human) IP-WB Antibody Pair
H00002288-PW1	FKBP4 (Human) IP-WB Antibody Pair	H00026035-PW1	GLCE (Human) IP-WB Antibody Pair
H00002288-PW2	FKBP4 (Human) IP-WB Antibody Pair	H00011146-PW1	GLMN (Human) IP-WB Antibody Pair
H00002328-PW1	FMO3 (Human) IP-WB Antibody Pair	H00011146-PW2	GLMN (Human) IP-WB Antibody Pair
H00064122-PW1	FN3K (Human) IP-WB Antibody Pair	H00002739-PW1	GLO1 (Human) IP-WB Antibody Pair
H00002354-PW1	FOSB (Human) IP-WB Antibody Pair	H00002746-PW1	GLUD1 (Human) IP-WB Antibody Pair
H00002355-PW2	FOSL2 (Human) IP-WB Antibody Pair	H00027232-PW1	GNMT (Human) IP-WB Antibody Pair
H00003170-PW1	FOXA2 (Human) IP-WB Antibody Pair	H00008443-PW1	GNPAT (Human) IP-WB Antibody Pair
H00003170-PW2	FOXA2 (Human) IP-WB Antibody Pair	H00051280-PW1	GOLM1 (Human) IP-WB Antibody Pair
H00003171-PW1	FOXA3 (Human) IP-WB Antibody Pair	H00009527-PW1	GOSR1 (Human) IP-WB Antibody Pair
H00050943-PW1	FOXP3 (Human) IP-WB Antibody Pair	H00002805-PW1	GOT1 (Human) IP-WB Antibody Pair
H00002358-PW1	FPR2 (Human) IP-WB Antibody Pair	H00002811-PW1	GP1BA (Human) IP-WB Antibody Pair
H00002483-PW1	FRG1 (Human) IP-WB Antibody Pair	H00002262-PW1	GPC5 (Human) IP-WB Antibody Pair
H00002483-PW2	FRG1 (Human) IP-WB Antibody Pair	H00002262-PW2	GPC5 (Human) IP-WB Antibody Pair
H00010818-PW1	FRS2 (Human) IP-WB Antibody Pair	H00002262-PW3	GPC5 (Human) IP-WB Antibody Pair
H00079187-PW1	FSD1 (Human) IP-WB Antibody Pair	H00023171-PW1	GPD1L (Human) IP-WB Antibody Pair
H00079187-PW2	FSD1 (Human) IP-WB Antibody Pair	H00002821-PW1	GPI (Human) IP-WB Antibody Pair
H00002512-PW2	FTL (Human) IP-WB Antibody Pair	H00011321-PW1	GPN1 (Human) IP-WB Antibody Pair
H00002512-PW3	FTL (Human) IP-WB Antibody Pair	H00131601-PW1	GPR175 (Human) IP-WB Antibody Pair
H00002523-PW1	FUT1 (Human) IP-WB Antibody Pair	H00084706-PW1	GPT2 (Human) IP-WB Antibody Pair
H00002524-PW1	FUT2 (Human) IP-WB Antibody Pair	H00002880-PW1	GPX5 (Human) IP-WB Antibody Pair
H00002395-PW1	FXN (Human) IP-WB Antibody Pair	H00009402-PW1	GRAP2 (Human) IP-WB Antibody Pair
H00002395-PW2	FXN (Human) IP-WB Antibody Pair	H00009402-PW2	GRAP2 (Human) IP-WB Antibody Pair
H00007855-PW1	FZD5 (Human) IP-WB Antibody Pair	H00002886-PW1	GRB7 (Human) IP-WB Antibody Pair
H00050486-PW1	G0S2 (Human) IP-WB Antibody Pair	H00002886-PW2	GRB7 (Human) IP-WB Antibody Pair
H00010146-PW1	G3BP1 (Human) IP-WB Antibody Pair	H00064388-PW1	GREM2 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00002905-PW1	GRIN2C (Human) IP-WB Antibody Pair	H00003109-PW4	HLA-DMB (Human) IP-WB Antibody Pair
H00002870-PW1	GRK6 (Human) IP-WB Antibody Pair	H00003112-PW1	HLA-DOB (Human) IP-WB Antibody Pair
H00145258-PW1	GSC (Human) IP-WB Antibody Pair	H00003119-PW1	HLA-DQB1 (Human) IP-WB Antibody Pair
H00145258-PW2	GSC (Human) IP-WB Antibody Pair	H00003127-PW1	HLA-DRB5 (Human) IP-WB Antibody Pair
H00002934-PW1	GSN (Human) IP-WB Antibody Pair	H00003133-PW1	HLA-E (Human) IP-WB Antibody Pair
H00002936-PW1	GSR (Human) IP-WB Antibody Pair	H00003141-PW1	HLCS (Human) IP-WB Antibody Pair
H00002936-PW2	GSR (Human) IP-WB Antibody Pair	H00003142-PW1	HLX (Human) IP-WB Antibody Pair
H00002938-PW1	GSTA1 (Human) IP-WB Antibody Pair	H00003142-PW2	HLX (Human) IP-WB Antibody Pair
H00009446-PW1	GSTO1 (Human) IP-WB Antibody Pair	H00010362-PW1	HMG20B (Human) IP-WB Antibody Pair
H00002950-PW1	GSTP1 (Human) IP-WB Antibody Pair	H00003146-PW1	HMGB1 (Human) IP-WB Antibody Pair
H00002953-PW1	GSTT2 (Human) IP-WB Antibody Pair	H00003155-PW1	HMGCL (Human) IP-WB Antibody Pair
H00002954-PW1	GSTZ1 (Human) IP-WB Antibody Pair	H00003172-PW1	HNF4A (Human) IP-WB Antibody Pair
H00011036-PW1	GTF2A1L (Human) IP-WB Antibody Pair	H00009456-PW1	HOMER1 (Human) IP-WB Antibody Pair
H00011036-PW2	GTF2A1L (Human) IP-WB Antibody Pair	H00009455-PW1	HOMER2 (Human) IP-WB Antibody Pair
H00011036-PW3	GTF2A1L (Human) IP-WB Antibody Pair	H00003198-PW1	HOXA1 (Human) IP-WB Antibody Pair
H00002960-PW1	GTF2E1 (Human) IP-WB Antibody Pair	H00003211-PW1	HOXB1 (Human) IP-WB Antibody Pair
H00002965-PW1	GTF2H1 (Human) IP-WB Antibody Pair	H00003211-PW2	HOXB1 (Human) IP-WB Antibody Pair
H00002966-PW1	GTF2H2 (Human) IP-WB Antibody Pair	H00003211-PW3	HOXB1 (Human) IP-WB Antibody Pair
H00002967-PW1	GTF2H3 (Human) IP-WB Antibody Pair	H00003219-PW1	HOXB9 (Human) IP-WB Antibody Pair
H00002968-PW1	GTF2H4 (Human) IP-WB Antibody Pair	H00003219-PW2	HOXB9 (Human) IP-WB Antibody Pair
H00002969-PW1	GTF2I (Human) IP-WB Antibody Pair	H00003219-PW3	HOXB9 (Human) IP-WB Antibody Pair
H00009569-PW1	GTF2IRD1 (Human) IP-WB Antibody Pair	H00003224-PW1	HOXC8 (Human) IP-WB Antibody Pair
H00051512-PW1	GTSE1 (Human) IP-WB Antibody Pair	H00003231-PW1	HOXD1 (Human) IP-WB Antibody Pair
H00009626-PW1	GUCA1C (Human) IP-WB Antibody Pair	H00003231-PW2	HOXD1 (Human) IP-WB Antibody Pair
H00009626-PW2	GUCA1C (Human) IP-WB Antibody Pair	H00003234-PW1	HOXD8 (Human) IP-WB Antibody Pair
H00002982-PW1	GUCY1A3 (Human) IP-WB Antibody Pair	H00003248-PW1	HPGD (Human) IP-WB Antibody Pair
H00002987-PW1	GUK1 (Human) IP-WB Antibody Pair	H00003263-PW1	HPX (Human) IP-WB Antibody Pair
H00002987-PW2	GUK1 (Human) IP-WB Antibody Pair	H00003273-PW1	HRG (Human) IP-WB Antibody Pair
H00002987-PW3	GUK1 (Human) IP-WB Antibody Pair	H00009653-PW1	HS2ST1 (Human) IP-WB Antibody Pair
H00003002-PW1	GZMB (Human) IP-WB Antibody Pair	H00003291-PW1	HSD11B2 (Human) IP-WB Antibody Pair
H00003026-PW1	HABP2 (Human) IP-WB Antibody Pair	H00003294-PW1	HSD17B2 (Human) IP-WB Antibody Pair
H00003033-PW1	HADH (Human) IP-WB Antibody Pair	H00051478-PW1	HSD17B7 (Human) IP-WB Antibody Pair
H00003033-PW2	HADH (Human) IP-WB Antibody Pair	H00007923-PW1	HSD17B8 (Human) IP-WB Antibody Pair
H00051179-PW1	HAO2 (Human) IP-WB Antibody Pair	H00003283-PW1	HSD3B1 (Human) IP-WB Antibody Pair
H00003035-PW1	HARS (Human) IP-WB Antibody Pair	H00003283-PW2	HSD3B1 (Human) IP-WB Antibody Pair
H00026762-PW1	HAVCR1 (Human) IP-WB Antibody Pair	H00003305-PW1	HSPA1L (Human) IP-WB Antibody Pair
H00010767-PW1	HBS1L (Human) IP-WB Antibody Pair	H00003306-PW1	HSPA2 (Human) IP-WB Antibody Pair
H00003052-PW1	HCCS (Human) IP-WB Antibody Pair	H00003315-PW1	HSPB1 (Human) IP-WB Antibody Pair
H00003055-PW1	HCK (Human) IP-WB Antibody Pair	H00023640-PW1	HSPBP1 (Human) IP-WB Antibody Pair
H00003055-PW2	HCK (Human) IP-WB Antibody Pair	H00051504-PW1	HSPC152 (Human) IP-WB Antibody Pair
H00003059-PW1	HCLS1 (Human) IP-WB Antibody Pair	H00003329-PW1	HSPD1 (Human) IP-WB Antibody Pair
H00003062-PW1	HCRT2R (Human) IP-WB Antibody Pair	H00003364-PW1	HUS1 (Human) IP-WB Antibody Pair
H00003065-PW1	HDAC1 (Human) IP-WB Antibody Pair	H00003364-PW2	HUS1 (Human) IP-WB Antibody Pair
H00079885-PW1	HDAC11 (Human) IP-WB Antibody Pair	H00003373-PW1	HYAL1 (Human) IP-WB Antibody Pair
H00008841-PW1	HDAC3 (Human) IP-WB Antibody Pair	H00008372-PW1	HYAL3 (Human) IP-WB Antibody Pair
H00008841-PW2	HDAC3 (Human) IP-WB Antibody Pair	H00003383-PW1	ICAM1 (Human) IP-WB Antibody Pair
H00051564-PW1	HDAC7 (Human) IP-WB Antibody Pair	H00003383-PW2	ICAM1 (Human) IP-WB Antibody Pair
H00055869-PW1	HDAC8 (Human) IP-WB Antibody Pair	H00003386-PW1	ICAM4 (Human) IP-WB Antibody Pair
H00055869-PW2	HDAC8 (Human) IP-WB Antibody Pair	H00003418-PW1	IDH2 (Human) IP-WB Antibody Pair
H00009709-PW1	HERPUD1 (Human) IP-WB Antibody Pair	H00003423-PW1	IDS (Human) IP-WB Antibody Pair
H00009709-PW2	HERPUD1 (Human) IP-WB Antibody Pair	H00003434-PW1	IFIT1 (Human) IP-WB Antibody Pair
H00003073-PW1	HEXA (Human) IP-WB Antibody Pair	H00003433-PW1	IFIT2 (Human) IP-WB Antibody Pair
H00003081-PW1	HGD (Human) IP-WB Antibody Pair	H00003437-PW1	IFIT3 (Human) IP-WB Antibody Pair
H00003081-PW2	HGD (Human) IP-WB Antibody Pair	H00010410-PW1	IFITM3 (Human) IP-WB Antibody Pair
H00064399-PW1	HHIP (Human) IP-WB Antibody Pair	H00010410-PW2	IFITM3 (Human) IP-WB Antibody Pair
H00003091-PW1	HIF1A (Human) IP-WB Antibody Pair	H00010410-PW3	IFITM3 (Human) IP-WB Antibody Pair
H00003092-PW1	HIP1 (Human) IP-WB Antibody Pair	H00003439-PW1	IFNA1 (Human) IP-WB Antibody Pair
H00147746-PW1	HIPK4 (Human) IP-WB Antibody Pair	H00003446-PW1	IFNA10 (Human) IP-WB Antibody Pair
H00003098-PW1	HK1 (Human) IP-WB Antibody Pair	H00003451-PW1	IFNA17 (Human) IP-WB Antibody Pair
H00003105-PW1	HLA-A (Human) IP-WB Antibody Pair	H00003452-PW1	IFNA21 (Human) IP-WB Antibody Pair
H00003107-PW1	HLA-C (Human) IP-WB Antibody Pair	H00003441-PW1	IFNA4 (Human) IP-WB Antibody Pair
H00003108-PW1	HLA-DMA (Human) IP-WB Antibody Pair	H00003442-PW1	IFNA5 (Human) IP-WB Antibody Pair
H00003109-PW1	HLA-DMB (Human) IP-WB Antibody Pair	H00003443-PW1	IFNA6 (Human) IP-WB Antibody Pair
H00003109-PW2	HLA-DMB (Human) IP-WB Antibody Pair	H00003445-PW1	IFNA8 (Human) IP-WB Antibody Pair
H00003109-PW3	HLA-DMB (Human) IP-WB Antibody Pair	H00003456-PW1	IFNB1 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00338376-PW1	IFNE (Human) IP-WB Antibody Pair	H00030818-PW2	KCNIP3 (Human) IP-WB Antibody Pair
H00338376-PW2	IFNE (Human) IP-WB Antibody Pair	H00386618-PW1	KCTD4 (Human) IP-WB Antibody Pair
H00003459-PW1	IFNGR1 (Human) IP-WB Antibody Pair	H00010657-PW1	KHDRBS1 (Human) IP-WB Antibody Pair
H00003460-PW1	IFNGR2 (Human) IP-WB Antibody Pair	H00003795-PW1	KHK (Human) IP-WB Antibody Pair
H00007866-PW1	IFRD2 (Human) IP-WB Antibody Pair	H00009933-PW1	KIAA0020 (Human) IP-WB Antibody Pair
H00010644-PW1	IGF2BP2 (Human) IP-WB Antibody Pair	H00023392-PW1	KIAA0368 (Human) IP-WB Antibody Pair
H00003489-PW1	IGFBP6 (Human) IP-WB Antibody Pair	H00057214-PW1	KIAA1199 (Human) IP-WB Antibody Pair
H00003500-PW1	IGHG1 (Human) IP-WB Antibody Pair	H00057214-PW2	KIAA1199 (Human) IP-WB Antibody Pair
H00003514-PW1	IGKC (Human) IP-WB Antibody Pair	H00026153-PW1	KIF26A (Human) IP-WB Antibody Pair
H00003543-PW1	IGLL1 (Human) IP-WB Antibody Pair	H00011004-PW1	KIF2C (Human) IP-WB Antibody Pair
H00064375-PW1	IKZF4 (Human) IP-WB Antibody Pair	H00084063-PW1	KIRREL2 (Human) IP-WB Antibody Pair
H00003588-PW1	IL10RB (Human) IP-WB Antibody Pair	H00028999-PW1	KLF15 (Human) IP-WB Antibody Pair
H00003594-PW1	IL12RB1 (Human) IP-WB Antibody Pair	H00011279-PW1	KLF8 (Human) IP-WB Antibody Pair
H00003598-PW1	IL13RA2 (Human) IP-WB Antibody Pair	H00054758-PW1	KLHDC4 (Human) IP-WB Antibody Pair
H00003600-PW1	IL15 (Human) IP-WB Antibody Pair	H00005655-PW1	KLK10 (Human) IP-WB Antibody Pair
H00027189-PW1	IL17C (Human) IP-WB Antibody Pair	H00003817-PW1	KLK2 (Human) IP-WB Antibody Pair
H00053342-PW1	IL17D (Human) IP-WB Antibody Pair	H00003817-PW2	KLK2 (Human) IP-WB Antibody Pair
H00112744-PW1	IL17F (Human) IP-WB Antibody Pair	H00003821-PW1	KLRC1 (Human) IP-WB Antibody Pair
H00008809-PW1	IL18R1 (Human) IP-WB Antibody Pair	H00003821-PW2	KLRC1 (Human) IP-WB Antibody Pair
H00003553-PW1	IL1B (Human) IP-WB Antibody Pair	H00010219-PW1	KLRG1 (Human) IP-WB Antibody Pair
H00003553-PW2	IL1B (Human) IP-WB Antibody Pair	H00003827-PW1	KNG1 (Human) IP-WB Antibody Pair
H00003556-PW1	IL1RAP (Human) IP-WB Antibody Pair	H00003827-PW2	KNG1 (Human) IP-WB Antibody Pair
H00003557-PW1	IL1RN (Human) IP-WB Antibody Pair	H00003866-PW1	KRT15 (Human) IP-WB Antibody Pair
H00050604-PW1	IL20 (Human) IP-WB Antibody Pair	H00003872-PW1	KRT17 (Human) IP-WB Antibody Pair
H00051561-PW1	IL23A (Human) IP-WB Antibody Pair	H00003875-PW1	KRT18 (Human) IP-WB Antibody Pair
H00011009-PW1	IL24 (Human) IP-WB Antibody Pair	H00054474-PW1	KRT20 (Human) IP-WB Antibody Pair
H00246778-PW1	IL27 (Human) IP-WB Antibody Pair	H00003885-PW1	KRT34 (Human) IP-WB Antibody Pair
H00003560-PW1	IL2RB (Human) IP-WB Antibody Pair	H00003886-PW1	KRT35 (Human) IP-WB Antibody Pair
H00009235-PW1	IL32 (Human) IP-WB Antibody Pair	H00003852-PW1	KRT5 (Human) IP-WB Antibody Pair
H00003576-PW1	IL8 (Human) IP-WB Antibody Pair	H00003856-PW1	KRT8 (Human) IP-WB Antibody Pair
H00003576-PW2	IL8 (Human) IP-WB Antibody Pair	H00003889-PW1	KRT83 (Human) IP-WB Antibody Pair
H00003576-PW3	IL8 (Human) IP-WB Antibody Pair	H00008942-PW1	KYNU (Human) IP-WB Antibody Pair
H00003576-PW4	IL8 (Human) IP-WB Antibody Pair	H00091133-PW1	L3MBTL4 (Human) IP-WB Antibody Pair
H00003576-PW5	IL8 (Human) IP-WB Antibody Pair	H00003903-PW1	LAIR1 (Human) IP-WB Antibody Pair
H00003608-PW1	ILF2 (Human) IP-WB Antibody Pair	H00051056-PW1	LAP3 (Human) IP-WB Antibody Pair
H00003609-PW1	ILF3 (Human) IP-WB Antibody Pair	H00055323-PW1	LARP6 (Human) IP-WB Antibody Pair
H00003614-PW1	IMPDH1 (Human) IP-WB Antibody Pair	H00003932-PW1	LCK (Human) IP-WB Antibody Pair
H00003615-PW1	IMPDH2 (Human) IP-WB Antibody Pair	H00051451-PW1	LCMT1 (Human) IP-WB Antibody Pair
H00003621-PW1	ING1 (Human) IP-WB Antibody Pair	H00003934-PW1	LCN2 (Human) IP-WB Antibody Pair
H00083729-PW1	INHBE (Human) IP-WB Antibody Pair	H00003939-PW1	LDHA (Human) IP-WB Antibody Pair
H00003628-PW1	INPP1 (Human) IP-WB Antibody Pair	H00003948-PW1	LDHC (Human) IP-WB Antibody Pair
H00003628-PW2	INPP1 (Human) IP-WB Antibody Pair	H00011061-PW1	LECT1 (Human) IP-WB Antibody Pair
H00003632-PW1	INPP5A (Human) IP-WB Antibody Pair	H00010637-PW1	LEFTY1 (Human) IP-WB Antibody Pair
H00128239-PW1	IQGAP3 (Human) IP-WB Antibody Pair	H00007044-PW1	LEFTY2 (Human) IP-WB Antibody Pair
H00003654-PW1	IRAK1 (Human) IP-WB Antibody Pair	H00007044-PW2	LEFTY2 (Human) IP-WB Antibody Pair
H00003654-PW2	IRAK1 (Human) IP-WB Antibody Pair	H00003954-PW1	LETM1 (Human) IP-WB Antibody Pair
H00003654-PW3	IRAK1 (Human) IP-WB Antibody Pair	H00003955-PW1	LFNG (Human) IP-WB Antibody Pair
H00003659-PW1	IRF1 (Human) IP-WB Antibody Pair	H00003956-PW1	LGALS1 (Human) IP-WB Antibody Pair
H00003661-PW1	IRF3 (Human) IP-WB Antibody Pair	H00003958-PW1	LGALS3 (Human) IP-WB Antibody Pair
H00003661-PW2	IRF3 (Human) IP-WB Antibody Pair	H00003959-PW1	LGALS3BP (Human) IP-WB Antibody Pair
H00003663-PW1	IRF5 (Human) IP-WB Antibody Pair	H00005641-PW1	LGMN (Human) IP-WB Antibody Pair
H00010379-PW1	IRF9 (Human) IP-WB Antibody Pair	H00005641-PW2	LGMN (Human) IP-WB Antibody Pair
H00009636-PW1	ISG15 (Human) IP-WB Antibody Pair	H00010859-PW1	LILRB1 (Human) IP-WB Antibody Pair
H00003669-PW1	ISG20 (Human) IP-WB Antibody Pair	H00079727-PW1	LIN28 (Human) IP-WB Antibody Pair
H00003669-PW2	ISG20 (Human) IP-WB Antibody Pair	H00286826-PW1	LIN9 (Human) IP-WB Antibody Pair
H00003696-PW1	ITGB8 (Human) IP-WB Antibody Pair	H00003988-PW1	LIPA (Human) IP-WB Antibody Pair
H00003697-PW1	ITIH1 (Human) IP-WB Antibody Pair	H00009388-PW1	LIPG (Human) IP-WB Antibody Pair
H00003702-PW1	ITK (Human) IP-WB Antibody Pair	H00009516-PW1	LITAF (Human) IP-WB Antibody Pair
H00003705-PW1	ITPK1 (Human) IP-WB Antibody Pair	H00004000-PW1	LMNA (Human) IP-WB Antibody Pair
H00008681-PW1	JMJD7-PLA2G4B (Human) IP-WB Antibody Pair	H00004001-PW1	LMNB1 (Human) IP-WB Antibody Pair
H00003725-PW1	JUN (Human) IP-WB Antibody Pair	H00004016-PW1	LOXL1 (Human) IP-WB Antibody Pair
H00003726-PW1	JUNB (Human) IP-WB Antibody Pair	H00004023-PW1	LPL (Human) IP-WB Antibody Pair
H00010524-PW1	KAT5 (Human) IP-WB Antibody Pair	H00079705-PW1	LRRK1 (Human) IP-WB Antibody Pair
H00003753-PW1	KCNE1 (Human) IP-WB Antibody Pair	H00090678-PW1	LRSAM1 (Human) IP-WB Antibody Pair
H00030818-PW1	KCNIP3 (Human) IP-WB Antibody Pair	H00004046-PW1	LSP1 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00054585-PW1	LZTFL1 (Human) IP-WB Antibody Pair	H00004312-PW1	MMP1 (Human) IP-WB Antibody Pair
H00054585-PW2	LZTFL1 (Human) IP-WB Antibody Pair	H00004319-PW1	MMP10 (Human) IP-WB Antibody Pair
H00009587-PW1	MAD2L1BP (Human) IP-WB Antibody Pair	H00004321-PW1	MMP12 (Human) IP-WB Antibody Pair
H00009587-PW2	MAD2L1BP (Human) IP-WB Antibody Pair	H00004321-PW2	MMP12 (Human) IP-WB Antibody Pair
H00010459-PW1	MAD2L2 (Human) IP-WB Antibody Pair	H00004322-PW1	MMP13 (Human) IP-WB Antibody Pair
H00004110-PW1	MAGEA11 (Human) IP-WB Antibody Pair	H00004322-PW2	MMP13 (Human) IP-WB Antibody Pair
H00004102-PW1	MAGEA3 (Human) IP-WB Antibody Pair	H00004314-PW1	MMP3 (Human) IP-WB Antibody Pair
H00004103-PW1	MAGEA4 (Human) IP-WB Antibody Pair	H00004314-PW2	MMP3 (Human) IP-WB Antibody Pair
H00004107-PW1	MAGEA8 (Human) IP-WB Antibody Pair	H00004316-PW1	MMP7 (Human) IP-WB Antibody Pair
H00028986-PW1	MAGEH1 (Human) IP-WB Antibody Pair	H00004317-PW1	MMP8 (Human) IP-WB Antibody Pair
H00028986-PW2	MAGEH1 (Human) IP-WB Antibody Pair	H00027304-PW1	MOCS3 (Human) IP-WB Antibody Pair
H00004128-PW1	MAOA (Human) IP-WB Antibody Pair	H00004340-PW1	MOG (Human) IP-WB Antibody Pair
H00004129-PW1	MAOB (Human) IP-WB Antibody Pair	H00004342-PW1	MOS (Human) IP-WB Antibody Pair
H00005605-PW1	MAP2K2 (Human) IP-WB Antibody Pair	H00004352-PW1	MPL (Human) IP-WB Antibody Pair
H00006416-PW1	MAP2K4 (Human) IP-WB Antibody Pair	H00004355-PW1	MPP2 (Human) IP-WB Antibody Pair
H00005607-PW1	MAP2K5 (Human) IP-WB Antibody Pair	H00004355-PW2	MPP2 (Human) IP-WB Antibody Pair
H00006885-PW1	MAP3K7 (Human) IP-WB Antibody Pair	H00004356-PW1	MPP3 (Human) IP-WB Antibody Pair
H00010454-PW1	MAP3K7IP1 (Human) IP-WB Antibody Pair	H00004357-PW1	MPST (Human) IP-WB Antibody Pair
H00023118-PW1	MAP3K7IP2 (Human) IP-WB Antibody Pair	H00010205-PW1	MPZL2 (Human) IP-WB Antibody Pair
H00005600-PW1	MAPK11 (Human) IP-WB Antibody Pair	H00055686-PW1	MREG (Human) IP-WB Antibody Pair
H00006300-PW1	MAPK12 (Human) IP-WB Antibody Pair	H00055686-PW2	MREG (Human) IP-WB Antibody Pair
H00005603-PW1	MAPK13 (Human) IP-WB Antibody Pair	H00065993-PW1	MRPS34 (Human) IP-WB Antibody Pair
H00005603-PW2	MAPK13 (Human) IP-WB Antibody Pair	H00000931-PW1	MS4A1 (Human) IP-WB Antibody Pair
H00005603-PW3	MAPK13 (Human) IP-WB Antibody Pair	H00004440-PW1	MS11 (Human) IP-WB Antibody Pair
H00005597-PW1	MAPK6 (Human) IP-WB Antibody Pair	H00124540-PW1	MSI2 (Human) IP-WB Antibody Pair
H00005599-PW1	MAPK8 (Human) IP-WB Antibody Pair	H00010232-PW1	MSLN (Human) IP-WB Antibody Pair
H00005599-PW2	MAPK8 (Human) IP-WB Antibody Pair	H00010232-PW2	MSLN (Human) IP-WB Antibody Pair
H00005599-PW3	MAPK8 (Human) IP-WB Antibody Pair	H00004487-PW1	MSX1 (Human) IP-WB Antibody Pair
H00005599-PW4	MAPK8 (Human) IP-WB Antibody Pair	H00004487-PW2	MSX1 (Human) IP-WB Antibody Pair
H00005599-PW5	MAPK8 (Human) IP-WB Antibody Pair	H00004487-PW3	MSX1 (Human) IP-WB Antibody Pair
H00005599-PW6	MAPK8 (Human) IP-WB Antibody Pair	H00004487-PW4	MSX1 (Human) IP-WB Antibody Pair
H00009261-PW1	MAPKAPK2 (Human) IP-WB Antibody Pair	H00004487-PW5	MSX1 (Human) IP-WB Antibody Pair
H00009261-PW2	MAPKAPK2 (Human) IP-WB Antibody Pair	H00004488-PW1	MSX2 (Human) IP-WB Antibody Pair
H00009261-PW3	MAPKAPK2 (Human) IP-WB Antibody Pair	H00092140-PW1	MTDH (Human) IP-WB Antibody Pair
H00008550-PW1	MAPKAPK5 (Human) IP-WB Antibody Pair	H00004524-PW1	MTHFR (Human) IP-WB Antibody Pair
H00010982-PW1	MAPRE2 (Human) IP-WB Antibody Pair	H00051537-PW1	MTP18 (Human) IP-WB Antibody Pair
H00004141-PW1	MARS (Human) IP-WB Antibody Pair	H00004595-PW1	MUTYH (Human) IP-WB Antibody Pair
H00053615-PW1	MBD3 (Human) IP-WB Antibody Pair	H00004597-PW1	MVD (Human) IP-WB Antibody Pair
H00027349-PW1	MCAT (Human) IP-WB Antibody Pair	H00004598-PW1	MVK (Human) IP-WB Antibody Pair
H00004171-PW1	MCM2 (Human) IP-WB Antibody Pair	H00004602-PW1	MYB (Human) IP-WB Antibody Pair
H00004171-PW2	MCM2 (Human) IP-WB Antibody Pair	H00010514-PW1	MYBBP1A (Human) IP-WB Antibody Pair
H00010445-PW1	MCRS1 (Human) IP-WB Antibody Pair	H00004609-PW1	MYC (Human) IP-WB Antibody Pair
H00004192-PW1	MDK (Human) IP-WB Antibody Pair	H00004615-PW2	MYD88 (Human) IP-WB Antibody Pair
H00004194-PW1	MDM4 (Human) IP-WB Antibody Pair	H00004633-PW1	MYL2 (Human) IP-WB Antibody Pair
H00051102-PW1	MECR (Human) IP-WB Antibody Pair	H00004653-PW1	MYOC (Human) IP-WB Antibody Pair
H00054797-PW1	MED18 (Human) IP-WB Antibody Pair	H00011143-PW1	MYST2 (Human) IP-WB Antibody Pair
H00004221-PW1	MEN1 (Human) IP-WB Antibody Pair	H00029104-PW1	N6AMT1 (Human) IP-WB Antibody Pair
H00004223-PW1	MEOX2 (Human) IP-WB Antibody Pair	H00027163-PW1	NAAA (Human) IP-WB Antibody Pair
H00004240-PW1	MFGE8 (Human) IP-WB Antibody Pair	H00004665-PW1	NAB2 (Human) IP-WB Antibody Pair
H00004240-PW2	MFGE8 (Human) IP-WB Antibody Pair	H00008883-PW1	NAE1 (Human) IP-WB Antibody Pair
H00004248-PW1	MGAT3 (Human) IP-WB Antibody Pair	H00004674-PW1	NAP1L2 (Human) IP-WB Antibody Pair
H00011282-PW1	MGAT4B (Human) IP-WB Antibody Pair	H00004675-PW1	NAP1L3 (Human) IP-WB Antibody Pair
H00010724-PW1	MGEA5 (Human) IP-WB Antibody Pair	H00093100-PW1	NAPRT1 (Human) IP-WB Antibody Pair
H00004255-PW1	MGMT (Human) IP-WB Antibody Pair	H00004677-PW1	NARS (Human) IP-WB Antibody Pair
H00004255-PW2	MGMT (Human) IP-WB Antibody Pair	H00000009-PW1	NAT1 (Human) IP-WB Antibody Pair
H00004258-PW1	MGST2 (Human) IP-WB Antibody Pair	H00000010-PW1	NAT2 (Human) IP-WB Antibody Pair
H00004258-PW2	MGST2 (Human) IP-WB Antibody Pair	H00004684-PW1	NCAM1 (Human) IP-WB Antibody Pair
H00004277-PW1	MICB (Human) IP-WB Antibody Pair	H00004692-PW1	NDN (Human) IP-WB Antibody Pair
H00004282-PW1	MIF (Human) IP-WB Antibody Pair	H00004692-PW2	NDN (Human) IP-WB Antibody Pair
H00004286-PW1	MITF (Human) IP-WB Antibody Pair	H00004692-PW3	NDN (Human) IP-WB Antibody Pair
H00008195-PW1	MKKS (Human) IP-WB Antibody Pair	H00051079-PW1	NDUFA13 (Human) IP-WB Antibody Pair
H00008195-PW2	MKKS (Human) IP-WB Antibody Pair	H00004716-PW1	NDUFB10 (Human) IP-WB Antibody Pair
H00023609-PW1	MKRN2 (Human) IP-WB Antibody Pair	H00004712-PW1	NDUFB6 (Human) IP-WB Antibody Pair
H00004291-PW1	MLF1 (Human) IP-WB Antibody Pair	H00004719-PW1	NDUFS1 (Human) IP-WB Antibody Pair
H00055904-PW1	MLL5 (Human) IP-WB Antibody Pair	H00004720-PW1	NDUFS2 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00004722-PW1	NDUFS3 (Human) IP-WB Antibody Pair	H00116448-PW4	OLIG1 (Human) IP-WB Antibody Pair
H00004722-PW2	NDUFS3 (Human) IP-WB Antibody Pair	H00010215-PW1	OLIG2 (Human) IP-WB Antibody Pair
H00004747-PW1	NEFL (Human) IP-WB Antibody Pair	H00010215-PW2	OLIG2 (Human) IP-WB Antibody Pair
H00252969-PW1	NEIL2 (Human) IP-WB Antibody Pair	H00010215-PW3	OLIG2 (Human) IP-WB Antibody Pair
H00004751-PW1	NEK2 (Human) IP-WB Antibody Pair	H00010215-PW4	OLIG2 (Human) IP-WB Antibody Pair
H00004751-PW2	NEK2 (Human) IP-WB Antibody Pair	H00004973-PW1	OLR1 (Human) IP-WB Antibody Pair
H00004752-PW1	NEK3 (Human) IP-WB Antibody Pair	H00004958-PW1	OMD (Human) IP-WB Antibody Pair
H00091754-PW1	NEK9 (Human) IP-WB Antibody Pair	H00010133-PW1	OPTN (Human) IP-WB Antibody Pair
H00004745-PW1	NELL1 (Human) IP-WB Antibody Pair	H00023594-PW1	ORC6L (Human) IP-WB Antibody Pair
H00004760-PW1	NEUROD1 (Human) IP-WB Antibody Pair	H00023594-PW2	ORC6L (Human) IP-WB Antibody Pair
H00004761-PW1	NEUROD2 (Human) IP-WB Antibody Pair	H00005004-PW1	ORM1 (Human) IP-WB Antibody Pair
H00050674-PW1	NEUROG3 (Human) IP-WB Antibody Pair	H00009180-PW1	OSMR (Human) IP-WB Antibody Pair
H00050674-PW2	NEUROG3 (Human) IP-WB Antibody Pair	H00078990-PW1	OTUB2 (Human) IP-WB Antibody Pair
H00004771-PW1	NF2 (Human) IP-WB Antibody Pair	H00005013-PW1	OTX1 (Human) IP-WB Antibody Pair
H00084901-PW1	NFATC2IP (Human) IP-WB Antibody Pair	H00005018-PW1	OXA1L (Human) IP-WB Antibody Pair
H00084901-PW2	NFATC2IP (Human) IP-WB Antibody Pair	H00005018-PW2	OXA1L (Human) IP-WB Antibody Pair
H00085409-PW1	NKD2 (Human) IP-WB Antibody Pair	H00009943-PW1	OXSRI (Human) IP-WB Antibody Pair
H00001482-PW1	NKX2-5 (Human) IP-WB Antibody Pair	H00009943-PW2	OXSRI (Human) IP-WB Antibody Pair
H00001482-PW2	NKX2-5 (Human) IP-WB Antibody Pair	H00009943-PW3	OXSRI (Human) IP-WB Antibody Pair
H00022829-PW1	NLGN4Y (Human) IP-WB Antibody Pair	H00009943-PW4	OXSRI (Human) IP-WB Antibody Pair
H00004830-PW1	NME1 (Human) IP-WB Antibody Pair	H00009943-PW5	OXSRI (Human) IP-WB Antibody Pair
H00004831-PW1	NME2 (Human) IP-WB Antibody Pair	H00005034-PW1	P4HB (Human) IP-WB Antibody Pair
H00008382-PW1	NME5 (Human) IP-WB Antibody Pair	H00005036-PW1	PA2G4 (Human) IP-WB Antibody Pair
H00008382-PW2	NME5 (Human) IP-WB Antibody Pair	H00005042-PW1	PABPC3 (Human) IP-WB Antibody Pair
H00009111-PW1	NMI (Human) IP-WB Antibody Pair	H00053354-PW1	PANK1 (Human) IP-WB Antibody Pair
H00004837-PW1	NNMT (Human) IP-WB Antibody Pair	H00009060-PW1	PAPSS2 (Human) IP-WB Antibody Pair
H00064127-PW1	NOD2 (Human) IP-WB Antibody Pair	H00050855-PW1	PARD6A (Human) IP-WB Antibody Pair
H00051491-PW1	NOP16 (Human) IP-WB Antibody Pair	H00005071-PW1	PARK2 (Human) IP-WB Antibody Pair
H00004856-PW1	NOV (Human) IP-WB Antibody Pair	H00005080-PW1	PAX6 (Human) IP-WB Antibody Pair
H00004856-PW2	NOV (Human) IP-WB Antibody Pair	H00005089-PW1	PBX2 (Human) IP-WB Antibody Pair
H00010811-PW1	NOXA1 (Human) IP-WB Antibody Pair	H00005094-PW1	PCBP2 (Human) IP-WB Antibody Pair
H00004862-PW1	NPAS2 (Human) IP-WB Antibody Pair	H00005096-PW1	PCCB (Human) IP-WB Antibody Pair
H00010577-PW1	NPC2 (Human) IP-WB Antibody Pair	H00005097-PW1	PCDH1 (Human) IP-WB Antibody Pair
H00004867-PW1	NPHP1 (Human) IP-WB Antibody Pair	H00005105-PW1	PCK1 (Human) IP-WB Antibody Pair
H00007827-PW1	NPHS2 (Human) IP-WB Antibody Pair	H00005105-PW2	PCK1 (Human) IP-WB Antibody Pair
H00007827-PW2	NPHS2 (Human) IP-WB Antibody Pair	H00005106-PW1	PCK2 (Human) IP-WB Antibody Pair
H00055666-PW1	NPLOC4 (Human) IP-WB Antibody Pair	H00005111-PW1	PCNA (Human) IP-WB Antibody Pair
H00004879-PW1	NPPB (Human) IP-WB Antibody Pair	H00005127-PW1	PCTK1 (Human) IP-WB Antibody Pair
H00004885-PW1	NPTX2 (Human) IP-WB Antibody Pair	H00005127-PW2	PCTK1 (Human) IP-WB Antibody Pair
H00001728-PW1	NQO1 (Human) IP-WB Antibody Pair	H00005134-PW1	PDCD2 (Human) IP-WB Antibody Pair
H00008856-PW1	NR1I2 (Human) IP-WB Antibody Pair	H00009659-PW1	PDE4DIP (Human) IP-WB Antibody Pair
H00009970-PW1	NR1I3 (Human) IP-WB Antibody Pair	H00009659-PW2	PDE4DIP (Human) IP-WB Antibody Pair
H00007026-PW1	NR2F2 (Human) IP-WB Antibody Pair	H00005152-PW1	PDE9A (Human) IP-WB Antibody Pair
H00002649-PW1	NR6A1 (Human) IP-WB Antibody Pair	H00005154-PW1	PDGFA (Human) IP-WB Antibody Pair
H00008829-PW1	NRP1 (Human) IP-WB Antibody Pair	H00080310-PW2	PDGFD (Human) IP-WB Antibody Pair
H00050814-PW1	NSDHL (Human) IP-WB Antibody Pair	H00080310-PW3	PDGFD (Human) IP-WB Antibody Pair
H00004909-PW1	NTF4 (Human) IP-WB Antibody Pair	H00005159-PW1	PDGFRB (Human) IP-WB Antibody Pair
H00004521-PW1	NUDT1 (Human) IP-WB Antibody Pair	H00008050-PW1	PDHX (Human) IP-WB Antibody Pair
H00000318-PW1	NUDT2 (Human) IP-WB Antibody Pair	H00002923-PW1	PDIA3 (Human) IP-WB Antibody Pair
H00011163-PW1	NUDT4 (Human) IP-WB Antibody Pair	H00005165-PW1	PDK3 (Human) IP-WB Antibody Pair
H00011162-PW1	NUDT6 (Human) IP-WB Antibody Pair	H00009124-PW1	PDLIM1 (Human) IP-WB Antibody Pair
H00004926-PW1	NUMA1 (Human) IP-WB Antibody Pair	H00064236-PW1	PDLIM2 (Human) IP-WB Antibody Pair
H00004928-PW1	NUP98 (Human) IP-WB Antibody Pair	H00009260-PW1	PDLIM7 (Human) IP-WB Antibody Pair
H00115861-PW1	NXNL1 (Human) IP-WB Antibody Pair	H00009260-PW2	PDLIM7 (Human) IP-WB Antibody Pair
H00004938-PW1	OAS1 (Human) IP-WB Antibody Pair	H00005170-PW1	PDPK1 (Human) IP-WB Antibody Pair
H00004950-PW1	OCLN (Human) IP-WB Antibody Pair	H00005173-PW1	PDYN (Human) IP-WB Antibody Pair
H00004950-PW2	OCLN (Human) IP-WB Antibody Pair	H00079849-PW1	PDZD3 (Human) IP-WB Antibody Pair
H00004953-PW1	ODC1 (Human) IP-WB Antibody Pair	H00005175-PW1	PECAM1 (Human) IP-WB Antibody Pair
H00004968-PW1	OGG1 (Human) IP-WB Antibody Pair	H00010455-PW1	PECI (Human) IP-WB Antibody Pair
H00004969-PW1	OGN (Human) IP-WB Antibody Pair	H00005184-PW1	PEPD (Human) IP-WB Antibody Pair
H00011339-PW1	OIP5 (Human) IP-WB Antibody Pair	H00005184-PW2	PEPD (Human) IP-WB Antibody Pair
H00011339-PW2	OIP5 (Human) IP-WB Antibody Pair	H00008863-PW1	PER3 (Human) IP-WB Antibody Pair
H00116448-PW1	OLIG1 (Human) IP-WB Antibody Pair	H00008863-PW2	PER3 (Human) IP-WB Antibody Pair
H00116448-PW2	OLIG1 (Human) IP-WB Antibody Pair	H00005192-PW1	PEX10 (Human) IP-WB Antibody Pair
H00116448-PW3	OLIG1 (Human) IP-WB Antibody Pair	H00005196-PW1	PF4 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00005208-PW1	PFKFB2 (Human) IP-WB Antibody Pair	H00005475-PW1	PPEF1 (Human) IP-WB Antibody Pair
H00005209-PW1	PFKFB3 (Human) IP-WB Antibody Pair	H00005481-PW1	PPIID (Human) IP-WB Antibody Pair
H00005209-PW2	PFKFB3 (Human) IP-WB Antibody Pair	H00009360-PW1	PPIG (Human) IP-WB Antibody Pair
H00005223-PW1	PGAM1 (Human) IP-WB Antibody Pair	H00009647-PW1	PPM1F (Human) IP-WB Antibody Pair
H00005225-PW1	PGC (Human) IP-WB Antibody Pair	H00009647-PW2	PPM1F (Human) IP-WB Antibody Pair
H00025796-PW1	PGLS (Human) IP-WB Antibody Pair	H00333926-PW1	PPM1J (Human) IP-WB Antibody Pair
H00005236-PW1	PGM1 (Human) IP-WB Antibody Pair	H00333926-PW2	PPM1J (Human) IP-WB Antibody Pair
H00005238-PW1	PGM3 (Human) IP-WB Antibody Pair	H00153743-PW1	PPP1R2P3 (Human) IP-WB Antibody Pair
H00001911-PW1	PHC1 (Human) IP-WB Antibody Pair	H00005507-PW1	PPP1R3C (Human) IP-WB Antibody Pair
H00055274-PW1	PHF10 (Human) IP-WB Antibody Pair	H00005519-PW1	PPP2R1B (Human) IP-WB Antibody Pair
H00026227-PW1	PHGDH (Human) IP-WB Antibody Pair	H00005521-PW1	PPP2R2B (Human) IP-WB Antibody Pair
H00005261-PW1	PHKG2 (Human) IP-WB Antibody Pair	H00005521-PW2	PPP2R2B (Human) IP-WB Antibody Pair
H00000401-PW1	PHOX2A (Human) IP-WB Antibody Pair	H00028227-PW1	PPP2R3B (Human) IP-WB Antibody Pair
H00055361-PW1	PI4K2A (Human) IP-WB Antibody Pair	H00028227-PW2	PPP2R3B (Human) IP-WB Antibody Pair
H00008554-PW1	PIAS1 (Human) IP-WB Antibody Pair	H00005535-PW1	PPP3R2 (Human) IP-WB Antibody Pair
H00009063-PW1	PIAS2 (Human) IP-WB Antibody Pair	H00023532-PW1	PRAME (Human) IP-WB Antibody Pair
H00009463-PW1	PICK1 (Human) IP-WB Antibody Pair	H00000639-PW1	PRDM1 (Human) IP-WB Antibody Pair
H00009463-PW2	PICK1 (Human) IP-WB Antibody Pair	H00010935-PW1	PRDX3 (Human) IP-WB Antibody Pair
H00005283-PW1	PIGH (Human) IP-WB Antibody Pair	H00010549-PW1	PRDX4 (Human) IP-WB Antibody Pair
H00055011-PW1	PIH1D1 (Human) IP-WB Antibody Pair	H00025824-PW1	PRDX5 (Human) IP-WB Antibody Pair
H00055011-PW2	PIH1D1 (Human) IP-WB Antibody Pair	H00005551-PW1	PRF1 (Human) IP-WB Antibody Pair
H00011040-PW1	PIM2 (Human) IP-WB Antibody Pair	H00005551-PW2	PRF1 (Human) IP-WB Antibody Pair
H00008396-PW1	PIP4K2B (Human) IP-WB Antibody Pair	H00051422-PW1	PRKAG2 (Human) IP-WB Antibody Pair
H00079837-PW1	PIP4K2C (Human) IP-WB Antibody Pair	H00005589-PW1	PRKCSH (Human) IP-WB Antibody Pair
H00023396-PW1	PIP5K1C (Human) IP-WB Antibody Pair	H00008575-PW1	PRKRA (Human) IP-WB Antibody Pair
H00051268-PW1	PIPOX (Human) IP-WB Antibody Pair	H00008575-PW2	PRKRA (Human) IP-WB Antibody Pair
H00005307-PW1	PITX1 (Human) IP-WB Antibody Pair	H00005613-PW1	PRKX (Human) IP-WB Antibody Pair
H00005313-PW1	PKLR (Human) IP-WB Antibody Pair	H00005613-PW2	PRKX (Human) IP-WB Antibody Pair
H00063876-PW1	PKNOX2 (Human) IP-WB Antibody Pair	H00055170-PW1	PRMT6 (Human) IP-WB Antibody Pair
H00051365-PW1	PLA1A (Human) IP-WB Antibody Pair	H00054496-PW1	PRMT7 (Human) IP-WB Antibody Pair
H00023659-PW1	PLA2G15 (Human) IP-WB Antibody Pair	H00005626-PW1	PROP1 (Human) IP-WB Antibody Pair
H00011145-PW1	PLA2G16 (Human) IP-WB Antibody Pair	H00005627-PW1	PROS1 (Human) IP-WB Antibody Pair
H00005327-PW1	PLAT (Human) IP-WB Antibody Pair	H00005639-PW1	PRRG2 (Human) IP-WB Antibody Pair
H00059338-PW1	PLEKHA1 (Human) IP-WB Antibody Pair	H00010942-PW1	PRSS21 (Human) IP-WB Antibody Pair
H00005347-PW1	PLK1 (Human) IP-WB Antibody Pair	H00005652-PW1	PRSS8 (Human) IP-WB Antibody Pair
H00005347-PW2	PLK1 (Human) IP-WB Antibody Pair	H00005660-PW1	PSAP (Human) IP-WB Antibody Pair
H00010733-PW2	PLK4 (Human) IP-WB Antibody Pair	H00005688-PW1	PSMA7 (Human) IP-WB Antibody Pair
H00005352-PW1	PLOD2 (Human) IP-WB Antibody Pair	H00005696-PW1	PSMB8 (Human) IP-WB Antibody Pair
H00005354-PW1	PLP1 (Human) IP-WB Antibody Pair	H00005696-PW2	PSMB8 (Human) IP-WB Antibody Pair
H00005359-PW1	PLSCR1 (Human) IP-WB Antibody Pair	H00005711-PW1	PSMD5 (Human) IP-WB Antibody Pair
H00005359-PW2	PLSCR1 (Human) IP-WB Antibody Pair	H00005713-PW1	PSMD7 (Human) IP-WB Antibody Pair
H00057048-PW1	PLSCR3 (Human) IP-WB Antibody Pair	H00005721-PW1	PSME2 (Human) IP-WB Antibody Pair
H00005373-PW1	PMM2 (Human) IP-WB Antibody Pair	H00005725-PW1	PTBP1 (Human) IP-WB Antibody Pair
H00010654-PW1	PMVK (Human) IP-WB Antibody Pair	H00005728-PW1	PTEN (Human) IP-WB Antibody Pair
H00010654-PW2	PMVK (Human) IP-WB Antibody Pair	H00009536-PW1	PTGES (Human) IP-WB Antibody Pair
H00010654-PW3	PMVK (Human) IP-WB Antibody Pair	H00009536-PW2	PTGES (Human) IP-WB Antibody Pair
H00009240-PW1	PNMA1 (Human) IP-WB Antibody Pair	H00080142-PW1	PTGES2 (Human) IP-WB Antibody Pair
H00005409-PW1	PNMT (Human) IP-WB Antibody Pair	H00005740-PW1	PTGIS (Human) IP-WB Antibody Pair
H00005409-PW2	PNMT (Human) IP-WB Antibody Pair	H00005742-PW1	PTGS1 (Human) IP-WB Antibody Pair
H00087178-PW1	PNPT1 (Human) IP-WB Antibody Pair	H00005743-PW1	PTGS2 (Human) IP-WB Antibody Pair
H00057645-PW1	POGK (Human) IP-WB Antibody Pair	H00005741-PW1	PTH (Human) IP-WB Antibody Pair
H00023649-PW1	POLA2 (Human) IP-WB Antibody Pair	H00005741-PW2	PTH (Human) IP-WB Antibody Pair
H00023649-PW2	POLA2 (Human) IP-WB Antibody Pair	H00005741-PW3	PTH (Human) IP-WB Antibody Pair
H00005423-PW1	POLB (Human) IP-WB Antibody Pair	H00005741-PW4	PTH (Human) IP-WB Antibody Pair
H00011232-PW1	POLG2 (Human) IP-WB Antibody Pair	H00005741-PW5	PTH (Human) IP-WB Antibody Pair
H00005429-PW1	POLH (Human) IP-WB Antibody Pair	H00005741-PW6	PTH (Human) IP-WB Antibody Pair
H00051426-PW1	POLK (Human) IP-WB Antibody Pair	H00005747-PW1	PTK2 (Human) IP-WB Antibody Pair
H00051082-PW1	POLR1D (Human) IP-WB Antibody Pair	H00005747-PW2	PTK2 (Human) IP-WB Antibody Pair
H00005434-PW1	POLR2E (Human) IP-WB Antibody Pair	H00005770-PW1	PTPN1 (Human) IP-WB Antibody Pair
H00005437-PW1	POLR2H (Human) IP-WB Antibody Pair	H00005771-PW1	PTPN2 (Human) IP-WB Antibody Pair
H00005443-PW1	POMC (Human) IP-WB Antibody Pair	H00084867-PW1	PTPN5 (Human) IP-WB Antibody Pair
H00005452-PW1	POU2F2 (Human) IP-WB Antibody Pair	H00005778-PW1	PTPN7 (Human) IP-WB Antibody Pair
H00005459-PW1	POU4F3 (Human) IP-WB Antibody Pair	H00284119-PW1	PTRF (Human) IP-WB Antibody Pair
H00005463-PW1	POU6F1 (Human) IP-WB Antibody Pair	H00005805-PW1	PTS (Human) IP-WB Antibody Pair
H00005471-PW1	PPAT (Human) IP-WB Antibody Pair	H00080324-PW1	PUS1 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00080324-PW2	PUS1 (Human) IP-WB Antibody Pair	H00006049-PW1	RNF6 (Human) IP-WB Antibody Pair
H00005818-PW1	PVRL1 (Human) IP-WB Antibody Pair	H00006050-PW1	RNH1 (Human) IP-WB Antibody Pair
H00005819-PW1	PVRL2 (Human) IP-WB Antibody Pair	H00006051-PW1	RNPEP (Human) IP-WB Antibody Pair
H00025945-PW1	PVRL3 (Human) IP-WB Antibody Pair	H00051765-PW1	RP6-213H19.1 (Human) IP-WB Antibody Pair
H00005822-PW1	PWP2 (Human) IP-WB Antibody Pair	H00006117-PW1	RPA1 (Human) IP-WB Antibody Pair
H00005831-PW1	PYCR1 (Human) IP-WB Antibody Pair	H00029935-PW1	RPA4 (Human) IP-WB Antibody Pair
H00065263-PW1	PYCR1 (Human) IP-WB Antibody Pair	H00079657-PW1	RPAP3 (Human) IP-WB Antibody Pair
H00005836-PW1	PYGL (Human) IP-WB Antibody Pair	H00009045-PW1	RPL14 (Human) IP-WB Antibody Pair
H00005859-PW1	QARS (Human) IP-WB Antibody Pair	H00009045-PW2	RPL14 (Human) IP-WB Antibody Pair
H00025797-PW1	QPCT (Human) IP-WB Antibody Pair	H00055197-PW1	RPRD1A (Human) IP-WB Antibody Pair
H00057403-PW1	RAB22A (Human) IP-WB Antibody Pair	H00006198-PW1	RPS6KB1 (Human) IP-WB Antibody Pair
H00010981-PW1	RAB32 (Human) IP-WB Antibody Pair	H00006199-PW1	RPS6KB2 (Human) IP-WB Antibody Pair
H00010981-PW2	RAB32 (Human) IP-WB Antibody Pair	H00057521-PW1	RPTOR (Human) IP-WB Antibody Pair
H00023682-PW1	RAB38 (Human) IP-WB Antibody Pair	H00010325-PW1	RRAGB (Human) IP-WB Antibody Pair
H00116442-PW1	RAB39B (Human) IP-WB Antibody Pair	H00006237-PW1	RRAS (Human) IP-WB Antibody Pair
H00005865-PW1	RAB3B (Human) IP-WB Antibody Pair	H00006237-PW2	RRAS (Human) IP-WB Antibody Pair
H00338382-PW1	RAB7B (Human) IP-WB Antibody Pair	H00006241-PW1	RRM2 (Human) IP-WB Antibody Pair
H00023637-PW1	RABGAP1 (Human) IP-WB Antibody Pair	H00006241-PW2	RRM2 (Human) IP-WB Antibody Pair
H00005888-PW1	RAD51 (Human) IP-WB Antibody Pair	H00023212-PW1	RRS1 (Human) IP-WB Antibody Pair
H00025788-PW1	RAD54B (Human) IP-WB Antibody Pair	H00065078-PW1	RTN4R (Human) IP-WB Antibody Pair
H00005897-PW1	RAG2 (Human) IP-WB Antibody Pair	H00000864-PW1	RUNX3 (Human) IP-WB Antibody Pair
H00005899-PW1	RALB (Human) IP-WB Antibody Pair	H00010856-PW1	RUVBL2 (Human) IP-WB Antibody Pair
H00010928-PW1	RALBP1 (Human) IP-WB Antibody Pair	H00006256-PW1	RXRA (Human) IP-WB Antibody Pair
H00055103-PW1	RALGPS2 (Human) IP-WB Antibody Pair	H00006256-PW2	RXRA (Human) IP-WB Antibody Pair
H00010411-PW1	RAPGEF3 (Human) IP-WB Antibody Pair	H00023429-PW1	RYBP (Human) IP-WB Antibody Pair
H00005919-PW1	RARRES2 (Human) IP-WB Antibody Pair	H00023429-PW2	RYBP (Human) IP-WB Antibody Pair
H00005917-PW1	RARS (Human) IP-WB Antibody Pair	H00023429-PW3	RYBP (Human) IP-WB Antibody Pair
H00005921-PW1	RASA1 (Human) IP-WB Antibody Pair	H00006278-PW1	S100A7 (Human) IP-WB Antibody Pair
H00023551-PW1	RASD2 (Human) IP-WB Antibody Pair	H00006278-PW2	S100A7 (Human) IP-WB Antibody Pair
H00025780-PW1	RASGRP3 (Human) IP-WB Antibody Pair	H00006280-PW1	S100A9 (Human) IP-WB Antibody Pair
H00011228-PW1	RASSF8 (Human) IP-WB Antibody Pair	H00006291-PW1	SAA4 (Human) IP-WB Antibody Pair
H00011228-PW2	RASSF8 (Human) IP-WB Antibody Pair	H00006291-PW2	SAA4 (Human) IP-WB Antibody Pair
H00010181-PW1	RBM5 (Human) IP-WB Antibody Pair	H00025939-PW1	SAMHD1 (Human) IP-WB Antibody Pair
H00023543-PW1	RBM9 (Human) IP-WB Antibody Pair	H00010284-PW1	SAP18 (Human) IP-WB Antibody Pair
H00005937-PW1	RBMS1 (Human) IP-WB Antibody Pair	H00054938-PW1	SARS2 (Human) IP-WB Antibody Pair
H00001104-PW1	RCC1 (Human) IP-WB Antibody Pair	H00060485-PW1	SAV1 (Human) IP-WB Antibody Pair
H00001104-PW2	RCC1 (Human) IP-WB Antibody Pair	H00060485-PW2	SAV1 (Human) IP-WB Antibody Pair
H00005957-PW1	RCVRN (Human) IP-WB Antibody Pair	H00010609-PW1	SC65 (Human) IP-WB Antibody Pair
H00005957-PW2	RCVRN (Human) IP-WB Antibody Pair	H00051540-PW1	SCLY (Human) IP-WB Antibody Pair
H00005957-PW3	RCVRN (Human) IP-WB Antibody Pair	H00006322-PW1	SCML1 (Human) IP-WB Antibody Pair
H00005957-PW4	RCVRN (Human) IP-WB Antibody Pair	H00006342-PW1	SCP2 (Human) IP-WB Antibody Pair
H00005962-PW1	RDX (Human) IP-WB Antibody Pair	H00057147-PW1	SCYL3 (Human) IP-WB Antibody Pair
H00005971-PW1	RELB (Human) IP-WB Antibody Pair	H00006385-PW1	SDC4 (Human) IP-WB Antibody Pair
H00005984-PW1	RFC4 (Human) IP-WB Antibody Pair	H00006386-PW1	SDCBP (Human) IP-WB Antibody Pair
H00005993-PW1	RFX5 (Human) IP-WB Antibody Pair	H00006386-PW2	SDCBP (Human) IP-WB Antibody Pair
H00008490-PW1	RGS5 (Human) IP-WB Antibody Pair	H00195814-PW1	SDR16C5 (Human) IP-WB Antibody Pair
H00006009-PW1	RHEB (Human) IP-WB Antibody Pair	H00113675-PW1	SDSL (Human) IP-WB Antibody Pair
H00253260-PW1	RICTOR (Human) IP-WB Antibody Pair	H00113675-PW2	SDSL (Human) IP-WB Antibody Pair
H00083732-PW1	RIOK1 (Human) IP-WB Antibody Pair	H00113675-PW3	SDSL (Human) IP-WB Antibody Pair
H00055781-PW1	RIOK2 (Human) IP-WB Antibody Pair	H00006396-PW1	SEC13 (Human) IP-WB Antibody Pair
H00008780-PW1	RIOK3 (Human) IP-WB Antibody Pair	H00055829-PW1	SELS (Human) IP-WB Antibody Pair
H00008780-PW2	RIOK3 (Human) IP-WB Antibody Pair	H00123228-PW1	SENP8 (Human) IP-WB Antibody Pair
H00006017-PW1	RLBP1 (Human) IP-WB Antibody Pair	H00001731-PW1	SEPT1 (Human) IP-WB Antibody Pair
H00006013-PW1	RLN1 (Human) IP-WB Antibody Pair	H00001731-PW2	SEPT1 (Human) IP-WB Antibody Pair
H00140432-PW1	RNF113B (Human) IP-WB Antibody Pair	H00151011-PW1	SEPT10 (Human) IP-WB Antibody Pair
H00140432-PW2	RNF113B (Human) IP-WB Antibody Pair	H00005413-PW1	SEPT5 (Human) IP-WB Antibody Pair
H00055905-PW1	RNF114 (Human) IP-WB Antibody Pair	H00051156-PW1	SERPINA10 (Human) IP-WB Antibody Pair
H00051132-PW1	RNF12 (Human) IP-WB Antibody Pair	H00145264-PW1	SERPINA12 (Human) IP-WB Antibody Pair
H00051132-PW2	RNF12 (Human) IP-WB Antibody Pair	H00145264-PW2	SERPINA12 (Human) IP-WB Antibody Pair
H00168433-PW1	RNF133 (Human) IP-WB Antibody Pair	H00000012-PW1	SERPINA3 (Human) IP-WB Antibody Pair
H00168433-PW2	RNF133 (Human) IP-WB Antibody Pair	H00000012-PW2	SERPINA3 (Human) IP-WB Antibody Pair
H00255488-PW1	RNF144B (Human) IP-WB Antibody Pair	H00001992-PW1	SERPINB1 (Human) IP-WB Antibody Pair
H00056254-PW1	RNF20 (Human) IP-WB Antibody Pair	H00005055-PW1	SERPINB2 (Human) IP-WB Antibody Pair
H00064320-PW1	RNF25 (Human) IP-WB Antibody Pair	H00006317-PW1	SERPINB3 (Human) IP-WB Antibody Pair
H00140545-PW1	RNF32 (Human) IP-WB Antibody Pair	H00006317-PW2	SERPINB3 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00006317-PW3	SERPINB3 (Human) IP-WB Antibody Pair	H00009655-PW1	SOCS5 (Human) IP-WB Antibody Pair
H00008710-PW1	SERPINB7 (Human) IP-WB Antibody Pair	H00009306-PW1	SOCS6 (Human) IP-WB Antibody Pair
H00003053-PW1	SERPIND1 (Human) IP-WB Antibody Pair	H00054937-PW1	SOHLH2 (Human) IP-WB Antibody Pair
H00005054-PW1	SERPINE1 (Human) IP-WB Antibody Pair	H00054937-PW2	SOHLH2 (Human) IP-WB Antibody Pair
H00005054-PW2	SERPINE1 (Human) IP-WB Antibody Pair	H00006652-PW1	SORD (Human) IP-WB Antibody Pair
H00005345-PW1	SERPINF2 (Human) IP-WB Antibody Pair	H00009580-PW1	SOX13 (Human) IP-WB Antibody Pair
H00000710-PW1	SERPING1 (Human) IP-WB Antibody Pair	H00011063-PW1	SOX30 (Human) IP-WB Antibody Pair
H00027244-PW1	SESN1 (Human) IP-WB Antibody Pair	H00006662-PW1	SOX9 (Human) IP-WB Antibody Pair
H00083667-PW1	SESN2 (Human) IP-WB Antibody Pair	H00006672-PW1	SP100 (Human) IP-WB Antibody Pair
H00002810-PW1	SFN (Human) IP-WB Antibody Pair	H00003431-PW1	SP110 (Human) IP-WB Antibody Pair
H00006439-PW1	SFTPFB (Human) IP-WB Antibody Pair	H00003431-PW2	SP110 (Human) IP-WB Antibody Pair
H00006444-PW1	SGCD (Human) IP-WB Antibody Pair	H00008877-PW1	SPHK1 (Human) IP-WB Antibody Pair
H00006445-PW1	SGCG (Human) IP-WB Antibody Pair	H00008877-PW2	SPHK1 (Human) IP-WB Antibody Pair
H00006445-PW2	SGCG (Human) IP-WB Antibody Pair	H00056848-PW1	SPHK2 (Human) IP-WB Antibody Pair
H00084197-PW1	SGK196 (Human) IP-WB Antibody Pair	H00006690-PW1	SPINK1 (Human) IP-WB Antibody Pair
H00010110-PW1	SGK2 (Human) IP-WB Antibody Pair	H00006693-PW1	SPN (Human) IP-WB Antibody Pair
H00010110-PW2	SGK2 (Human) IP-WB Antibody Pair	H00006697-PW1	SPR (Human) IP-WB Antibody Pair
H00151648-PW1	SGOL1 (Human) IP-WB Antibody Pair	H00010252-PW1	SPRY1 (Human) IP-WB Antibody Pair
H00151648-PW2	SGOL1 (Human) IP-WB Antibody Pair	H00006728-PW1	SRP19 (Human) IP-WB Antibody Pair
H00151648-PW3	SGOL1 (Human) IP-WB Antibody Pair	H00006732-PW1	SRPK1 (Human) IP-WB Antibody Pair
H00006448-PW1	SGSH (Human) IP-WB Antibody Pair	H00063826-PW1	SRR (Human) IP-WB Antibody Pair
H00006449-PW1	SGTA (Human) IP-WB Antibody Pair	H00085464-PW1	SSH2 (Human) IP-WB Antibody Pair
H00006449-PW2	SGTA (Human) IP-WB Antibody Pair	H00117178-PW1	SSX2IP (Human) IP-WB Antibody Pair
H00117157-PW1	SH2D1B (Human) IP-WB Antibody Pair	H00006487-PW1	ST3GAL3 (Human) IP-WB Antibody Pair
H00009047-PW1	SH2D2A (Human) IP-WB Antibody Pair	H00006480-PW1	ST6GAL1 (Human) IP-WB Antibody Pair
H00063898-PW1	SH2D4A (Human) IP-WB Antibody Pair	H00006775-PW1	STAT4 (Human) IP-WB Antibody Pair
H00053358-PW1	SHC3 (Human) IP-WB Antibody Pair	H00006776-PW1	STAT5A (Human) IP-WB Antibody Pair
H00399694-PW1	SHC4 (Human) IP-WB Antibody Pair	H00006780-PW1	STAU1 (Human) IP-WB Antibody Pair
H00006474-PW1	SHOX2 (Human) IP-WB Antibody Pair	H00008987-PW1	STBD1 (Human) IP-WB Antibody Pair
H00059307-PW1	SIGIRR (Human) IP-WB Antibody Pair	H00010963-PW1	STIP1 (Human) IP-WB Antibody Pair
H00000946-PW1	SIGLEC6 (Human) IP-WB Antibody Pair	H00008576-PW1	STK16 (Human) IP-WB Antibody Pair
H00000946-PW2	SIGLEC6 (Human) IP-WB Antibody Pair	H00009263-PW1	STK17A (Human) IP-WB Antibody Pair
H00006490-PW1	SILV (Human) IP-WB Antibody Pair	H00008859-PW1	STK19 (Human) IP-WB Antibody Pair
H000051547-PW1	SIRT7 (Human) IP-WB Antibody Pair	H00008859-PW2	STK19 (Human) IP-WB Antibody Pair
H00006502-PW1	SKP2 (Human) IP-WB Antibody Pair	H00202374-PW1	STK32A (Human) IP-WB Antibody Pair
H00006563-PW1	SLC14A1 (Human) IP-WB Antibody Pair	H00065975-PW1	STK33 (Human) IP-WB Antibody Pair
H00000292-PW1	SLC25A5 (Human) IP-WB Antibody Pair	H00065975-PW2	STK33 (Human) IP-WB Antibody Pair
H00376497-PW1	SLC27A1 (Human) IP-WB Antibody Pair	H00065975-PW3	STK33 (Human) IP-WB Antibody Pair
H00006523-PW1	SLC5A1 (Human) IP-WB Antibody Pair	H00065975-PW4	STK33 (Human) IP-WB Antibody Pair
H00057152-PW1	SLURP1 (Human) IP-WB Antibody Pair	H00065975-PW5	STK33 (Human) IP-WB Antibody Pair
H00004091-PW1	SMAD6 (Human) IP-WB Antibody Pair	H00011329-PW1	STK38 (Human) IP-WB Antibody Pair
H00056916-PW1	SMARCAD1 (Human) IP-WB Antibody Pair	H00083931-PW1	STK40 (Human) IP-WB Antibody Pair
H00006598-PW1	SMARCB1 (Human) IP-WB Antibody Pair	H00083931-PW2	STK40 (Human) IP-WB Antibody Pair
H00006605-PW1	SMARCE1 (Human) IP-WB Antibody Pair	H00003925-PW1	STMN1 (Human) IP-WB Antibody Pair
H00009126-PW1	SMC3 (Human) IP-WB Antibody Pair	H00030968-PW1	STOML2 (Human) IP-WB Antibody Pair
H00006607-PW1	SMN2 (Human) IP-WB Antibody Pair	H00092335-PW1	STRADA (Human) IP-WB Antibody Pair
H00055512-PW1	SMPD3 (Human) IP-WB Antibody Pair	H00000412-PW1	STS (Human) IP-WB Antibody Pair
H00027293-PW1	SMPDL3B (Human) IP-WB Antibody Pair	H00008676-PW1	STX11 (Human) IP-WB Antibody Pair
H00057154-PW1	SMURF1 (Human) IP-WB Antibody Pair	H00023673-PW1	STX12 (Human) IP-WB Antibody Pair
H00056950-PW1	SMYD2 (Human) IP-WB Antibody Pair	H00023673-PW2	STX12 (Human) IP-WB Antibody Pair
H00006615-PW1	SNAI1 (Human) IP-WB Antibody Pair	H00002054-PW1	STX2 (Human) IP-WB Antibody Pair
H00006591-PW1	SNAI2 (Human) IP-WB Antibody Pair	H00006810-PW1	STX4 (Human) IP-WB Antibody Pair
H00006591-PW2	SNAI2 (Human) IP-WB Antibody Pair	H00008417-PW1	STX7 (Human) IP-WB Antibody Pair
H00008773-PW1	SNAP23 (Human) IP-WB Antibody Pair	H00055359-PW1	STYK1 (Human) IP-WB Antibody Pair
H00006622-PW1	SNCA (Human) IP-WB Antibody Pair	H00006817-PW1	SULT1A1 (Human) IP-WB Antibody Pair
H00079753-PW1	SNIP1 (Human) IP-WB Antibody Pair	H00027284-PW1	SULT1B1 (Human) IP-WB Antibody Pair
H00006626-PW1	SNRPA (Human) IP-WB Antibody Pair	H00006783-PW1	SULT1E1 (Human) IP-WB Antibody Pair
H00010073-PW1	SNUPN (Human) IP-WB Antibody Pair	H00006820-PW1	SULT2B1 (Human) IP-WB Antibody Pair
H00006642-PW1	SNX1 (Human) IP-WB Antibody Pair	H00006821-PW1	SUOX (Human) IP-WB Antibody Pair
H00029916-PW1	SNX11 (Human) IP-WB Antibody Pair	H00023075-PW1	SWAP70 (Human) IP-WB Antibody Pair
H00029916-PW2	SNX11 (Human) IP-WB Antibody Pair	H00008224-PW1	SYN3 (Human) IP-WB Antibody Pair
H00257364-PW1	SNX33 (Human) IP-WB Antibody Pair	H00006881-PW1	TAF10 (Human) IP-WB Antibody Pair
H00051429-PW1	SNX9 (Human) IP-WB Antibody Pair	H00008148-PW1	TAF15 (Human) IP-WB Antibody Pair
H00051429-PW2	SNX9 (Human) IP-WB Antibody Pair	H00009014-PW1	TAF1B (Human) IP-WB Antibody Pair
H00006646-PW1	SOAT1 (Human) IP-WB Antibody Pair	H00010629-PW1	TAF6L (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00010629-PW2	TAF6L (Human) IP-WB Antibody Pair	H00057761-PW1	TRIB3 (Human) IP-WB Antibody Pair
H00006888-PW1	TALDO1 (Human) IP-WB Antibody Pair	H00010626-PW1	TRIM16 (Human) IP-WB Antibody Pair
H00051347-PW1	TAOK3 (Human) IP-WB Antibody Pair	H00010346-PW1	TRIM22 (Human) IP-WB Antibody Pair
H00006895-PW1	TARBP2 (Human) IP-WB Antibody Pair	H00008805-PW1	TRIM24 (Human) IP-WB Antibody Pair
H00006895-PW2	TARBP2 (Human) IP-WB Antibody Pair	H00007706-PW1	TRIM25 (Human) IP-WB Antibody Pair
H00023435-PW1	TARDBP (Human) IP-WB Antibody Pair	H00022954-PW1	TRIM32 (Human) IP-WB Antibody Pair
H00006905-PW1	TBCE (Human) IP-WB Antibody Pair	H00084675-PW1	TRIM55 (Human) IP-WB Antibody Pair
H00006907-PW1	TBL1X (Human) IP-WB Antibody Pair	H00055223-PW1	TRIM62 (Human) IP-WB Antibody Pair
H00026608-PW1	TBL2 (Human) IP-WB Antibody Pair	H00084676-PW2	TRIM63 (Human) IP-WB Antibody Pair
H00006919-PW1	TCEA2 (Human) IP-WB Antibody Pair	H00493829-PW1	TRIM72 (Human) IP-WB Antibody Pair
H00006920-PW1	TCEA3 (Human) IP-WB Antibody Pair	H00493829-PW2	TRIM72 (Human) IP-WB Antibody Pair
H00009338-PW1	TCEAL1 (Human) IP-WB Antibody Pair	H00493829-PW3	TRIM72 (Human) IP-WB Antibody Pair
H00009338-PW2	TCEAL1 (Human) IP-WB Antibody Pair	H00493829-PW4	TRIM72 (Human) IP-WB Antibody Pair
H00008115-PW1	TCL1A (Human) IP-WB Antibody Pair	H00009322-PW1	TRIP10 (Human) IP-WB Antibody Pair
H00008115-PW2	TCL1A (Human) IP-WB Antibody Pair	H00007205-PW1	TRIP6 (Human) IP-WB Antibody Pair
H00006947-PW1	TCN1 (Human) IP-WB Antibody Pair	H00008848-PW1	TSC22D1 (Human) IP-WB Antibody Pair
H00006947-PW2	TCN1 (Human) IP-WB Antibody Pair	H00001831-PW1	TSC22D3 (Human) IP-WB Antibody Pair
H00006949-PW1	TCOF1 (Human) IP-WB Antibody Pair	H00116461-PW1	TSEN15 (Human) IP-WB Antibody Pair
H00055775-PW1	TDP1 (Human) IP-WB Antibody Pair	H00007251-PW1	TSG101 (Human) IP-WB Antibody Pair
H00007004-PW1	TEAD4 (Human) IP-WB Antibody Pair	H00081629-PW1	TSSK3 (Human) IP-WB Antibody Pair
H00056157-PW1	TEX13A (Human) IP-WB Antibody Pair	H00081629-PW2	TSSK3 (Human) IP-WB Antibody Pair
H00051368-PW1	TEX264 (Human) IP-WB Antibody Pair	H00081629-PW3	TSSK3 (Human) IP-WB Antibody Pair
H00007024-PW1	TFCP2 (Human) IP-WB Antibody Pair	H00083983-PW1	TSSK6 (Human) IP-WB Antibody Pair
H00007045-PW1	TGFB1 (Human) IP-WB Antibody Pair	H00007265-PW1	TTC1 (Human) IP-WB Antibody Pair
H00007052-PW1	TGM2 (Human) IP-WB Antibody Pair	H00123016-PW1	TTC8 (Human) IP-WB Antibody Pair
H00055145-PW1	THAP1 (Human) IP-WB Antibody Pair	H00007272-PW1	TTK (Human) IP-WB Antibody Pair
H00007072-PW1	TIA1 (Human) IP-WB Antibody Pair	H00079861-PW1	TUBAL3 (Human) IP-WB Antibody Pair
H00026230-PW1	TIAM2 (Human) IP-WB Antibody Pair	H00081027-PW1	TUBB1 (Human) IP-WB Antibody Pair
H00007076-PW1	TIMP1 (Human) IP-WB Antibody Pair	H00056604-PW1	TUBB4Q (Human) IP-WB Antibody Pair
H00114609-PW1	TIRAP (Human) IP-WB Antibody Pair	H00051174-PW1	TUBD1 (Human) IP-WB Antibody Pair
H00007083-PW1	TK1 (Human) IP-WB Antibody Pair	H00051174-PW2	TUBD1 (Human) IP-WB Antibody Pair
H00007086-PW1	TKT (Human) IP-WB Antibody Pair	H00010844-PW1	TUBGCP2 (Human) IP-WB Antibody Pair
H00030012-PW1	TLX3 (Human) IP-WB Antibody Pair	H00027229-PW1	TUBGCP4 (Human) IP-WB Antibody Pair
H00011070-PW1	TMEM115 (Human) IP-WB Antibody Pair	H00007284-PW1	TUFM (Human) IP-WB Antibody Pair
H00055863-PW1	TMEM126B (Human) IP-WB Antibody Pair	H00007284-PW2	TUFM (Human) IP-WB Antibody Pair
H00029767-PW1	TMOD2 (Human) IP-WB Antibody Pair	H00007288-PW1	TULP2 (Human) IP-WB Antibody Pair
H00008795-PW1	TNFRSF10B (Human) IP-WB Antibody Pair	H00005756-PW1	TWF1 (Human) IP-WB Antibody Pair
H00008795-PW2	TNFRSF10B (Human) IP-WB Antibody Pair	H00005756-PW2	TWF1 (Human) IP-WB Antibody Pair
H00004982-PW1	TNFRSF11B (Human) IP-WB Antibody Pair	H00011344-PW1	TWF2 (Human) IP-WB Antibody Pair
H00000608-PW1	TNFRSF17 (Human) IP-WB Antibody Pair	H00007291-PW1	TWIST1 (Human) IP-WB Antibody Pair
H00010673-PW1	TNFSF13B (Human) IP-WB Antibody Pair	H00007291-PW2	TWIST1 (Human) IP-WB Antibody Pair
H00010673-PW2	TNFSF13B (Human) IP-WB Antibody Pair	H00007291-PW3	TWIST1 (Human) IP-WB Antibody Pair
H00008740-PW1	TNFSF14 (Human) IP-WB Antibody Pair	H00007291-PW4	TWIST1 (Human) IP-WB Antibody Pair
H00008995-PW1	TNFSF18 (Human) IP-WB Antibody Pair	H00117581-PW1	TWIST2 (Human) IP-WB Antibody Pair
H00008995-PW2	TNFSF18 (Human) IP-WB Antibody Pair	H00117581-PW2	TWIST2 (Human) IP-WB Antibody Pair
H00010188-PW1	TNK2 (Human) IP-WB Antibody Pair	H00007295-PW1	TXN (Human) IP-WB Antibody Pair
H00010188-PW2	TNK2 (Human) IP-WB Antibody Pair	H00023071-PW1	TXNDC4 (Human) IP-WB Antibody Pair
H00007139-PW1	TNNT2 (Human) IP-WB Antibody Pair	H00010190-PW1	TXNDC9 (Human) IP-WB Antibody Pair
H00007148-PW1	TNXB (Human) IP-WB Antibody Pair	H00007297-PW1	TYK2 (Human) IP-WB Antibody Pair
H00054472-PW1	TOLLIP (Human) IP-WB Antibody Pair	H00001890-PW1	TYMP (Human) IP-WB Antibody Pair
H00010953-PW1	TOMM34 (Human) IP-WB Antibody Pair	H00007298-PW1	TYMS (Human) IP-WB Antibody Pair
H00010953-PW2	TOMM34 (Human) IP-WB Antibody Pair	H00007298-PW2	TYMS (Human) IP-WB Antibody Pair
H00007157-PW1	TP53 (Human) IP-WB Antibody Pair	H00007299-PW1	TYR (Human) IP-WB Antibody Pair
H00007157-PW2	TP53 (Human) IP-WB Antibody Pair	H00007306-PW1	TYRP1 (Human) IP-WB Antibody Pair
H00009540-PW1	TP53I3 (Human) IP-WB Antibody Pair	H00053347-PW1	UBASH3A (Human) IP-WB Antibody Pair
H00001200-PW1	TPP1 (Human) IP-WB Antibody Pair	H00009246-PW1	UBE2L6 (Human) IP-WB Antibody Pair
H00008717-PW1	TRADD (Human) IP-WB Antibody Pair	H00027338-PW1	UBE2S (Human) IP-WB Antibody Pair
H00007185-PW1	TRAF1 (Human) IP-WB Antibody Pair	H00029089-PW1	UBE2T (Human) IP-WB Antibody Pair
H00007185-PW2	TRAF1 (Human) IP-WB Antibody Pair	H00029089-PW2	UBE2T (Human) IP-WB Antibody Pair
H00010758-PW1	TRAF3IP2 (Human) IP-WB Antibody Pair	H00029089-PW3	UBE2T (Human) IP-WB Antibody Pair
H00007188-PW1	TRAF5 (Human) IP-WB Antibody Pair	H00023304-PW1	UBR2 (Human) IP-WB Antibody Pair
H00007189-PW1	TRAF6 (Human) IP-WB Antibody Pair	H00007371-PW1	UCK2 (Human) IP-WB Antibody Pair
H00010131-PW1	TRAP1 (Human) IP-WB Antibody Pair	H00054963-PW1	UCKL1 (Human) IP-WB Antibody Pair
H00001787-PW1	TRDMT1 (Human) IP-WB Antibody Pair	H00054600-PW1	UGT1A9 (Human) IP-WB Antibody Pair
H00028951-PW1	TRIB2 (Human) IP-WB Antibody Pair	H00054600-PW2	UGT1A9 (Human) IP-WB Antibody Pair

Antibody Pairs for IP-WB

Cat. Num.	Product Name	Cat. Num.	Product Name
H00051720-PW1	UIMC1 (Human) IP-WB Antibody Pair	H00081029-PW1	WNT5B (Human) IP-WB Antibody Pair
H00007369-PW1	UMOD (Human) IP-WB Antibody Pair	H00007478-PW1	WNT8A (Human) IP-WB Antibody Pair
H00007369-PW2	UMOD (Human) IP-WB Antibody Pair	H00011059-PW1	WWP1 (Human) IP-WB Antibody Pair
H00011274-PW1	USP18 (Human) IP-WB Antibody Pair	H00011060-PW1	WWP2 (Human) IP-WB Antibody Pair
H00009099-PW1	USP2 (Human) IP-WB Antibody Pair	H00000331-PW1	XIAP (Human) IP-WB Antibody Pair
H00010868-PW1	USP20 (Human) IP-WB Antibody Pair	H00007507-PW1	XPA (Human) IP-WB Antibody Pair
H00023032-PW1	USP33 (Human) IP-WB Antibody Pair	H00057510-PW1	XPO5 (Human) IP-WB Antibody Pair
H00081839-PW1	VANGL1 (Human) IP-WB Antibody Pair	H00007517-PW1	XRCC3 (Human) IP-WB Antibody Pair
H00009218-PW1	VAPA (Human) IP-WB Antibody Pair	H00007518-PW1	XRCC4 (Human) IP-WB Antibody Pair
H00114990-PW1	VASN (Human) IP-WB Antibody Pair	H00002547-PW1	XRCC6 (Human) IP-WB Antibody Pair
H00007411-PW1	VBP1 (Human) IP-WB Antibody Pair	H00010413-PW1	YAP1 (Human) IP-WB Antibody Pair
H00007421-PW1	VDR (Human) IP-WB Antibody Pair	H00051087-PW1	YBX2 (Human) IP-WB Antibody Pair
H00007423-PW1	VEGFB (Human) IP-WB Antibody Pair	H00007525-PW1	YES1 (Human) IP-WB Antibody Pair
H00007423-PW2	VEGFB (Human) IP-WB Antibody Pair	H00009189-PW1	ZBED1 (Human) IP-WB Antibody Pair
H00389136-PW1	VGLL3 (Human) IP-WB Antibody Pair	H00084936-PW1	ZFYVE19 (Human) IP-WB Antibody Pair
H00007429-PW1	VIL1 (Human) IP-WB Antibody Pair	H00084936-PW2	ZFYVE19 (Human) IP-WB Antibody Pair
H00007431-PW1	VIM (Human) IP-WB Antibody Pair	H00009372-PW1	ZFYVE9 (Human) IP-WB Antibody Pair
H00007433-PW1	VIPR1 (Human) IP-WB Antibody Pair	H00064393-PW1	ZMAT3 (Human) IP-WB Antibody Pair
H00008876-PW1	VNN1 (Human) IP-WB Antibody Pair	H00079698-PW1	ZMAT4 (Human) IP-WB Antibody Pair
H00008876-PW2	VNN1 (Human) IP-WB Antibody Pair	H00079698-PW2	ZMAT4 (Human) IP-WB Antibody Pair
H00057617-PW1	VPS18 (Human) IP-WB Antibody Pair	H00007690-PW1	ZNF131 (Human) IP-WB Antibody Pair
H00023339-PW1	VPS39 (Human) IP-WB Antibody Pair	H00007566-PW1	ZNF18 (Human) IP-WB Antibody Pair
H00007443-PW1	VRK1 (Human) IP-WB Antibody Pair	H00007566-PW2	ZNF18 (Human) IP-WB Antibody Pair
H00128434-PW1	VSTM2L (Human) IP-WB Antibody Pair	H00007569-PW1	ZNF182 (Human) IP-WB Antibody Pair
H00128434-PW2	VSTM2L (Human) IP-WB Antibody Pair	H00007762-PW1	ZNF215 (Human) IP-WB Antibody Pair
H00007448-PW1	VTN (Human) IP-WB Antibody Pair	H00007775-PW1	ZNF232 (Human) IP-WB Antibody Pair
H00023063-PW1	WAPAL (Human) IP-WB Antibody Pair	H00010127-PW1	ZNF263 (Human) IP-WB Antibody Pair
H00007453-PW1	WARS (Human) IP-WB Antibody Pair	H00010781-PW1	ZNF266 (Human) IP-WB Antibody Pair
H00007453-PW2	WARS (Human) IP-WB Antibody Pair	H00023567-PW1	ZNF346 (Human) IP-WB Antibody Pair
H00010352-PW1	WARS2 (Human) IP-WB Antibody Pair	H00059348-PW1	ZNF350 (Human) IP-WB Antibody Pair
H00010810-PW1	WASF3 (Human) IP-WB Antibody Pair	H00252884-PW1	ZNF396 (Human) IP-WB Antibody Pair
H00114049-PW1	WBSCR22 (Human) IP-WB Antibody Pair	H00084307-PW1	ZNF397 (Human) IP-WB Antibody Pair
H00057728-PW1	WDR19 (Human) IP-WB Antibody Pair	H00084838-PW1	ZNF496 (Human) IP-WB Antibody Pair
H00091833-PW1	WDR20 (Human) IP-WB Antibody Pair	H00090441-PW1	ZNF622 (Human) IP-WB Antibody Pair
H00091833-PW2	WDR20 (Human) IP-WB Antibody Pair	H00007579-PW1	ZSCAN20 (Human) IP-WB Antibody Pair
H00010785-PW1	WDR4 (Human) IP-WB Antibody Pair	H00007589-PW1	ZSCAN21 (Human) IP-WB Antibody Pair
H00008839-PW1	WISP2 (Human) IP-WB Antibody Pair	H00007589-PW2	ZSCAN21 (Human) IP-WB Antibody Pair
H00008839-PW2	WISP2 (Human) IP-WB Antibody Pair	H00011130-PW1	ZWINT (Human) IP-WB Antibody Pair
H00089780-PW1	WNT3A (Human) IP-WB Antibody Pair	H00007791-PW1	ZYX (Human) IP-WB Antibody Pair

IP&ChIP Kits

Cat. Num.	Product Name	Cat. Num.	Product Name
KA0669	Acetyl-Histone H3 ChIP Kit	KA0569	Protein G Agarose Immunoprecipitation Kit
KA0670	Acetyl-Histone H4 ChIP Kit	KA0671	Tissue Acetyl-Histone H3 ChIP Kit
KA0685	Chromatin Immunoprecipitation (ChIP) Kit	KA0672	Tissue Acetyl-Histone H4 ChIP Kit
KA0692	Methylated DNA Immunoprecipitation Kit	KA0686	Tissue Chromatin Immunoprecipitation (ChIP) Kit
KA0694	Methyl-CpG Binding domain Protein 2 ChIP Kit	KA0693	Tissue Methylated DNA Immunoprecipitation Kit
KA0696	Methyl-Histone H3-K27 ChIP Kit	KA0695	Tissue Methyl-CpG Binding domain Protein 2 ChIP Kit
KA0697	Methyl-Histone H3-K4 ChIP Kit	KA0699	Tissue Methyl-Histone H3-K27 ChIP Kit
KA0698	Methyl-Histone H3-K9 ChIP Kit	KA0700	Tissue Methyl-Histone H3-K4 ChIP Kit
KA1514	One-Step ChIP Kit	KA0701	Tissue Methyl-Histone H3-K9 ChIP Kit
KA0568	Protein A Agarose Immunoprecipitation Kit	KA0702	Tri-Methyl-Histone H3-K9 ChIP Kit

Precipitor™ 16

High Throughput Immunoprecipitation and Purification

Abnova's Precipitor™ 16 system is an automated magnetic bead platform for high throughput precipitation and purification of proteins. Combining 96 deep well plate with affinity conjugated magnetic beads, Precipitor™ 16 easily handles 16 different assays simultaneously by transferring beads from one well to the next for mixing, binding, washing, and elution reactions via the robotic action of parallel magnetic rods. It simplifies the routine yet labor-intensive process and addresses the needs of rigorous proteomic screening and biomarker discovery applications such as immunoprecipitation (IP, ChIP, RIP), recombinant protein purification, and protein-protein interaction. Precipitor™ 16 delivers reproducible and consistent results by obviating the drawbacks of manual operation. Its integrated on-screen display allows easy change of parameters tailored to your experiment. Moreover, you can select from a large scope of available antibody reagents to accelerate your research!

Applications

- Immunoprecipitation
- Chromatin-Immunoprecipitation
- RNA-Immunoprecipitation
- Recombinant Protein Purification
- Protein-Protein Interaction


Advantages

- Rapid, multi-step automation
- Compatible with 96 well plate format
- Reproducibility and consistency of results
- Multiplex analysis of up to 16 analytes
- No external device or PC required



Affinity Magnetic Precipitation Principle

● Magnetic protein A bead ● Antibody ● Target ● Non-specific




- (A) Magnetic rods
(B) Cover strip rack
(C) 96 deep well plate rack
(D) 8-channel cover strip
(E) 96 deep well plate


Step 1
Mix antibody with protein A beads



Step 4
Wash with PBS/T




Step 2
Blocking to remove non-specific binding




Step 5
Elute




Step 3
React with sample lysate



Step 6
Release



Demonstration of Precipitor™



Accessory Reagents

- Naked Bead
- Indirect Antibody Conjugation
 - Protein A Bead
 - Protein G Bead
- Direct Antibody Conjugation
 - IP and ChIP Ab Bead
 - Goat Anti-Mouse Bead
- Protein Purification
 - Anti-His Bead
 - Anti-GST Bead
 - Anti-FC Bead
 - Anti-Myc Bead

See Abnova website under "Systems" for more accessory reagents.

Precipitor™ Plus

Large-Scale Protein and Nucleic Acid Purification

Precipitor™ Plus is a fully automated magnetic bead platform for large-scale protein and nucleic acid purification. Precipitor™ Plus can process up to 5 ml of sample volume and handle 1 to 24 independent samples from a wide range of sources including DNA, RNA, protein, cell/tissue lysates, whole blood, serum, and plasma. The robotic arm consisting of parallel magnetic rods transfers the sample from one well to the next via the action of magnetic beads for mixing, binding, washing and elution reactions. The system contains pre-programmed reaction protocols which the user can select directly from the on-screen display. It also allows for modification and creation of new protocols via an external PC. The whole reaction process takes only 30-120 minutes depending on the specific experiment. Precipitor™ Plus simplifies the labor intensive operation, shortens cycle time, and delivers faster results for your research and clinical studies.

Applications


- Immunoprecipitation
- Chromatin-Immunoprecipitation
- RNA-Immunoprecipitation
- Recombinant Protein Purification
- Nucleic Acid Purification

Advantages

- Process up to 5 ml of sample volume
- Handle 1 to 24 independent samples
- Deliver reproducible and consistent results
- Allow walk-away operation
- Ease of use and flexible platform




(A) Magnetic rods
(B) Cover strip rack
(C) 24 deep well plate rack
(D) 8-channel cover strip
(E) 24 deep well plate




Principle

- Magnetic protein A bead
- Antibody
- Target
- Non-specific


Step 1 Mix antibody with protein A beads




Step 2 Blocking to remove non-specific binding




Step 3 React with sample lysate




Step 4 Wash with PBS/T



Step 5 Elute



Step 6 Release



Accessory

Magnetic Beads

- Naked Bead
- Indirect Antibody Conjugation
 - Protein A Bead
 - Protein G Bead
- Direct Antibody Conjugation
 - IP and ChIP Ab Bead
 - Goat Anti-Mouse Bead
- Protein Purification
 - Anti-His Bead
 - Anti-GST Bead
 - Anti-FC Bead
 - Anti-Myc Bead



DNA/RNA Purification Kits

- Blood DNA Extraction Kit (Pre-filled)
- Blood RNA Extraction Kit (Pre-filled)
- Tissue DNA Extraction Kit (Pre-filled)
- Tissue RNA Extraction Kit (Pre-filled)
- Rice DNA Extraction Kit (Pre-filled)
- Plant DNA Extraction Kit (Pre-filled)
- Plant RNA Extraction Kit (Pre-filled)

See Abnova website under "Systems" for more accessory reagents.

